

Hematology Devices (for in vitro diagnostic use)

ABX VET Pack

04/02/09
A95A00249BEN

REF 0604052

REAGENT 1 0,5L

REAGENT 2 0,3L

REAGENT 3 3,4L

HORIBA ABX SAS
BP 7290
34184 Montpellier
cedex 4 - France

Exclusive use:

ABX Micros ABC Vet
SCIL ABC Vet
ABX Micros ESV60
SCIL Vet ABC Plus

1. Functions^a

R1: enzymatic solution with proteolytic action for the cleaning of blood cell counters.

R2: erythrocyte lysing agent for white blood cells counting and differentiation and hemoglobin determination.

R3: buffered isotonic solution for the determination and the differentiation of blood cells, and the measurement of hematocrit.

Measurement procedure to be followed in using the device:

principle of the method, specific analytical performance characteristics, analytical sensitivity, diagnostic sensitivity, analytical specificity, diagnostic specificity, accuracy, repeatability, reproducibility (including control of known relevant interference), limits of detection, limitations of the method and information about the use of available reference measurement procedures and materials by the user: see «Section: Specifications» in the instrument User Manual.

2. Conservation & Expiration

Storage conditions: stored at 5°C (41°F) to 25°C (77°F).

Expiration date: refer to «expiration date» reagent packaging label.

3. Measurements, Principles & Results

Directions for use: see «Section: Reagent Location and connection» in the instrument User Manual.

Measuring Principles : see «Section: Technology» in the instrument User Manual.

Results: refer to the instrument User Manual.

Performance data: see «Section: Specifications» in the instrument User Manual.

Note: if performance changes, call your HORIBA Medical representative.

a.Modification from A to B : ABX Micros ESV60, SCIL Vet ABC Plus

4. Composition & Handling precautions

R1:	Organic Buffer < 20%
	Proteolytic enzyme < 1%
	Detergent..... < 5%
	Preservative < 25%
	pH: 9,6 +/- 0,4 (T = 20°C)
	Resistivity: 72 +/- 2 Ω (T = 20°C)
	Description: Colorless liquid.
R2:	Preservative <25%
	Quaternary Ammonium salt < 20%
	Sodium Carbonate.. 0,2%
	Surfactant..... 0,2%
	pH: 10,5 +/- 0,5 (T = 20°C)
	Resistivity: 55 +/- 1 Ω (T = 20°C)
	Description: Colorless liquid.
R3:	Sodium fluoride..... < 3%
	Sodium Azide..... < 0,1%
	Sodium Hydroxyde < 1%
	Dimethylolurea..... 0,1%
	pH: 7 +/- 0,1 (T = 20°C)
	Resistivity: 60 +/- 6 Ω (T = 20°C)
	Description: Limpid and odorless aqueous solution.

Handling Precautions: avoid contact with eyes, skin and clothing. Wear laboratory gloves when handling the product. The product may be harmful if ingested or inhaled. Keep the pack closed when not in use. Please refer to the MSDS associated with these reagents.



Use established, good laboratory working practices when handling the used Pack.

Specimen Collection and Mixing: see «Section: Specimen collection and Mixing» in the instrument User Manual.

5. Limitations & Waste disposal

Limitations: see «Section: Specifications» in the instrument User Manual.

Safe Waste Disposal: follow your laboratory's protocol when neutralizing and disposing of waste. Please refer to the MSDS associated with these reagents.

Follow your local and/or national guidelines for biohazard waste disposal.