
**LAtrol™ Abnormal
Control Plasma**

REF 816A

IVD

**LAtrol™ Normal
Control Plasma**

REF 816N

EC **REP**

Obelis s.a
Bd. Général Wahis 53, 1030 Brussels, BELGIUM

REF 816A: 10 vials each containing 0.5 mL of LATrol™ Abnormal Control
10 fläschchen mit je 0.5 mL LATrol™ Abnormale Kontrolle
10 flacons contenant chacun 0,5 mL de LATrol™ Contrôle Anormal
10 viales con 0,5 mL de Control Anormal LATrol™ cada uno
10 flaconi contenenti ciascuno 0,5 mL di LATrol™ Controllo Anormale

REF 816N: 10 vials each containing 1.0 mL LATrol™ Normal Control
10 fläschchen mit je 1.0 mL LATrol™ Normal Kontrolle
10 flacons contenant chacun 1,0 mL de LATrol™ Contrôle Normal
10 viales con 1,0 mL de Control Normal LATrol™ cada uno
10 flaconi contenenti ciascuno 1,0 mL di LATrol™ Controllo Normale

ENGLISH

INTENDED USE

The LATrol™ Abnormal Control (REF 816A) and LATrol™ Normal Control (REF 816N) plasmas have been developed for use as part of daily quality control procedures for Lupus Anticoagulant (LA) testing. These control plasmas are designed to be used with ACTICLOT® dPT™ (REF 824), DVVtest® (REF 810/825), and DVVconfirm® (REF 815/815L).

DESCRIPTION

LATrol Abnormal Control is a lyophilised preparation of a Lupus Anticoagulant plasma with added buffer, which has been determined to be positive for LA in accordance with the revised criteria of the Scientific and Standardization Committee's (SSC) Subcommittee on Lupus Anticoagulant detection.^{1,2}

LATrol Normal Control is a lyophilised preparation of a multi-donor normal plasma pool with added buffer. This normal plasma may also be used in mixing tests to determine the presence of inhibitors or factor deficiencies in patient plasmas that test positive using ACTICLOT dPT and DVVtest.

Source plasmas for these products are processed in a manner consistent with established procedures to ensure the plasma is platelet-poor.²

WARNING

The source material for these plasmas has been found to be non-reactive for Hepatitis B Surface Antigen (HBsAg), Hepatitis C Virus (HCV) and Human Immunodeficiency Virus Type 1 and Type 2 (HIV-1, HIV-2) using US Food and Drug Administration cleared methods. As no known test method can provide complete assurance that products derived from human blood will not transmit HBsAg, HCV, HIV-1, HIV-2 or other blood-borne pathogens, these plasmas should be handled as recommended for any potentially infectious human specimen.

REAGENT PREPARATION AND STABILITY

These lyophilised plasmas are stable until the date indicated on the vial label when stored at 2°C to 8°C. Reconstitute these plasmas with purified, deionized or distilled water as follows:

REF	816A	816N
Volume H ₂ O	0.5 mL	1.0 mL

After reconstitution, re-stopper the vial, swirl the contents gently and allow the vials to stand at room temperature for at least 15 minutes to ensure complete dissolution. Once reconstituted, these plasmas will remain stable for:

	8 hours
	2° - 8°C

Mix gently prior to use. Refer to the ACTICLOT dPT, DVVtest and DVVconfirm package inserts for complete assay procedures.

EXPECTED RESULTS/PERFORMANCE CHARACTERISTICS

BioMedica Diagnostics offers Instrument Applications for several coagulation analysers. These Instrument Applications may contain platform specific programming and performance data which differ from that provided in this Instructions for Use. In these cases, the information contained in the Instrument Application supersedes the information in this Instruction for Use. Please consult the specific manufacturer's instrument manual for complete operating instructions.

LAtrol Abnormal Control and LAtrol Normal Control plasmas have been used in a quality control protocol for ACTICLOT dPT, DVVtest and DVVconfirm reagents. Assay values were generated using the ACL[®] 300R, BCT[®], CA-7000, MLA[®] 900C, ST4, STA Compact[®] and Yumizen G800 coagulation analyzers.

The following clotting times are for example only. Actual clotting times may vary among laboratories and are dependent upon the method and instrumentation employed for clotting end-point detection. The protocol included multiple testing over several days. The overall coefficients of variation for the normal control and abnormal control are listed. The limits of quality control acceptability for these control plasmas should be established by each laboratory as part of a complete quality control protocol.

Table 1. LATrol™ Precision Data with ACTICLOT® dPT™

Coagulation Analyzer	Control	dPT Screening (sec)	Intra-Assay CV (%)	Inter-Assay CV (%)	dPT Confirmatory (sec)	Intra-Assay CV (%)	Inter-Assay CV (%)
ACL® 300R	816N	32.8	2.5	5.1	30.2	3.8	5.3
	816A	63.8	1.9	7.1	36.9	3.2	3.8
BCT®	816N	47.5	0.5	3.2	51.6	1.7	4.5
	816A	89.2	0.6	5.2	61.9	1.2	3.7
CA 7000	816N	44.3	0.9	ND	37.9	2.0	Nd
	816A	84.1	5.4	ND	48.6	3.3	Nd
MLA® 900C	816N	27.9	2.5	3.7	27.2	2.8	4.1
	816A	51.6	2.4	8.6	30.5	1.5	3.7
ST4	816N	38.1	0.7	ND	37.3	3.8	ND
	816A	60.6	6.5	ND	44.0	6.7	ND
STA Compact®	816N	40.2	0.8	3.4	39.7	0.9	4.3
	816A	77.9	1.1	7.2	46.0	1.0	4.8

ND – Not Determined

Table 2. LATrol™ Precision Data with DVVtest® and DVVconfirm®

Coagulation Analyzer	Control	DVVtest (sec)	Intra-Assay CV (%)	Inter-Assay CV (%)	DVVconfirm (sec)	Intra-Assay CV (%)	Inter-Assay CV (%)
ACL® 300R	816N	30.4	1.4	ND	31.3	0.6	ND
	816A	64.2	3.2	ND	35.1	1.4	ND
BCT®	816N	31.8	0.3	2.4	33.2	0.5	3.7
	816A	63.8	0.5	2.7	43.0	0.6	4.8
MLA® 900C	816N	31.0	0.9	2.9	31.7	1.3	2.3
	816A	63.6	1.4	3.9	34.5	0.9	3.2
ST4	816N	36.6	0.9	ND	34.4	1.0	ND
	816A	80.7	1.2	ND	41.1	1.2	ND
STA Compact®	816N	36.2	0.7	3.3	34.4	0.7	3.6
	816A	71.9	0.9	3.8	38.9	1.3	5.6
Yumizen G800	816N	40.3	2.0	2.2	31.9	1.8	2.2
	816A	103.4	3.3	4.3	44.1	3.6	4.5

ND – Not Determined

REFERENCES

1. CLSI. *Laboratory Testing for the Lupus Anticoagulant; Approved Guideline*. CLSI document H60-A. Wayne, PA: Clinical and Laboratory Standards Institute, 2014.

2. Pengo V, Tripodi A, Reber G, Rand JH, Ortel TL, Galli M, de Groot PG. Update of the guidelines for lupus anticoagulant detection. *J Thromb Haemost* 2009, 7: 1737–1740.

ACL and MLA are registered trademarks of Instrumentation Laboratory, SpA.

BCT is a registered trademark of Dade Behring Inc.

STA Compact is a registered trademark of Diagnostica Stago SA

Yumizen is a mark of HORIBA ABX SAS, Montpellier, France

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Yumizen G800 es una marca registrada de HORIBA ABX SAS, Montpellier, France