Hematology Devices (for in vitro diagnostic use) ABX CRP Std

A01A00071B / I1002840000A

www.horiba-abx.com

Explore the future

HORIBA GROUP

Exclusive use:

Micros CRP 200
Micros CRP
Micros 45
Pentra 60
Pentra 60 C+
Pentra 80
Pentra XL 80
Pentra 120
Pentra 120 Retic
Pentra DX 120

REF

0501016

CAL

2x1mL

IVD

CE

HORIBA ABX
BP 7290 - 34184 Montpellier
cedex 4 - France

Date: 03/01/06

FORM 0867 - rev 1

S.A.S au capital de 41.700.000 Euros - RCS Montpellier 328 031 042 - SIRET 328 031 042 000 42 - APE 332 B

EN

1. Functions

ABX CRP Std is used for calibration of CRP concentration measurement on ABX MICROS CRP 200 instruments, by immuno-turbidimetry technique.^a

2. Summary

ABX CRP Std is a human serum - based Calibrator, used on instruments requiring calibration on a periodic basis. Shall be used in conjunction with ABX CRP REA reagent (see ABX CRP REA leaflet)^b

The CRP values are defined by standard IFCC method, by using international reference CRM 470.

3. Principles of measurement

ABX CRP Std allows calibration of instrument, previously to measurements taken on whole blood sample (see «Instructions for use» chap 5). The CRP in whole blood sample is measured by immuno-turbidimetry using ABX CRP REA.

Refer to ABX CRP REA reagent leaflet for principles of measurement.

4. Warning and precautions

Potentially biohazardous material.
For in vitro diagnostic use.

Each human blood donor unit used in the compounding of the serum base for the above referenced products has been tested by an

FDA approved method and found to be non reactive for HIV 1/2 (HIV) Antibody, HCV Antibody, Hepatitis B Surface Antigen (HBsAg). The above referenced products were only formed with these donors.

Because no known test method can provide total assurance that products derived from human blood will not transmit infectious diseases, products containing materials from human sources should be handled as if potentially infectious.

Use safe laboratory procedures as outlined in Biosafety in Microbiological and Biomedical Laboratories. Contains < 0,1% Sodium Azide as preservative. Avoid swallowing and contact with skin or mucous membranes.

5. Instruction for use

See section «Quality Assurance» in the instrument User Manual: «CRP Calibration»

After use seal tightly the bottle. Pay attention not to mistake the cap. The color of the ABX CRP Std cap is blue.^c

6. Storage conditions

Storage temperature: 2 - 10°C (do not freeze)^d

Stability before opening: refer to «Expiry date» packaging label.

Stability after opening: 3 months (or until the «Expiry date» whatever comes first)

7. Assigned values

Refer to ABX CRP Std «Value» packaging label.^e

8. References

Biosafety in Microbiological and Biomedical Laboratories. DC/NIH manual, 1993.

a.Modification from index A to B : immuno-turbidimetry technique
b.Reagent to used in conjunction with abx crp Std
c.sealing information
d.Do not freeze
e.Assigned values addition

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