

ACTICHROME[®] Heparin (Anti-FXa)

REF 832

Obelis s.a
Boulevard Général Wahis 53, 1030 Brussels, BELGIUM

INTENDED USE

ACTICHROME[®] Heparin (Anti-FXa) is a chromogenic assay intended for the quantitative measurement of therapeutic heparin, by measurement of factor Xa activity, in plasma from patients undergoing anticoagulant therapy. The assay is for use by laboratory professionals and may be performed manually, or by using semi-automated or automated coagulation analysers. The assay is for in vitro diagnostic use.

EXPLANATION OF THE TEST

The inhibitory effect of antithrombin III (AT-III) on thrombin, factor Xa and other coagulation serine proteases in plasma is increased several thousand-fold by heparin. This inhibition accounts for the anticoagulant effect of heparin. Low molecular weight heparin (LMWH) preparations appear to catalyze the reaction between factor Xa and AT-III more readily than the reaction between thrombin (FIIa) and AT-III. Unfractionated heparin (UFH) catalyzes both reactions equally. This makes the factor Xa inhibition test the most useful assay as it is applicable to the widest variety of heparin preparations.

PRINCIPLE OF THE METHOD

ACTICHROME Heparin (Anti-FXa) is a multi-step, chromogenic assay for measuring heparin activity by measuring residual factor Xa activity. In the first step, AT-III is added to the patient plasma and incubated at 37°C. In the second step, bovine factor Xa is added to the mixture and again incubated at 37°C. Heparin present in the patient sample will inhibit factor Xa activity. Last, a chromogenic substrate, highly specific for factor Xa, is added to the mixture followed by a final incubation at 37°C. Factor Xa activity present hydrolyses the substrate, releasing a p-nitroaniline (pNA) chromophore which turns the mixture yellow. The assay reaction is stopped by adding acetic acid, which turns the mixture blue. The absorbance of the pNA in the reaction solution at 405 nm is measured and compared to those values obtained from a standard curve generated using known heparin activity levels. The rate of factor Xa inhibition is directly proportional to the heparin concentration since both factor Xa and AT-III are in excess. The residual factor Xa activity is inversely proportional to the heparin concentration.^{1,2}

REAGENTS

The kit contains reagents sufficient to perform 200 tests using an automated coagulation analyzer, 100 tests if a manual end-point method is used.

R1 Bovine Factor Xa Reagent: 4 vials (lyophilized).

R2 Human Antithrombin III Reagent: 4 vials (lyophilized)

R3 SPECTROZYME[®] FXa: 4 vials each containing 4 µmoles of substrate (lyophilized).

WARNING

This product contains human source material that has been found to be non-reactive for Hepatitis B Surface Antigen (HBsAg), Hepatitis C Virus (HCV) and Human Immunodeficiency Virus Type 1 and Type 2 (HIV-1, HIV-2) using registered methods. As no known test method can provide complete assurance that products derived from human specimens will not transmit HBsAg, HCV, HIV-1, HIV-2 or other blood-borne pathogens, this reagent should be handled as recommended for any potentially infectious human specimen.

This product contains animal source material. As no known test method can provide complete assurance that products derived from animal specimens will not transmit blood-borne pathogens, this reagent should be handled as recommended for any potentially infectious human specimen.

Dispose of discarded material and packing in accordance with all applicable local regulations.

If any serious incident occurs when using ACTICHROME Heparin (Anti-FXa), contact BioMedica Diagnostics or your local distributor, and report to your local competent authority.

Bovine Factor Xa	Warning		CONT Tris-(hydroxymethyl)aminomethane H315, H319, H335, P261, P264, P280, P302 + P352, P305 + P351 + P338, P337 + P313
Human Antithrombin III	Warning		H315, H319, H335, P261, P264, P280, P302 + P352, P305 + P351 + P338, P337 + P313
SPECTROZYME® FXa	Warning		CONT N-Methoxycarbonyl-D-cyclohexyl-glycyl-glycyl-arginine-para-nitroanilide acetate H315, H319, H335, P261, P264, P280, P302 + P352, P305 + P351 + P338, P337 + P313

Hazard Statements	H315 Causes skin irritation. H319 Causes serious eye irritation. H335 May cause respiratory irritation.
Precautionary Statements:	P261 Avoid breathing dust/fume/gas/mist/vapors/spray. P264 Wash thoroughly after handling. P280 Wear protective gloves/ protective clothing/ eye protection/ face protection. P302 + P352 IF ON SKIN: Wash with plenty of water. P305 + P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. P337 + P313 If eye irritation persists: Get medical advice/attention.

REAGENT PREPARATION AND STORAGE

Lyophilized reagents are stable until the expiration date indicated on the label when stored as instructed.

- R1 Bovine Factor Xa Reagent:** Reconstitute with 5 mL deionized water. Reconstituted reagent is stable for 4 weeks at –20°C.
- R2 Human Antithrombin III Reagent:** Reconstitute with 5 mL deionized water. Reconstituted reagent is stable for 4 weeks at –20°C.
- R3 SPECTROZYME® FXa:** Reconstitute with 5 mL deionized water. Reconstituted reagent is stable for 4 weeks at –20°C.

SPECIMEN COLLECTION AND PREPARATION

Only citrate collected platelet poor plasma may be used for this assay. Do Not Use EDTA collected plasma. See "Collection, Transport and Processing of Blood Specimens for Testing Plasma-based Coagulation Assays and Molecular Hemostasis Assays; Approved Guidelines-Fifth Edition", CLSI Document H21-A5, Vol. 28, No. 5, January 2008. Plasma collection should be performed as follows:

- Collect 9 parts of blood into 1 part of 3.2% (0.109 M) trisodium citrate anticoagulant solution.
- Centrifuge the blood sample at 3,000 x g for 15 minutes.
- Plasma should be stored at 2° - 8°C and assayed within 2 hours. Alternatively, plasma may be stored at –70°C for up to 6 months
- Frozen plasma should be thawed rapidly at 37°C. Thawed plasmas should be stored at 2°-8°C and assayed within 2 hours.

PROCEDURE

Materials Provided – See Reagents

Materials Required But Not Provided

Automated coagulation analyser set at 405 nm
Spectrophotometer set at 405 nm and laboratory time deionized water
glacial acetic acid, 0.9% NaCl
25 µL, 200 µL, 500 µL and 2500 µL pipettes
37°C water or dry bath
plastic test tubes or cuvettes
commercially available normal plasma
Heparin – See 'Assay Calibration'

Assay Calibration

Pooled normal human plasma (at least 10 normal donors) which has been collected in the same way as plasmas to be tested, may be used in the preparation of the heparin standards. Pooled normal human plasma (from at least 20 normal donors) that has been collected in the same manner as plasmas to be tested or commercially available normal plasma may be used for preparation of the heparin standards. The specific heparin preparation used during therapy must be added to the pooled normal plasma to prepare the standards.

Prepare an 8.0 USP unit/mL solution of heparin in saline (0.9% NaCl). When starting with a 1000 unit/mL heparin stock solution, pipette 40 µL of the heparin stock into 4.96 mL of the saline. Then prepare the heparin standards using the 8.0 USP unit/mL solution as shown in Table 1.

Table 1 – Preparation of Heparin Standards

Heparin Standard	Pooled Normal Plasma	Volume of Heparin Standard
0.8 unit/mL	900 µL	100 µL of 8.0 USP unit/mL
0.4 unit/mL	500 µL	500 µL of 0.8 USP unit/mL
0.2 unit/mL	500 µL	500 µL of 0.4 USP unit/mL
0.0 unit/mL	500 µL	None

Assay Procedure

ACTICHROME Heparin (Anti-FXa) may be performed manually, or by using semi-automated or automated coagulation analysers.

Automated Instrumentation Method

BioMedica Diagnostics offers instrument applications for performing ACTICHROME Heparin (Anti-FXa) on several coagulation analysers. These instrument applications may contain platform specific programming and performance data which differ from that provided in this Instructions for Use. In these cases, the information contained in the instrument application supersedes the information in this Instruction for Use. Please consult the specific manufacturer's instrument manual for complete operating instructions. The performance of this assay was determined on an MLA ELECTRA 900C™ optical analyzer.

Manual End-point Method

1. Add 200 µL of AT-III to a plastic test tube.
2. Add 25 µL of plasma sample or heparin standard. Mix and incubate at 37°C for 2 minutes.
3. Add 200 µL of Bovine Factor Xa. Mix and incubate at 37°C for exactly 1 minute.
4. Add 200 µL of SPECTROZYME FXa. Mix and incubate at 37°C for exactly 5 minutes.
5. Add 200 µL of glacial acetic acid. Mix.
6. Add 200 µL of deionized water (optional).*
7. Read absorbance at 405 nm in a 1 cm cuvette against a blank prepared in the following order:
 - 200 µL acetic acid
 - 200 µL AT-III
 - 25 µL pooled normal plasma
 - 200 µL factor Xa
 - 200 µL SPECTROZYME FXa
 - * 200 µL water (optional)

* Some spectrophotometers require a minimum of 1 mL volume in cuvette.


QUALITY CONTROL

Quality control samples should be run according to Laboratory Guidelines. Control samples should be run on freshly reconstituted reagent vials and each time the assay is run. Commercial heparin control plasma may be used for quality control of the assay. Alternatively, pooled normal plasma containing 0.5 and 0.25 units/mL heparin can be made and used as in-house controls. If quality control samples fail to yield heparin concentrations within the correct ranges in any given test, the test should be repeated with freshly reconstituted reagents.

RESULTS

Construct a standard curve by plotting the mean absorbance at a wavelength of 405 nm value for each heparin standard versus its corresponding concentration and draw the best-fit calibration curve. Interpolate the heparin in the unknown samples directly from the standard curve. A standard curve should be generated each time the assay is performed. The following curve is for demonstration purposes only.

Representative Standard Curve

LIMITATIONS OF THE PROCEDURE

Different heparin preparations may elicit different slopes in the standard curve. The same heparin administered to the patient must be used to generate the standard curve.

Platelet contamination in plasma samples may result in the release of platelet factor 4 (PF4), a potent heparin neutralizer. PF4 released into plasma may result in an underestimation of the heparin concentration. Freezing and thawing of platelets is known to release PF4 from the platelets. If plasma is to be frozen and thawed prior to use in the assay, then the plasma must be prepared as described in the sample preparation section so that it is platelet poor.

INTERFERENCES

Icteric, lipemic and hemolyzed samples may interfere with assay results. Studies have shown that these samples result in an underestimation of the heparin concentration³.

PERFORMANCE CHARACTERISTICS

Accuracy

A study performed using an automated instrument to evaluate ACTICHROME Heparin (Anti-FXa) against another commercially available chromogenic heparin assay (using unfractionated heparin) gave a correlation coefficient of 0.913 for 92 samples tested, with a regression equation of $y = 0.84x + 0.003$. The Standard Error of Estimation = 0.07.

A study performed using an automated instrument to evaluate ACTICHROME Heparin (Anti-FXa) against another commercially available chromogenic heparin assay (using low molecular weight heparin) gave a correlation coefficient of 0.895 for 26 samples tested, with a regression equation of $y = 1.36x - 0.139$. The Standard Error of Estimation = 0.14.

Precision

The following estimates of coefficient of variation (CV) were observed using an automated spectrophotometer.

Table 2 – Intra-Assay and Inter-Assay Variations

UFH Concentration	Intra-Assay Variation	Inter-Assay Variation
0.44 units/mL	3.3 %	4.1 %
0.25 units/mL	5.6 %	4.8 %

LMW Heparin Concentration	Intra-Assay Variation	Inter-Assay Variation
0.48 units/mL	2.9 %	6.2 %
0.24 units/mL	5.1 %	9.3 %

Sensitivity

ACTICHROME Heparin (Anti-FXa) is designed to give a linear standard curve for heparin levels between 0 and 0.8 units/mL.

Specificity

Specificity of the assay has been ensured by the use of purified bovine factor Xa and purified human AT-III. In addition, the chromogenic substrate used to measure residual factor Xa activity is highly specific for factor Xa.

SUMMARY OF SAFETY AND PERFORMANCE

The Summary of Safety and Performance (SSP) is available in the European Database on Medical Devices (Eudamed), at <https://ec.europa.eu/tools/eudamed> under the Basic UDI-DI (GMN) 0062794400265Q, or upon request.

REFERENCES

1. Teien, A. N., Lie, M. and Abildgaard, U. Assay of heparin in plasma using a chromogenic substrate. *Thrombosis Research* **8**: 413-416, 1976.
2. Teien, A. N. and Lie, M. Evaluation of an amidolytic heparin assay method: Increased sensitivity by adding purified antithrombin III. *Thrombosis Research* **10**: 399-410, 1977.
3. Data on file at BioMedica Diagnostics.

MLA is a registered trademark of Instrumentation Laboratory, SpA

Changes from the previous revision are noted by red dotted lines in the left margin.

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MLA ist ein ist ein geschützter Markenname von Instrumentation Laboratory, SpA