ABX Pentra Amylase CP

REF	A11A01628
REAGENT 1	26 mL
REAGENT 2	6,5 mL

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Reference

- Lorentz K. α -Amylase. In: Thomas L, editor. Clinical laboratory diagnostics. 1st ed. Frankfurt: TH-Book Verlagsgesellschaft (1998): 192-202.
- Moss DW, Henderson AR. Digestive enzymes of pancreatic origin. In: Burtis CA, Ashwood ER, editors. Tietz Textbook of Clinical Chemistry. 3rd ed. Philadelphia: WB Saunders Company (1999): 689-698.
- Kruse-Jarres JD, Kaiser C, Hafkenscheid JC, Hohenwallner W, Stein W, Bohner J et al. Evaluation of a new alpha-amylase assay using 4,6-ethylidene-(G7)-1-4-nitrophenyl-(G1)-alpha-D maltoheptaoside as substrate. J. Clin. Chem. Biochem. (1989) **27**: 103-113.
- Lorentz K. Approved recommendation on IFCC methods for the measurement of catalytic concentration of enzymes; Part 9. IFCC method for α -Amylase (1,4- α -D-Glucan-4-glucanohydrolase, EC 3.2.1.1). Clin. Chem. Lab. Med. (1998) **38**: 185-203.
- Guder WG, Zawta B, the Quality of Diagnostics Samples. Samples: From the Patient to the Laboratory. 1st ed. Guder WG, Narayanan S, Zawta B. (WHILEY-VCH, Darmstadt, Germany) (2001): 16-17.
- Guder WG, Zawta B, the Quality of Diagnostics Samples. Samples: From the Patient to the Laboratory. 1st ed. Guder WG, Narayanan S, Zawta B. (WHILEY-VCH, Darmstadt, Germany) (2001): 50-51.
- Roberts WL, McMillin GA, Burtis CA, Bruns DE. Reference Information for the Clinical Laboratory, Tietz Textbook of Clinical Chemistry and Molecular Diagnostics. 4th Ed., Burtis C.A., Ashwood E.R., Bruns D.E. (Elsevier Saunders eds. St louis, USA) (2006): 2255.
- Junge W, Waldenström J, Bouman A et al. Evaluation of the Assays for Total and Pancreatic alpha-Amylase based on 100% cleavage of Et-G7-PNP at 6 European Clinical Centres (Poster Medlab 97); Bâle, Suisse: 12th IFCC European Congress of Clinical Chemistry (1997): 17-22.
- Council Directive (2000/54/EC). Official Journal of the European Communities. No. L262 from October 17, 2000: 21-45.
- Evaluation of detection capability for clinical laboratory measurement procedures. Approved Guideline, 2nd ed., CLSI (NCCLS) document EP17-A2 (2012) **32** (8).
- Vassault A, Grafmeyer D, Naudin C et al. Protocole de validation de techniques (document B). Ann. Biol. Clin. (1986) **44**: 686-745.
- Evaluation of Precision Performance of Quantitative Measurement Method. Approved Guideline, CLSI (NCCLS) document EP5-A2 (2004) **24** (25).
- Evaluation of Linearity of Quantitative Measurement Procedures. 2nd Edition, CLSI (NCCLS) guideline EP06-Ed2 (2020) **40** (16).
- Measurement Procedure Comparison and Bias Estimation Using Patient Samples. Approved Guideline, 3rd ed., CLSI (NCCLS) document EP09c (2018) **38** (12).
- Passing H, Bablok W. A new biometrical procedure for testing the equality of measurements from two different analytical methods. J. Clin. Chem. Clin. Biochem. (1983) **21**: 709-720.
- Young DS. Effects of Drugs on Clinical Laboratory Tests. 4th Edition, Washington, DC, AACC Press (1997) **3**: 143-163.
- Young DS. Effects of Preanalytical Variables on Clinical Laboratory Tests. 2nd Edition, Washington, DC, AACC Press (1997) **3**: 120-132.