Utilisation

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Paramètres de test

Test:	Chol	Chemistry:	Cholesterol
Chemistry No.:	210	Print Name:	Cholesterol
Reaction Type:	Endpoint	Reaction Direction:	Positive
Pri. Wave:	510 nm	Sec. Wave:	670 nm
Decimal.:	0	Samp. Type:	Serum
Blank Time:		Reaction Time:	18 20
Unit:	mg/dL	Incubation Time:	0

	Sample Vol.	Aspirated	Diluent	Reagent Vol.	Diluent
Standard;	2	uL	uL	R1: 200	uL
Decreased;		uL	uL		uL
Increased;		uL	uL		uL

Linearity Range (Standard);	0-500	Linearity Limit:	
Linearity Range (Decreased);		Substrate Depletion:	
Linearity Range (Increased);		Mixed Blank Abs.:	-40000 40000
R1 Blank Abs.:	-40000 40000	On-board Stability:	30 Day (s)
Blank Response	-40000 40000	Reagent Alarm Limit:	5
Twin Chemistry:			

Prozone Check:		
Q1:	Q2:	Q3:
Q4:	PC:	ABS:

Use Qualitative Result:	Range:	Flag:
-------------------------	--------	-------

Slope Offset:	Slope	Offset	Unit
	1	0	mg/dL

Pretreatment:	
Pretreat Sample Vol.:	uL
Pretreat Reagent Vol.:	uL

Ref. Range:	Sample Type:	Gender:	Age Range:	Ref. Range:	Critical Range:	Unit:
-------------	--------------	---------	------------	-------------	-----------------	-------

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Chem:	Chol	Calibrator	Conc.	Pos	Lot No.
Calibration Setting		Water	0.0	W	
Math Model:	Two-Point Linear	Chem Cal	*	*	
Factor:	Replicates: 2				
Acceptance Limits					
Cal Time:	336 hr.				
Slope Diff:	SD:				
Sensitivity:	Repeatability:				* User Defined
Deter Coeff:					
Auto Calib:					
	<input type="checkbox"/> Cal Time				

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Calculs (Exemple)

Abs. = Absorbance

$$\frac{\text{Abs. (Patient)}}{\text{Abs. (Standard)}} \times \text{Concentration of Std.} = \text{Cholesterol (mg/dl)}$$

Exemple: Abs. (Patient) = 0.40, Abs. (Standard) = 0.32, Concentration of Standard = 200 mg/dl

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Symboles

Use by (YYYY-MM-DD)	Lot and batch code
Catalog number	Manufacturer
In vitro diagnostic medical device	Temperature limitation
Consult instructions for use	Rx Only: Prescription Use Only
CE mark	Authorized representative in the European Community

12-C7510-160 Manufactured by HORIBA Instruments Incorporated - Pointe Brand 5449 Research Drive Canton, MI 48188

Manufactured by HORIBA Instruments Incorporated – Pointe Brand
5449 Research Drive, Canton, MI 48188

European Authorized Representative:

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