

## Intended Use

For the quantitative determination of total protein concentration in serum using the Yumizen C230 and Yumizen C240 analyzers. **Rx Only.**

## Method History

The color reaction of protein molecules with cupric ions, known as the Biuret color reaction, has been known since 1878. Since the Riegler<sup>1</sup> publications of 1914, several attempts have been made to stabilize the cupric ions in the alkaline reagent. Kingsley,<sup>2,3</sup> modified the procedure in 1939 and 1942 to include the use of sodium potassium tartrate as a complexing agent. This procedure was later modified by Weichselbaum<sup>4</sup> and Gornall.<sup>5</sup> The present method is based on these modifications.

## Principle



Protein in serum forms a violet colored complex when reacted with cupric ions in an alkaline solution. The intensity of the violet color is proportional to the amount of protein present when compared to a solution with known protein concentration.

## Reagent Content

Sodium Hydroxide 600mM, Copper Sulfate 12mM, Sodium Potassium Tartrate 32mM, Potassium Iodide 30mM, Non-reactive ingredients.

## Precautions

1. This reagent is for *in vitro* diagnostic use only.
2. Avoid ingestion. DO NOT PIPETTE BY MOUTH. In case of ingestion drink large amounts of water and seek medical attention quickly.
3. Avoid contact with skin and eyes. The reagent contains sodium hydroxide which is corrosive. In case of contact with skin, flush with water. For eyes, seek medical attention.

## Reagent Preparation

Reagent comes in a ready to use form.

## Reagent Storage

Store reagent at room temperature(15-30°C). The reagent is stable until the expiration date appearing on the label when stored as directed.

## Reagent Deterioration

The reagent should be a clear, pale blue solution. Turbidity or the presence of a black precipitate indicates reagent deterioration and should not be used.

## Specimen Collection and Storage

1. Unhemolyzed serum is the specimen of choice.
2. Gross hemolysis will cause elevated results because of the released hemoglobin as well as the increase in background color.
3. Lipemic sera cause elevated results. A serum Blank should be performed.
4. Samples with bromosulfophthalein (BSP) will result in falsely elevated results.<sup>8</sup>
5. Protein in serum is stable for one week at room temperature (18-25°C) and for at least one month refrigerated (2-8°C) when guarded against evaporation.<sup>6</sup>

## Interferences

Young, et al.<sup>7</sup> has reviewed a number of drugs and substances that may affect protein concentrations.

## Materials Provided

Total Protein reagent

## Materials Required but not Provided

1. Yumizen C230 / Yumizen C240 Analyzer
2. Yumizen C230 / Yumizen C240 Operation Manual
3. Pointe Chemistry Calibrator, catalog number C7506-50
4. Pointe Chemistry control, catalog number C7592-100

## Test Parameters

Test:	TP	Chemistry:	Total Protein
Chemistry No.:	229	Print Name:	Total Protein
Reaction Type:	Endpoint	Reaction Direction:	Positive
Pri. Wave:	546 nm	Sec. Wave:	670 nm
Decimal.:	0.1	Samp. Type:	Serum
Blank Time:		Reaction Time:	10 11
Unit:	g/dL	Incubation Time:	0

	Sample Vol.	Aspirated	Diluent	Reagent Vol.	Diluent
Standard;	4	uL	uL	180	uL uL
Decreased;		uL	uL	uL	
Increased;		uL	uL	uL	

Linearity Range (Standard);	1.0-15.0	Linearity Limit:	
Linearity Range (Decreased);		Substrate Depletion:	
Linearity Range (Increased);		Mixed Blank Abs.:	- 40000 40000
R1 Blank Abs.:	- 40000 40000	On-board Stability:	30 Day (s)
Blank Response	- 40000 40000	Reagent Alarm Limit:	5
Twin Chemistry:			

Prozone Check:		
Q1:	Q2:	Q3:
Q4:	PC:	ABS:

Use Qualitative Result:	
Range:	Flag:

Slope Offset:			
Slope	Offset	Unit	
1	0	g/dL	

Pretreatment:			
Pretreat Sample Vol.:	uL	Pretreat Reagent Vol.:	uL

Ref. Range:			
Sample Type:	Gender:	Age Range:	Ref. Range: Critical Range: Unit:

# Pointe Total Protein (Biuret) Reagent Set

## Calibration Setup Parameters

Chem: TP																	
Calibration Setting	<table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <th>Calibrator</th> <th>Conc.</th> <th>Pos</th> <th>Lot No.</th> </tr> <tr> <td>Water</td> <td>0.0</td> <td>W</td> <td></td> </tr> <tr> <td>Chem Cal</td> <td>*</td> <td>*</td> <td></td> </tr> <tr> <td> </td> <td> </td> <td> </td> <td> </td> </tr> </table>	Calibrator	Conc.	Pos	Lot No.	Water	0.0	W		Chem Cal	*	*					
Calibrator	Conc.	Pos	Lot No.														
Water	0.0	W															
Chem Cal	*	*															
Math Model: Two-Point Linear																	
Factor: Replicates: 2																	
Acceptance Limits																	
Cal Time: 336 hr.																	
Slope Diff: SD:																	
Sensitivity: Repeatability: * User Defined																	
Deter Coeff:																	
Auto Calib.																	
<input type="checkbox"/> Cal Time																	

## Calibration

Use an NIST-traceable serum calibrator. The procedure should be calibrated according to the instrument manufacturer's calibration instructions. If control results are found to be out of range, the procedure should be re-calibrated.

## Quality Control

- Use control sera with known total protein concentrations to monitor the integrity of the reaction.
- Quality control requirements should be performed in conformance with local, state, and/or Federal regulations or accreditation requirements.

## Calculation (Example)

Abs. = Absorbance

$$\frac{\text{Abs. of Unknown}}{\text{Abs. of Standard}} \times \text{Conc. of standard} = \text{Total Protein (g/dl)}$$

Example: Abs. of Unknown = 0.350, Abs. of Standard = 0.400  
Concentration of Standard = 8 g/dl

$$\text{Then: } \frac{0.350}{0.400} \times 8 = 7.00 \text{ g/dl}$$

## Limitations

- Samples with values above 15.0 g/dl should be diluted 1:1 with 0.9% saline, re-run and result multiplied by two.
- The Biuret procedure is not sensitive at low ranges (<1 g/dl). Do not use for urine or spinal fluid.

## Expected Values<sup>8</sup>

6.2 – 8.5 g/dl

- The effect of posture, when blood is drawn, varies with the individual but recumbent values are usually lower than ambulatory. Differences may be as much as 1.2 g/dl.
- It is strongly recommended that each laboratory establish its own range.

## Performance










- Linearity: 1.0 – 15.0 g/dl
- Comparison: A study was performed between the Yumizen 200 series analyzers and a similar analyzer and method, resulting in a correlation coefficient of 0.996 with a regression equation of  $y=0.919x + 0.36$ .
- Precision: Precision studies were performed using the Yumizen 200 series analyzers following a modification of the guidelines which are contained in NCCLS document EP5-T2.<sup>9</sup>





Within Run			Day to Day		
Mean	S.D.	C.V.%	Mean	S.D.	C.V.%
4.77	0.13	2.7	4.78	0.07	1.5
7.46	0.16	2.1	7.63	0.14	1.8

## References

- Riegler, E., Anal. Chem. 53:242 (1914).
- Kingsley, G.R., J. Biol. Chem. 131:197 (1939).
- Kingsley, G.R., J. Lab. Clin. Med. 27:840 (1942).
- Weichselbaum, T., Amer. J. Clin. Path. 16:40 (1946).
- Gornall, A., et al, J. Biol. Chem. 177:752 (1949).
- Henry, R.J., et al, Clinical Chemistry: Principles and Technics, Harper & Row, New York, p.415 (1974).
- Young, D.S., et al, Clin. Chem. 21:1D (1975).
- Tietz, N.W., Fundamentals of Clinical Chemistry Philadelphia, W.B. Saunders, pp. 299, (1976).
- NCCLS document "Evaluation of Precision Performance of Clinical Chemistry Devices", 2<sup>nd</sup> Ed. (1992).

## Symbol Key

 Use by (YYYY-MM-DD)	 Lot and batch code
 Catalog number	 Manufacturer
 In vitro diagnostic medical device	 Temperature limitation
 Consult instructions for use	<b>Rx Only:</b> Prescription Use Only
 CE mark	 Authorized representative in the European Community

 12-T7528-120
  Manufactured by  
 HORIBA Instruments Incorporated - Pointe Brand  
 5449 Research Drive Canton, MI 48188
 


Manufactured by HORIBA Instruments Incorporated – Pointe Brand  
5449 Research Drive, Canton, MI 48188

European Authorized Representative:  
Obelis s.a.

Boulevard Général Wahis 53  
1030 Brussels, BELGIUM

Tel: (32)2.732.59.54 Fax:(32)2.732.60.03 email: mail@obelis.net

## Certified to Perform Reagents

The Pointe reagents are certified to be manufactured according to specified parameters. Any Pointe reagent product not meeting specifications through its listed expiration date will be remedied immediately without charge.