## Utilisation

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Test:	TP	Chemistry:	Total Protein
Chemistry No.:	229	Print Name:	Total Protein
Reaction Type:	Endpoint	Reaction Direction:	Positive
Pri. Wave:	546 nm	Sec. Wave:	670 nm
Decimal.:	0.1	Samp. Type:	Serum
Blank Time:		Reaction Time:	10 11
Unit:	g/dL	Incubation Time:	0

	Sample Vol.	Aspirated	Diluent	Reagent Vol.	Diluent
Standard;	4	uL	uL	180	uL uL
Decreased;		uL	uL	uL	
Increased;		uL	uL	uL	

Linearity Range (Standard):	1.0-15.0	Linearity Limit:	
Linearity Range (Decreased):		Substrate Depletion:	
Linearity Range (Increased):		Mixed Blank Abs.:	- 40000 40000
R1 Blank Abs.:	- 40000 40000	On-board Stability:	30 Day (s)
Blank Response	- 40000 40000	Reagent Alarm Limit:	5
Twin Chemistry:			

Prozone Check:		
Q1:	Q2:	Q3:
Q4:	PC:	ABS:

Use Qualitative Result:	
Range:	Flag:

Slope Offset:			
Slope	Offset	Unit	
1	0	g/dL	

Pretreatment:	
Pretreat Sample Vol.:	uL Pretreat Reagent Vol.: uL

Ref. Range:	
Sample Type:	Gender: Age Range: Ref. Range: Critical Range: Unit:

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Chem:	TP				
Calibration Setting		Calibrator	Conc.	Pos	Lot No.
Math Model:	Two-Point Linear	Water	0.0	W	
Factor:	Replicates: 2	Chem Cal	*	*	
Acceptance Limits					
Cal Time:	336 hr.				
Slope Diff:	SD:				
Sensitivity:	Repeatability:				* User Defined
Deter Coeff:					
Auto Calib.					
	<input type="checkbox"/> Cal Time				

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## Calcul (Exemple)

Abs. = Absorbance

$$\frac{\text{Abs. of Unknown}}{\text{Abs. of Standard}} \times \text{Conc. of standard} = \text{Total Protein (g/dl)}$$

Example: Abs. of Unknown = 0.350, Abs. of Standard = 0.400  
Concentration of Standard = 8 g/dl

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## Symbol Key

Use by (YYYY-MM-DD)	<b>LOT</b> Lot and batch code
<b>REF</b> Catalog number	Manufacturer
<b>IVD</b> In vitro diagnostic medical device	Temperature limitation
Consult instructions for use	<b>Rx Only:</b> Prescription Use Only
CE mark	<b>EC REP</b> Authorized representative in the European Community

**REF** 12-T7528-120 Manufactured by HORIBA Instruments Incorporated - Pointe Brand 5449 Research Drive Canton, MI 48188 **IVD**

Manufactured by HORIBA Instruments Incorporated – Pointe Brand  
5449 Research Drive, Canton, MI 48188

European Authorized Representative:

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