Utilisation

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Paramètres de test

Test:	URIC	Chemistry:	Uric Acid
Chemistry No.:	231	Print Name:	Uric Acid
Reaction Type:	Endpoint	Reaction Direction:	Positive
Pri. Wave:	510 nm	Sec. Wave:	670 nm
Decimal.:	0.1	Samp. Type:	Serum
Blank Time:		Reaction Time:	35 37
Unit:	mg/dL	Incubation Time:	0

	Sample Vol.	Aspirated	Diluent	Reagent Vol.	Diluent
Standard:	4	uL	uL	180	uL uL
Decreased:		uL	uL	uL	
Increased:		uL	uL	uL	

Linearity Range (Standard) 0-20	Linearity Limit:
Linearity Range (Decreased):	Substrate Depletion:
Linearity Range (Increased):	Mixed Blank Abs.: - 40000 40000
R1 Blank Abs.: 40000 40000	On-board Stability: 30 Day (s)
BlankResponse 40000 40000	Reagent Alarm Limit: 5
Twin Chemistry:	

Prozone Check:		
Q1:	Q2:	Q3:
Q4:	PC:	ABS:

Use Qualitative Result:	
Range:	Flag:

Slope Offset:			
Slope	Offset	Unit	
1	0	mg/dL	

Pretreatment:	
Pretreat Sample Vol.:	Pretreat Reagent Vol.:
uL	uL

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Ref. Range:
Sample Type: Gender: Age Range: Ref. Range: Critical Range: Unit:

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Chem: UA																					
Calibration Setting	<table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th>Calibrator</th> <th>Conc.</th> <th>Pos</th> <th>Lot No.</th> </tr> </thead> <tbody> <tr> <td>Water</td> <td>0.0</td> <td>W</td> <td></td> </tr> <tr> <td>Chem Cal</td> <td>*</td> <td>*</td> <td></td> </tr> <tr> <td> </td> <td> </td> <td> </td> <td> </td> </tr> <tr> <td> </td> <td> </td> <td> </td> <td> </td> </tr> </tbody> </table>	Calibrator	Conc.	Pos	Lot No.	Water	0.0	W		Chem Cal	*	*									
Calibrator	Conc.	Pos	Lot No.																		
Water	0.0	W																			
Chem Cal	*	*																			
Math Model: Two-Point Linear																					
Factor: Replicates: 2																					
Acceptance Limits																					
Cal Time: 336 hr.																					
Slope Diff: SD:																					
Sensitivity: Repeatability: * User Defined																					
Deter Coeff:																					
Auto Calib.																					
<input type="checkbox"/> Cal Time																					

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Calculs (Exemple)

A = Absorbance

$$\frac{A(\text{Unk})}{A(\text{Std})} \times \text{Conc. of Std (mg/dl)} = \text{Uric Acid (mg/dl)}$$

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Symboles

Use by (YYYY-MM-DD)	LOT Lot and batch code
REF Catalog number	Manufacturer
IVD In vitro diagnostic medical device	Temperature limitation
Consult instructions for use	Rx Only: Prescription Use Only
CE mark	EC REP Authorized representative in the European Community

REF 12-U7581-120	Manufactured by HORIBA Instruments Incorporated - Pointe Brand 5449 Research Drive Canton, MI 48188	IVD
Manufactured by HORIBA Instruments Incorporated – Pointe Brand 5449 Research Drive, Canton, MI 48188		
European Authorized Representative: Obelis s.a. Boulevard Général Wahis 53 1030 Brussels, BELGIUM Tel: (32)2.732.59.54 Fax:(32)2.732.60.03 email: mail@obelis.net		

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 **CE** CE mark

 **Manufacturer**

 **EC REP** Authorized representative in the European Community

 **IVD** *In vitro* diagnostic medical device