

Intended Use

For the quantitative determination of Calcium in serum or heparinized plasma using the Yumizen C560 analyzer. **Rx Only.**

Clinical Significance^{1,2}

Increased serum calcium may be observed in hyperparathyroidism, vitamin D intoxication, multiple myeloma and some neoplastic diseases of bone. Decreased serum calcium may be observed in hypoparathyroidism, vitamin D deficiency, steatorrhea, nephrosis, and nephritis.


Method History

Various methodologies have been developed for the determination of calcium including flame photometry, fluorescent, gravimetric and titrimetric procedures, ion selective electrodes, and atomic absorption. Atomic absorption has been recommended as the reference method but it requires expensive instrumentation.³

Specific dye binding methodologies have become popular for calcium determination because they are rapid, convenient and inexpensive. Procedures using the dyes alizarin 3-sulfonate and methylthymol blue have been described.^{4,5} A method using *o*-cresolphthalein complexone as the chromagen was developed in 1966 by Connerty and Biggs, and modified by Gitelman in 1967 and Baginski, et al, in 1973.^{6,7,8} *o*-Cresolphthalein complexone procedures have been widely used for the determination of calcium.

The present procedure uses Arsenazo III and has been modified to provide a highly sensitive and stable reagent system. Magnesium interference is prevented by the inclusion of 8-hydroxyquinoline sulfonate. The reagent is provided as a convenient ready to use liquid.

Principle



Calcium reacts with Arsenazo III in a slightly alkaline medium to form a purple-colored complex which absorbs at 650 nm. The intensity of the color is proportional to the calcium concentration.

Reagents

Calcium reagent: Arsenazo III \geq 0.15mM, 8-Hydroxyquinoline Sulfonate 5.0mM, Buffer, Surfactant.

Reagent Preparation

Reagent is ready to use.

Reagent Storage and Stability

Store reagent at room temperature (15-30°C). The reagent is stable until the expiration date appearing on the label when stored as directed. Manufacturer studies have shown reagent is stable for 30 days once placed in the refrigerated reagent carousel (2-10°C), however reagent stability may vary based on individual laboratory conditions.

Reagent Deterioration

Do not use if the reagent has become noticeably turbid.

Precautions and Hazards

1. This reagent is for *in vitro* diagnostic use only.
2. Reagent may be irritating to the skin. Avoid contact. Flush with water if contact occurs.

Hazards:

Hazard Classifications: Reproductive Toxicity (Category 2)

Hazard Statements: H361: Suspected of damaging fertility or the unborn child

Precautionary Statements: **Prevention:** P202 Do not handle until all safety precautions have been read and understood.

P280 Wear protective gloves/protective clothing/eye protection/face protection. **Response:** P308 + P313 IF exposed or concerned: Get medical advice/attention. **Storage:** P404 Store in a closed container. **Disposal:** P501: Dispose of contents to an approved waste disposal plant. **Refer to the Safety Data Sheet for this product (SDS-C7529) available by calling 1-734-487-8300.**



Signal Word: Warning

Specimen Collection and Storage

1. Fresh, unhemolyzed serum is the preferred specimen.
2. Heparinized plasma may also be used.
3. Anticoagulants other than heparin should not be used.⁹
4. Remove serum from clot as soon as possible since red cells can absorb calcium.¹⁰
5. Older serum specimens containing visible precipitate should not be used.^{11,12}
6. Serum calcium is stable for 24 hours at room temperature, one week at 2-8°C, and up to five months frozen (-15 to -25°C) and protected from evaporation.¹³ Specimens should not be thawed and refrozen.

Interferences

1. Substances that contain or complex with calcium cause inaccurate results.¹⁴
2. Glass tubes often are coated with a residue containing calcium. They should be acid-washed before use. Alternatively, plastic tubes may be used.
3. Bilirubin up to 20 mg/dl and hemoglobin to 500 mg/dl do not interfere.
4. Severe lipemia may cause elevated results. A serum blank should be run for greatest accuracy.
5. For a comprehensive review of interferences see Young, et al.¹⁵

Pointe Calcium (Arsenazo) Reagent Set

Materials Provided

Calcium Reagent

Materials Required but not Provided

1. Yumizen C560 Analyzer
2. Yumizen C560 Operation manual
3. Chemistry Calibrator, catalog number C7506-50
4. Chemistry control, catalog number C7592-100

Limitations

Samples with calcium values exceeding 15mg/dl¹⁶ should be diluted with an equal volume of saline, the assay repeated, and the result multiplied by two. Severely lipemic samples should be run with a serum blank for greatest accuracy.

Calibration

Use an NIST-traceable serum calibrator. The procedure should be calibrated according to the instrument manufacturer's calibration instructions. If control results are found to be out of range, the test may need to be re-calibrated. Under typical operating conditions manufacturer calibration stability studies have shown the calibration curve will be stable for at least 14 days.

Quality Control

The integrity of the reaction should be monitored by use of normal and abnormal control sera with known calcium concentrations. These controls should be run at least with every working shift in which calcium assays are performed. It is recommended that each laboratory establish its own frequency of control determination. Quality control requirements should be performed in conformance with local, state, and/or Federal regulations or accreditation requirements.

Expected Value

Adults: 8.5-10.4 mg/dl¹⁷

Newborns: 7.8-11.2 mg/dl¹⁸

It is strongly recommended that each laboratory establish its own reference range.

Performance

1. Assay Range: 0.1-15 mg/dl¹⁶
2. Comparison: A study was performed between the Yumizen C560 and a similar analyzer using this method, resulting in the following:

Method	Calcium
N	139
Mean Calcium (mg/dL)	9.32
Range (mg/dL)	0.6-14.5
Standard Deviation	3.89
Regression Analysis	$y = 1.093x - 1.02$
Correlation Coefficient	0.9808

3. Precision: Precision studies were performed using the Yumizen C560 analyzer following a modification of the guidelines which are contained in NCCLS document EP5-T2.¹⁹

Sample	Within Day		
	LOW	MID	HIGH
N	20	20	20
Mean	5.04	9.32	12.98
Standard Deviation	0.07	0.13	0.16
Coefficient of Variation (%)	1.3%	1.4%	1.2%

Sample	Total		
	LOW	MID	HIGH
N	40	40	40
Mean	5.21	9.35	12.95
Standard Deviation	0.17	0.15	0.24
Coefficient of Variation (%)	3.3%	1.6%	1.9%

4. Sensitivity: 2SD limit of detection (95% Conf) = 0.1 mg/dL

References

1. Tietz, N.W., Fundamentals of Clinical Chemistry, Philadelphia, W.B. Saunders, p. 149 (1984).
2. Henry, J.B., Clinical Diagnosis and Management by Laboratory Methods, Philadelphia, W.B. Saunders, p. 149 (1984).
3. Cali, J.P., et al, N.B.S., Sp. Publication 260:36 (1972).
4. Connerty, H.V. and Biggs, A.R., Am. J. Clin. Chem. 11:716 (1965).
5. Gindler, E.M. and King, J.D., Am. J. Clin. Path. 58:376 (1972).
6. Connerty, H.V. and Biggs, A.R., Am. J. Clin. Path. 45:290 (1966).
7. Gitelman, H.J., Anal. Biochem. 18:521 (1967).
8. Baginski, E.S., et al, Clin. Chem. Acta 46:49 (1973).
9. Richerich, R., Clinical Chemistry: Theory and Practice, New York, Academic Press, p. 304 (1969).
10. Peters, J.P., Van Slyke, D.D., Quantitative Clinical Chemistry – Vol. 2, Baltimore, Williams and Wilkins, (1932).
11. Chen, P.S., et al, Anal. Chem. 26:1967 (1954).
12. Tayeau, F., et al, Bull. Soc. Pharm. Bordeaux, 95:206 (1956).
13. Henry, R.J., et al, Clinical Chemistry: Principles and Technics, Hagerstown (MD), Harper and Row, p. 669 (1974).
14. Tietz, N.W., Fundamentals of Clinical Chemistry, Philadelphia, W.B. Saunders, pp. 903-904 (1984).
15. Young, D.S., et al, Clin. Chem. 21:1D (1975).
16. HORIBA Medical records.
17. Tietz, N.W., Fundamentals of Clinical Chemistry, Philadelphia, W.B. Saunders, p. 1208 (1984).
18. Meites, Samuel, Pediatric Clinical Chemistry, Washington DC, AACC Press, p. 81 (1989).
19. NCCLS document "Evaluation of Precision Performance of Clinical Chemistry Devices", 2nd Ed. (1992).

CHEMISTRY PARAMETERS

Chem:	CA	No.:	209	Sample Type:	Serum
Chemistry:	Calcium (Arsenazo)			Print Name:	CA
Reaction Type:	End Point			Reaction Direction:	Positive
Pri Wave:	660			Sec Wave:	
Unit:	mg/dL			Decimal:	0.1
Blank Time:	10 12			Reaction Time:	27 29
	Sample Vol.	Aspirated	Diluent	Reagent Vol.	Diluent
Standard:	1.5 ul	-- ul	-- ul	R1: 125 ul	-- ul
Decreased:	-- ul	-- ul	-- ul	R2: -- ul	-- ul
Increased:	-- ul	-- ul	-- ul	R3: -- ul	-- ul
	<input type="checkbox"/> Sample Blank	<input checked="" type="checkbox"/> Auto Rerun		R4: -- ul	-- ul
 <u>Slope/Offset Adjustment</u>					
	Slope: 1	Offset: 0			

Linearity Range (Standard)	0.1	15	Linearity Limit:	
Linearity Range (Decreased)	---	---	Substrate Depletion:	
Linearity Range (Increased)	---	---	Mixed Blank Abs:	
R1 Blank Abs:	---	---	Uncapping Time	
Blank Response:	---	---	Reagent Alarm Limit:	
Twin Chemistry:			<input type="checkbox"/> Enzyme Linear Extension	
<input type="checkbox"/> Prozone Check		<input type="radio"/> Rate Check	<input type="radio"/> Antigen Addition	
Q1:	Q2:	Q3:	Q4:	
PC:	ABS:			

Pointe Calcium (Arsenazo) Reagent Set

CALIBRATION PARAMETERS

Calibrator Definition						
Calibrator:	*	Lot No.:			*	
Exp Date:	*					
Carousel		Pos				
Sample Carousel 1	*					
Sample Carousel 2						
Sample Carousel 3						
Reagent/Calibration						
<u>Calibrator</u>	<u>Pos</u>	<u>Lot No</u>	<u>Exp Date</u>	<u>Chem</u>	<u>Conc</u>	<u>Unit</u>
Water	W	*	*	CA	0	mg/dL
Chemistry Calibrator	*	*	*	CA	*	mg/dL
Calibration Setup						
Chem:	CA					
<u>Calibration Settings</u>						
Math Model:	Two-Point Linear					
Factor:	Replicates:		2			
<u>Acceptance Limits</u>						
Cal Time:	336	Hour				
Slope Diff:	---	SD:	---			
Sensitivity :	---	Repeatability:	---			
Deter Coeff:	---					
<u>Auto Calib.</u>						
<input type="checkbox"/> Bottle Changed	<input type="checkbox"/> Lot Changed	<input type="checkbox"/> Cal Time				

It is recommended that two levels of control material be assayed daily.
* Indicates user defined parameter.

REF 14-C7529-360



Manufactured by
HORIBA Instruments Incorporated-Pointe Brand
5449 Research Drive Canton, MI 48188



Certified to Perform Reagents

The Pointe reagents are certified to be manufactured according to specified parameters. Any Pointe reagent product not meeting specifications through its listed expiration date will be remedied immediately without charge.

Manufactured by HORIBA Instruments Incorporated – Pointe Brand
5449 Research Drive, Canton, MI 48188

European Authorized Representative:
Obelis s.a.

Boulevard Général Wahis 53
1030 Brussels, BELGIUM

Tel: (32)2.732.59.54 Fax:(32)2.732.60.03 email: mail@obelis.net



Symbol Key



Use by (YYYY-MM-DD)



Lot and batch code



Catalog number



Manufacturer



Temperature limitation



Consult instructions for use



In vitro diagnostic medical device **Rx Only:** Prescription Use Only