

Intended Use


For the quantitative determination of Carbon Dioxide in serum using the Yumizen C560 analyzer. For *in vitro* diagnostic use only. **Rx Only.**

Method History

Early methods for the determination of carbon dioxide were based on either volumetric or manometric determination of the CO₂ released from a sample by acid treatment. These methods used the instruments of Van Slyke^{1,2} until they were replaced by the Natelson microgasometer,³ which still uses manometric determination of total CO₂.

Methods have been developed for Auto Analyzers⁴ but these suffer from baseline drift⁵ and require equipment which many laboratories do not have. Enzymatic methods for CO₂ have been introduced by Wilson,⁶ Menson⁷ and Norris⁸ using phosphoenolpyruvate carboxylase. The present procedure is an enzymatic assay utilizing Phosphoenolpyruvate Carboxylase (PEPC) and a NADH analog.

Principle



Carbon Dioxide (in the form of bicarbonate ions) reacts with phosphoenolpyruvate (PEP), in the presence of phosphoenolpyruvate carboxylase (PEPC), to form oxaloacetate. The cofactor then in the presence of malate dehydrogenase (MDH) is oxidized by the oxaloacetate. The decrease in absorbance monitored between 405 and 415 nm resulting is proportional to the amount of CO₂ in the sample.

Clinical Significance⁵

The measurement of Carbon Dioxide is useful in the assessment of acid-base balance disturbances. Elevated CO₂ is observed in metabolic alkalosis and compensated respiratory acidosis. Low CO₂ is observed in compensated respiratory alkalosis and metabolic acidosis. Differentiation between the metabolic and respiratory conditions is only possible through additional laboratory determinations.

Reagents

CO₂ reagent: PEP 6mM, Magnesium Ions 10mM, NADH analog, MDH (porcine) ≥ 1200U/L, PEPC (microbial) ≥ 200U/L, Buffer, pH 7.4 ± 0.1 non-reactive stabilizers with surfactants and preservative.

Reagent Preparation

Reagent provided as a ready to use liquid.

Reagent Storage and Stability

Reagent is stable until expiration date indicated on vial label when stored tightly capped at 2-8°C. (15 months from date of manufacture).

Reagent Deterioration

1. Reagent should appear clear and pale yellow in color.
2. Do not use if reagent appears to be turbid, this would indicate deterioration.

Precautions and Hazards

1. Reagents are for *in vitro* diagnostic use only.
2. Do not ingest. Toxicity has not been established.
3. Do not pipet by mouth to avoid CO₂ contamination from the expired air.

Hazards:

Hazard Classifications: Not a hazardous substance or mixture.

Pictogram: Not required.

Signal Word: Not required.

Hazard Statements: Not a hazardous substance or mixture.

Precautionary Statements: Not a hazardous substance or mixture

Refer to the Safety Data Sheet for this product (SDS-C7502) available by calling 1-734-487-8300

Specimen Collection and Storage

1. Fresh, unhemolyzed serum collected under anaerobic conditions is the recommended specimen.
2. The sample may be stored in ice water under anaerobic conditions for up to one hour.⁹

Interferences

1. Interferences were evaluated for this carbon dioxide method on a Yumizen C560 analyzer. No interference was observed by bilirubin up to 20.0 mg/dl, hemoglobin up to 219.9 mg/dl and lipemia (intralipid) up to 1000 mg/dl. (Using a criteria of >10% variance from control. CO₂ level was 26, 23, 25 mmol/L respectively)
2. CO₂ from air or the breath of the analyst is a major interference in this assay. Reagent handling, specimen collection, and all storage instructions must be strictly followed to minimize this interference.
3. A number of conditions and substances have been reported to affect serum Carbon Dioxide levels.^{10,11,12}

Pointe Carbon Dioxide Reagent Set

Materials Provided

Carbon Dioxide Reagent

Materials Required but not Provided

1. Yumizen C560 Analyzer
2. Yumizen C560 Operation manual
3. Chemistry Calibrator, catalog number C7506-50
4. Chemistry control, catalog number C7592-100

Limitations

1. Samples exceeding 40 mmol/L must be diluted 1:1 with saline, re-assayed, and the result multiplied by two.
2. Carbon Dioxide contamination must be avoided. Keep reagent tightly capped when not in use.

Calibration

Use an NIST-traceable serum calibrator. The procedure should be calibrated according to the instrument manufacturer's calibration instructions. If control results are found to be out of range, the test may need to be re-calibrated. Under typical operating conditions manufacturer calibration stability studies have shown the calibration curve will be stable for at least 1 day.

Quality Control

To monitor the reliability of results, two levels of control sera with known Carbon Dioxide values should be run with patient samples. Quality control requirements should be performed in conformance with local, state, and/or Federal regulations or accreditation requirements.

Expected Values ⁹

23-34 mmol/L

It is strongly recommended that each laboratory determine its own reference range.

Performance

1. Assay Range: 2 - 40 mmol/L
2. Comparison: A study was performed between the Yumizen C560 and a similar analyzer using this method, resulting in the following:

Method	Carbon Dioxide
N	97
Mean CO ₂ (mmol/L)	22.5
Range (mmol/L)	4-39
Standard Deviation	6.4
Regression Analysis	$y = 0.962x - 2.2$
Correlation Coefficient	0.9543

3. Precision: Within Day precision was investigated by running two samples in replicates of 20 on the same day. Day to Day results were obtained by performing one run per day over a span of 20 days. Precision studies were performed using the Yumizen C560 analyzer following a modification of the guidelines which are contained in NCCLS document EP5-T2.¹³

Sample	Within Day		
	LOW	MID	HIGH
N	20	20	20
Mean	7.8	23.9	30.9
Standard Deviation	0.4	0.5	0.6
Coefficient of Variation (%)	5.7%	2.1%	2.0%

Sample	Total		
	LOW	MID	HIGH
N	40	40	40
Mean	10.0	22.6	29.4
Standard Deviation	1.1	2.0	1.0
Coefficient of Variation (%)	10.5%	8.7%	3.3%

4. Sensitivity: 2SD limit of detection (95% Conf) = 1 mmol/L

References

1. Van Slyke, D.D. and Stadie, W.C., J. Biol. Chem. 49:1 (1921).
2. Van Slyke, D.D. and Neil, J.M., J. Biol. Chem. 61:523 (1924).
3. Natelson, S., Microtechniques of Clinical Chemistry, C. Thomas, Springfield, IL. P.147 (1961).
4. Skeggs, L.T. Jr., Am. J. Clin. Path. 33:181 (1960).
5. Tietz, N.W., Fundamentals of Clinical Chemistry, W.B. Saunders, Philadelphia, PA., pp 884-887 (1982).
6. Wilson, W., et al, Clin. Chem. 19:640 (1973).
7. Menson, R.C., et al, Clin. Chem. 20:872 (1974).
8. Norris, K.A., et al, Clin. Chem. 21:1093 (1975).
9. Henry, R.J., Clinical Chemistry: Principles and Technics, Harper & Row, New York, NY, p784 (1974).
10. Young, D.S., et al, Clin. Chem. 21:1D (1975).
11. Martin, E.W., In Hazard of Medication (Alexander, S.F., Farage, D.J., and Hassan, W.E., Jr. eds.), J.B. Lippincott Co., Philadelphia, PA., and Toronto, Canada, p. 169 (1971).
12. Constantino, N.V., and Kabat, H.F., Am. J. Hosp. Pharm. 30:24 (1973).
13. NCCLS document "Evaluation of Precision Performance of Clinical Chemistry Devices", 2nd Ed. (1992)

CHEMISTRY PARAMETERS

Chem:	CO2	No.:	208	Sample Type:	Serum
Chemistry:	Carbon Dioxide			Print Name:	CO2
Reaction Type:	Fixed Time			Reaction Direction:	Negative
Pri Wave:	412			Sec Wave:	505
Unit:	mmol/L			Decimal:	0
Blank Time:	10	12		Reaction Time:	18 41
	Sample Vol.	Aspirated	Diluent	Reagent Vol.	Diluent
Standard:	1.5 ul	-- ul	-- ul	R1: 150 ul	-- ul
Decreased:	-- ul	-- ul	-- ul	R2: -- ul	-- ul
Increased:	-- ul	-- ul	-- ul	R3: -- ul	-- ul
	<input type="checkbox"/> Sample Blank	<input checked="" type="checkbox"/> Auto Rerun		R4: -- ul	-- ul
<u>Slope/Offset Adjustment</u>					
Slope: 1		Offset: 0			

Linearity Range (Standard)	2	40	Linearity Limit:
Linearity Range (Decreased)	---	---	Substrate Depletion:
Linearity Range (Increased)	---	---	Mixed Blank Abs:
R1 Blank Abs:	---	---	Uncapping Time
Blank Response:	---	---	Reagent Alarm Limit:
Twin Chemistry:			<input type="checkbox"/> Enzyme Linear Extension
<input type="checkbox"/> Prozone Check		<input type="radio"/> Rate Check	<input type="radio"/> Antigen Addition
Q1:	Q2:	Q3:	Q4:
PC:	ABS:		

Pointe Carbon Dioxide Reagent Set

CALIBRATION PARAMETERS

Calibrator Definition						
Calibrator:	*	Lot No.:	*			
Exp Date:	*					
Carousel		Pos				
Sample Carousel 1	*					
Sample Carousel 2						
Sample Carousel 3						
Reagent/Calibration						
<u>Calibrator</u>	<u>Pos</u>	<u>Lot No</u>	<u>Exp Date</u>	<u>Chem</u>	<u>Conc</u>	<u>Unit</u>
Water	W	*	*	CO2	0	mmol/L
Chemistry Calibrator	*	*	*	CO2	*	mmol/L
Calibration Setup						
Chem:	CO2					
Calibration Settings						
Math Model:	Two-Point Linear					
Factor:		Replicates:	2			
Acceptance Limits						
Cal Time:	24	Hour				
Slope Diff:	---	SD:	---			
Sensitivity :	---	Repeatability:	---			
Deter Coeff:	---					
Auto Calib.						
<input type="checkbox"/> Bottle Changed	<input type="checkbox"/> Lot Changed	<input type="checkbox"/> Cal Time				

It is recommended that two levels of control material be assayed daily.
* Indicates user defined parameter.

REF 14-C7502-160



Manufactured by
HORIBA Instruments Incorporated-Pointe Brand
5449 Research Drive Canton, MI 48188



Certified to Perform Reagents

The Pointe reagents are certified to be manufactured according to specified parameters. Any Pointe reagent product not meeting specifications through its listed expiration date will be remedied immediately without charge.

Manufactured by HORIBA Instruments Incorporated – Pointe Brand
5449 Research Drive, Canton, MI 48188

European Authorized Representative:
Obelis s.a.

Boulevard Général Wahis 53
1030 Brussels, BELGIUM

Tel: (32)2.732.59.54 Fax:(32)2.732.60.03 email: mail@obelis.net



Symbol Key



Use by (YYYY-MM-DD)



Lot and batch code



Catalog number



Manufacturer



Temperature limitation



Consult instructions for use



In vitro diagnostic medical device **Rx Only:** Prescription Use Only