## Utilisation

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8. NCCLS document, "Protection of Laboratory Workers from Infectious Disease Transmitted by Blood, Body Fluids, and Tissue", 2<sup>nd</sup> Ed. (1991).
9. NCCLS document, "Procedures for the Collection of Diagnostic Blood Specimens by Skin Puncture", 3<sup>rd</sup> Ed. (1991).
10. NCCLS document, "National Evaluation Protocols for Interference Testing", Evaluation Protocol Number 7, Vol. 4, No. 8, (June 1984).
11. Young, D.S., et al, *Clin Chem* 21:1D, 1975.
12. NCCLS document, "Evaluation of Precision Performance of Clinical Chemistry Devices", 2<sup>nd</sup> Ed. (1992)
13. Liuzzo, G., et al, *N Eng J Med*, 331:417-424, 1994.
14. U.S. Patent nos. 6,248,597; 6,828, 158.

**PARAMETRAGE CHIMIE**

Chem:	CRP	No.:	214	Sample Type:	Serum
Chemistry:	CRP			Print Name:	CRP
Reaction Type:	End Point			Reaction Direction:	Positive
Pri Wave:	570			Sec Wave:	
Unit:	mg/L			Decimal	0.1
Blank Time:	47      49			Reaction Time:	80      82
	Sample Vol.	Aspirated	Diluent	Reagent Vol.	Diluent
Standard:	4.7 UI	--- ul	--- ul	R1: 120 ul	--- ul
Decreased:	4.7 UI	30 ul	120 ul	R2: 120 ul	-- ul
Increased:	--- UI	--- ul	--- ul	R3: --- ul	-- ul
	<input type="checkbox"/> Sample Blank	<input checked="" type="checkbox"/> Auto Rerun		R4: --- ul	--- ul
<b><u>Slope/Offset Adjustment</u></b>					
Slope: 1		Offset: 0			

Linearity Range (Standard)	0.1	320	Linearity Limit:
Linearity Range (Decreased)	---	---	Substrate Depletion:
Linearity Range (Increased)	---	---	Mixed Blank Abs:
R1 Blank Abs:	---	---	Uncapping Time
Blank Response:	---	---	Reagent Alarm Limit:
Twin Chemistry:			<input type="checkbox"/> Enzyme Linear Extension
<input type="checkbox"/> Prozone Check		<input type="radio"/> Rate Check	<input type="radio"/> Antigen Addition
Q1:	Q2:	Q3:	Q4:
PC:	ABS:		

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## PARAMETRAGE D'ETALONNAGE

<b>Calibrator Definition</b>						
Calibrator:	*	Lot No.:			*	
Exp Date:	*					
<b>Carousel</b>	<b>Pos</b>					
Sample Carousel 1	*					
Sample Carousel 2						
Sample Carousel 3						
<b>Reagent/Calibration</b>						
<u>Calibrator</u>	<u>Pos</u>	<u>Lot No</u>	<u>Exp Date</u>	<u>Chem</u>	<u>Conc</u>	<u>Unit</u>
Water	W	*	*	*	*	mg/L
CRP Cal 1	*	*	*	*	*	mg/L
CRP Cal 2	*	*	*	*	*	mg/L
CRP Cal 3	*	*	*	*	*	mg/L
CRP Cal 4	*	*	*	*	*	mg/L
CRP Cal 5	*	*	*	*	*	mg/L
<b>Calibration Setup</b>						
Chem:	CRP					
<u>Calibration Settings</u>						
Math Model:	Spline					
Factor:	Replicates:		2			
<u>Acceptance Limits</u>						
Cal Time:	336	Hour				
Slope Diff:	---	SD:	---			
Sensitivity :	---	Repeatability:	---			
Deter Coeff:	---					
<u>Auto Calib.</u>						
<input type="checkbox"/> Bottle Changed		<input type="checkbox"/> Lot Changed		<input type="checkbox"/> Cal Time		

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**REF** 14-C7568-480

Manufactured by  
HORIBA Instruments Incorporated-Pointe Brand  
5449 Research Drive Canton, MI 48188

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Manufactured by HORIBA Instruments Incorporated – Pointe Brand  
5449 Research Drive, Canton, MI 48188

European Authorized Representative:  
Obelis s.a.

Boulevard Général Wahis 53  
1030 Brussels, BELGIUM

Tel: (32)2.732.59.54 Fax:(32)2.732.60.03 email: mail@obelis.net

### Symbol Key

Use by (YYYY-MM-DD)	<b>LOT</b> Lot and batch code	<b>REF</b> Catalog number
Manufacturer	Temperature limitation	Consult instructions for use
<b>IVD</b> In vitro diagnostic medical device <b>Rx Only:</b> Prescription Use Only		