Utilisation

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Calibration

The autoHDL/LDL™ Cholesterol Calibrator is required for calibration. The value of the autoHDL/LDL™ calibrator was assigned by procedures traceable to the National Reference System for Cholesterol (NRS/CHOL). Refer to autoHDL/LDL™ Cholesterol Calibrator kit package insert for instructions. If control results are found to be out of range, the test may need to be re-calibrated. Under typical operating conditions manufacturer calibration stability studies have shown the calibration curve will be stable for at least 14 days.

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Sample	Within Day			Total		
	LOW	MID	HIGH	LOW	MID	HIGH
N	20	20	20	40	40	40
Moyenne	32.6	67.7	135.6	32.8	67.4	132.6
Écart-type	0.6	0.7	0.8	6.6	0.6	1.2
Coefficient de variation (%)	1.8%	1.1%	0.6%	1.9%	1.8%	5.0%

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Références

- Gotto A.M., Lipoprotein metabolism and the etiology of hyperlipidemia, Hospital practice, 23:Suppl.1,4 (1988).
- Crouse J.R. et al, Studies of low density lipoprotein molecular weight in human beings with coronary artery disease, J. Lipid Res., 26:566 (1985).
- Badimon J.J., Badimon L., Fuester V., Regression of Atherosclerotic Lesions by High-Density Lipoprotein Plasma Fraction in the Cholesterol-Fed Rabbit, Journal of Clinical Investigation, 1990; 85:1234-41.
- Castelli, W.P., et al, Cholesterol and other lipids in coronary heart disease, Circulation, 55:767 (1977).
- Barr, D.P., Russ E.M., Eder H.A., Protein-lipid relationships in human plasma, Am. J. Med. 11:480 (1951).
- Gordon, T., et al, High density lipoprotein as a protective factor against coronary heart disease, Am. J. Med., 62:707 (1977).
- Williams, P., Robinson D., Baily A., High density lipoprotein and coronary risk factor, Lancet, 1:72 (1979).
- Kannel W.B., Castelli W.P., Gordon, T., Cholesterol in the prediction of atherosclerotic disease; New perspectives based on the Framingham study, Am. Intern. Med., 90:85 (1979).
- National Institute on Health publication No. 93-3095, September 1993.
- Warnick G. Russell, Wood Peter D., National Cholesterol Education Program Recommendations for Measurement of High-Density Lipoprotein Cholesterol; Executive Summary Clinical Chemistry, Vol.41, No. 10 (1995).
- Grundy, S.M., et al, Summary of the Second Report of the National Cholesterol Education Program (NCEP) Expert Panel on Detection, Evaluation, and Treatment of High Blood Cholesterol in Adults (Adult Treatment Panel II), JAMA 1993, 269: 23, 3015-3023.
- National Committee for Clinical Laboratory Standards, National Evaluation Protocols for Interference Testing, Evaluation Protocol Number 7, Vol.4, No. 8, June 1984.
- Young, D.S., Effects of Drugs on Clinical Laboratory Tests, 3rd. Ed., AACC Press, Washington DC, 1990, 3-104 thru 3-106.
- Tietz, N.W., Clinical Guide to Laboratory Tests, W.B. Saunders Co., Philadelphia, 1986, p.256.
- NCCLS document "Evaluation of Precision Performance of Clinical Chemistry Devices" 2nd Ed. 1992.

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CHEMISTRY PARAMETERS

Chem:	HDL	No.:	218	Sample Type:	Serum	
Chemistry:	autoHDL Cholesterol	Print Name:	HDL	Reaction Direction:	Positive	
Reaction Type:	Fixed Time	Pri Wave:	605	Sec Wave:	660	
Unit:	mg/dL	Decimal:	0	Reaction Time:	51 82	
Blank Time:	47 49	Sample Vol.	Aspirated	Diluent	Reagent Vol.	Diluent
Standard:	2.0 ul	-- ul	-- ul	-- ul	R1: 120 ul	-- ul
Decreased:	-- ul	-- ul	-- ul	-- ul	R2: 40 ul	-- ul
Increased:	-- ul	-- ul	-- ul	-- ul	R3: -- ul	-- ul
	<input type="checkbox"/> Sample Blank	<input checked="" type="checkbox"/> Auto Rerun			R4: -- ul	-- ul
<u>Slope/Offset Adjustment</u>						
Slope: 1		Offset: 0				

Linearity Range (Standard)	2	150	Linearity Limit:
Linearity Range (Decreased)	---	---	Substrate Depletion:
Linearity Range (Increased)	---	---	Mixed Blank Abs:
R1 Blank Abs:	---	---	Uncapping Time
Blank Response:	---	---	Reagent Alarm Limit:
Twin Chemistry:			<input type="checkbox"/> Enzyme Linear Extension
<input type="checkbox"/> Prozone Check		<input type="radio"/> Rate Check	<input type="radio"/> Antigen Addition
Q1:	Q2:	Q3:	Q4:
PC:	ABS:		

CALIBRATION PARAMETERS

Calibrator Definition						
Calibrator:	*	Lot No.:	*			
Exp Date:	*					
Carousel		Pos				
Sample Carousel 1	*					
Sample Carousel 2						
Sample Carousel 3						
Reagent/Calibration						
<u>Calibrator</u>	<u>Pos</u>	<u>Lot No</u>	<u>Exp Date</u>	<u>Chem</u>	<u>Conc</u>	<u>Unit</u>
Water	W	*	*	HDL	0	mg/dL
autoHDL/LDL Calibrator	*	*	*	HDL	*	mg/dL
Calibration Setup						
Chem:	HDL					
Calibration Settings						
Math Model:	Two-Point Linear					
Factor:		Replicates:	2			
Acceptance Limits						
Cal Time:	336	Hour				
Slope Diff:	---	SD:	---			
Sensitivity :	---	Repeatability:	---			
Deter Coeff:	---					
Auto Calib.						
<input type="checkbox"/> Bottle Changed	<input type="checkbox"/> Lot Changed	<input type="checkbox"/> Cal Time				

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REF 14-H7545-300

Manufactured for HORIBA
Instruments Incorporated – Pointe Brand
5449 Research Drive Canton, MI 48188

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Manufactured by HORIBA Instruments Incorporated – Pointe Brand
5449 Research Drive, Canton, MI 48188

European Authorized Representative:

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