

Intended Use

The Liquid Stable 2-Part Homocysteine Reagent is intended for *in vitro* quantitative determination of total homocysteine in human serum and plasma on the Yumizen C560 analyzer. The device can assist in the diagnosis and treatment of patients suspected of having hyperhomocysteinemia and homocystinuria. **Rx Only.**

Clinical Significance

Homocysteine (HCY) is a thiol-containing amino acid produced by the intracellular demethylation of methionine. Homocysteine is exported into plasma where it circulates, mostly in its oxidized form, bound to plasma proteins as a protein-HCY mixed disulfide with albumin (protein-SS-HCY).^{1,5} Smaller amounts of reduced homocysteine and the disulfide homocysteine (HCY-SS-HCY) are present. Total homocysteine (tHCY) represents the sum of all the HCY species found in serum or plasma (free plus protein bound). Homocysteine is metabolized to either cysteine or methionine. In the vitamin B6 trans-sulphuration pathway, homocysteine is irreversibly catabolized to cysteine. A major part of homocysteine is remethylated to methionine, mainly by the folate and cobalamin-dependent enzyme methionine synthase. Homocysteine accumulates and is excreted into blood when these reactions are impaired.^{3,5} Severely elevated concentrations of total homocysteine are found in subjects with homocystinuria, a rare genetic disorder of the enzymes involved in the metabolism of homocysteine. Patients with homocystinuria exhibit mental retardation, early arteriosclerosis and arterial and venous thromboembolism.^{2,6} Other less severe genetic defects which lead to moderately elevated levels of total homocysteine are also found.⁷⁻⁹

Epidemiological studies have investigated the relationship between elevated homocysteine levels and cardiovascular disease (CVD). A meta-analysis of 27 of these studies, including more than 4000 patients, estimated that a 5 µmol/L increase in total homocysteine was associated with an odds ratio for coronary artery disease (CAD) of 1.6 (95% confidence interval [CI], 1.4 to 1.7 for men and 1.8 (95% CI 1.3 to 1.9) for women; the odds ratio for cerebrovascular disease was 1.5 (95% CI 1.3 to 1.9). The risk associated with a 5 µmol/L increase in total homocysteine was the same as that associated with 0.5 mmol/L (20 mg/dL) increase in cholesterol. Peripheral arterial disease also showed a strong association.¹⁰

Hyperhomocysteinemia, elevated levels of homocysteine, can be associated with an increased risk of CVD. There have also been many published reports of prospective studies on the relationship between hyperhomocysteinemia and risk of CVD in men and women who were initially healthy. End points were based on a cardiovascular event such as acute myocardial infarction, stroke, CAD, or mortality. The results of eleven of these nested case-control studies reviewed by Cattaneo¹¹ were equivocal where five of the studies support the association with risk and six do not. More recently homocysteine levels were determined in a prospective study of post-menopausal women who participated in the Women's Health Study. Specimens from 122 women, who subsequently developed cardiovascular events, were tested for homocysteine and compared to a control group of 244 women who were matched for age and smoking status. The women in the control group remained free of disease during the three year follow-up period. The results demonstrated that post-menopausal women who developed cardiovascular events had significantly higher baseline homocysteine levels. Those with levels in the highest quartile had a two-fold increase in risk of any cardiovascular event. Elevated baseline homocysteine levels were shown to be an independent risk factor.¹² Also, homocysteine levels were determined in 1933 elderly men and women for the Framingham Heart Study cohort and demonstrated that elevated levels of homocysteine are independently associated with increased rates of all-cause and CVD mortality.¹³

Patients with chronic renal disease experience an excess morbidity and mortality due to arteriosclerotic CVD. Elevated concentration of homocysteine is a frequently observed finding in the blood of these patients.

Although such patients lack some of the vitamins involved in the metabolism of homocysteine, the elevated HCY levels are mainly due to impaired HCY removal from the blood by the Kidneys.^{14,15}


Recent evidence has also implicated elevated blood levels of homocysteine in miscarriages and birth defects.¹⁶

Drugs such as methotrexate, carbamazepine, phenytoin, nitrous oxide, and 6-azauridine triacetate interfere with HCY metabolism and may give elevated levels of HCY.¹⁷

Test Summary and Principle

Bound or dimerized homocysteine (oxidized form) is reduced to free homocysteine, which then reacts with serine catalyzed by cystathionine beta-synthase (CBS) to form cystathionine. Cystathionine in turn is broken down by cystathionine beta-lyase (CBL) to form homocysteine, pyruvate and ammonia. Pyruvate is then converted by lactate dehydrogenase (LDH) to lactate with nicotinamide adenine dinucleotide (NADH) as coenzyme. The rate of NADH conversion to NAD⁺ is directly proportional to the concentration of homocysteine (Δ A340 nm).

Reduction: Dimerized homocysteine, mixed disulfide, and protein-bound forms of HCY in the sample are reduced to form free HCY by the use of tris [2-carboxyethyl] phosphine (TCEP).



Enzymatic Conversion: Free HCY is converted to cystathionine by the use of cystathionine beta-synthase and excess serine. The cystathionine is then broken down to homocysteine, pyruvate and ammonia. Pyruvate is converted to lactate via lactate dehydrogenase with NADH as coenzyme. The rate of NADH conversion to NAD⁺ (Δ A340 nm) is directly proportional to the concentration of homocysteine.

Reagents

R1 reagent: NADH (0.47 mM), LDH (38 KU/L), Serine (0.76 mM), Trizma Base 1-10%, Trizma Hydrochloride 1-10%, Sodium Azide < 1%. Reductant (TCEP:2.9 mM)

R2 reagent: Cycling Enzymes CBS (0.748 KU/L) and CBL (16.4 KU/L) Sodium Azide < 1%.

Calibrator 1: Aqueous homocysteine blank (0 µmol/L).

Calibrator 2: Aqueous homocysteine solution (28 µmol/L).

Reagent Preparation


R1 and R2 are packaged ready to use. The reagents are stable until the expiration date specified on the label.

Indications of Deterioration

The reagents should be clear of particulate material. They should be discarded if they become turbid.

Calibrator Preparation and Use

The calibrators are prepared gravimetrically and are traceable to Standard Reference Material NIST SRM 1955, confirmed by a designated measurement procedure (HPLC). The calibrators are supplied in the kit and are provided ready to use. Values are printed on the labels. Calibration stability studies have shown the calibration curve will be stable for at least 14 days.



Pointe Homocysteine Reagent

Precautions and Hazards

1. Adhere strictly to the instructions in this insert, particularly for handling and storage conditions.
2. Reagent 1 and Reagent 2 contain sodium azide which can react with lead or copper plumbing to form highly explosive metal azides. On disposal, flush with large quantities of water to prevent azide build-up.

Hazards:

R1: Hazard Classifications: Acute Toxicity, Oral (Category 4)

Hazard Statements: H302: Harmful if swallowed

Precautionary Statements: **Prevention:** P264: Wash skin thoroughly after handling. P270: Do not eat, drink or smoke when using this product. **Response:** P330: Rinse Mouth. P301 + P312: IF SWALLOWED: call a POISON CENTER or doctor/physician IF you feel unwell. **Disposal:** P501: Dispose of contents into sewer system after diluting with large volumes of water, if in accordance with local regulations.

R2: Hazard Classifications: Acute Toxicity, Oral (Category 4)

Hazard Statements: H302: Harmful if swallowed

Precautionary Statements: **Prevention:** P264: Wash skin thoroughly after handling. P270: Do not eat, drink or smoke when using this product. **Response:** P330: Rinse Mouth. P301 + P312: IF SWALLOWED: call a POISON CENTER or doctor/physician IF you feel unwell. **Disposal:** P501: Dispose of contents into sewer system after diluting with large volumes of water, if in accordance with local regulations.

Cal 1 and Cal 2: Hazard Classifications: Not a hazardous substance or mixture.

Pictogram and Signal Word: Not required.

Hazard Statements: Not a hazardous substance or mixture.

Precautionary Statements: Not a hazardous substance or mixture. **Refer to the Safety Data Sheet for this product (SDS-H7575) available by calling 1-734-487-8300**



R22: Harmful if swallowed.

R32: Contact with acids liberates very toxic gas.

S36/37/39: Wear suitable protective clothing, gloves and eye/face protection.

S29/35: Do not empty into drains; dispose of this material and its container in a safe way.

S46: If swallowed, seek medical advice immediately and show this container or label.

Reagent Storage

1. Store kit components at 2-8°C and use until the expiry date on the labels. Do not use expired reagents. Once placed on board reagent is stable for 30 days.
2. Reagents may be used on multiple occasions until the expiry date on the labels. Reagents **must** be returned to 2-8°C storage between use.
3. Do not mix different reagent kit lot numbers.
4. **DO NOT FREEZE REAGENTS.**
5. Do not expose Reagent 1 and Reagent 2 to light during on-board use.
6. Avoid contamination of reagents. Use a new disposable pipette tip for each reagent or sample manipulation.

Specimen Collection and Handling

1. Serum (collected in serum or serum separator tubes) and plasma (collected in potassium EDTA or lithium heparin tubes) may be used for the measurement of homocysteine. However, it is not recommended to use individual patient results from serum, heparinized plasma and EDTA plasma interchangeably.²⁷ Additionally matrix differences between serum and serum separator tubes and plasma tubes have been reported.¹⁹
To minimize increases in homocysteine concentration from synthesis by red blood cells, process specimens as follows:
 - Place all specimens (serum and plasma) on ice after collection and prior to processing. Serum may clot more slowly and the volume may be reduced.¹⁷
 - All specimens may be kept on ice for up to 6 hours prior to separation by centrifugation.¹⁷
 - Separate red blood cells from serum or plasma by centrifugation and transfer to a sample cup or other clean container.**Note:** Specimens not placed on ice immediately may exhibit a 10-20% increase in homocysteine concentration.¹⁸
2. If the assay will be performed within 2 weeks after collection, the specimen should be stored at 2-8°C. If the testing will be delayed more than 2 weeks, the specimen should be stored frozen at -20°C or colder. Specimens have been shown to be stable at -20°C for 8 months.^{17,19}
3. It is the responsibility of the operator to verify the correct specimen type(s) is (are) used in the liquid stable 2-Part Homocysteine Reagent.
4. Inspect all samples (specimens, calibrators and controls) for bubbles. Remove bubbles prior to analysis.
5. Specimens containing particulate matter (fibrin, red blood cells, or other matter) and visibly lipemic specimens should not be used with the assay. Results from these specimens may be inaccurate.
6. Mix specimens **thoroughly** after thawing by low speed vortexing or by gentle inversion to ensure consistency in results. Avoid repeated freezing and thawing. Specimens showing particulate matter, erythrocytes, or turbidity should be centrifuged before testing.
7. On-board instrument storage. EDTA plasma samples can be stored for 3 hours on-board the AU400. The other recommended sample tubes for use on the assay have not been tested.

Materials Provided

Homocysteine Reagent R1 and R2, Calibrators

Materials Required But Not Provided

1. Yumizen C560 Analyzer.
2. Yumizen C560 Operation manual.
3. Pipettes



Signal Word: Warning



Signal Word: Warning

Limitations

- The linear range of the liquid stable 2-Part Homocysteine Reagent when run as directed is 1-46 µmol/L. Specimens > 46 µmol/L should be diluted 1 part specimen to 2 parts Cal 0 µmol/L or 1 part specimen to 9 parts Cal 0 µmol/L as appropriate.
- The Reagents should be clear. Discard if turbid.
- Cystathionine is measured with homocysteine, but in the general population the cystathionine level (0.065 to 0.3 µmol/L) has a negligible effect. In very rare cases, end stage renal disease and patients with severe metabolic disturbances, cystathionine levels may rise dramatically and in severe cases cause greater than 20% interference.^{25,26}
- Hydroxylamine, present in several iron reagents may carryover (reagent probe or reaction cuvette) and cause falsely low results. Routine rinsing procedures are not adequate to eliminate this problem in most cases. Possible solutions would include special washing protocols, changing to an iron assay that used ascorbic acid as reductant or running iron and homocysteine assays on separate instruments.
- Carbamazepine, methotrexate, phenytoin, nitrous oxide, or 6-azauridine triacetate may affect the homocysteine concentration.¹⁷
- Samples with raised protein levels show > 10% difference compared to results obtained from normal samples and should be avoided.
- Note: Specimens from patients who are on drug therapy involving S-adenosyl-methionine may show falsely elevated levels of homocysteine. Patients who are taking methotrexate, carbamazepine, phenytoin, nitrous oxide, anticonvulsants, or 6-azauridine triacetate, may have elevated levels of homocysteine due to their effect on the pathway.
- Specimens containing particulate matter (fibrin, red blood cells, or other matter) and visibly lipemic specimens should not be used with the assay. Results from these specimens may be inaccurate.

Quality Control

Ensure that adequate maintenance and calibration is performed according to the manufacturer's instructions.

Assayed control materials with values for homocysteine in both the normal and abnormal ranges should be tested to validate reagent performance. Users should ensure that they are fully acquainted with the instructions for the assay, particularly the Precautions and the Reagent Storage sections. Users should demonstrate that they obtain performance specifications for precision and reportable range of test results comparable to those established by the manufacturer before reporting patient test results.

A Pointe Homocysteine Control Kit (H7575-CTL) containing low, medium and high controls is also available from HORIBA Medical for use with the liquid stable 2-Part Homocysteine Reagent.

Expected Values

Reference Range: The reference range should be determined by each laboratory to confirm the characteristics of the population being tested. As a point of reference the following data may be used until the laboratory has analyzed a sufficient number of specimens to determine its own reference range. The HCY concentration in plasma or serum of healthy individuals varies with age, gender, geographical area and genetic factors. Scientific literature reports reference values for adult male and females between 5 and 15 µmol/L, men having higher values than women, and post menopausal woman having higher homocysteine values than pre-menopausal women.^{17,20,21} HCY values will normally increase with age, giving a reference range among an elderly population (> 60 years) of 5-20 µmol/L.²² In countries with folic acid fortification programs, reduced levels of HCY may be observed.^{23,24}

Performance Characteristics

- Assay Range: 1-46 µmol/L.
- Correlation: A study was performed between the Yumizen C560 and a similar analyzer using this method, resulting in the following:

Method	Homocysteine
N	86
Mean Homocysteine (µmol/L)	11.88
Range (µmol/L)	1.5-43.8
Standard Deviation	7.31
Regression Analysis	$y = 0.967x + 0.50$
Correlation Coefficient	0.9969

- Precision: Precision studies were performed following a modification of the guidelines contained in the NCCLS document EP5-T2.²⁹

Sample	Within Day			Total		
	LOW	MID	HIGH	LOW	MID	HIGH
N	20	20	20	40	40	40
Mean	7.18	13.07	24.90	5.88	11.05	23.08
Standard Deviation	0.28	0.44	0.34	0.29	0.62	1.18
Coefficient of Variation (%)	4.0%	3.3%	1.4%	4.9%	5.6%	5.1%

- Sensitivity: 2SD limit of detection (95% Conf) = 0.2 µmol/L

Established on the AU400®:

- Dilution Linearity: The dilution linearity of the liquid stable 2-Part Homocysteine Reagent gives a % recovery range of 91–104% for all samples across the range of the assay (1–46 µmol/L) on the OLYMPUS AU400.
Samples > 46 µmol/L exhibit mean recovery of 100% ± 11% of the expected result when diluted into the assay range.

Pointe Homocysteine Reagent

6. Analytical Specificity: The specificity of the liquid stable 2-Part Homocysteine Reagent was assessed according to guidance in the CLSI Document EP7-A2³¹ for the interfering substances listed in the following table:

Interfering Substance	Interfering Substance Concentration	% Interference
Bilirubin	20 mg/dL	≤ ±10
Hemoglobin	500 mg/dL	≤ ±10
Red Blood Cell	0.4%	≤ ±10
Triglyceride (Intralipid solution)	500 mg/dL	≤ ±10
Glutathione	1000 µmol/L	≤ ±10
Methionine	800 µmol/L	≤ ±10
Cysteine	200 µmol/L	≤ ±10
Pyruvate	1250 µmol/L	≤ ±10

None of these substances interfered significantly in the assay.

Samples with raised protein levels show > 10% difference compared to results obtained from normal samples and should be avoided.

Refer to References section of this product insert (ref 17) for possible interferences caused by drugs, disease or preanalytical variables.

References

- McCully KS. Vascular Pathology of Homocysteinemia: Implications for the Pathogenesis of Arteriosclerosis. *Am J Pathol* 1969;56:111-122
- Malinow MR. Plasma Homocyst(e)ine and Arterial Occlusive Diseases: A Mini-Review. *Clin Chem* 1995;41:173-176
- Ueland PM. Homocysteine Species as Components of Plasma Redox Thiol Status. *Clin Chem* 1995;41:340-342
- Perry IJ, Refsum H, Morris RW, et al. Prospective Study of Serum Total Homocysteine Concentration and Risk of Stroke in Middle-aged British Men. *The Lancet* 1995;346:1395-1398
- Finkelstein JD. Methionine Metabolism in Mammals. *J Nutr Biochem* 1990;1:228-237
- Mudd SH, Levy HL, Skovby F. Disorders of Transsulfuration. In: Scriver CR, Beaudet AL, Sly WS, et al., eds *The Metabolic and Molecular Basis of Inherited Disease*. New York: McGraw-Hill, 1995;1279-1327
- Clarke R, Daly L, Robinson K, et al. Hyperhomocysteinemia: An Independent Risk Factor for Vascular Disease. *N Engl J Med* 1991;324:1149-1155
- Deloughery TG, Evans A, Sadeghi A, et al. Common Mutation in Methylenetetrahydrofolate Reductase. *Circulation* 1996;94:3074-3078
- Schmitz C, Lindpaintner K, Verhoef P, et al. Genetic Polymorphism of Methylenetetrahydrofolate Reductase and Myocardial Infarction. *Circulation* 1996;94:1812-1814
- Boushey CJ, Beresford SAA, Omenn GS, et al. A Quantitative Assessment of Plasma Homocysteine as a Risk Factor for Vascular Disease. *JAMA* 1995;274:1049-1057
- Cattaneo M. Hyperhomocysteinemia, Artherosclerosis and Thrombosis. *Thromb Haemost* 1999;81:165-176
- Ridker PM, Manson JE, Buring JE, et al. Homocysteine and Risk of Cardiovascular Disease Among Postmenopausal Women. *JAMA* 1999;281:1817-1821
- Bostom AG, Silbershatz H, Rosenberg IH, et al. Nonfasting Plasma Total Homocysteine Levels and All-Cause and Cardiovascular Disease Mortality in Elderly Framingham Men and Women. *Arch Intern Med* 1999;159:1077-1080
- Guttormsen AB, Svarstad E, Ueland PM, et al. Elimination of Homocysteine from Plasma in Subjects with Endstage Renal Failure. *Irish J Med Sci* 1995;164 (Suppl. 15):8-9
- Bostom AG, Lathrop L. Hyperhomocysteinemia in End-stage Renal Disease: Prevalence, Etiology, and Potential Relationship to Arteriosclerotic Outcomes. *Kidney Int* 1997;52:10-20
- Rosenquist TH, Ratashak SA, Selhub J. Homocysteine Induces Congenital Defects of the Heart and Neural Tube: Effect of Folic Acid. *Proc Natl Acad Sci USA* 1996;93:15227-15232
- Ueland PM, Refsum H, Stabler SP, et al. Total Homocysteine in Plasma or Serum: Methods and Clinical Applications. *Clin Chem* 1993;39:1764-1779
- Ueland PM, Refsum H. Plasma Homocysteine, A Risk Factor for Vascular Disease: Plasma Levels in Health, Disease, and Drug Therapy. *J Lab Clin Med* 1989;114:473-501
- Fiskerstrand T, Refsum H, Kvalheim G, et al. Homocysteine and Other Thiols in Plasma and Urine: Automated Determination and Sample Stability. *Clin Chem* 1993;39:263-271
- Nehler MR, Taylor LM Jr, Porter JM. Homocysteinemia as a Risk Factor for Atherosclerosis: A Review. *Cardiovascular Pathol* 1997;6:1-9
- Lussier-Cacan S, Xhignesse M, Piolot A, et al. Plasma Total Homocysteine in Healthy Subjects: Sex-Specific Relation with Biological Traits. *Am J Clin Nutr* 1996;64:587-593
- Clarke R, Woodhouse P, Ulvik A, et al. Variability and Determinants of Total Homocysteine Concentrations in Plasma in an Elderly Population. *Clin Chem* 1998;44:102-107
- Jacques PF, Selhub J, Bostom AG, et al. The Effect of Folic Acid Fortification on Plasma Folate and Total Homocysteine Concentrations. *N Engl J Med* 1999;340:1449-1454
- Lawrence JM, Petitti DB, Watkins M and Umekubo MA. Trends in Serum Folate after Food Fortification. *The Lancet* 1999;354:915-916
- Herrmann W, Schorr H, Obeid R, et al. Disturbed Homocysteine and Methionine Cycle Intermediates S-adenosylhomocysteine and S-adenosylmethionine are Related to Degree of Renal Insufficiency in Type 2 Diabetes. *Clin Chem* 2005;51:1-7
- Obeid R, Kuhlmann MK, Kohler H, et al. Response of Homocysteine, Cystathionine, and Methylmalonic Acid to Vitamin Treatment in Dialysis Patients. *Clin Chem* 2005;51:196-201
- Refsum H, Smith AD, Ueland PM, et al. Facts and recommendations about total homocysteine determinations: an expert opinion. *Clin Chem* 2004;50(1):3-32
- National Committee for Clinical Laboratory Standards. *Method Comparison and Bias Estimation using Patient Samples; Approved Guideline-Second Edition*. NCCLS document EP9-A2. Wayne, PA: NCCLS, 2002
- National Committee for Clinical Laboratory Standards. *Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline—Second Edition*. NCCLS Document EP5-A2, Wayne, PA: NCCLS, 2004
- National Committee for Clinical Laboratory Standards. *Protocols for the Determination of Limits of Detection and Limits of Quantitation; Approved Guideline*. NCCLS Document EP17-A. Wayne, PA: NCCLS, 2004.
- Clinical Laboratory Standards Institute. *Interference Testing in Clinical Chemistry; Approved Guideline—Second Edition*. CLSI Document EP7-A2. Wayne, PA: CLSI, 2005.

CHEMISTRY PARAMETERS

Chem:	HCY			No.:	220	Sample Type:	Serum
Chemistry:	Homocysteine					Print Name:	HCY
Reaction Type:	End Point					Reaction Direction:	Negative
Pri Wave:	340					Sec Wave:	380
Unit:	μmol/L					Decimal:	0.1
Blank Time:	47	49				Reaction Time:	80 82
	Sample Vol.	Aspirated	Diluent			Reagent Vol.	Diluent
Standard:	7.7 ul	--- ul	--- ul			R1: 120 ul	--- ul
Decreased:	--- ul	--- ul	--- ul			R2: 12 ul	-- ul
Increased:	--- ul	--- ul	--- ul			R3: --- ul	-- ul
	<input type="checkbox"/> Sample Blank	<input checked="" type="checkbox"/> Auto Rerun				R4: --- ul	--- ul
<u>Slope/Offset Adjustment</u>							
	Slope: 1		Offset: 0				

Linearity Range (Standard)	1	46	Linearity Limit:
Linearity Range (Decreased)	---	---	Substrate Depletion:
Linearity Range (Increased)	---	---	Mixed Blank Abs:
R1 Blank Abs:	---	---	Uncapping Time
Blank Response:	---	---	Reagent Alarm Limit:
Twin Chemistry:			<input type="checkbox"/> Enzyme Linear Extension
<input type="checkbox"/> Prozone Check		<input type="radio"/> Rate Check	<input type="radio"/> Antigen Addition
Q1:	Q2:	Q3:	Q4:
PC:	ABS:		

Pointe Homocysteine Reagent

CALIBRATION PARAMETERS

Calibrator Definition						
Calibrator:	*	Lot No.:	*			
Exp Date:	*					
Carousel		Pos				
Sample Carousel 1	*					
Sample Carousel 2						
Sample Carousel 3						
Reagent/Calibration						
<u>Calibrator</u>	<u>Pos</u>	<u>Lot No</u>	<u>Exp Date</u>	<u>Chem</u>	<u>Conc</u>	<u>Unit</u>
Homocysteine Cal 1	*	*	*	HCY	*	µmol/L
Homocysteine Cal 2	*	*	*	HCY	*	µmol/L
Calibration Setup						
Chem:	HCY					
<u>Calibration Settings</u>						
Math Model:	Two-Point Linear					
Factor:	Replicates:		2			
<u>Acceptance Limits</u>						
Cal Time:	336	Hour				
Slope Diff:	---	SD:	---			
Sensitivity :	---	Repeatability:	---			
Deter Coeff:	---					
<u>Auto Calib.</u>						
<input type="checkbox"/> Bottle Changed	<input type="checkbox"/> Lot Changed	<input type="checkbox"/> Cal Time				

It is recommended that two levels of control material be assayed daily.
* Indicates user defined parameter.

REF 14-H7575-144



Manufactured for HORIBA
Instruments Incorporated – Pointe Brand
5449 Research Drive Canton, MI 48188



Certified to Perform Reagents

The Pointe reagents are certified to be manufactured according to specified parameters. Any Pointe reagent product not meeting specifications through its listed expiration date will be remedied immediately without charge.

Manufactured by HORIBA Instruments Incorporated – Pointe Brand
5449 Research Drive, Canton, MI 48188

European Authorized Representative:

Obelis s.a.

Boulevard Général Wahis 53

1030 Brussels, BELGIUM

Tel: (32)2.732.59.54 Fax:(32)2.732.60.03 email: mail@obelis.net



Symbol Key



Use by (YYYY-MM-DD)



Lot and batch code



Catalog number



Manufacturer



Temperature limitation



Consult instructions for use



In vitro diagnostic medical device **Rx Only:** Prescription Use Only