

Intended Use

For the quantitative determination of Inorganic Phosphorus in serum using the Yumizen C560 analyzer. **Rx Only.**

Method History

The measurement of inorganic phosphorus in serum is usually accomplished by forming a phosphomolybdate complex and in turn reducing it to a molybdenum blue color complex. Methods differ as to the choice of reducing agents: stannous chloride¹, phenylhydrazine², aminonaphtholsulfonic acid³, ascorbic acid⁴, p-methylaminophenolsulfate⁵, N-phenyl-p-phenylenediamine⁶ and ferrous sulfate.⁷ These methods suffered from color instability, deproteinization steps and complexity of performance⁸. The addition of a surfactant eliminated the need to prepare a protein-free filtrate, accelerated color production, stabilized the color and simplified the procedure. Many of the components in these reagents were unstable and had to be stored separately. The quantitative measurement of unreduced phosphomolybdate complexes was first reported by Simonsen in 1946.⁹ Daly and Ertingshausen¹⁰ adapted that technique for the determination of inorganic phosphorus in 1972. Amador and Urban¹¹ modified this procedure further the same year. The present method is a modification of the above procedure using a single, stable reagent performing in the UV range.

Principle

Inorganic Phosphorus + H₂SO₄ + Ammonium Molybdate \longrightarrow Unreduced Phosphomolybdate Complex

Inorganic phosphorus reacts with ammonium molybdate in an acid medium to form a phosphomolybdate complex that absorbs light at 340nm. The absorbance at this wavelength is directly proportional to the amount of inorganic phosphorus present in the sample.

Reagents

Ammonium Molybdate 0.48 mM, Sulfuric Acid 220 mM with surfactant

Precautions and Hazards

This reagent is for *in vitro* diagnostic use only.

Hazards:

Hazard Classifications: Skin Corrosion/Irritation (Category 1), Serious eye damage/eye irritation (Category 1), Specific target organ toxicity, single exposure; Respiratory System (Category 1), Specific target organ toxicity, repeat exposure; Respiratory System (Category 1)

Hazard Statements: H314: Causes severe skin burns and eye damage, H319: Causes serious eye damage, H370: Causes damage to organs, H372: Causes damage to organs through prolonged or repeated exposure

Precautionary Statements

Prevention : P260: Do not breathe dust/fume/gas/mist/vapors/spray, P264: Wash hands thoroughly after handling, P270: Do not eat, drink or smoke when using this product, P280: Wear protective gloves/protective clothing/eye protection/face protection

Response : P314: Get medical advice/attention if you feel unwell, P363: Wash contaminated clothing before reuse, P301 + P330 + P331: IF SWALLOWED: Rinse mouth. Do NOT induce vomiting, P303 + P361 + P353: IF ON SKIN (or hair): remove/Take off immediately all contaminated clothing. Rinse SKIN with water/shower, P304 + P340: IF INHALED: Remove victim to fresh air and Keep at rest in a position comfortable for breathing, P305 + P351 + P338: IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing, P307 + P311: IF exposed: call a POISON CENTER or doctor/physician.

Storage : P404 : Store in a closed container.

Disposal : P501 : Dispose of contents into sewer system after diluting with large volumes of water, if in accordance with local regulations. **Refer to the Safety Data Sheet for this product (SDS-P7516) available by calling 1-734-487-8300.**



Signal Word: Danger

Reagent Preparation

Reagent comes in a ready to use form.

Reagent Storage and Stability

Store reagent at refrigerator temperature (2-8°C). The reagent is stable until the expiration date appearing on the label when stored as directed. Manufacturer studies have shown reagent is stable for 30 days once placed in the refrigerated reagent carousel (2-10°C), however reagent stability may vary based on individual laboratory conditions.

Reagent Deterioration

Do not use reagent if:

1. Reagent read against water has an absorbance greater than 0.500 at 340 nm.
2. The reagent fails to recover stated control values.

Specimen Collection and Storage

1. Unhemolyzed serum is specimen of choice.
2. Plasma should not be used since anticoagulants may produce falsely low values.¹²
3. Hemolyzed sample may give falsely high values.
4. Serum should be removed from the red cell clot as soon as possible.¹³
5. Serum inorganic phosphorus is stable for one week refrigerated and for three weeks frozen.^{13,14}

Interferences

For a comprehensive list of substances that interfere with the measurement of Inorganic Phosphorus see Young, et al.¹⁵

Pointe Inorganic Phosphorus Reagent Set (UV)

Materials Provided

Inorganic Phosphorus Reagent

Materials Required but not Provided

1. Yumizen C560 Analyzer
2. Yumizen C560 operation manual
3. Chemistry Calibrator, catalog number C7506-50
4. Chemistry control, catalog number C7592-100

Calibration

Use Pointe Chemistry Calibrator (Catalog Number C7506-50). The procedure should be calibrated according to the instrument manufacturer's calibration instructions. If control results are found to be out of range, the test may need to be re-calibrated. Under typical operating conditions manufacturer calibration stability studies have shown the calibration curve will be stable for at least 14 days.

Quality Control

The integrity of the reaction should be monitored by use of normal and abnormal control sera with known concentrations of inorganic phosphorus. Quality control requirements should be performed in conformance with local, state, and/or Federal regulations or accreditation requirements.

SI Units

To obtain results in SI Units (mmol/L), multiply the results in mg/dl by the factor 0.323.

Example: $3.4 \text{ mg/dl} \times 0.323 = 1.09 \text{ mmol/L}$.

Limitations

Detergents containing phosphate should not be used for cleaning glassware used in this procedure. Lipemic and icteric samples require a serum blank. For maximum accuracy a serum blank should be run with each sample.

Expected Values

Adults: 2.5-4.8mg/dl¹⁶
Children: 4.0-7.0mg/dl¹⁷

Values are decreased during menstrual period and after meals.¹⁷ It is strongly recommended that each laboratory establish its own normal values.

Performance

1. Assay Range: 0.0-12.0 mg/dL.
2. Correlation: A study was performed between the Yumizen C560 and a similar analyzer using this method, resulting in the following:

Method	Phosphorus
N	96
Mean Phosphorus (mg/dL)	4.30
Range (mg/dL)	0.5-9.9
Standard Deviation	1.83
Regression Analysis	$y = 0.936x + 0.25$
Correlation Coefficient	0.9724

3. Precision: Precision studies were performed following a modification of the guidelines contained in the NCCLS document EP5-T2.¹⁸

Sample	Within Day			Total		
	LOW	MID	HIGH	LOW	MID	HIGH
N	20	20	20	40	40	40
Mean	3.78	7.31	10.89	3.81	7.38	11.06
Standard Deviation	0.04	0.04	0.04	0.13	0.31	0.31
Coefficient of Variation (%)	1.1%	0.6%	0.4%	3.5%	4.2%	2.8%

4. Sensitivity: 2SD limit of detection (95% Conf) = 0.0 mg/dL

References

1. Osmond, M.F., Bull. Soc. Chim. 47:745 (1887).
2. Taylor, A.E., Miller, C.W., J. Biol., Chem 18:215 (1914).
3. Fiske, C.H., Subbarow, Y., J. Biol. Chem. 66:275 (1925).
4. Lowry, O.H., Lopez, J.A., J. Biol. Chem. 162:421 (1946).
5. Power, M.H., Standard Methods of Clinical Chemistry New York, Academic Press, (1953).
6. Dryer, R.L., et al, J. Biol. Chem. 225:177 (1957).
7. Tausky, H.H., Shorr, E., J. Biol. Chem. 202:675 (1953).
8. Martinek, R.G., J. Am. Med. Tech. 32:337 (1970).
9. Simonsen, D.G., et al, J. Biol. Chem. 166:747 (1946).
10. Daly, J.A., Ertingshausen, G., Clin. Chem. 18:263 (1972).
11. Amador, E., Urban, J., Clin. Chem. 18:601 (1972).
12. Goldenberg, H. Fernandez, A. Clin. Chem. 12:871 (1966).
13. Henry, R.J., et al, Clinical Chemistry: Principles and Technics, New York, Harper & Row, pp.122:143 (1964).
14. Hansk, A., Kao, J., Clin. Chem. 14:58 (1968).
15. Young, D.S., et al, Clin. Chem., 21:1D, (1975).

16. Henry, R.J., et al, Clinical Chemistry: Principles and Technics, 2nd Ed., Hagerstown (MD), Harper & Row, p.728 (1974).
17. Tietz, N.W., Fundamentals of Clinical Chemistry, Philadelphia, W.B. Saunders, p.917 (1976).
18. NCCLS document "Evaluation of Precision Performance of Clinical Chemistry Devices", 2nd Ed. (1992).

CHEMISTRY PARAMETERS

Chem:	PHOS	No.:	227	Sample Type:	Serum
Chemistry:	Inorganic Phosphorus			Print Name:	PHOS
Reaction Type:	End Point			Reaction Direction:	Positive
Pri Wave:	340			Sec Wave:	
Unit:	mg/dL			Decimal	0.1
Blank Time:	10	12		Reaction Time:	27 29
	Sample Vol.	Aspirated	Diluent	Reagent Vol.	Diluent
Standard:	1.5 ul	-- ul	-- ul	R1: 150 ul	-- ul
Decreased:	-- ul	-- ul	-- ul	R2: -- ul	-- ul
Increased:	-- ul	-- ul	-- ul	R3: -- ul	-- ul
	<input type="checkbox"/> Sample Blank	<input checked="" type="checkbox"/> Auto Rerun		R4: -- ul	-- ul
 <u>Slope/Offset Adjustment</u>					
Slope: 1		Offset: 0			

Linearity Range (Standard)	0	12	Linearity Limit:
Linearity Range (Decreased)	---	---	Substrate Depletion:
Linearity Range (Increased)	---	---	Mixed Blank Abs:
R1 Blank Abs:	---	---	Uncapping Time
Blank Response:	---	---	Reagent Alarm Limit:
Twin Chemistry:			<input type="checkbox"/> Enzyme Linear Extension
<input type="checkbox"/> Prozone Check		<input type="radio"/> Rate Check	<input type="radio"/> Antigen Addition
Q1:	Q2:	Q3:	Q4:
PC:	ABS:		

Pointe Inorganic Phosphorus Reagent Set (UV)

CALIBRATION PARAMETERS

Calibrator Definition						
Calibrator:	*	Lot No.:			*	
Exp Date:	*					
Carousel		Pos				
Sample Carousel 1	*					
Sample Carousel 2						
Sample Carousel 3						
Reagent/Calibration						
<u>Calibrator</u>	<u>Pos</u>	<u>Lot No</u>	<u>Exp Date</u>	<u>Chem</u>	<u>Conc</u>	<u>Unit</u>
Water	W	*	*	PHOS	0	mg/dL
Chemistry Calibrator	*	*	*	PHOS	*	mg/dL
Calibration Setup						
Chem:	PHOS					
<u>Calibration Settings</u>						
Math Model:	Two-Point Linear					
Factor:		Replicates:	2			
<u>Acceptance Limits</u>						
Cal Time:	336	Hour				
Slope Diff:	---	SD:	---			
Sensitivity :	---	Repeatability:	---			
Deter Coeff:	---					
<u>Auto Calib.</u>						
<input type="checkbox"/> Bottle Changed	<input type="checkbox"/> Lot Changed	<input type="checkbox"/> Cal Time				

It is recommended that two levels of control material be assayed daily.
* Indicates user defined parameter.

REF 14-P7516-174



Manufactured by
HORIBA Instruments Incorporated-Pointe Brand
5449 Research Drive Canton, MI 48188



Certified to Perform Reagents

The Pointe reagents are certified to be manufactured according to specified parameters. Any Pointe reagent product not meeting specifications through its listed expiration date will be remedied immediately without charge.

Manufactured by HORIBA Instruments Incorporated – Pointe Brand
5449 Research Drive, Canton, MI 48188

European Authorized Representative:
Obelis s.a.

Boulevard Général Wahis 53
1030 Brussels, BELGIUM

Tel: (32)2.732.59.54 Fax:(32)2.732.60.03 email: mail@obelis.net



Symbol Key

Use by (YYYY-MM-DD)	LOT Lot and batch code	REF Catalog number
Manufacturer	Temperature limitation	Consult instructions for use
IVD <i>In vitro</i> diagnostic medical device Rx Only: Prescription Use Only		