

## Intended Use

For the quantitative determination of total bilirubin in serum using the Yumizen C560 analyzer. **Rx Only.**

## Method History


Since the introduction of the diazo method for bilirubin determination by Ehrlich in 1883<sup>1</sup>, several modifications have been proposed to enhance the reaction. The Malloy and Evelyn method<sup>2</sup> employs methanol to catalyze the azo-coupling reaction of the indirect bilirubin, as well as to keep the azobilirubin in solution. A serious disadvantage of this method lies in the fact that protein may be precipitated by the methanol solution to yield falsely lowered results.

In 1938, Jendrassik and Grof.<sup>3</sup> presented an assay that gave reliable results. The method is, however, cumbersome and involves several pipetting steps.

The method presented here was developed by Wahlefeld et al.<sup>4</sup> A detergent is used to accelerate the reaction and to avoid protein precipitation. The diazo reagent is 2,5-dichlorophenyldiazonium tetrafluoroborate (DPD) that reacts very rapidly in coupling with bilirubin under acidic conditions. The resulting procedure is simple, yet exhibits good correlation when compared with the method of Jendrassik and Grof.

## Principle

Total bilirubin is coupled with a diazonium salt (DPD) in a strongly acid medium (pH 1 – 2).



The intensity of the color of the azobilirubin produced is proportional to the total bilirubin concentration and can be measured photometrically.

## Reagents

Total bilirubin R1 reagent: acid buffer 50 mmol/L, Surfactant. Total bilirubin R2 reagent: acid buffer >30 mmol/L, >2.0 mmol/L DPD and stabilizers.

## Reagent Preparation

Reagents provided as ready to use liquids.

## Reagent Storage and Stability

1. Packaged reagents are stored at 2-8°C. The reagents are stable until the expiration date appearing on the label when stored as directed.
2. Manufacturer studies have shown reagent is stable for 30 days once placed in the refrigerated reagent carousel (2-10°C), however reagent stability may vary based on individual laboratory conditions.
3. Do not freeze reagents.
4. Avoid exposure to direct sunlight.

## Reagent Deterioration

1. Do not use if reagents show evidence of contamination (turbidity)
2. The R2 may develop very slight precipitation that does not affect performance and will re-dissolve if the R2 is warmed gently.
3. R2 reagent containing a precipitate that does not re-dissolve and results in product discoloration should not be used.
4. Do not use if reagent fails to achieve assigned assay values of fresh control sera.

## Precautions

1. Reagents are toxic and corrosive. Do not pipette by mouth. Avoid contact with skin and clothing.
2. This reagent is for *in vitro* diagnostic use only.

## Hazards:

**R1 and R2:** Hazard Classifications: Skin Corrosion/Irritation (Category 1), Serious eye damage/eye irritation (Category 1)

Hazard Statements: H314: Causes severe skin burns and eye damage, H318: Causes serious eye irritation

Precautionary Statements: **Prevention:** P260: Do not breathe dust/fume/gas/mist/vapors/spray. P264: Wash skin thoroughly after handling. P280: Wear protective gloves/protective clothing/eye protection/face protection. **Response:** P310: Immediately call a POISON CENTER or doctor/physician. P363: Wash contaminated clothing before reuse. P301+P330+P331 : If SWALLOWED : Rinse mouth. Do NOT induce vomiting. P303+P361+P353 : IF ON SKIN (or hair) : Remove/Take off immediately all contaminated clothing. Rinse SKIN with water/shower. P304+P340: IF INHALED: Remove victim to fresh air and keep at rest in a position comfortable for breathing. P305+P351+P338: IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do.

Continue rinsing. **Storage:** P404: Store in a closed container. **Disposal:** P501: Dispose of contents into sewer system after diluting with large volumes of water, if in accordance with local regulations. **Refer to the Safety Data Sheet for this product (SDS-HB979) available by calling 1-734-487-8300.**



**Signal Word: Danger**

## Specimen Collection and Storage

1. Fresh, unhemolyzed serum is recommended.
2. Samples should be analyzed within two hours of collection if kept at room temperature in the dark and within twelve hours if kept refrigerated (2-8°C) and protected from light.<sup>5</sup>
3. Bilirubin in serum is stable for three months when stored frozen (-20°C) and protected from light.<sup>5</sup>
4. Direct sunlight may cause up to a 50% decrease in bilirubin within one hour.<sup>6</sup>
5. Specimen collection should be carried out in accordance with NCCLS M29-T2. No method can offer complete assurance that human blood samples will not transmit infection. Therefore, all blood samples should be considered potentially infectious.

## Interferences

1. All interference studies were performed according to the procedures recommended in NCCLS guideline No. EP7-P for interference testing in clinical chemistry.<sup>7</sup>
2. Serum hemoglobin has been shown to interfere with results. Serum Triglycerides up to 1000 mg/dl do not interfere with results.
3. A number of drugs and substances affect bilirubin results. See Young, et al.<sup>8</sup>

# Pointe Total Bilirubin Reagent Set

## Materials Provided

Total Bilirubin R1 reagent, Total Bilirubin R2 reagent

## Materials Required but not Provided

1. Yumizen C560 Analyzer
2. Yumizen C560 Operation manual
3. Chemistry Calibrator, catalog number C7506-50
4. Chemistry control, catalog number C7592-100

## Calibration

Use Pointe Chemistry Calibrator (Catalog Number C7506-50). Follow instrument application instructions for calibration. Refer to instrument manual instructions for calibration procedures and frequency. It is recommended that each laboratory determine its own frequency of calibration. If control results are found to be out of range, the test may need to be re-calibrated. Under typical operating conditions manufacturer calibration stability studies have shown the calibration curve will be stable for at least 7 days.

## Quality Control

The validity of the reaction should be monitored by use of the control sera with known normal and abnormal total bilirubin values. These controls should be run at least with every working shift in which total bilirubin assays are performed. It is recommended that each laboratory establish its own frequency of control determination. Quality control requirements should be performed in conformance with local, state, and/or Federal regulations or accreditation requirements.

## Limitations

1. Samples with values above 30 mg/dl must be diluted 1:1 with isotonic saline, re-assayed and the final answer multiplied by two.
2. Serum triglycerides up to 1000 mg/dl do not interfere with results.
3. Serum hemoglobin levels have been shown to interfere with results.

## Performance

1. Assay Range: 0.0-30.0 mg/dl
2. Comparison: A study was performed between the Yumizen C560 and a similar analyzer using this method, resulting in the following:

Method	Total Bilirubin
N	149
Mean T Bilirubin (mg/dL)	4.38
Range (mg/dL)	0.0-29.7
Standard Deviation	7.18
Regression Analysis	$y = 0.988x + 0.12$
Correlation Coefficient	0.9940

3. Precision: Precision studies were performed using the Yumizen C560 analyzer following a modification of the guidelines which are contained in NCCLS document EP5-T2.<sup>9</sup>

Sample	Within Day			Total		
	LOW	MID	HIGH	LOW	MID	HIGH
N	20	20	20	40	40	40
Mean	0.70	4.33	21.38	0.87	4.34	22.04
Standard Deviation	0.00	0.09	0.20	0.07	0.16	0.71
Coefficient of Variation (%)	0.0%	2.1%	0.9%	8.3%	3.8%	3.2%

4. Sensitivity: 2SD Limit of Detection (95% Conf) = 0.0 mg/dL

## Expected Values<sup>10</sup>

Total: Adults and infants older than 1 month: 0.2 –1.0 mg/dl

Infants: Full Term Newborn

Up to 24hrs: 2.0-6.0 mg/dl

Up to 48hrs: 6.0-10.0 mg/dl

Days 3-5: 4.0-8.0 mg/dl

## References

1. Ehrlich, P., Charite Ann. 8:140 (1883).
2. Malloy, H.T., Evelyn, K.A., J. Biol. Chem. 119:481 (1937).
3. Jendrassik, L., Grof, P., Biochem. Zeitschr. 297:81 (1938).
4. Wahlefeld AW, et al. Scand J Clin Lab Invest. 29 Supplement 126(1972).
5. Martinek, R.G., Clin. Chim. Acta 13:161 (1966).
6. Tietz, N.W., Fundamentals of Clinical Chemistry, Philadelphia, W.B. Saunders, p.1028 (1976).
7. NCCLS document, "National Evaluation Protocols for Interference Testing", Evaluation Protocol Number 7, Vol. 4, No. 8, (June 1984).
8. Young, D.S., Effects of Preanalytical Variables on Clinical Laboratory Tests, Washington DC, AACC Press, (1997)
9. NCCLS document "Evaluation of Precision Performance of Clinical Chemistry Devices", 2<sup>nd</sup> Ed. (1992).
10. Tietz, Textbook of Clinical Chemistry, Philadelphia, W.B. Saunders, 3<sup>rd</sup> Ed., p. 1170 (1999)

### CHEMISTRY PARAMETERS

Chem:	TBIL	No.:	207	Sample Type:	Serum
Chemistry:	Total Bilirubin	Print Name:	TBIL	Reaction Direction:	Positive
Reaction Type:	End Point	Sec Wave:	605	Decimal	0.1
Pri Wave:	546	Reaction Time:	80	82	
Unit:	mg/dL	Reagent Vol.	Diluent		
Blank Time:	47	49			
	Sample Vol.	Aspirated	Diluent		
Standard:	2.0 ul	--- ul	--- ul	R1:	120 ul --- ul
Decreased:	--- ul	--- ul	--- ul	R2:	31 ul -- ul
Increased:	--- ul	--- ul	--- ul	R3:	--- ul -- ul
	<input type="checkbox"/> Sample Blank	<input checked="" type="checkbox"/> Auto Rerun		R4:	--- ul --- ul
<b>Slope/Offset Adjustment</b>					
Slope: 1		Offset: 0			

Linearity Range (Standard)	0	30	Linearity Limit:
Linearity Range (Decreased)	---	---	Substrate Depletion:
Linearity Range (Increased)	---	---	Mixed Blank Abs:
R1 Blank Abs:	---	---	Uncapping Time
Blank Response:	---	---	Reagent Alarm Limit:
Twin Chemistry:			<input type="checkbox"/> Enzyme Linear Extension
<input type="checkbox"/> Prozone Check		<input type="radio"/> Rate Check	<input type="radio"/> Antigen Addition
Q1:	Q2:	Q3:	Q4:
PC:	ABS:		

# Pointe Total Bilirubin Reagent Set

## CALIBRATION PARAMETERS

<b>Calibrator Definition</b>						
Calibrator:	*	Lot No.:			*	
Exp Date:	*					
<b>Carousel</b>	<b>Pos</b>					
Sample Carousel 1	168					
Sample Carousel 2						
Sample Carousel 3						
<b>Reagent/Calibration</b>						
<u>Calibrator</u>	<u>Pos</u>	<u>Lot No</u>	<u>Exp Date</u>	<u>Chem</u>	<u>Conc</u>	<u>Unit</u>
Water	W	*	*	TBIL	0	mg/dL
Chemistry Calibrator	*	*	*	TBIL	*	mg/dL
<b>Calibration Setup</b>						
Chem:	TBIL					
<u>Calibration Settings</u>						
Math Model:	Two-Point Linear					
Factor:		Replicates:	2			
<u>Acceptance Limits</u>						
Cal Time:	168	Hour				
Slope Diff:	---	SD:	---			
Sensitivity :	---	Repeatability:	---			
Deter Coeff:	---					
<u>Auto Calib.</u>						
<input type="checkbox"/> Bottle Changed	<input type="checkbox"/> Lot Changed	<input type="checkbox"/> Cal Time				

It is recommended that two levels of control material be assayed daily. \* Indicates user defined parameter.

**REF** 14-HB979-365



Manufactured by  
HORIBA Instruments Incorporated-Pointe Brand  
5449 Research Drive Canton, MI 48188



### Certified to Perform Reagents

The Pointe reagents are certified to be manufactured according to specified parameters. Any Pointe reagent product not meeting specifications through its listed expiration date will be remedied immediately without charge.

Manufactured by HORIBA Instruments Incorporated – Pointe Brand  
5449 Research Drive, Canton, MI 48188

European Authorized Representative:  
Obelis s.a.

Boulevard Général Wahis 53  
1030 Brussels, BELGIUM

Tel: (32)2.732.59.54 Fax:(32)2.732.60.03 email: mail@obelis.net



### Symbol Key



Use by (YYYY-MM-DD)



Lot and batch code



Catalog number



Manufacturer



Temperature limitation



Consult instructions for use



In vitro diagnostic medical device **Rx Only:** Prescription Use Only