

Procedure for External Quality Control tubes

Important

This document is applicable to the following versions of analysers:



- Yumizen H500 since v2.2
- Yumizen H550 since v3
- HELO 1 (Yumizen H1500 / H2500 v1.4)
- HELO 2 (Yumizen H1500 / H2500 v1.7)



To run an EQC tube on any other versions, please refer to your HORIBA customer representative.

I. To run an EQC tube on Yumizen H500/H550

Follow the instructions of the user manual section « To Run an External Control Blood Sample ».

The **Analysis** and the **Analysis mode** of your EQC tube can be find in section Annex I.

II. To run an EQC tube in rack mode on Yumizen H1500/H2500

If not already done, create a « predefined test » for your EQC tube (see section V)

Follow the instructions of the user manual section « To Run an External Control Blood Sample ».

If you have a T6000 track and if your rack has not been automatically routed, repeat the above instructions on all analysers.

If for any reason you are unable to pass your EQC in one or several analysers in rack mode, then follow the instructions of section III for these devices.

III. To run an EQC tube in STAT mode (only Yumizen H1500/H2500)

Follow the instructions of the user manual section « To Run a Control Blood Sample in Stat Mode ».

The specific **Test**, **Mode** and **Sample Type** for your EQC tube can be find in Annex II.

If you don't find on screen the required **Sample Type**, follow instructions of section VI.

IV. To recover the results

Follow the instructions of the user manual section « To Export Results ».

V. To create a predefined test for you EQC tube

Follow the instructions of the user manual section « To Create an External Control Blood Sample ».

The specific **Test**, **Mode** and **Sample Type** for your EQC tube can be find in Annex II.

If you don't find on screen the required **Sample Type**, follow instructions of section VI.

VI. To create a custom Analysis Profile (Sample Type)

Retrieve the information related to your EQC program in Annex III :

- The **Name** of the new profile to create.
- The initial **Sample Type** from which to start.
- The list of **Modifications** to be applied.

Follow the instructions of the user manual section « To Create an Analysis Profile », applying the elements above recovered.

Annex I. EQC settings for Yumizen H500/H550

Table 1 - UK NEQAS

EQC tube	Analysis	Analysis Mode
FB	CBC	EQC Control
DL	DIF	EQC Control

Table 2 – CAP

EQC tube	Analysis	Analysis Mode
FH16	DIF	EQC Control

Table 3 - EQAS

EQC tube	Analysis	Analysis Mode	
		Yumizen H500 v2.2	Yumizen H500 / H550 since v3
BC90	CBC	EQC Patient	EQC Patient mode 1

Table 4 - RIQAS

EQC tube	Analysis	Analysis Mode	
		Yumizen H500 v2.2	Yumizen H500 / H550 since v3
RQ9118	CBC	EQC Patient	EQC Patient mode 1

Table 5 - Biologie-Propective

EQC tube	Analysis	Analysis Mode	
		Yumizen H500 v2.2	Yumizen H500 / H550 since v3
NUMSANG	DIR + PLTo	EQC Patient	EQC Patient mode 1

Table 6 - ProBioQual

EQC tube	Analysis	Analysis Mode
HD	CBC	EQC Control
HDF	DIF	EQC Control

Annex II. EQC settings for Yumizen H1500/H2500

Table 7 - UK NEQAS

EQC tube	Test	Mode	Sample Type (also called Analysis Profile)
FB	CBC	CONTROL	UK NEQAS ¹
DL	DIF	CONTROL	UK NEQAS ¹
RE	RET	BLOOD	STANDARD

Table 8 - CAP

EQC tube	Test	Mode	Sample Type (also called Analysis Profile)
FH16	DIF	CONTROL	CAP ²
RTQ	RET	CONTROL	DIFFTROL

Table 9 - EQAS

EQC tube	Test	Mode	Sample Type (also called Analysis Profile)
BC90	CBC	BLOOD	STANDARD

Table 10 - RIQAS

EQC tube	Test	Mode	Sample Type (also called Analysis Profile)
RQ9118	CBC	BLOOD	STANDARD

Table 11 - INSTAND

EQC tube	Test	Mode	Sample Type (also called Analysis Profile)
211	CBC	CONTROL	DIFFTROL
216	RET	CONTROL	MINOTROL RETIC

Table 12 - Biologie-Prospective

EQC tube	Test	Mode	Sample Type (also called Analysis Profile)
NUMSANG	DIR, PLTo	BLOOD	STANDARD

¹ See Table 14 Annex III

² See Table 15 Annex III

Table 13 - ProBioQual

EQC tube	Test	Mode	Sample Type (also called Analysis Profile)
HD	CBC, PLTo	CONTROL	DIFFTROL
HDF	DIF, PLTo	CONTROL	DIFFTROL
HJ	RET	BLOOD	STANDARD

Annex III. Specific analysis profile for Yumizen H1500/H2500

Table 14 – UK NEQAS

Name	Initial Sample Type (also called Analysis Profile)	Modifications	
		HELO 1 (Yumizen H1500/H2500 v1.4)	HELO 2 (Yumizen H1500/H2500 v1.7)
UK NEQAS	DIFFTROL	<ul style="list-style-type: none"> LMNE matrix -LGauto: False LMNE matrix -rBFL: 10* LMNE matrix -rLG: 10 	<ul style="list-style-type: none"> LMNE matrix -LGauto: False LMNE matrix -rBFL: 10 LMNE matrix -rLGI: 10 LMNE matrix -rLGh: 10 LMNE matrix -rLeftGaussLgMax: 35 TNC Distribution -rTNC1 : 4 nrbcLympProp: 0

Table 15 – CAP

Name	Initial Sample Type (also called Analysis Profile)	Modifications	
		HELO 1 (Yumizen H1500/H2500 v1.4)	HELO 2 (Yumizen H1500/H2500 v1.7)
CAP	DIFFTROL	<ul style="list-style-type: none"> BASO Distribution - calibpBaso: 0.01 WBC Low count - WbvSynthesis.basoLowWbcRatio: 0.1 WBC Noise - BASO.maxEventsBG: 12 000 RBC Abnormal indices - RbcPltSynthesis.MchcMax: 500 LMNE Matrix - rLeftGaussLgMax: 35 LMNE matrix - rLG : 15 LMNE matrix - rBFL : 15 LMNE matrix - LGauto : Desactivated 	TBD