

REF A11A01754

REAGENT 29 mL

IVD **CE** Rx Only

HORIBA ABX SAS
Parc Euromédecine
Rue du Caducée
BP 7290
34184 Montpellier Cedex 4
FRANCE



ABX Pentra Deproteinizer CP

- Pentra C200
- Pentra C400
- ABX Pentra 400

Cleaning solution used on HORIBA Medical clinical chemistry analyzers.

Intended Use

ABX Pentra Deproteinizer CP is a cleaning solution used on HORIBA Medical clinical chemistry analyzers. Clinical laboratories use.

Reagents

ABX Pentra Deproteinizer CP is ready-to-use. This reagent contains sodium hypochlorite (29% v/v at 15% of active chlorine) as well as sodium hydroxide (< 0.5%).

The pH of the solution is greater than 11.

ABX Pentra Deproteinizer CP should be used according to this notice. The manufacturer cannot guarantee its performance if used otherwise.

Handling

Remove the cap of the cassette, place the cassette in the refrigerated reagent compartment of the instrument.

For further details, please refer to the maintenance procedures outlined in the User Manual for Pentra C200, Pentra C400 or ABX Pentra 400 clinical chemistry analyzer.

Materials required but not provided

- Automated clinical chemistry analyzer.
- Standard laboratory equipment.

Storage and Stability

Stability before opening:

Stable up to the expiry date on the label if stored at 2-8°C.

Stability after opening:

Stable for 1 month if installed on the refrigerated compartment of the instrument.

Waste Management

Please refer to local legal requirements.

General Precautions ^a

- This reagent is for professional *in vitro* diagnostic use only.
For laboratory use.
- For prescription use only.
- This reagent is classified as hazardous in compliance with regulation (EC) N°.1272/2008.

^aModification: general precautions modification.

ABX Pentra Deproteinizer CP

- **Danger:** Due to sodium hypochlorite presence.
 - H290:** May be corrosive to metals.
 - H318:** Causes serious eye damage.
 - H400:** Very toxic to aquatic life.
 - H411:** Toxic to aquatic life with long lasting effects.
 - P280:** Wear eye or face protection.
 - P234:** Keep only in original packaging.
 - P273:** Avoid release to the environment.
 - P305 + P351 + P338:** IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
 - P310:** Immediately call a POISON CENTER or doctor/physician.
 - P390:** Absorb spillage to prevent material damage.
 - P391:** Collect spillage.
- Do not swallow. Avoid contact with skin and mucous membranes.
- Observe the standard laboratory precautions for use.
- The reagent cassettes are disposable and should be disposed of in accordance with the local legal requirements.
- Please refer to the SDS associated with the reagent.
- Do not use the product if there is visible evidence of biological, chemical or physical deterioration.
- Do not use the product if the recommended storage conditions, including temperature, are not followed.
- User must be trained by a HORIBA Medical representative before attempting to operate the device.
- It is the user's responsibility to verify that this document is applicable to the reagent used.
- For technical assistance, you can call +33 (0)4 67 14 15 16.
- Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the country in which the user and/or the patient is established.