

REF 1300141429

REAGENT 1 2 x 25 mL

REAGENT 2 2 x 8 mL

IVD  2797

HORIBA ABX SAS
Parc Euromédecine
Rue du Caducée
BP 7290
34184 Montpellier Cedex 4
FRANCE

Yumizen C Creatinine Jaffé

- Yumizen C230
- Yumizen C240

Diagnostic reagent for quantitative *in vitro* determination of creatinine in serum, plasma and urine by colorimetry.

Intended Use

Yumizen C Creatinine Jaffé reagent is a diagnostic reagent for quantitative *in vitro* determination of Creatinine in human serum, plasma and urine based on a kinetic method using alkaline picrate (Jaffé method). Clinical laboratories use.

Creatinine measurements are used in the diagnosis and treatment of renal diseases, and as a calculation basis for measuring other urine analytes.

Assessing physiologic and pathologic variations of activity of creatinine in human serum, plasma and urine is useful for screening or follow-up of these diseases.

Clinical Interest

Creatinine measurements are used in the diagnosis and treatment of renal diseases and as a calculation basis for measuring other urine analytes.

Method

In 1886, Jaffe developed an assay for creatinine based upon the reaction between creatinine and sodium picrate (1). In 1904, Folin (2) used this reaction for the quantitative determination of creatinine in urine. Kinetic procedures based on the observed reaction rates various substances, including creatinine, with alkaline picrate have been proposed by Fabing (3) and Soldin (4). This improved Jaffe chemistry is a kinetic procedure which does not require deproteinization of the sample and is formulated to reduce the interference in serum proteins.

Creatinine + alkaline picrate \longrightarrow creatinine-picrate complex

At an alkaline pH, creatinine reacts with picrate to form Janovsky complex.

The rate of increase in absorbance at 510 nm due to the formation of creatinine-picrate complex is directly

proportional to the creatinine concentration present in the sample.

Reagents

Yumizen C Creatinine Jaffé is ready-to-use.

Reagent 1 (R1):

Sodium hydroxide 0.25 mol/L
Surfactants

Reagent 2 (R2):

Picric acid 31.4 mmol/L

Yumizen C Creatinine Jaffé should be used according to this notice. The manufacturer cannot guarantee its performance if used otherwise.

Handling

1. Remove the caps of the cassettes.
2. If present, remove foam by using a plastic pipette.
3. Place the cassettes into the refrigerated reagent compartment.

Calibrator

For calibration, use:

ABX Pentra Multical (A11A01652) (not included)
10 x 3 mL (lyophilisate)

Yumizen C Creatinine Jaffé

Control

For internal quality control, use:

- **ABX Pentra N MultiControl** (1300054414) (not included)
10 x 5 mL (lyophilisate)
- **ABX Pentra P MultiControl** (1300054415) (not included)
10 x 5 mL (lyophilisate)
- **Yumizen C Urine Level 1 Control** (1300023946) (not included)
6 x 5 mL
- **Yumizen C Urine Level 2 Control** (1300023947) (not included)
6 x 5 mL

Each control should be assayed daily and/or after a calibration.

The frequency of controls and the confidence intervals should correspond to laboratory guidelines and country-specific directives. You should follow federal, state and local guidelines for testing quality control materials. The results must be within the range of the defined confidence limits. Each laboratory should establish a procedure to follow if the results exceed these confidence limits.

Materials Required but not Provided

- Automated clinical chemistry analyzer: Yumizen C230/C240
- Calibrator: **ABX Pentra Multical** (A11A01652)
- Controls:
 - ABX Pentra N MultiControl** (1300054414)
 - ABX Pentra P MultiControl** (1300054415)
 - Yumizen C Urine Level 1 Control** (1300023946)
 - Yumizen C Urine Level 2 Control** (1300023947)
- Standard laboratory equipment.

Specimen

This device intended testing population is general population.

Specimen types

- Fresh, clear serum.
- Plasma in lithium heparin.
- Fresh centrifuged urine.

Anticoagulants other than those listed have not been tested by HORIBA and are therefore not recommended for use with this assay.

24h urine has to be collected without additive.

Stability

Serum, plasma (5)

- At 20-25°C: 7 days
- At 4-8°C: 7 days
- At -20°C: 3 months

Urine (6)

- At 20-25°C: 2 days
- At 4-8°C: 6 days
- At -20°C: 6 months

Reference Range

Each laboratory should establish its own reference ranges. The values given here are used as guidelines only.

Serum, plasma (7)

Men	Women
8 - 13 mg/L	6 - 12 mg/L
0.8 - 1.3 mg/dL	0.6 - 1.2 mg/dL
71 - 115 µmol/L	53 - 106 µmol/L

Urine (24 hours) (8)

Men	Women
14 - 26 mg/kg/day	11 - 20 mg/kg/day
124 - 230 µmol/kg/day	97 - 177 µmol/kg/day

Clinical sensitivity and specificity, positive predictive value and negative predictive value are not commonly reported for this analyte. This is largely attributed to the fact that this analyte is not sole indicator for the intended purpose and patient treatment decision making. To arrive at a diagnosis and a course of treatment, results from others routine clinical chemistry tests should be used in conjunction with other diagnostic information and the attending health-care professional's evaluation of the patient's condition.

Storage and Stability

Stability before opening:

Stable up to the expiry date on the label if stored at 18-26°C.

Yumizen C Creatinine Jaffé

Stability after opening:

Refer to the paragraph "Performance on Yumizen C230/C240".

Waste Management

Please refer to local legal requirements.

General Precautions

- This reagent is for professional *in vitro* diagnostic use only.
For laboratory use.
- For prescription use only.
- This reagent is classified as hazardous in compliance with regulation (EC) N°.1272/2008.
- **Reagent 1 and 2 (R1 and R2):**
Warning: This reagent is obtained from substances of animal origin. Consequently, it should be treated as potentially infectious and handled with the appropriate cautions in accordance with good laboratory practices (9).
- **Reagent 1 (R1):**
Warning
H290: May be corrosive to metals.
H315: Causes skin irritation.
H319: Cause serious eye irritation.
P234: Keep only in original container.
P264: Wash hands thoroughly after handling.
P280: Wear protective gloves/protective clothing/eye protection/face protection.
P302 + P352: IF ON SKIN: Wash with plenty of soap and water.
P305 + P351 + P338: IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
P321: Specific treatment (see [***] on this label).
P332 + P313: If skin irritation occurs: Get medical advice/attention.
P337 + P313: If eye irritation persists: Get medical advice/attention.
P362 + P364: Take off contaminated clothing and wash it before reuse.
P390: Absorb spillage to prevent material damage.
P406: Store in corrosive resistant container with a resistant inner liner.
- Observe the standard laboratory precautions for use.
- The reagent cassettes are disposable and should be disposed of in accordance with the local legal requirements.
- Please refer to the SDS associated with the reagent.
- Do not use the product if there is visible evidence of biological, chemical or physical deterioration.

- Do not use the product if the recommended storage conditions, including temperature, are not followed.
- User must be trained by a HORIBA representative before attempting to operate the device.
- It is the user's responsibility to verify that this document is applicable to the reagent used.
- For technical assistance, you can call +33 (0)4 67 14 15 16.
- Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the country in which the user and/or the patient is established.

Performance on Yumizen C230/C240

Serum, plasma

The performance data listed below have been obtained on the Yumizen C230/C240 analyzer.

Number of tests: approximately 2 x 155 tests

On Board Reagent Stability

Once opened, the reagent cassette placed in the refrigerated Yumizen C230/C240 compartment is stable for 14 days.

Sample volume: 13 µL/test

Lowest Detectable Level

The lowest detectable level represents the lowest measurable level of analyte that can be distinguished from zero. It is calculated as the absolute mean plus three standard deviations of 20 replicates of an analyte free sample. The lowest detectable level is estimated at 4.29 µmol/L (0.05 mg/dL).

Limit of Quantitation

The limit of quantitation is determined according to CLSI (NCCLS), EP17-A2 protocol (10) and equals 10.00 µmol/L (0.11 mg/dL).

Accuracy and Precision

Repeatability (*within-run precision*)

Repeatability according to the recommendations found in the CLSI (NCCLS), EP05-A3 protocol (11) with samples tested 20 times:

- 2 controls
- 3 specimens (low / medium / high levels)

Yumizen C Creatinine Jaffé

	Mean value µmol/L	Mean value mg/dL	CV %
Control specimen 1	109.57	1.24	1.2
Control specimen 2	343.89	3.89	1.9
Specimen 1	46.84	0.53	2.9
Specimen 2	158.32	1.79	1.8
Specimen 3	586.81	6.63	1.2

Reproducibility (total precision)

Reproducibility according to the recommendations found in the CLSI (NCCLS), EP05-A3 protocol (11) with samples tested in duplicate for 20 days (2 series per day):

- 2 controls
- 3 specimens (low / medium / high levels)

	Mean value µmol/L	Mean value mg/dL	CV %
Control specimen 1	110.75	1.25	3.9
Control specimen 2	342.46	3.87	2.5
Specimen 1	43.42	0.49	4.1
Specimen 2	146.43	1.65	3.2
Specimen 3	555.10	6.27	2.3

Measuring Range

The assay confirmed a measuring range from 10 µmol/L (0.11 mg/dL) to 1600 µmol/L (18.08 mg/dL).

The measuring range is extended up to 6400 µmol/L (72.32 mg/dL) with the automatic post-dilution.

The reagent linearity has been assessed up to 1600 µmol/L (18.08 mg/dL) according to the recommendations found in the CLSI (NCCLS), EP06-Ed2 protocol (12).

Correlation

Patient samples: Serum

Number of patient samples: 100

Specimens are correlated with a commercial reagent taken as reference according to the recommendations found in the CLSI (NCCLS), EP09c protocol (13).

Values ranged from 35.94 µmol/L (0.41 mg/dL) to 1172.76 µmol/L (13.25 mg/dL).

The equation for the allometric line obtained using Passing-Bablok regression procedure (14) is:

$$Y = 0.9076 X + 10.970 \text{ (µmol/L)}$$

$$Y = 0.9076 X + 0.124 \text{ (mg/dL)}$$

with a correlation coefficient $r^2 = 0.999$.

Interferences

Haemoglobin: No significant influence is observed up to 579 µmol/L (1000 mg/dL).

Triglycerides: No significant influence is observed up to a triglyceride concentration of 11.13 mmol/L (973.88 mg/dL).

Total Bilirubin: No significant influence is observed up to 331.10 µmol/L (19.37 mg/dL).

Direct Bilirubin: No significant influence is observed up to 281.95 µmol/L (16.49 mg/dL).

Other limitations are given by Young as a list of drugs and preanalytical variables known to affect this methodology (15, 16).

Calibration Stability

The reagent is calibrated on Day 0. The calibration stability is checked by testing 2 control specimens.

The calibration stability is 1 day.

Note: A recalibration is recommended when reagent lots change, and when quality control results fall outside the range established.

Conversion Factor

$$\mu\text{mol/L} \times 0.0113 = \text{mg/dL}$$

Urine

The performance data listed below have been obtained on the Yumizen C230/C240 analyzer.

Number of tests: approximately 2x 155 tests

On Board Reagent Stability

Once opened, the reagent cassette placed in the refrigerated Yumizen C230/C240 compartment is stable for 14 days.

Sample volume: 13 µL/test

Lowest Detectable Level

The lowest detectable level represents the lowest measurable level of analyte that can be distinguished from zero. It is calculated as the absolute mean plus three standard deviations of 20 replicates of an analyte free sample. The lowest detectable level is estimated at 39.76 µmol/L (0.45 mg/dL).

Limit of Quantitation

The limit of quantitation is determined according to CLSI (NCCLS), EP17-A2 protocol (10) and equals 100 µmol/L (1.13 mg/dL).

Yumizen C Creatinine Jaffé

Accuracy and Precision

Repeatability (within-run precision)

Repeatability according to the recommendations found in the CLSI (NCCLS), EP05-A3 protocol (11) with samples tested 20 times:

- 2 controls
- 3 specimens (low / medium / high levels)

	Mean value μmol/L	Mean value mg/dL	CV %
Control specimen 1	5791.55	65.44	1.2
Control specimen 2	11428.05	129.14	0.7
Specimen 1	1140.26	12.88	1.0
Specimen 2	8256.13	93.29	0.7
Specimen 3	19887.20	224.73	0.7

Reproducibility (total precision)

Reproducibility according to the recommendations found in the CLSI (NCCLS), EP05-A3 protocol (11) with samples tested in duplicate for 20 days (2 series per day):

- 2 controls
- 3 specimens (low / medium / high levels)

	Mean value μmol/L	Mean value mg/dL	CV %
Control specimen 1	5737.92	64.84	2.6
Control specimen 2	11277.55	127.44	2.4
Specimen 1	1055.37	11.93	2.9
Specimen 2	8025.72	90.69	2.0
Specimen 3	19916.96	225.06	2.0

Measuring Range

The assay confirmed a measuring range from 100 μmol/L (1.13 mg/dL) to 25000 μmol/L (282.5 mg/dL).

The measuring range is extended up to 100000 μmol/L (1130 mg/dL) with the automatic post-dilution.

The reagent linearity has been assessed up to 25000 μmol/L (282.5 mg/dL) according to the recommendations found in the CLSI (NCCLS), EP06-Ed2 protocol (12).

Correlation

Patient samples: urine

Number of patient samples: 99

Specimens are correlated with a commercial reagent taken as reference according to the recommendations found in the CLSI (NCCLS), EP09c protocol (13).

Values ranged from 2142.82 μmol/L (24.21 mg/dL) to 18687.14 μmol/L (211.16 mg/dL).

The equation for the allometric line obtained using Passing-Bablok regression procedure (14) is:

$$Y = 0.8707 X + 301.1 \text{ (}\mu\text{mol/L)}$$

$$Y = 0.8707 X + 3.40 \text{ (mg/dL)}$$

with a correlation coefficient $r^2 = 0.996$.

Interferences

Haemoglobin: No significant influence is observed up to 579 μmol/L (1000 mg/dL).

Triglycerides: No significant influence is observed up to a triglyceride concentration of 11.39 mmol/L (996.63 mg/dL).

Direct Bilirubin: No significant influence is observed up to 522.10 μmol/L (30.54 mg/dL).

Other limitations are given by Young as a list of drugs and preanalytical variables known to affect this methodology (15, 16).

Calibration Stability

The reagent is calibrated on Day 0. The calibration stability is checked by testing 2 control specimens.

The calibration stability is 1 day.

Note: A recalibration is recommended when reagent lots change, and when quality control results fall outside the range established.

Conversion Factor:

$$\mu\text{mol/L} \times 0.0113 = \text{mg/dL}$$

Reference

1. Jaffe M. Hoppe Selyer's. Z. Physiol. Chem. (1886) **10**: 391-400.
2. Folin O. Beitrag zur Chemie des Kreatinins und Kreatins im Harn. Z. Physiol. Chem. (1904) **41**: 223-242.
3. Fabing DL, Ertingshausen G. Automated Reaction Rate Method for the Determination of Serum Creatinine with the Centrifichem. Clin. Chem. (1971) **17**: 391.
4. Soldin S, Henderson L, Hill G. The Effect of Bilirubin and Ketones on Reaction Rate Methods for the Measurement of Creatinine. Clin. BioChem. (1978): 82-86.
5. Use of anticoagulants in diagnostic laboratory investigations. WHO publication WHO/DIL/LAB/99.1 Rev. 2 (2002): 28.

Yumizen C Creatinine Jaffé

6. Use of anticoagulants in diagnostic laboratory investigations. WHO publication WHO/DIL/LAB/99.1 Rev. 2 (2002): 46.
7. Tietz NW. Clinical guide to laboratory tests, 3rd Ed, (WB. Saunders eds. Philadelphia USA), (1995): 186.
8. Roberts WL, McMillin GA, Burtis CA, Bruns DE. Reference Information for the Clinical Laboratory, Tietz Textbook of Clinical Chemistry and Molecular Diagnostics. 4th Ed; Burtis CA, Ashwood ER, Bruns DE, (Elsevier Saunders eds. St Louis, USA), (2006): 2264.
9. Council Directive (2000/54/EC). Official Journal of the European Communities. No. L262 from October 17, 2000: 21-45.
10. Evaluation of detection capability for clinical laboratory measurement procedures. Approved Guideline, 2nd ed., CLSI (NCCLS) document EP17-A2 (2012) **32** (8).
11. Evaluation of Precision of Quantitative Measurement Procedures. Approved Guideline, CLSI (NCCLS) document EP05-A3 (2014) **24** (25).
12. Evaluation of Linearity of Quantitative Measurement Procedures. 2nd Edition, CLSI (NCCLS) guideline EP06-Ed2 (2020) **40** (16).
13. Measurement Procedure Comparison and Bias Estimation Using Patient Samples. Approved Guideline, 3rd ed., CLSI (NCCLS) document EP09c (2018) **38** (12).
14. Passing H, Bablok W. A new biometrical procedure for testing the equality of measurements from two different analytical methods. J. Clin. Chem. Clin. Biochem. (1983) **21**: 709-720.
15. Young DS. Effects of Drugs on Clinical Laboratory Tests. 5th Edition, Washington, DC, AACC Press (2000).
16. Young DS. Effects of Preanalytical Variables on Clinical Laboratory Tests. 2nd Edition, Washington, DC, AACC Press (1997) **3**: 120-132.