BioMedica Diagnostics Inc.
P.O. Box 1030, 94 Wentworth Road
Windsor, NS B0N 2T0 Canada
1.902.798.5105
www.biomedicadiagnostics.com

DVVtest[®] 10

DVVtest[®] 25

REF 810/825

DVVconfirm[®] 5

DVVconfirm[®] 10

REF 815/815L

Obelis s.a
Boulevard Général Wahis 53, 1030 Brussels, BELGIUM

ENGLISH

INTENDED USE

DVVtest® is a single-step dilute Russell's Viper Venom Time (dRVVT) test intended for the determination of the presence of Lupus Anticoagulants (LA) in human plasma. DVVtest may be performed manually, or by using semi-automated or automated coagulation analysers. The reagent is for in vitro diagnostic use.

DVVconfirm® is a single-step dilute Russell's Viper Venom Time (dRVVT) test intended for the determination of the presence of Lupus Anticoagulants (LA) in human plasma. DVVconfirm may be performed manually, or by using semi-automated or automated coagulation analysers. The reagent is for in vitro diagnostic use.

EXPLANATION OF THE TEST

The identification of Lupus Anticoagulants (LA) in plasma is a diagnostic hallmark of the Antiphospholipid Syndrome (APS) which is characterized clinically by arterial/venous thrombosis, recurrent spontaneous fetal loss, thrombocytopenia and/or neurological disorders.¹ LA may also be induced as a result of administration of chlorpromazine, procainamide, thiorazine and certain antibiotics. In 1952, Conley and Hartmann² first described the presence of LA in patients with Systemic Lupus Erythematosus (SLE). It is now recognized that LA are more commonly found in patients who do not have SLE.³ LA autoantibodies are specifically directed against a variety of phospholipid binding proteins including β 2-glycoprotein I (β 2GPI), prothrombin, and annexin V which are complexed to different anionic phospholipids (e.g. cardiolipin, phosphatidylinositol and phosphatidylserine⁴). Lupus Anticoagulants are immunoglobulins of the IgG, IgM and IgA isotypes that prolong one or more of the *in vitro* phospholipid-dependent coagulation tests (e.g. activated partial thromboplastin time (aPTT), dilute prothrombin time (dPT), Textarin time or dRVVT). The International Society for Thrombosis and Hemostasis,

Scientific Subcommittee Criteria on Lupus Anticoagulants and Phospholipid-dependent Antibodies has recommended that LA be diagnosed using coagulation-based screening tests and a confirmatory test containing a high phospholipid concentration.⁵ BioMedica Diagnostics has developed DVVtest as a primary screening diagnostic test for LA and DVVconfirm as the companion high phospholipid-containing test to confirm the diagnosis of LA. DVVtest and DVVconfirm are specially formulated to reduce the effect of heparin to maximize the sensitivity and specificity of diagnosing LA.¹⁴

PRINCIPLE OF THE METHOD

DVVtest reagent contains activator isolated from Russell's Viper Venom that directly activates Factor X to Factor Xa in the presence of phospholipids and calcium. Factor Xa cleaves prothrombin to thrombin, which converts fibrinogen to fibrin leading to detectable clot formation in plasma. This direct activation of Factor X bypasses the Contact and Intrinsic pathways in the coagulation cascade, thereby excluding interference from deficiencies of Factors VIII, IX, XI and XII, or their respective inhibitors. A positive DVVtest is indicated by a significant prolongation of the phospholipid-dependent clotting time (> 2 standard deviations above the laboratory's normal reference mean). If a patient is suspected to have LA, the DVVtest may also be performed on samples with a normal APTT because the dilution and type of phospholipids contained in the DVVtest reagent increases the test's sensitivity and specificity to LA.⁶

DVVconfirm is the high phospholipid containing coagulation assay used in conjunction with DVVtest to confirm for the presence of LA in plasma.⁷ DVVconfirm is formulated with activator isolated from Russell's Viper Venom and a high concentration of phospholipids. The clotting time of a plasma containing LA should be significantly shorter with DVVconfirm as compared to DVVtest. The presence of LA in plasma samples is confirmed when the ratio of the DVVtest clotting time to the DVVconfirm clotting time is greater than the range of the laboratory's normal reference DVVtest/DVVconfirm Ratio (mean normal ratio \pm 2 S.D.).

REAGENTS

The DVVtest and DVVconfirm reagents are supplied lyophilized and consist of proprietary mixtures of activator isolated from Russell's Viper Venom, calcium, and phospholipids, as well as inert additives and preservatives. Unopened vials stored at 2°- 8°C are stable until the expiration date printed on the label.

The DVVtest reagents are packaged in two vial sizes: REF 810, comprises 10 vials, with each vial containing reagent for 20 tests as performed by most automated methods; REF 825, comprises 10 vials, with each vial containing reagent for 50 tests as performed by most automated methods.

The DVVconfirm reagents are packaged in two vial sizes: REF 815, comprises 10 vials, with each vial containing reagent for 10 tests as performed by most automated methods; REF 815L, comprises 10 vials, with each vial containing reagent for 20 tests as performed by most automated methods.

WARNING

DVVtest®	Warning		H319, H412, P264, P273, P280, P305 + P351 + P338, P337+ P313		
DVVconfirm®	Danger		<table border="1"><tr><td>CONT</td><td>Imidazole</td></tr></table> H315, H318, H360, H412, P202, P280, P273, P302 + P352, P332 + P313, P305 + P351 + P338	CONT	Imidazole
CONT	Imidazole				

H315 Causes skin irritation.

H318 Causes serious eye damage.

H319 Causes serious eye irritation.

P202 Do not handle until all safety precautions have been read and understood.

H360 May damage fertility or the unborn child.

H412 Harmful to aquatic life with long lasting effects.

P264 Wash thoroughly after handling.

P273 Avoid release to the environment.

P280 Wear protective gloves/protective clothing/eye protection/face protection.

P302 + P352 IF ON SKIN: Wash with plenty of water.

P305 + P351 + P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

P332 + P313 If skin irritation occurs: Get medical advice/attention.

P337 + P313 If eye irritation persists: Get medical advice/attention.

SAMPLE COLLECTION AND PREPARATION

Blood samples may be drawn using syringes or evacuated siliconized tubes. Mix nine parts of blood with one part of 3.2 or 3.8% (109 or 129 mM) of the dihydrate form of trisodium citrate. If the patient's haematocrit is greater than 55%, the DVVtest and DVVconfirm results may be inaccurate, and require an adjustment of the blood to anticoagulant ratio.⁸ After mixing, centrifuge the samples at a minimum of 5,000-x g for 10 minutes. It is imperative that all plasmas tested using this reagent be platelet poor, having fewer than 10^4 platelets/ μL , especially prior to their freezing, otherwise some shortening of the test results may occur which could mask the presence of LA⁵. This can be achieved by means of the above centrifugation step or by filtering the plasma through a 0.22 micron filter. Collected plasma may then be stored at 2°- 8°C but should be tested within four hours. Alternatively, plasma may be rapidly frozen and stored at -70°C for up to six months. Frozen plasmas must be thawed rapidly at 37°C before use, and tested immediately, or stored for no more than two hours at 2°- 8°C.⁸

MATERIALS REQUIRED BUT NOT PROVIDED

LAtrol™ Abnormal Control Plasma (REF 816A), LAtrol™ Normal Control Plasma (REF 816N)¹⁰ or equivalent
Pipettes that deliver volumes from 100 µL up to 5.0 mL

Purified deionized or distilled water

REAGENT PREPARATION AND STABILITY

NOTE: The DVVconfirm lyophilizate has a powdery appearance; however, the DVVtest lyophilizate may appear crystalline.

Reconstitute each vial of the reagent with deionized water according to the instructions on the label.

REF	810	825	815	815L
Volume H ₂ O	2.0 mL	5.0 mL	1.0 mL	2.0 mL

Mix the reagent well and allow to stand at room temperature for at least 15 minutes to ensure complete dissolution. Once reconstituted, the reagent will remain stable for:

	24 hours	5 days	1 month
	20°- 25°C	2°- 8°C	-20°C

Reconstituted reagent may be recovered from instrument reagent reservoirs provided the 24-hour room temperature stability time has not passed.

Reconstituted reagent may be frozen in plastic storage tubes to ensure its stability.

To use frozen reagent, thaw the reagent for at least 10 minutes at 37°C and vortex well before use.

TEST PROCEDURES

DVVtest and DVVconfirm may be performed manually, or by using semi-automated or automated coagulation analysers.

Equal volumes of reagent and sample are added (i.e., 100 µL DVVtest or DVVconfirm and 100 µL sample plasma). DVVtest or DVVconfirm results are reported in time (seconds). Reconstituted reagent should be mixed well before use. Both reagent and sample plasma should be incubated at 37°C for at least two minutes prior to testing. Use dedicated reagent tubing and place magnetic stir bars in reagent reservoirs. Perform a rinsing and/or cleaning procedure after using the reagent, especially on instruments that do not have dedicated reagent tubing, before performing other clotting or chromogenic assays.

BioMedica Diagnostics offers Instrument Applications for several coagulation analysers. These Instrument Applications may contain platform specific programming and performance data which differ from that provided in this Instructions for Use. In these cases, the information contained in the Instrument Application supersedes the information in this Instruction for Use. Please consult the specific manufacturer's instrument manual for complete operating instructions.

QUALITY CONTROL


Both normal and abnormal LA control plasmas should be tested with each group of assays, with a change in personnel or work shift, and with every 40 samples.¹¹ Control plasmas must be platelet poor containing fewer than 10^4 platelets/ μ L. LATrol Abnormal and Normal Plasmas (REF 816A and 816N)¹⁰ may be used for quality control. Values for both the normal and abnormal LA control plasmas should fall within the laboratory's established control ranges before patient samples are tested.⁸

NORMAL REFERENCE RANGE

Each laboratory should establish its own DVVtest and DVVconfirm normal reference range (mean value \pm 2 S.D.). A minimum of 20 plasmas from healthy blood donors representative of the patient population, including both men and women and spanning the adult age range, should be used (see EXPECTED VALUES). When establishing the normal reference range, handle the normal plasma samples in the same manner as test samples. If frozen samples are tested exclusively, the normal ranges should be established on frozen normal samples only. It is not recommended to test mixed plasma populations of fresh and frozen samples, either for the reference ranges or for routine testing. A new normal reference range MUST be established whenever there is a change in reagent lot number or instrument or at least once per year (whichever comes first). Obtain data over a period of several days to determine day-to-day variation.

INTERPRETATION AND REPORTING RESULTS

LA Decision Flow Chart

DVVtest Results

If the DVVtest time of the sample plasma is within the established normal reference range for the laboratory (mean \pm 2 S.D.), the test is negative for LA.

If the DVVtest time of the sample plasma is greater than 2 S.D.'s above the mean of the normal reference, the test is positive for LA.

NOTE: A positive DVVtest indicates either LA are present or the plasma may contain a Factor II, V or X deficiency. When the DVVtest is positive, the presence of LA should be confirmed by retesting using DVVconfirm to obtain a DVVtest/DVVconfirm Ratio (see below). The DVVconfirm test should be performed on the same patient sample and on the same day that the DVVtest has been performed. Mixing studies should be performed to determine the presence of a factor deficiency.

DVVconfirm Results

If a sample tested positive for DVVtest, then the DVVconfirm test should be performed. The results from DVVconfirm testing can be reported using DVVtest/DVVconfirm Ratio or Normalized Ratio.

1) The DVVtest/DVVconfirm Ratio is calculated by dividing the DVVtest time (seconds) by the DVVconfirm time (seconds) as follows:

$$\text{DVVtest/DVVconfirm Ratio} = \text{DVVtest (sec)} \div \text{DVVconfirm (sec)}$$

The DVVtest/DVVconfirm Ratio must be determined for normal reference plasmas. The ratio of the test plasma when greater than 2 standard deviations above the mean of the laboratory's normal reference DVVtest/DVVconfirm Ratio confirms a positive test for LA. If the DVVtest/DVVconfirm Ratio result is within the normal reference ratio range, then the test is negative for LA, or other abnormalities are present.

2) The Normalized Ratio is calculated as follows:

$$\text{Normalized Ratio} = \frac{\text{DVVtest (sample)} \div \text{DVVtest (mean normal reference)}}{\text{DVVconfirm (sample)} \div \text{DVVconfirm (mean normal reference)}}$$

Example Only: Patient DVVtest: 58.9 seconds
 Patient DVVconfirm: 32.2 seconds
 Mean DVVtest: 35.8 seconds
 Mean DVVconfirm: 32.8 seconds

$$\text{Normalized Ratio} = \frac{58.9 \div 35.8}{32.2 \div 32.8} = \frac{1.64}{0.98} = 1.67$$

Interpretation of Normalized Ratio Results

Normalized Ratio	LA Status
> 2.0	LA is strongly present
1.5 to 2.0	LA is moderately present
1.2 to 1.5	LA is weakly present

Test results may be reported either as positive/negative for LA, or by reporting the DVVtest/DVVconfirm Ratio or Normalized Ratio. When reporting the Ratio results, the value for the reference range must also be reported.

EXPECTED VALUES

Normal Reference Ranges (± 2 S.D.)¹²

INSTRUMENT:	MLA 1000C	ST4	ACL 300+
DVVtest Time (sec)	28 - 45	29 - 51	28 - 47
DVVconfirm Time (sec)	30 - 44	21 - 46	28 - 40
DVVtest/DVVconfirm Ratio	0.77 - 1.21	0.95 - 1.47	0.83 - 1.35
n	21	25	20

Reference ranges were generated by BioMedica Diagnostics in accordance with CLSI Document H21-A5.⁸ Donors included both men and women, spanning the adult age range and having normal PT and APTT values. These results should be used only as a guide. Results will vary among laboratories, and each testing site must establish its own normal reference range for all instruments and methods employed.

INTERFERING SUBSTANCES

Samples suspected of containing inhibitors of Factors II, V or X may give false positive results. If DVVtest is performed on plasmas that contain high levels of Factor VIII (greater than 200%), the reagent may not identify the presence of LA, yielding a false negative.¹²

Oral anticoagulants and other Vitamin K antagonists may prolong the DVVtest and DVVconfirm times. In 2006, an update to the classification criteria for antiphospholipid syndrome (APS) was issued.¹³ The update recommends that if the patient is on an oral anticoagulant, the patient sample should be diluted 1:2 with normal

plasma (1 part patient plasma + 1 part normal plasma) before performing the test, providing that the INR (international normalized ratio) is < 3.5. When the INR > 3.5, LA testing should not be attempted.

If the laboratory does not know that the patient is on an oral anticoagulant, the test results for the DVVconfirm must be looked at closely. If the DVVconfirm does not show significant correction (a shorter clot time) relative to the DVVtest and is prolonged relative to the DVVconfirm time for normal plasma, then the presence of an oral anticoagulant is suspected. In this case the patient sample should be mixed 1:2 with normal plasma and the DVVtest performed again. If the DVVtest clot time is still prolonged, the DVVconfirm should also be performed on the 1:2 mix to confirm if the DVVtest results are phospholipid-dependent.

DVVtest and DVVconfirm contain agents that neutralize standard unfractionated heparin up to and including 1.0 U/mL. Plasmas containing heparin levels greater than 1.0 U/mL may give elevated results with these tests. Some low molecular weight heparins (LMWH) may interfere with DVVtest. Studies have shown that plasma spiked with 2.0 units/mL of Fragmin® (dalteparin sodium injection) did not interfere with assay results. Plasma spiked with Lovenox® (enoxaparin sodium injection) at concentrations greater than 0.25 units/mL interfered with assay results.¹²

LIMITATIONS OF THE PROCEDURE

Do not test hemolyzed, lipemic, or icteric samples. Samples suspected to have Factor II, V or X deficiencies should be tested using mixing studies.

Samples should be tested using at least one other coagulation-based LA test as no single test can absolutely identify the presence of Lupus Anticoagulants in every plasma.⁵

PERFORMANCE CHARACTERISTICS

Precision

The DVVtest and DVVconfirm assays have been performed in a quality control program by BioMedica Diagnostics and two independent laboratories. The overall coefficient of variation (CV) for normal plasmas was less than 4.0% for DVVtest and less than 5.0% for DVVconfirm. For abnormal plasmas the overall CV was less than 6.5% for DVVtest and less than 5.0% for DVVconfirm.¹²

Specificity

The DVVtest and DVVconfirm reagents were tested¹² using both normal and abnormal plasmas as described below.

Plasma Sample	Positive DVVtest/DVVconfirm
Lupus anticoagulant plasma	100% (17/17)
normal plasma	2.1% (2/96)
heparinized plasma	0% (0/2)
factor deficient plasma	0% (0/8)
factor VIII inhibitor plasma	0% (0/2)

SUMMARY OF SAFETY AND CLINICAL PERFORMANCE

The Summary of Safety and Clinical Performance (SCCP) is available in the European Database on Medical Devices (Eudamed), at <https://ec.europa.eu/tools/eudamed> under the following Basic UDI-DI (GMN):

REF 810 0062794400205C

REF 815 0062794400215E

REF 825 0062794400255N

REF 815L 0062794400315H

REFERENCES

1. Harris, E. N., Asherson, R. A., and Hughes, G. R. V. Antiphospholipid antibodies - autoantibodies with a difference. *Ann Revs Med* 1988, **39**: 261-271.
2. Conley, C. L. and Hartmann, R. C. A hemorrhagic disorder caused by a circulating anticoagulant in patients with disseminated lupus erythematosus. *J Clin Invest* 1952, **31**: 621-622.
3. Triplett, D. A., Brandt, J. T. and Maas, R. L. The laboratory heterogeneity of lupus -anticoagulants. *Arch Path Lab Med* 1985; **109**: 946-951.
4. Thiagarajan, P., Shapiro, S. S and DeMarco, L. Monoclonal immunoglobulin M lambda coagulation inhibitor with phospholipid specificity: Mechanism of a lupus anticoagulant. *J Clin Invest* 1980, **66**: 397-405.
5. Pengo V, Tripodi A, Reber G, Rand JH, Ortel TL, Galli M, de Groot PG. Update of the guidelines for lupus anticoagulant detection. *J Thromb Haemost* 2009; **7**: 1737-40.
6. Exner, T., Papadopoulos, G. and Koutts, J. Use of a simplified dilute Russell viper venom time (dRVVT) confirms heterogeneity among 'Lupus Anticoagulants.' *Blood Coag Fibrinol* 1990, **1**: 259-266.
7. Thiagarajan, P., Pengo, V. and Shapiro, S. S. The use of the dilute Russell viper venom time for the diagnosis of lupus anticoagulants. *Blood* 1986, **68**: 869-874.
8. CLSI. *Collection, transport and processing of blood specimens for testing plasma-based coagulation assays and molecular hemostasis assays – 5th Edition*. CLSI document H21-A5. Wayne, PA: Clinical and Laboratory Standards Institute; 2008.
9. Brandt, J. T. "Q & A". *CAP Today*, August 1992; Page 61.
10. Available from BioMedica Diagnostics, or your local distributor
11. Passey, R. B. Walking the straight and narrow on quality control. *MLO* 1993; **2**: 39-43.
12. Data on file, BioMedica Diagnostics.

13. Miyakis, S., *et al.* International consensus statement on an update of the classification criteria for definite antiphospholipid syndrome (APS). *Journal of Thrombosis and Haemostasis* 2006, **4**: 295-306.
14. CLSI. *Laboratory Testing for the Lupus Anticoagulant; Approved Guideline*. CLSI document H60-A. Wayne, PA: Clinical and Laboratory Standards Institute; 2014.

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LITERATUR

1. Harris, E. N., Asherson, R. A., and Hughes, G. R. V. Antiphospholipid antibodies - autoantibodies with a difference. *Ann Revs Med* 1988, **39**: 261-271.
2. Conley, C. L and Hartmann, R. C. A hemorrhagic disorder caused by a circulating anticoagulant in patients with disseminated lupus erythematosus. *J Clin Invest* 1952, **31**: 621-622.

3. Triplett, D. A., Brandt, J. T. and Maas, R. L. The laboratory heterogeneity of lupus anticoagulants. *Arch Path Lab Med* 1985, **109**: 946-951.
4. Thiagarajan, P., Shapiro, S. S. and DeMarco L. Monoclonal immunoglobulin M lambda coagulation inhibitor with phospholipid specificity: Mechanism of a lupus anticoagulant. *J Clin Invest* 1980, **66**: 397-405.
5. Pengo V, Tripodi A, Reber G, Rand JH, Ortel TL, Galli M, de Groot PG. Update of the guidelines for lupus anticoagulant detection. *J Thromb Haemost* 2009; **7**: 1737–40.
6. Exner, T., Papadopoulos, G. and Koutts, J. Use of a simplified dilute Russell viper venom time (dRVVT) confirms heterogeneity among 'Lupus Anticoagulants.' *Blood Coag Fibrinol* 1990, **1**: 259-266.
7. Thiagarajan, P., Pengo, V. and Shapiro, S. S. The use of the dilute Russell viper venom time for the diagnosis of lupus anticoagulants. *Blood* 1986, **68**: 869-874.
8. CLSI. *Collection, transport and processing of blood specimens for testing plasma-based coagulation assays and molecular hemostasis assays – 5th Edition*. CLSI document H21-A5. Wayne, PA: Clinical and Laboratory Standards Institute; 2008.
9. Brandt, J. T. "Q & A". *CAP Today*, August 1992; Page 61.
10. Erhältlich über BioMedica Diagnostics, or your local distributor.
11. Passey, R. B. Walking the straight and narrow on quality control. *MLO* 1993, **2**: 39-43.
12. Daten verfügbar bei BioMedica Diagnostics
13. Miyakis, S., *et al.* International consensus statement on an update of the classification criteria for definite antiphospholipid syndrome (APS). *Journal of Thrombosis and Haemostasis* 2006, **4**: 295-306.
14. CLSI. *Laboratory Testing for the Lupus Anticoagulant; Approved Guideline*. CLSI document H60-A. Wayne, PA: Clinical and Laboratory Standards Institute; 2014.

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4. Thiagarajan, P., Shapiro, S. S. and DeMarco, L. Monoclonal immunoglobulin M lambda coagulation inhibitor with phospholipid specificity: Mechanism of a lupus anticoagulant. *J Clin Invest* 1980, **66**: 397-405.
5. Pengo V, Tripodi A, Reber G, Rand JH, Ortel TL, Galli M, de Groot PG. Update of the guidelines for lupus anticoagulant detection. *J Thromb Haemost* 2009; **7**: 1737–40.
6. Exner, T., Papadopoulos, G. and Koutts, J. Use of a simplified dilute Russell viper venom time (dRVVT) confirms heterogeneity among 'Lupus Anticoagulants.' *Blood Coag Fibrinol* 1990, **1**: 259-266.
7. Thiagarajan, P., Pengo, V. and Shapiro, S. S. The use of the dilute Russell viper venom time for the diagnosis of lupus anticoagulants. *Blood* 1986, **68**: 869-874.
8. CLSI. *Collection, transport and processing of blood specimens for testing plasma-based coagulation assays and molecular hemostasis assays – 5th Edition*. CLSI document H21-A5. Wayne, PA: Clinical and Laboratory Standards Institute; 2008.
9. Brandt, J. T. "Q & A". *CAP Today*, August 1992; Page 61.
10. Disponible chez BioMedica Diagnostics, ou votre distributeur local.
11. Passey, R. B. Walking the straight and narrow on quality control. *MLO* 1993, **2**: 39-43.
12. Documentation interne, BioMedica Diagnostics.
13. Miyakis, S., *et al.* International consensus statement on an update of the classification criteria for definite antiphospholipid syndrome (APS). *Journal of Thrombosis and Haemostasis* 2006, **4**: 295-306.
14. CLSI. *Laboratory Testing for the Lupus Anticoagulant; Approved Guideline*. CLSI document H60-A. Wayne, PA: Clinical and Laboratory Standards Institute; 2014.

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2. Conley, C. L and Hartmann, R. C. A hemorrhagic disorder caused by a circulating anticoagulant in patients with disseminated lupus erythematosus. *J Clin Invest* 1952, **31**: 621-622.
3. Triplett, D. A., Brandt, J. T. and Maas, R. L. The laboratory heterogeneity of lupus anticoagulants. *Arch Path Lab Med* 1985, **109**: 946-951.
4. Thiagarajan, P., Shapiro, S. S. and DeMarco, L. Monoclonal immunoglobulin M lambda coagulation inhibitor with phospholipid specificity: Mechanism of a lupus anticoagulant. *J Clin Invest* 1980, **66**: 397-405.
5. Pengo V, Tripodi A, Reber G, Rand JH, Ortel TL, Galli M, de Groot PG. Update of the guidelines for lupus anticoagulant detection. *J Thromb Haemost* 2009; **7**: 1737-40.
6. Exner, T., Papadopoulos, G. and Koutts, J. Use of a simplified dilute Russell viper venom time (dRVVT) confirms heterogeneity among 'Lupus Anticoagulants.' *Blood Coag Fibrinol* 1990, **1**: 259-266.
7. Thiagarajan, P., Pengo, V. and Shapiro, S. S. The use of the dilute Russell viper venom time for the diagnosis of lupus anticoagulants. *Blood* 1986, **68**: 869-874.
8. CLSI. *Collection, transport and processing of blood specimens for testing plasma-based coagulation assays and molecular hemostasis assays – 5th Edition*. CLSI document H21-A5. Wayne, PA: Clinical and Laboratory Standards Institute; 2008.
9. Brandt, J. T. "Q & A". *CAP Today*, August 1992; Page 61.
10. Disponibili presso l'BioMedica Diagnostics, o presso il vostro distributore locale.
11. Passey, R. B. Walking the straight and narrow on quality control. *MLO* 1993, **2**: 39-43.
12. Dati disponibili su file presso la BioMedica Diagnostics.
13. Miyakis, S., *et al.* International consensus statement on an update of the classification criteria for definite antiphospholipid syndrome (APS). *Journal of Thrombosis and Haemostasis* 2006, **4**: 295-306.
14. CLSI. *Laboratory Testing for the Lupus Anticoagulant; Approved Guideline*. CLSI document H60-A. Wayne, PA: Clinical and Laboratory Standards Institute; 2014.

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REF 815L 0062794400315H

REFERENCES

1. Harris, E. N., Asherson, R. A. and Hughes, G. R. V. Antiphospholipid antibodies - autoantibodies with a difference. *Ann Revs Med* 1988, **39**: 261-271.
2. Conley, C. L and Hartmann, R. C. A hemorrhagic disorder caused by a circulating anticoagulant in patients with disseminated lupus erythematosus. *J Clin Invest* 1952, **31**: 621-622.
3. Triplett, D. A., Brandt, J. T. and Maas, R. L. The laboratory heterogeneity of lupus -anticoagulants. *Arch Path Lab Med* 1985, **109**: 946-951.
4. Thiagarajan, P., Shapiro, S. S. and DeMarco, L. Monoclonal immunoglobulin M lambda coagulation inhibitor with phospholipid specificity: Mechanism of a lupus anticoagulant. *J Clin Invest* 1980, **66**: 397-405.
5. Pengo V, Tripodi A, Reber G, Rand JH, Ortel TL, Galli M, de Groot PG. Update of the guidelines for lupus anticoagulant detection. *J Thromb Haemost* 2009; **7**: 1737-40.
6. Exner, T., Papadopoulos, G. and Koutts, J. Use of a simplified dilute Russell viper venom time (dRVVT) confirms heterogeneity among 'Lupus Anticoagulants.' *Blood Coag Fibrinol* 1990, **1**: 259-266.
7. Thiagarajan, P., Pengo, V. and Shapiro, S. S. The use of the dilute Russell viper venom time for the diagnosis of lupus anticoagulants. *Blood* 1986, **68**: 869-874.
8. CLSI. *Collection, transport and processing of blood specimens for testing plasma-based coagulation assays and molecular hemostasis assays – 5th Edition*. CLSI document H21-A5. Wayne, PA: Clinical and Laboratory Standards Institute; 2008.
9. Brandt, J. T. "Q & A". *CAP Today*, August 1992, Page 61.
10. Available from BioMedica Diagnostics, or your local Distributor
11. Passey, R. B. Walking the straight and narrow on quality control. *MLO* 1993, **2**: 39-43.
12. Data on file, BioMedica Diagnostics.

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