



**Hematology Analyzer - HELO 2.X**

## **Output Format for Host Connection**

Ref: RAA090BEN



# Output Format for Host Connection

---



**HORIBA ABX SAS**  
Parc Euromédecine  
Rue du Caducée  
BP 7290  
34184 Montpellier Cedex 4  
FRANCE



# Contents

---

<b>1. Foreword</b> .....	<b>1</b>
1.1. Document Update.....	1
<b>2. Connection</b> .....	<b>2</b>
2.1. Ethernet Connection.....	2
<b>3. ASTM Format</b> .....	<b>4</b>
3.1. Introduction.....	4
3.2. Connection Specifications (LIS01-A2).....	5
3.3. Records General Format Specifications (LIS2-A2).....	7
3.4. Special Characteristics for HORIBA Data.....	19
3.5. Curves and Matrix Transmission.....	32
3.6. Example of Data Frame.....	37



# 1. Foreword

## 1.1. Document Update

---

### 1.1.1. Revisions

Internal Reference	Software Version used for Documentation	Document Date Issued
RAA090AEN	1.7.0	11/2023
RAA090AAEN	1.7.1	11/2024
RAA090BEN	1.7.3	08/2025

This document applies to the latest software version listed and higher versions.

### 1.1.2. What's New?

Update	Chapter
Updates	<a href="#">Output Data Characteristics</a> <a href="#">Order Record</a> <a href="#">Manufacturer Record</a> <a href="#">Parameters</a> <a href="#">Units</a>
New chapter added.	<a href="#">Universal Test ID Definition</a>

## 2. Connection

### 2.1. Ethernet Connection

---

#### 2.1.1. Ethernet Connection Overview

---



The network connection mode is available for the ASTM format on the Yumizen H1500/H2500.


---

The implementation of network-based communication is based on the Windows Socket standard.

The data transmitted between the client and the server takes the form of ASTM high level packets.

This connection is made through the RJ45 connector on the instrument.

The keepalive feature is activated on the Yumizen H1500/H2500. The instrument sends a keepalive signal to the Host every 15 minutes in order to preserve the connection.




#### 2.1.2. To Configure the Instrument Network

Access: **Home** > **Settings** > **System** > **Interfaces**

Make sure that the RJ45 cable is connected from the instrument to the laboratory network.

- The instrument network settings must be configured only if a fixed IP address is required.
- Ask the laboratory IT for network settings (free IP address, subnet and gateway).

1. Press **Edit**.
2. Scroll down to see the **Edit** button of the **Network** area.
3. In the **Network** area, select the line with **True** in the **Link** column.  
The information in the **Link** column shows you which network card is connected.  
**True** means connected and **False** not connected.



4. Press **Edit**.
5. Choose **Manual configuration** as the IP address must be fixed.
6. Type the configuration defined by the laboratory IT:
  - **IP address**
  - **Subnet mask**
  - **Default gateway**
7. Press **Confirm** to validate or **Cancel** to go back to the screen without saving changes.



To go back to the default values in all screen areas:

1. Press **Reinitialize**.
2. Press **Confirm**.

### 2.1.3. To Configure the Connection to the Host

Access: **Home > Settings > System > Data management**

The IP address provided by your laboratory must be fixed and set during the first installation for each instrument.

1. Press **Edit**.
2. In the **Data management** field, select **Yes** to allow the connection to the Host.
3. In the **DM IP address** field, enter the fixed IP address of the Host.
4. In the **TCP port** field, enter the fixed TCP port number of the Host.
5. In the **Connection mode** drop-down list, choose **Query**.
6. In the **Protocol** drop-down list, choose **P8000** or **ASTM E1394**.  
If the protocol **ASTM E1394** is selected, M (Manufacturer Information) records are not managed and the following information is not transmitted: curves and matrices, traceability, statistics, remote command.
7. In the **Parameters units**:
  - Choose **Current** if you want the instrument to send your results to the Host in the units you chose in the **Home > Settings > System > General** screen.
  - Or choose **International System (SI)** if you want the instrument to send your results to the Host always in international system.
8. Press **Confirm** to validate or **Cancel** to go back to the screen without saving changes.
9. Test the connection by pressing **Test** in the **Connection test** area.  
If the connection is failed, make sure you have the right fixed IP address and repeat the procedure again.



To go back to the default values in all screen areas:

1. Press **Reinitialize**.
2. Press **Confirm**.

## 3. ASTM Format

The HORIBA analyzers format corresponds to the ASTM specifications LIS01-A2 & LIS2-A2:

- LIS01-A2: Standard specification for low level protocol to transfer messages between clinical and laboratory instruments and computer systems.
- LIS2-A2: Standard specification for transferring information between clinical and laboratory instruments and computer systems.

### 3.1. Introduction

A connection between a computer (host) and a HORIBA instrument can be performed when the protocol, the format description and the connection mode have been properly setup.

Term	Definition
<ACK>	Acknowledgment (ASCII decimal 6)
[C1]	The most significant character of checksum
[C2]	The least significant character of checksum
[DATA]	The data contents of the record
<ENQ>	Inquire (ASCII decimal 5)
<ETB>	End of transmission block (ASCII decimal 23). For use only when a single record is too large to fit into one frame.
<ETX>	End of text (ASCII decimal 3). Required at the end of each record.
[frame number]	Single digit frame number "0" to "7", starts with "1".
<LF>	Line feed (ASCII decimal 10).
<NAK>	Negative acknowledgment (ASCII Decimal 21).
<STX>	Start of frame (ASCII decimal 2).

Term	Definition
Communications packet	All framing required for transmission of data. This framing includes: <STX>[frame number][DATA] [<ETB> or <ETX>][C1][C2] <LF>
Component Field	One of several related pieces of information within a field.
Field	A specific location within a record for a piece of information, indicated by a field delimiter and position.
Frame	A complete communications packet.
LIS	Laboratory Information System
Message	A collection of related information; a group of records that begins with a "Header" record and ends with a "Terminator" record. A single record could theoretically constitute a message, but within this context, a message always contains multiple records.
<EOT>	End of transmission (ASCII decimal 4)
<CR>	Carriage return (ASCII decimal 13)
Record	In reference to the low level protocol, a record is the message data (shown as [DATA]) as described within the communications packet. If the data is longer than 240 characters, then it must be split in two (or more) parts and sent in two (or more) communications packets. The intermediate packet uses the <ETB> character, and the ending packet uses the <ETX> character. No single communications packet contains more than one record. In reference to the message layer, a record can be one of the following codes: H (header), P (patient), O (order), R (result), L (terminator), C (comment).
Session	A total unit of communication activity used in this standard to indicate the events starting with the establishment phase and ending with the termination phase.
Test	A determination of a single analyte or a combination of values from other determinations or observations from which a variable or gradable result is derived.

## 3.2. Connection Specifications (LIS01-A2)

### 3.2.1. Hardware and Software Characteristics

The hardware settings of the interface is an ethernet connection via an RJ45 cable.

### 3.2.2. Output Data Characteristics

- Characters: ASCII
- Maximum message length: 247 characters
- The numerical values must use the character "." as decimal separator

### 3.2.3. Communication Protocol

#### Standard control characters

Control String	Hexadecimal value
<ENQ>	\$05
<ACK>	\$06
<NAK>	\$15
<STX>	\$02
<ETX>	\$03
<ETB>	\$17
<CR>	\$0D

Control String	Hexadecimal value
<LF>	\$0A
<EOT>	\$04

#### Typical discussion between the instrument and the host

Instrument	<>	Host
<ENQ>	>	
	<	<ACK>
<STX>1...Data...<CR><ETX>xx <CR><LF>	>	
	<	<ACK>
<STX>2...Data...<CR><ETX>xx <CR><LF>	>	
	<	<ACK>
<EOT>	>	

#### Typical discussion between the host and the instrument

Instrument	<>	Host
	<	<ENQ>
<ACK>	>	
	<	<STX>1...Data...<CR><ETX>xx <CR><LF>
<ACK>	>	
	<	<STX>2...Data...<CR><ETX>xx <CR><LF>
<ACK>	>	
	<	<EOT>

#### Discussion with conflict between the instrument and the host

No answer from Host for an <ENQ>

- Timeout: 15 seconds
- In case of conflict: 1 second before a new transmission, up to 3 transmissions. Host timeout: 20 seconds
- In case of negative answer <NAK>: No time before a new transmission, up to 6 transmissions

Instrument	<>	Host
<ENQ>	>	
	<	<ENQ>
Wait 1 second		Wait 20 seconds
<ENQ>	>	
	<	
	...	
<EOT>	>	

#### Defect packet during discussion between instrument and host

Instrument	<>	Host
<ENQ>	>	
	<	<ACK>
<STX>1...Data...<CR><ETX>xx <CR><LF>	>	
	<	<NAK>
<STX>1...Data...<CR><ETX>xx <CR><LF>	>	
	<	<ACK>
<STX>2...Data...<CR><ETX>xx <CR><LF>	>	
	<	<ACK>
<EOT>	>	

### 3.2.4. ASTM Data Frame Format

A sequential number located after the <STX> character is inserted into each data frame. The frame number is set to 1 when the transfer phase is initialized and is incremented by 1 for each frame up to 7 and then returns to 0.

The frame number allows the receiver to distinguish new and re-transmitted frames. In case of re-transmitted frame (after a <NAK> response from the host), the frame number is not incremented: <STX>1...Data...<CR><ETX>xx<CR><LF>.

#### Frame format

ASTM field	Definition	Transmitted data	# of bytes	Comments
0	STX	\$02	1	
1	Frame number	1 to 7, 0, ...	1	Frame number is set to 1, incremented by 1 for each frame up to 7, and then returns to 0
2	Data message		240 max.	Header, Patient, Order, Result and Comment messages
3	End of data message ETX if end frame		1	
4	Checksum		2	
5	CRLF	\$0D \$0A	2	

#### Frame checksum

According to LIS01-A2, the frame checksum (<STX>1...Data...<CR><ETX>xx<CR><LF>) is defined as modulo 256 of ASCII values sum between <STX> not included and <ETX> included characters: 1...Data...<CR><ETX>.

### 3.3. Records General Format Specifications (LIS2-A2)

Data frames encapsulate records defined by the LIS-A2 norm, records themselves encapsulate ASTM fields.

Record ID	ASTM Definition
H	Header
P	Patient
O	Order
R	Result
C	Comment
Q	Query (Request information order)
M	Manufacturer information
L	Terminator record

#### 3.3.1. Structure of Records

##### Structure of records for order transmission

- H (Header)
- P (Patient)
- .. C (Patient Comments) optional
- .. O (Order)
- ... C (Order Comments) optional
- L (Terminator)

The transmission of an order without patient record, even if not allowed in LIS2-A2, is accepted and managed in the Yumizen H1500/H2500.

##### Instrument patient file modification by host

- H (Header)
- P (Patient)
- .. C (Patient Comments) optional
- L (Terminator)

##### Structure of records for result transmission

- H (Header)
- P (Patient)
- .. C (Patient Comments) optional
- .. O (Order)
- ... C (Order Comments) optional
- ... C (Run Alarms) optional
- ... M (Curves and Matrix points)
- ... M (Curves and Matrix points)
- ... M (Traceability)
- ... R (Result)
- .... C (Flag Result) optional
- ... R (Result)
- .... C (Flag Result) optional
- .....
- .....
- ... R (Result)
- .... C (Flag Result) optional
- L (Terminator)

##### Structure of records for statistics

- H (Header)
- .. M (Statistics)
- .. M (Statistics)
- L (Terminator)

##### Structure of records for remote command

- H (Header)
- .. M (Remote command)
- L (Terminator)

### 3.3.2. Description of Records

Only fields described with their specified length are used by HORIBA instruments.

The length of a field can be less than the maximum value but must not be more.

Delimiters must be used even if a field is empty.

Field inside records are separated by "|" (ASCII \$7C).

Component inside fields are separated by "^" (ASCII \$5E).

Repeated fields inside records are separated by "\".

#### 3.3.2.1. Alphanumeric Data

UTF-8 encoding is used for alphanumeric fields.

When alphanumeric data is sent, all the characters below 0x20 are replaced by an escape sequence with the following format: &Xhhh&.

"hhh" is the hexadecimal value of ASCII character completed with zero on 4 digits.

For example, <ETB> should be replaced by: <&X0017&>.

When alphanumeric data is received, the escape sequence &Xhhhh& is converted to the corresponding characters.

When alphanumeric data is transmitted, all delimiters characters they can contain must be replaced by their corresponding escape sequence as below:

Delimiter	Escape sequence
Field delimiter	&F&
Component delimiter	&S&
Repeat delimiter	&R&
Escape delimiter	&E&

#### 3.3.2.2. Records to Send

Fields that are not used are sent empty.

When sending records, the instrument sends only non-empty components, i.e. without component delimiters for the last empty components of the field.

#### 3.3.2.3. Received Records

If a field value, length, delimiter of a received record does not correspond to the required input type, the instrument generates an error log, and can ignore the record and its following ones (depending on the error and the message).

For more information about errors and messages, refer to the *Appendices > Error Management* chapter.

#### 3.3.2.4. Header Record

ASTM field	Definition	Transmitted data	Field max. length	Input type	Repeat delimiter	Mandatory
6.1	Record Type	H	1	Fixed	No	Yes
		Like standard:				
		<ul style="list-style-type: none"> <li>■   Field delimiter</li> <li>■ \ Repeat delimiter</li> <li>■ ^ Component delimiter</li> <li>■ &amp; Escape delimiter</li> </ul>				
6.2	Delimiters definition		4	Text	No	Yes
6.3	Message Control ID					
6.4	Access Password					
6.5	Sender Name (from instrument to host)	MHR1^SerialNumber^Software version	39 (15^12^10)	Fixed^Alphanumeric^Fixed	No	Yes
6.5	Sender Name (from host to instrument)	Host name	32	Alphanumeric	No	No

ASTM field	Definition	Transmitted data	Field max. length	Input type	Repeat delimiter	Mandatory
6.6	Sender Address					
6.7	Reserved					
6.8	Sender Telephone Nb					
6.9	Characteristics of Sender					
6.10	Receiver ID					
6.11	Comments or Special Instructions					
6.12	Processing ID	P: Patient message Q: Quality control message D: Technician	1	Fixed list	No	Yes
6.13	ASTM Version Nb	LIS2-A2	9	Fixed	No	Yes
6.14	Date and Time of message					

There should not be the field delimiter between fields 6.1 and 6.2 (as it is in the value of field 6.2).

In case of a response to a request (query, ...), the field 6.5 should be an exact copy from the field 6.10 sent in the request.

### 3.3.2.5. Patient Record

ASTM field	Definition	Transmitted data	Field max. length	Input type	Repeat delimiter	Mandatory
7.1	Record Type	P	1	Fixed	No	Yes
7.2	Sequence Number	1, 2, ...	2	Numeric	No	Yes

ASTM field	Definition	Transmitted data	Field max. length	Input type	Repeat delimiter	Mandatory
7.3	Practice Assigned Patient ID					
7.4	Laboratory Assigned Patient ID	Patient Id	25	Alphanumeric	No	No
7.5	Patient ID No 3					
7.6	Patient Name	Name^First name	41 (20^20)	Alphanumeric	No	No
7.7	Mother's Maiden Name	YYYYMMDD^A GE^U				
7.8	Birth date	<ul style="list-style-type: none"> <li>■ YYYYMMDD: Date of birth</li> <li>■ AGE: Patient age</li> <li>■ U: Unit of age</li> </ul>	14 (8^3^1)	Date^Numeric^ Fixed list	No	No
7.9	Patient Sex	M = Male F = Female U = Unknown	1	Fixed list	No	No
7.10	Patient Race-Ethnic Origin	Ethnic origin	20	Alphanumeric	No	No
7.11	Patient Address					
7.12	Reserved					
7.13	Patient Telephone Nb					
7.14	Attending Physician ID	PhysicianID^PhysicianName	20^30	Alphanumeric^ Alphanumeric	No	No^No
7.15	Special Field 1					
7.16	Special Field 2					
7.17	Patient Height					
7.18	Patient Weight					

ASTM field	Definition	Transmitted data	Field max. length	Input type	Repeat delimiter	Mandatory
7.19	Patient's Known or Suspected Diagnosis					
7.20	Patient Active Medication					
7.21	Patient's Diet					
7.22	Practice Field 1					
7.23	Practice Field 2					
7.24	Admission and Discharge Dates					
7.25	Admission Status					
7.26	Location	Service name	20	Alphanumeric	No	No
7.27	Nature of Alternative Diagnostic Code and Classifiers					
7.28	Nature of Alternative Diagnostic Code and Classifiers					
7.29	Patient Religion					
7.30	Martial status					
7.31	Isolation Status					
7.32	Language					
7.33	Hospital Service					
7.34	Hospital Institution					
7.35	Dosage Category	Patient type	20	Alphanumeric	No	No

The following units can be used for the patient age transmitted in the Birth date field (7.8).

- Y (year)
- M (month)
- W (week)
- D (day)
- H (hour)

A patient age transmitted in hours (from host) is converted into days when received on the instrument.

### 3.3.2.6. Order Record

ASTM field	Definition	Transmitted data	Field max. length	Input type	Repeat delimiter	Mandatory
8.1	Record Type	O	1	Fixed	No	Yes
8.2	Sequence Number	1, 2, ...	2	Numeric	No	Yes
8.3	Sample ID	SampleID^RackLoadingNb^RackBarcodeID^RackPosition <ul style="list-style-type: none"> <li>■ Sample ID</li> <li>■ Number of runs of the rack</li> <li>■ Rack ID (barcode)</li> <li>■ Position of the sample on the rack</li> </ul>	31 (16^2^8^2)	Alphanumeric^Numeric^Alphanumeric^Numeric	No	Yes (sample ID or rack information)
8.4	Instrument Specimen ID					

ASTM field	Definition	Transmitted data	Field max. length	Input type	Repeat delimiter	Mandatory
8.5	Universal Test ID	^^^TestName (CBC or CBR or DIF or DIR or DIF_LV or RET or CBF or RBC_PLTO or SLIDE_EC) <i>For more detail, refer to the Universal Test ID Definition chapter</i>	13 (^^^10)	Fixed list	Yes	From instrument to host: Yes From host to instrument: No
8.6	Priority	R: routine S: STAT	1	Fixed list	No	No
8.7	Requested/Ordered Date and Time					
8.8	Specimen Collection Date and Time	YYYYMMDDH HMMSS	14	Date and time	No	No
8.9	Collection End Time					
8.10	Collection Volume					
8.11	Collector ID					
8.12	Action Code (from instrument to host)	X: test already in progress	1	Fixed list	No	No
8.12	Action Code (from host to instrument)	N: new order C: cancel order	1	Fixed list	No	Yes
8.13	Danger Code					
8.14	Relevant Clinical Information					
8.15	Date/Time Specimen Received					
8.16	Specimen Descriptor	SpecimenType ^^SpecimenLiq uid	42 (20^^20)	Alphanumeric	No	No
8.17	Ordering Physician					

ASTM field	Definition	Transmitted data	Field max. length	Input type	Repeat delimiter	Mandatory
8.18	Physician Tel Nb					
8.19	User Field 1					
8.20	User Field 2					
8.21	Laboratory Field 1 (from instrument to host)	Prefix^Value Prefix: FT (free text), AP (analysis profile), SP (slide profile)	41 (10^30)	Fixed list^Alphanumeric	Yes	No
8.21	Laboratory Field 1 (from host to instrument)					
8.22	Laboratory Field 2					
8.23	Date and Time Results reported or last modified					
8.24	Instrument Charge to Computer System					
8.25	Instrument Section ID					
8.26	Report Types (from instrument to host)	F: final P: preliminary result X: order cannot be done I: in instrument pending	1	Fixed list	No	Yes

ASTM field	Definition	Transmitted data	Field max. length	Input type	Repeat delimiter	Mandatory
8.26	Report Types (from host to instrument)	Q: response to request information Z: no record for this patient Y: no test for this record O: response to result or download	1	Fixed list	No	Yes
8.27	Reserved					
8.28	Location or Ward of Specimen Collection					
8.29	Nosocomial Infection Flag					
8.30	Specimen Service					
8.31	Specimen institution					

The order must comply with the following conditions, otherwise the received order is ignored:

- Sample ID data of field 8.3 matches Sample ID data of the pending query.
- At least one of the Universal Test ID fields (8.5) contain a known TestName data.
- In case of several Universal Test ID fields (8.5), the order can contain one blood test, one slide test and one or several unknown tests. The order cannot contain several blood tests.
- The Action code field (8.12) from host to instrument is N (new order).

If a received order contains Universal Test ID fields (8.5) with only unknown TestName data, the Yumizen H1500/H2500 sends back the order with the Report Types field (8.26) set to X and the received order is ignored.

### Remaining Tests

Remaining tests must be sent to the instrument as long as they are not completed.

For example, the instrument sends a blood test result and a slide reflex is requested by the host. The host must send an order which contains a Universal Test ID field (8.5) with the SLIDE\_EC TestName.

When the instrument sends the slide test result, the host must send an order with the Action code field (8.12) set to N (new order) and the Report Types field (8.26) set to Y (no test for this record). This means there are no more tests to run.

### Option: Yumizen T6000

When the instrument is used with a Yumizen T6000, all tests must be sent to the instrument, including tests which cannot be done on this instrument.

For example, if a slide is requested, the order must contain a Universal Test ID field (8.5) with the SLIDE\_EC TestName even if the instrument is not connected to the Yumizen SPS.

### 3.3.2.7. Result Record

ASTM field	Definition	Transmitted data	Field max. length	Input type	Repeat delimiter	Mandatory
9.1	Record Type	R	1	Fixed	No	Yes
9.2	Sequence Number	1, 2, ...	2	Numeric	No	Yes
9.3	Universal Test ID (from instrument to host)	^^^English result name^LOINC Code associated with the result frame, if available	24 (^^^10^10)	^^^Open list^Open list	No	^^^Yes^No
9.3	Universal Test ID (from host to instrument)					

ASTM field	Definition	Transmitted data	Field max. length	Input type	Repeat delimiter	Mandatory
9.4	Data or Measurement Value	Test result or --- (invalid) or +++ (above visibility range) or slide status Slide status: <ul style="list-style-type: none"> <li>■ BLOOD_DROP PED</li> <li>■ SLIDE_IN_RA CK</li> <li>■ SLIDE_EJECT ED</li> </ul>	16	Alphanumeric	No	No
9.5	Unit or Set of units	Unit text (ISO 2955 or specific) or - (if no unit)	10	Open list	No	Yes
9.6	Reference Ranges					
9.7	Result Abnormal Flag (from instrument to host)	L: below reference limits LL: below critical limits H: above reference limits HH: above critical limits <: below limit of quantification >: above linearity range >>: above visibility range X: invalid A: reject N: normal	2	Fixed list	No	Yes
9.7	Result Abnormal Flag (from host to instrument)					

ASTM field	Definition	Transmitted data	Field max. length	Input type	Repeat delimiter	Mandatory
9.8	Nature of Abnormality Testing					
9.9	Result Status (from instrument to host)	F: final result or below limit of quantification or above linearity range W: warning (suspicion on validity) X: above visibility range or invalid or reject or process flag	1	Fixed list	No	Yes
9.9	Result Status (from host to instrument)					
9.10	Date of Change in Normative Values or Units					
9.11	Operator Identification (from instrument to host)	Login^^User profile or LastName FirstName^^User profile User profile: <ul style="list-style-type: none"> <li>■ TECHNICIAN</li> <li>■ MANAGER</li> <li>■ USER</li> </ul>	63 (41^^20 )	Alphanumeric	No	Yes
9.11	Operator Identification (from host to instrument)					
9.12	Date/Time Test Starting	YYYYMMDDHHM MSS	14	Date	No	Yes
9.13	Date/Time Test Completed					

ASTM field	Definition	Transmitted data	Field max. length	Input type	Repeat delimiter	Mandatory
9.14	Device Identification (from instrument to host)	9380BDED579C	15	Alphanumeric	No	Yes
9.14	Device Identification (from host to instrument)					

The test result and its unit are transmitted in the Data or Measurement Value field (9.4) and in the Unit or Set of units field (9.5).

- Decimal separator: "."
- Symbol "10^" (e.g. 10^9/L) is sent as "1E" (e.g. 1E09/L)
- Symbol "mm^3" (e.g. mm^3) is sent as "3" (e.g. mm3)
- Symbol "μ" (e.g. μL) is sent as "u" (e.g. uL)

### 3.3.2.8. Comment Record

ASTM field	Definition	Transmitted data	Field max. length	Input type	Repeat delimiter	Mandatory
10.1	Record Type	C	1	Fixed	No	Yes
10.2	Sequence Number	1, 2, ...	2	Numeric	No	Yes
10.3	Comment Source	P: Practice I: Clinical instrument system L: Computer system (Filler)	1	Fixed list	No	Yes

ASTM field	Definition	Transmitted data	Field max. length	Input type	Repeat delimiter	Mandatory
		For result comment (after R frame): alarm	200	Open list	Yes	Yes
10.4	Comment Text (from instrument to host)	For order comment (after O frame): AlarmType^MeasurementType^AlarmMain^AlarmDetail Refer to the table below	103 (20^20^30^30)	Open list	Yes	Yes^No^Yes^No
		For patient comment or sample comment (after P or O frame): free text	200	Alphanumeric	No	Yes
10.4	Comment Text (from host to instrument)	Comments	200	Alphanumeric	No	Yes
10.5	Comment Type	I: Instrument flag comment	1	Fixed list	No	Yes

Alarm type	Measurement type	Alarm main	Alarm detail
P (process)	Not applicable	Process flag category: <ul style="list-style-type: none"> <li>■ BAD_IDENTIFICATION</li> <li>■ BAD_EXTRACTION</li> <li>■ BAD_SAMPLING</li> <li>■ BAD_CONDITIONS</li> <li>■ BAD_CALCULATION</li> <li>■ STOP_REQUESTED</li> <li>■ SLIDE_ERROR</li> </ul>	Not applicable
D (device)	Measurement type	Alarm main message	Alarm detail message

Alarm type	Measurement type	Alarm main	Alarm detail
S (sample)	Measurement type	Alarm main message	Alarm detail message
C (calculation)	Considered channel	Alarm technical name	Not applicable

Analytical alarms are transmitted in the Comment Text field (10.4) as follows:

- A first comment with the alarm type D (if the analytical alarm is linked to an instrument technical problem) or S (if the analytical alarm is linked to the blood sample), the measurement type and the alarm message (main + detail).
- A second comment with the alarm type C, the considered channel and the alarm technical name.

Refer to the [Analytical Alarm Messages \(Device Alarm\)](#) and [Analytical Alarm Messages \(Sample Alarm\)](#) chapters for the list of analytical alarms and associated transmitted data.

**Example of a message sent by Yumizen H1500/H2500:**

4C|1|||S^DIFF^WBC\_ABN\_MAT^SEP\_NEU\_EOS|

5C|2||C^LMNE^NeuEosSep|

**3.3.2.9. Request Information Record**

ASTM field	Definition	Transmitted data	Field max. length	Input type	Repeat delimiter	Mandatory
11.1	Record Type	Q	1	Fixed	No	Yes
11.2	Sequence Number	1..99	2	Numeric	No	Yes

ASTM field	Definition	Transmitted data	Field max. length	Input type	Repeat delimiter	Mandatory
11.3	Starting Range ID Number	^SampleID^RackLoadingNb^RackBarcodeID^RackPosition ■ Sample ID ■ Number of runs of the rack ■ Rack ID (barcode) ■ Position of the sample on the rack	32 (^16^2^8^2)	^Alphanumeric^Numeric^Alphanumeric^Numeric	No	Yes (sample ID or rack information)
11.4	End Of identifier List					
11.5	Universal Test ID	ALL	3	Fixed	No	Yes
11.6	Time limits					
11.7	Time Max limits					
11.8	Time Min limits					
11.9	Physician Name					
11.10	Telephone Number					
11.11	Reserved for User					
11.12	Reserved for User					
11.13	Request Information Status Codes	O: request for test information	1	Fixed	No	Yes

### 3.3.2.10. Manufacturer Record

#### Curves and Matrices

ASTM field	Definition	Transmitted data	Field max. length	Input type	Repeat delimiter	Mandatory
14.1	Record Type	M	1	Fixed	No	Yes
14.2	Sequence Number	1..99	2	Numeric	No	Yes
14.3	Message Type	MATRIX, HISTOGRAM	10	Fixed	No	Yes
14.4	Measurement Type	For MATRIX: DIFF, PLTO, RET, BFDIFF For HISTOGRAM: DIFF, PLT, WBC, RBC, BFRBC, BFDIFF	20	Alphanumeric	No	Yes
14.5	Name	Graphic name	20	Alphanumeric	No	Yes
14.6	Thresholds	Threshold Encoding Type^Threshold Data FLOATLE-stream/deflate:base64		Alphanumeric^ Alphanumeric	Yes	No^Yes
14.7	Points	Graphic Encoding Type^Graphic Data FLOATLE-stream/deflate:base64		Alphanumeric^ Alphanumeric	Yes	No^Yes

The graphic names for matrices and histograms are transmitted in the Name field (14.5) as follows:

- For the matrices: LMNERESABS, RBCPLTRESABS, RETRESFLUO, MPRESABS
- For the histograms: BASOALONGGRES, EOSALONGABS, EOSALONGGRES, LYMALONGABS, LYMALONGGRES, MONALONGABS, MONALONGGRES, NEUALONGABS, NEUALONGGRES, PLTALONGGRES, TNCALONGGRES, RBCALONGGRES, RBCPLTALONGGRES, MNALONGABS, MNALONGGRES, PNALONGABS, PNALONGGRES

#### Traceability for Reagents

ASTM field	Definition	Transmitted data	Field max. length	Input type	Repeat delimiter	Mandatory
14.1	Record Type	M	1	Fixed	No	Yes
14.2	Sequence Number	1..99	2	Numeric	No	Yes
14.3	Message Type	REAGENT	10	Fixed	No	Yes
14.4	Traceability Name	Reagent name	20	Open list	Yes	Yes
14.5	Traceability Information	Lot number^Opening Date^Expiration Date	39 (15^14^8)	Alphanumeric^ Date and time^Date	Yes	Yes^No^Yes s

The reagent names are transmitted in the Traceability Name field (14.4) as follows:

- For the instrument: BASOLYSE, CLEANER, DILUENT, FLUOCYTE, NUCEDIFF, LYSEBIO
- For the slide system: BUFFER, METHANOL, STAIN1, STAIN2, DILUENT

#### Traceability for Quality Control

ASTM field	Definition	Transmitted data	Field max. length	Input type	Repeat delimiter	Mandatory
14.1	Record Type	M	1	Fixed	No	Yes
14.2	Sequence Number	1..99	2	Numeric	No	Yes

ASTM field	Definition	Transmitted data	Field max. length	Input type	Repeat delimiter	Mandatory
14.3	Message Type	QC	10	Fixed	No	Yes
14.4	Traceability Name	Quality control indicator name EXTQC, WESTGARD, TIME2QC	20	Open list	Yes	Yes
14.5	Traceability Information	INDETERMINE D, CORRECT, INCORRECT	20	Alphanumeric	Yes	Yes

### Traceability for Settings

ASTM field	Definition	Transmitted data	Field max. length	Input type	Repeat delimiter	Mandatory
14.1	Record Type	M	1	Fixed	No	Yes
14.2	Sequence Number	1..99	2	Numeric	No	Yes
14.3	Message Type	SETTING	10	Fixed	No	Yes
14.4	Traceability Name	Setting name RUO, WBCDIFF	20	Open list	Yes	Yes
14.5	Traceability Information	For RUO: TRUE, FALSE For WBCDIFF: 2, 5, 6	20	Alphanumeric	Yes	Yes

### Statistics

ASTM field	Definition	Transmitted data	Field max. length	Input type	Repeat delimiter	Mandatory
14.1	Record Type	M	1	Fixed	No	Yes
14.2	Sequence Number	1..99	2	Numeric	No	Yes

ASTM field	Definition	Transmitted data	Field max. length	Input type	Repeat delimiter	Mandatory
14.3	Message Type	STATS	10	Fixed	No	Yes
14.4	Statistics Information	Statistic Type^Start Date Time^End Date Time^Session Session: SESSION, UPTIME	81 (30^14^ 14^20)	Open list^Date and time^Date and time^Open list	No	Yes
14.5	Statistics Item	Name^Info^Val ue	72 (30^20^ 20)	Open list^Open list^Alphanu meric	Yes	Yes

The statistic types are transmitted in the Statistics Information field (14.4) as follows:

- SAMPLE\_ORDER
- RACK
- PATIENT
- CONTROL
- RERUN\_REFLEX
- REPEAT\_CALI
- INVALID\_RUN
- FAILURE
- REAGENTS
- TECHNICAL\_CYCLES

This table gives details of the names and information available for each type of statistic.

Statistic Type	Statistic Name	Statistic Information
SAMPLE_ORDER	<ul style="list-style-type: none"> <li>■ SamplingNumber</li> <li>■ AnalyticalSamplingStatistics</li> <li>■ SlideSamplingStatistics</li> </ul>	
RACK	<ul style="list-style-type: none"> <li>■ RackLoadedFromLoader</li> <li>■ RackEjectedOnTray</li> <li>■ RackLoadedFromConveyor</li> <li>■ RackEjectedOnConveyor</li> </ul>	

Statistic Type	Statistic Name	Statistic Information	
PATIENT	■ RunReportPatientCBC		
	■ RunReportPatientDIF		
	■ RunReportPatientRET		
	■ RunReportPatientCBR		
	■ RunReportPatientDIR		
	■ RunReportPatientCBF		
	■ RunReportPatientDIF_LV		
	■ RunReportPatientRBC_PLTO		
	■ RunReportPatientSLIDE_E		
	■ RunReportPatientSLIDE_EC		
	■ RunReportPatientSLIDE_C		
	CONTROL	■ RunReportControlCBC	
		■ RunReportControlDIF	
■ RunReportControlRET			
■ RunReportControlCBR			
■ RunReportControlDIR			
■ RunReportControlCBF			
■ RunReportControlDIF_LV			
■ RunReportControlRBC_PLTO			
RERUN_REFLEX	■ RerunPatient		
	■ ReflexPatient		
	■ RerunControl		
	■ ReflexControl		
REPEAT_CALI	■ RunRepeatabilityDIF		
	■ RunRepeatabilityDIR		
	■ RunCalibrationDIF		
INVALID_RUN	■ BadIdentification		
	■ BadExtraction		
	■ BadSampling		
	■ BadCondition		
	■ BadCalculation		
	■ StopRequested		
	■ BadSmearStain		
FAILURE	■ FailureNumber	■ CIM ■ INST ■ SPS	

Statistic Type	Statistic Name	Statistic Information	
REAGENTS	■ Volume	■ BASOLYSE ■ BUFFER ■ CLEANER ■ DILUENT ■ FLUOCYTE ■ NUCEDIFF ■ LYSEBIO ■ METHANOL ■ STAIN1 ■ STAIN2 ■ RIBBON	
	TECHNICAL_CYCLES	■ AnalyzerStartup	
		■ AnalyzerShutdown	
		■ AnalyzerRecovery	
		■ AnalyzerBgNoiseCheck	
		■ AnalyzerCleaner	
		■ AnalyzerMinoclair	
		■ CIMStartup	
		■ CIMShutdown	
		■ CIMRecovery	
		■ SPSSStartupDaily	
■ SPSSStartupWeekly			
■ SPSSShutdownDaily			
■ SPSSShutdownWeekly			
■ SPSRecovery			

### Remote Command

ASTM field	Definition	Transmitted data	Field max. length	Input type	Repeat delimiter	Mandatory
14.1	Record Type	M	1	Fixed	No	Yes
14.2	Sequence Number	1..99	2	Numeric	No	Yes
14.3	Message Type	EXECUTE	10	Fixed	No	Yes

ASTM field	Definition	Transmitted data	Field max. length	Input type	Repeat delimiter	Mandatory
14.4	Command	Command Type^Command Name^Command ID Command Type: QC Command Name: EXTQC_INCORRECT, EXTQC_CORRECT	42 (10^20^10)	Open list^Open list^Alphanumeric	No	Yes
14.5	Parameter					
14.6	Status (from instrument to host)	Comment status: DONE, FAILED, IGNORED, ACCEPTED, UNKNOWN_COMMAND, UNKNOWN_PARAMETER	20	Open list	No	Yes
14.6	Status (from host to instrument)					

The host sends a remote command to indicate the external quality control status (correct or incorrect). The Command field (14.4) is: EXTQC\_CORRECT or EXTQC\_INCORRECT.

The instrument sends a first answer to the host to indicate if the external quality control status update is accepted or not. The Status field (14.6) is:

- ACCEPTED
- IGNORED: the command is ignored when the Processing ID field (6.12) in the header record is not P
- UNKNOWN\_COMMAND: unknown Command Type
- UNKNOWN\_PARAM: unknown Command Name

The instrument updates the external quality control status and sends a second answer to the host to indicate if the update has been successful or not. The Status field (14.6) is:

- DONE
- FAILED

### 3.3.2.11. Terminal Record

ASTM field	Definition	Transmitted data	Field max. length	Input type	Repeat delimiter	Mandatory
12.1	Record Type	L	1	Fixed	No	Yes
12.2	Sequence Number	1	1	Fixed	No	Yes
12.3	Termination code (from instrument to host)	N: Normal	1	Fixed	No	Yes
12.3	Termination code (from host to instrument)					

## 3.4. Special Characteristics for HORIBA Data

### 3.4.1. Data Presentation

The parameters correspond to the Universal Test ID field (9.3) and the units correspond to the Unit or Set of units field (9.5).

### 3.4.1.1. Parameters

LOINC Code: Logical Observation Identifiers Names & Codes

CBC Codes	LOINC Code	Definition
<b>RBC</b>	789-8	Red Blood Cells
<b>HGB</b>	718-7	Hemoglobin Concentration
<b>HCT</b>	4544-3	Hematocrit
<b>MCV</b>	787-2	Mean Corpuscular Volume
<b>MCH</b>	785-6	Mean Corpuscular Hemoglobin
<b>MCHC</b>	786-4	Mean Corpuscular Hemoglobin Concentration
<b>RDW-SD</b>	21000-5	Red Distribution Width Standard Deviation
<b>RDW-CV</b>	788-0	Red Distribution Width
<b>MIC</b>	X-MIC	Microcytic Red Blood Cells percentage (versus RBC)
<b>MAC</b>	X-MAC	Macrocytic Red Blood Cells percentage (versus RBC)
<b>PLT</b>	777-3	Platelets
<b>PCT</b>	51637-7	Plateletcrit
<b>PDW</b>	51631-0	Platelets Distribution Width
<b>MPV</b>	32623-1	Mean Platelet Volume
<b>P-LCC</b>	96354-6	Platelets - Large Cell Count
<b>P-LCR</b>	48386-7	Platelets - Large Cell Ratio
<b>WBC</b>	6690-2	White Blood Cells
<b>NRBC#</b>	771-6	Nucleated Red Blood Cells absolute value
<b>NRBC%</b>	58413-6	Nucleated Red Blood Cells percentage
<b>TNC</b>	50774-9	Total Nucleated cells

PLT opt Codes	LOINC Code	Definition
<b>PLT-Ox</b>	97995-5	Platelets from optical channel
<b>LPF</b>	97994-8	Large Platelet Fraction

DIF Codes	LOINC Code	Definition
<b>LYM#</b>	731-0	Lymphocytes absolute value
<b>LYM%</b>	736-9	Lymphocytes percentage
<b>MON#</b>	742-7	Monocytes absolute value
<b>MON%</b>	5905-5	Monocytes percentage
<b>NEU#</b>	751-8	Neutrophils absolute value
<b>NEU%</b>	770-8	Neutrophils percentage
<b>EOS#</b>	711-2	Eosinophils absolute value
<b>EOS%</b>	713-8	Eosinophils percentage
<b>BAS#</b>	704-7	Basophils absolute value
<b>BAS%</b>	706-2	Basophils percentage
<b>IMG#</b>	53115-2	Immature Granulocytic cells absolute value
<b>IMG%</b>	71695-1	Immature Granulocytic cells percentage
<b>IMM#</b>	X-IMM#	Immature Monocytic cells absolute value
<b>IMM%</b>	X-IMM%	Immature Monocytic cells percentage
<b>IML#</b>	X-IML#	Immature Lymphocytic cells absolute value
<b>IML%</b>	X-IML%	Immature Lymphocytic cells percentage
<b>ALY#</b>	43743-4	Atypical Lymphocytes absolute value
<b>ALY%</b>	42250-1	Atypical Lymphocytes percentage
<b>LIC#</b>	55432-9	Large Immature Cells absolute value
<b>LIC%</b>	55433-7	Large Immature Cells percentage

RET Codes	LOINC Code	Definition
<b>RET#</b>	14196-0	Reticulocytes absolute value
<b>RET%</b>	17849-1	Reticulocytes percentage
<b>RET-L</b>	X-RET-L	Reticulocytes with a low RNA content
<b>RET-M</b>	X-RET-M	Reticulocytes with a medium RNA content
<b>RET-H</b>	X-RET-H	Reticulocytes with a high RNA content
<b>CRC</b>	X-CRC	Corrected Reticulocyte Count
<b>MRV</b>	48706-6	Mean Reticulocyte Volume

RET Codes	LOINC Code	Definition
<b>RHCC</b>	X-RHCC	Reticulocyte Hemoglobin Cellular Content
<b>IRF</b>	X-IRF	Immature Reticulocyte Fraction

CBF Codes	LOINC Code	Definition
<b>BFRBC</b>	X-BFRBC	Red Blood Cells absolute value
<b>BFWBC</b>	X-BFWBC	White Blood Cells absolute value
<b>BFMN#</b>	X-BFMN#	Mononuclear absolute value
<b>BFMN%</b>	X-BFMN%	Mononuclear percentage
<b>BFPN#</b>	X-BFPN#	Polymorphonuclear absolute value
<b>BFPN%</b>	X-BFPN%	Polymorphonuclear percentage

### Technical Parameters

These parameters are not used for screening or diagnostic purposes. They may be used by technicians for troubleshooting.

CBC Codes	LOINC Code	Definition
<b>TNCHGB</b>	X-TNCHGB	WBC count in Chamber 2 TNC/HGB
<b>TNCDIF</b>	X-TNCDIF	WBC count in Chamber 4 LMNE/NRBC
<b>TNCBAS</b>	X-TNCBAS	WBC count in Chamber 3 BASO/TNC2

RET Codes	LOINC Code	Definition
<b>MFI</b>	X-MFI	Mean Fluorescence Index
<b>PIC</b>	X-PIC	Peak Index Channel (Fluorescence range)

CBF Codes	LOINC Code	Definition
<b>BFTNCHGB</b>	X-BFTNCHGB	WBC count in Chamber 2 TNC/HGB
<b>BFTNCDIF</b>	X-BFTNCDIF	WBC count in Chamber 4 LMNE/NRBC

### 3.4.1.2. Universal Test ID Definition

Universal Test ID	Sample type	Parameter list
CBC	Whole blood	RBC / HGB / HCT / MCV / MCH / MCHC / RDW-SD / RDW-CV / MIC / MAC / PLT / PCT / PDW / MPV / P-LCC / P-LCR / WBC / NRBC# / NRBC% / TNC / TNCHGB / TNCDIF / TNCBAS
DIF	Whole blood	CBC + LYM# / LYM% / MON# / MON% / NEU# / NEU% / EOS# / EOS% / BAS# / BAS% / IMG# / IMG% / IMM# / IMM% / IML# / IML% / ALY# / ALY% / LIC# / LIC%
DIF_LV	Whole blood	DIF - Low value
RET	Whole blood	RBC / RET# / RET% / RET-L / RET-M / RET-H / CRC / MRV / MFI / PIC / RHCC / IRF
CBR	Whole blood	CBC + RET
DIR	Whole blood	DIF + RET
RBC_PLTO	Whole blood	RBC / HGB / HCT / MCV / MCH / MCHC / RDW-SD / RDW-CV / MIC / MAC / PLT / PCT / PDW / MPV / P-LCC / P-LCR / PLT-Ox / LPF
CBF	Body fluids	BFRBC / BFWBC / BFTNCHGB / BFTNCDIF / BFMN# / BFMN% / BFPN# / BFPN%
SLIDE_EC	Other	Slide

### 3.4.1.3. Units

Default unit system: SI (international)

Codes: CBC	SI (international)	Conventional	mmol/L	Japan	China
<b>RBC</b>	10 <sup>12</sup> /L	10 <sup>6</sup> /mm <sup>3</sup>	10 <sup>12</sup> /L	10 <sup>4</sup> /μL	10 <sup>12</sup> /L
<b>HGB</b>	g/L	g/dL	mmol/L	g/dL	g/L
<b>HCT</b>	L/L	%	L/L	%	%
<b>MCV</b>	fL	fL	fL	fL	fL
<b>MCH</b>	pg	pg	fmol	pg	pg
<b>MCHC</b>	g/L	g/dL	mmol/L	g/dL	g/L

Codes: CBC	SI (international)	Conventional	mmol/L	Japan	China
RDW-SD	fL	fL	fL	fL	fL
RDW-CV	%	%	%	%	%
MIC	%	%	%	%	%
MAC	%	%	%	%	%
PLT	10 <sup>9</sup> /L	10 <sup>3</sup> /mm <sup>3</sup>	10 <sup>9</sup> /L	10 <sup>4</sup> /μL	10 <sup>9</sup> /L
PCT	%	%	%	%	%
PDW	fL	fL	fL	fL	fL
MPV	fL	fL	fL	fL	fL
P-LCC	10 <sup>9</sup> /L	10 <sup>3</sup> /mm <sup>3</sup>	10 <sup>9</sup> /L	10 <sup>4</sup> /μL	10 <sup>9</sup> /L
P-LCR	%	%	%	%	%
WBC	10 <sup>9</sup> /L	10 <sup>3</sup> /mm <sup>3</sup>	10 <sup>9</sup> /L	10 <sup>2</sup> /μL	10 <sup>9</sup> /L
NRBC#	10 <sup>9</sup> /L	10 <sup>3</sup> /mm <sup>3</sup>	10 <sup>9</sup> /L	10 <sup>2</sup> /μL	10 <sup>9</sup> /L
NRBC%	%	%	%	%	%
TNC	10 <sup>9</sup> /L	10 <sup>3</sup> /mm <sup>3</sup>	10 <sup>9</sup> /L	10 <sup>2</sup> /μL	10 <sup>9</sup> /L

Codes: PLT opt	SI (international)	Conventional	mmol/L	Japan	China
PLT-Ox	10 <sup>9</sup> /L	10 <sup>3</sup> /mm <sup>3</sup>	10 <sup>9</sup> /L	10 <sup>4</sup> /μL	10 <sup>9</sup> /L
LPF	%	%	%	%	%

Codes: DIF	SI (international)	Conventional	mmol/L	Japan	China
LYM#	10 <sup>9</sup> /L	10 <sup>3</sup> /mm <sup>3</sup>	10 <sup>9</sup> /L	10 <sup>2</sup> /μL	10 <sup>9</sup> /L
LYM%	%	%	%	%	%
MON#	10 <sup>9</sup> /L	10 <sup>3</sup> /mm <sup>3</sup>	10 <sup>9</sup> /L	10 <sup>2</sup> /μL	10 <sup>9</sup> /L
MON%	%	%	%	%	%
NEU#	10 <sup>9</sup> /L	10 <sup>3</sup> /mm <sup>3</sup>	10 <sup>9</sup> /L	10 <sup>2</sup> /μL	10 <sup>9</sup> /L
NEU%	%	%	%	%	%
EOS#	10 <sup>9</sup> /L	10 <sup>3</sup> /mm <sup>3</sup>	10 <sup>9</sup> /L	10 <sup>2</sup> /μL	10 <sup>9</sup> /L
EOS%	%	%	%	%	%
BAS#	10 <sup>9</sup> /L	10 <sup>3</sup> /mm <sup>3</sup>	10 <sup>9</sup> /L	10 <sup>2</sup> /μL	10 <sup>9</sup> /L

Codes: DIF	SI (international)	Conventional	mmol/L	Japan	China
BAS%	%	%	%	%	%
IMG#	10 <sup>9</sup> /L	10 <sup>3</sup> /mm <sup>3</sup>	10 <sup>9</sup> /L	10 <sup>2</sup> /μL	10 <sup>9</sup> /L
IMG%	%	%	%	%	%
IMM#	10 <sup>9</sup> /L	10 <sup>3</sup> /mm <sup>3</sup>	10 <sup>9</sup> /L	10 <sup>2</sup> /μL	10 <sup>9</sup> /L
IMM%	%	%	%	%	%
IML#	10 <sup>9</sup> /L	10 <sup>3</sup> /mm <sup>3</sup>	10 <sup>9</sup> /L	10 <sup>2</sup> /μL	10 <sup>9</sup> /L
IML%	%	%	%	%	%
ALY#	10 <sup>9</sup> /L	10 <sup>3</sup> /mm <sup>3</sup>	10 <sup>9</sup> /L	10 <sup>2</sup> /μL	10 <sup>9</sup> /L
ALY%	%	%	%	%	%
LIC#	10 <sup>9</sup> /L	10 <sup>3</sup> /mm <sup>3</sup>	10 <sup>9</sup> /L	10 <sup>2</sup> /μL	10 <sup>9</sup> /L
LIC%	%	%	%	%	%

Codes: RET	SI (international)	Conventional	mmol/L	Japan	China
RET#	10 <sup>9</sup> /L	10 <sup>6</sup> /mm <sup>3</sup>	10 <sup>9</sup> /L	10 <sup>4</sup> /μL	10 <sup>9</sup> /L
RET%	%	%	%	%	%
RET-L	%	%	%	%	%
RET-M	%	%	%	%	%
RET-H	%	%	%	%	%
CRC	%	%	%	%	%
MRV	fL	fL	fL	fL	fL
RHCC	pg	pg	fmol	pg	pg
IRF	ratio	ratio	ratio	ratio	ratio

Codes: CBF	SI (international)	Conventional	mmol/L	Japan	China
BFRBC	10 <sup>6</sup> /L	10 <sup>6</sup> /L	10 <sup>6</sup> /L	10 <sup>4</sup> /μL	10 <sup>6</sup> /L
BFWBC	10 <sup>6</sup> /L	10 <sup>6</sup> /L	10 <sup>6</sup> /L	10 <sup>2</sup> /μL	10 <sup>6</sup> /L
BFMN#	10 <sup>6</sup> /L	10 <sup>6</sup> /L	10 <sup>6</sup> /L	10 <sup>2</sup> /μL	10 <sup>6</sup> /L
BFMN%	%	%	%	%	%

Codes: CBF	SI (international)	Conventional	mmol/L	Japan	China
BFPN#	10 <sup>6</sup> /L	10 <sup>6</sup> /L	10 <sup>6</sup> /L	10 <sup>2</sup> /μL	10 <sup>6</sup> /L
BFPN%	%	%	%	%	%

### Technical Parameters

Codes: CBC	SI (international)	Conventional	mmol/L	Japan	China
TNCHGB	10 <sup>9</sup> /L	10 <sup>3</sup> /mm <sup>3</sup>	10 <sup>9</sup> /L	10 <sup>2</sup> /μL	10 <sup>9</sup> /L
TNCDIF	10 <sup>9</sup> /L	10 <sup>3</sup> /mm <sup>3</sup>	10 <sup>9</sup> /L	10 <sup>2</sup> /μL	10 <sup>9</sup> /L
TNCBAS	10 <sup>9</sup> /L	10 <sup>3</sup> /mm <sup>3</sup>	10 <sup>9</sup> /L	10 <sup>2</sup> /μL	10 <sup>9</sup> /L

Codes: RET	SI (international)	Conventional	mmol/L	Japan	China
MFI	ratio	ratio	ratio	ratio	ratio
PIC	-	-	-	-	-

Codes: CBF	SI (international)	Conventional	mmol/L	Japan	China
BFTNCHGB	10 <sup>6</sup> /L	10 <sup>6</sup> /L	10 <sup>6</sup> /L	10 <sup>2</sup> /μL	10 <sup>6</sup> /L
BFTNCDIF	10 <sup>6</sup> /L	10 <sup>6</sup> /L	10 <sup>6</sup> /L	10 <sup>2</sup> /μL	10 <sup>6</sup> /L

#### 3.4.1.4. Slide Status

The slide status is transmitted in the Data or Measurement Value field (9.4).

- BLOOD\_DROPPED
- SLIDE\_IN\_RACK
- SLIDE\_EJECTED

## 3.4.2. Alarms and Pathologies

### 3.4.2.1. Suspicion and Reject

When a result is suspected of being abnormal or false, it is not reliable and the instrument returns a flag in field 9.9.

Refer to [Description of Records](#).

### 3.4.2.2. Reference and Critical Ranges

Flags when result exceeds reference or critical ranges are transmitted through field 9.7, they should be compared, to obtain a full result information, to the ranges set by the user.

Refer to [Description of Records](#).

### 3.4.2.3. Analytical Alarm Messages (Device Alarm)

Analytical alarms are transmitted in the Comment Text field (10.4) as follows:

- A first comment with the alarm type D (if the analytical alarm is linked to an instrument technical problem) or S (if the analytical alarm is linked to the blood sample), the measurement type and the alarm message (main + detail).
- A second comment with the alarm type C, the considered channel and the alarm technical name.

This chapter lists device alarms in separate tables for each measurement type.

- The first column indicates the transmitted data in the first comment.
- The second column indicates the transmitted data in the second comment.



Some non-listed alarms type C may only be transmitted in the second comment and can be ignored.

**WBC**

Transmitted data Alarm type, measurement type, alarm message	Transmitted data Alarm type, considered channel, alarm technical name	Description
D^WBC^ANA_ERR^UNST_TNC	C^SYNTHESIS_WBC^UnstbResCountTnc1	WBC Analytical error TNC channel unstable measure
D^WBC^ANA_ERR^UNST_TNC	C^SYNTHESIS_WBC^UnstbWidthCumulTnc1	WBC Analytical error TNC channel unstable measure
D^WBC^ANA_ERR^CLOG_TNC	C^TNC^NoEvents	WBC Analytical error TNC channel clog?
D^WBC^ANA_ERR^CLOG_TNC	C^TNC^LowEvents	WBC Analytical error TNC channel clog?
D^WBC^ANA_ERR^UNST_BASO	C^BASO^UnstbResCount	WBC Analytical error BASO channel unstable measure
D^WBC^ANA_ERR^UNST_BASO	C^BASO^UnstbWidthCumul	WBC Analytical error BASO channel unstable measure
D^WBC^ANA_ERR^CLOG_BASO	C^BASO^NoEvents	WBC Analytical error BASO channel clog?
D^WBC^ANA_ERR^CLOG_BASO	C^BASO^LowEvents	WBC Analytical error BASO channel clog?
D^WBC^ANA_ERR^UNST_DIFF	C^LMNE^UnstbResCount	WBC Analytical error DIFF channel unstable measure
D^WBC^ANA_ERR^UNST_DIFF	C^SYNTHESIS_WBC^UnstbResCountLmne1	WBC Analytical error DIFF channel unstable measure
D^WBC^ANA_ERR^LIGHT_BEAM_DIFF	C^LMNE^LightShift	WBC Analytical error DIFF channel light beam error
D^WBC^ANA_ERR^CLOG_DIFF	C^LMNE^NoEvents	WBC Analytical error DIFF channel clog?
D^WBC^ANA_ERR^CLOG_DIFF	C^LMNE^LowEvents	WBC Analytical error DIFF channel clog?

Transmitted data Alarm type, measurement type, alarm message	Transmitted data Alarm type, considered channel, alarm technical name	Description
D^WBC^ANA_ERR^UNST_DIFF	C^LMNE^LowCorrel	WBC Analytical error DIFF channel unstable measure
D^WBC^ANA_ERR^UNST_DIFF	C^LMNE^LowBackgndNoiseExcess	WBC Analytical error DIFF channel unstable measure
D^WBC^ANA_ERR^BUBBLE_DIFF	C^SYNTHESIS_WBC^LowAbsCorrel	WBC Analytical error DIFF channel bubble?
D^WBC^ANA_ERR^UNST_DIFF	C^SYNTHESIS_WBC^LmneNoWbc	WBC Analytical error DIFF channel unstable measure
D^WBC^ANA_ERR^UNST_TNC	C^SYNTHESIS_WBC^TncNoWbc	WBC Analytical error TNC channel unstable measure
D^WBC^ANA_ERR^UNBAL_TNC_DIFF_BASO	C^SYNTHESIS_WBC^BasoLmneBalance2	WBC Analytical error TNC/DIFF/BASO channels unbalanced
D^WBC^ANA_ERR^UNBAL_TNC	C^SYNTHESIS_WBC^TncBasoBalance2	WBC Analytical error TNC channel unbalanced
D^WBC^ANA_ERR^UNBAL_DIFF	C^SYNTHESIS_WBC^TncLmneBalance2	WBC Analytical error DIFF channel unbalanced

**DIFF**

Transmitted data Alarm type, measurement type, alarm message	Transmitted data Alarm type, considered channel, alarm technical name	Description
D^DIFF^ANA_ERR^UNST_BASO	C^SYNTHESIS_WBC^BasoNoWbc	WBC Analytical error BASO channel unstable measure
D^DIFF^ANA_ERR^BAL_TNC_BASO	C^SYNTHESIS_WBC^BasoExcess	WBC Analytical error TNC/BASO channel balance

## RBC

Transmitted data Alarm type, measurement type, alarm message	Transmitted data Alarm type, considered channel, alarm technical name	Description
D^RBC^ANA_ERR^UNST_RBC	C^RBC^UnstbResCount	RBC Analytical error RBC channel unstable measure
D^RBC^ANA_ERR^CLOG_RBC	C^RBC^NoEvents	RBC Analytical error RBC channel clog?
D^RBC^ANA_ERR^CLOG_RBC	C^RBC^LowEvents	RBC Analytical error RBC channel clog?
D^RBC^ANA_ERR^UNST_RBC	C^RBC^Noise	RBC Analytical error RBC channel unstable measure
D^RBC^ANA_ERR^UNST_RBC	C^RBC^AbnormWidth	RBC Analytical error RBC channel unstable measure
D^RBC^ANA_ERR^CLOG_RBC	C^RBC^RbcMiss	RBC Analytical error RBC channel clog?
D^RBC^ANA_ERR^BAL_RBC_H GB	C^SYNTHESIS_RBC_PLT^Abno rmMch	RBC Analytical error RBC/HGB channel balance
D^RBC^ANA_ERR^BAL_RBC_H GB	C^SYNTHESIS_RBC_PLT^Abno rmMchc	RBC Analytical error RBC/HGB channel balance

## RBC\_OTH

Transmitted data Alarm type, measurement type, alarm message	Transmitted data Alarm type, considered channel, alarm technical name	Description
D^RBC_OTH^ANA_ERR^INTER F_HGB	C^SYNTHESIS_WBC^TncLmne Balance3	HGB Analytical error HGB channel interference
D^RBC_OTH^ANA_ERR^UNST_ HGB	C^HGB^UnstbMeasure	HGB Analytical error HGB channel unstable measure

Transmitted data Alarm type, measurement type, alarm message	Transmitted data Alarm type, considered channel, alarm technical name	Description
D^RBC_OTH^ANA_ERR^UNST_ HGB	C^HGB^UnstbBlank	HGB Analytical error HGB channel unstable measure
D^RBC_OTH^ANA_ERR^UNST_ HGB	C^HGB^BlankShift	HGB Analytical error HGB channel unstable measure

## PLT

Transmitted data Alarm type, measurement type, alarm message	Transmitted data Alarm type, considered channel, alarm technical name	Description
D^PLT^ANA_ERR^UNST_PLT	C^PLT^UnstbResCount	PLT Analytical error PLT channel unstable measure
D^PLT^ANA_ERR^CLOG_PLT	C^PLT^NoEvents	PLT Analytical error PLT channel clog?
D^PLT^ANA_ERR^CLOG_PLT	C^PLT^LowEvents	PLT Analytical error PLT channel clog?
D^PLT^ANA_ERR^UNST_PLT	C^PLT^Noise	PLT Analytical error PLT channel unstable measure
D^PLT^ANA_ERR^UNST_PLT	C^PLT^AbnormWidth	PLT Analytical error PLT channel unstable measure
D^PLT^ANA_ERR^UNST_PLT	C^PLT^NoiseEvents	PLT Analytical error PLT channel unstable measure
D^PLT^ANA_ERR^CLOG_PLT	C^PLT^PitMiss	PLT Analytical error PLT channel clog?

## PLTO

Transmitted data Alarm type, measurement type, alarm message	Transmitted data Alarm type, considered channel, alarm technical name	Description
D^PLTO^ANA_ERR^UNST_RBC	C^RBC^UnstbResCount	RBC Analytical error RBC channel unstable measure
D^PLTO^ANA_ERR^CLOG_RB C	C^RBC^NoEvents	RBC Analytical error RBC channel clog?
D^PLTO^ANA_ERR^CLOG_RB C	C^RBC^LowEvents	RBC Analytical error RBC channel clog?
D^PLTO^ANA_ERR^UNST_RBC	C^RBC^Noise	RBC Analytical error RBC channel unstable measure
D^PLTO^ANA_ERR^UNST_RBC	C^RBC^AbnormWidth	RBC Analytical error RBC channel unstable measure
D^PLTO^ANA_ERR^CLOG_RB C	C^RBC^RbcMiss	RBC Analytical error RBC channel clog?
D^PLTO^ANA_ERR^BAL_RBC_ HGB	C^SYNTHESIS_RBC_PLT^Abno rmMch	RBC Analytical error RBC/HGB channel balance
D^PLTO^ANA_ERR^BAL_RBC_ HGB	C^SYNTHESIS_RBC_PLT^Abno rmMchc	RBC Analytical error RBC/HGB channel balance
D^PLTO^ANA_ERR^UNST_PLT O	C^RBC_PLTO^UnstbResCount	PLTO Analytical error PLTO channel unstable measure
D^PLTO^ANA_ERR^LIGHT_BE AM_PLTO	C^RBC_PLTO^LightShift	PLTO Analytical error PLTO channel light beam error
D^PLTO^ANA_ERR^CLOG_PLT O	C^RBC_PLTO^NoEvents	PLTO Analytical error PLTO channel clog?
D^PLTO^ANA_ERR^CLOG_PLT O	C^RBC_PLTO^LowEvents	PLTO Analytical error PLTO channel clog?
D^PLTO^ANA_ERR^UNST_PLT O	C^RBC_PLTO^LowCorrel	PLTO Analytical error PLTO channel unstable measure

Transmitted data Alarm type, measurement type, alarm message	Transmitted data Alarm type, considered channel, alarm technical name	Description
D^PLTO^ANA_ERR^UNST_PLT O	C^RBC_PLTO^NoiseEvents	PLTO Analytical error PLTO channel unstable measure
D^PLTO^ANA_ERR^UNST_PLT O	C^RBC_PLTO^UnstbPltoCount	PLTO Analytical error PLTO channel unstable measure
D^PLTO^ANA_ERR^BAL_PLTO _RBC	C^SYNTHESIS_RBC_PLT^RbcR bcoBalance	PLTO Analytical error PLTO/RBC channel balance

## RET

Transmitted data Alarm type, measurement type, alarm message	Transmitted data Alarm type, considered channel, alarm technical name	Description
D^RET^ANA_ERR^UNST_RET	C^RET^UnstbResCount	RET Analytical error RET channel unstable measure
D^RET^ANA_ERR^LIGHT_BE M_RET	C^RET^LightShift	RET Analytical error RET channel light beam error
D^RET^ANA_ERR^CLOG_RET	C^RET^NoEvents	RET Analytical error RET channel clog?
D^RET^ANA_ERR^CLOG_RET	C^RET^LowEvents	RET Analytical error RET channel clog?
D^RET^ANA_ERR^UNST_RET	C^RET^LowCorrel	RET Analytical error RET channel unstable measure
D^RET^ANA_ERR^CLOG_RET	C^RET^RbcMiss	RET Analytical error RET channel clog?
D^RET^ANA_ERR^CLOG_RET	C^RET^RetMiss	RET Analytical error RET channel clog?

Transmitted data Alarm type, measurement type, alarm message	Transmitted data Alarm type, considered channel, alarm technical name	Description
D^RET^ANA_ERR^UNST_RET	C^SYNTHESIS_RBC_PLT^NcRetInterf	RET Analytical error RET channel unstable measure
D^RET^ANA_ERR^BAL_RBC_RET	C^SYNTHESIS_RBC_PLT^RbcRetRbcBalance	RET Analytical error RBC/RET channel balance

### BFWBC

Transmitted data Alarm type, measurement type, alarm message	Transmitted data Alarm type, considered channel, alarm technical name	Description
D^BFWBC^ANA_ERR^UNST_TNC	C^BF_TNC^UnstbResCount	WBC Analytical error TNC channel unstable measure
D^BFWBC^ANA_ERR^UNST_TNC	C^BF_TNC^UnstbWidthCumul	WBC Analytical error TNC channel unstable measure
D^BFWBC^ANA_ERR^CLOG_TNC	C^BF_TNC^NoEvents	WBC Analytical error TNC channel clog?
D^BFWBC^ANA_ERR^CLOG_TNC	C^BF_TNC^LowEvents	WBC Analytical error TNC channel clog?
D^BFWBC^ANA_ERR^UNST_DIFF	C^BF_LMNE^LowBackgndNoiseExcess	WBC Analytical error DIFF channel unstable measure
D^BFWBC^ANA_ERR^BAL_DIFF_TNC	C^SYNTHESIS_BF^BfTncLmneBalance	WBC Analytical error DIFF/TNC channel balance

### BFDIFF

Transmitted data Alarm type, measurement type, alarm message	Transmitted data Alarm type, considered channel, alarm technical name	Description
D^BFDIFF^ANA_ERR^UNST_DIFF	C^BF_LMNE^UnstbResCount	WBC Analytical error DIFF channel unstable measure
D^BFDIFF^ANA_ERR^LIGHT_BEAM_DIFF	C^BF_LMNE^LightShift	WBC Analytical error DIFF channel light beam error
D^BFDIFF^ANA_ERR^CLOG_DIFF	C^BF_LMNE^NoEvents	WBC Analytical error DIFF channel clog?
D^BFDIFF^ANA_ERR^CLOG_DIFF	C^BF_LMNE^LowEvents	WBC Analytical error DIFF channel clog?
D^BFDIFF^ANA_ERR^UNST_DIFF	C^BF_LMNE^LowCorrel	WBC Analytical error DIFF channel unstable measure

### BFRBC

Transmitted data Alarm type, measurement type, alarm message	Transmitted data Alarm type, considered channel, alarm technical name	Description
D^BFRBC^ANA_ERR^UNST_RBC	C^BF_RBC^UnstbResCount	RBC Analytical error RBC channel unstable measure
D^BFRBC^ANA_ERR^UNST_RBC	C^SYNTHESIS_BF^NoiseBfRbc	RBC Analytical error RBC channel unstable measure
D^BFRBC^ANA_ERR^UNST_RBC	C^SYNTHESIS_BF^AbnormWidthBfRbc	RBC Analytical error RBC channel unstable measure

### 3.4.2.4. Analytical Alarm Messages (Sample Alarm)

Analytical alarms are transmitted in the Comment Text field (10.4) as follows:

- A first comment with the alarm type D (if the analytical alarm is linked to an instrument technical problem) or S (if the analytical alarm is linked to the blood sample), the measurement type and the alarm message (main + detail).
- A second comment with the alarm type C, the considered channel and the alarm technical name.

This chapter lists device alarms in separate tables for each measurement type.

- The first column indicates the transmitted data in the first comment.
- The second column indicates the transmitted data in the second comment.



Some non-listed alarms type C may only be transmitted in the second comment and can be ignored.

#### WBC

Transmitted data Alarm type, measurement type, alarm message	Transmitted data Alarm type, considered channel, alarm technical name	Description
S^WBC^WBC_ABN_MAT^SEP_LYM_NRBC	C^LMNE^NrbcWbcInterf	WBC abn. matrix LYM/NRBC
S^WBC^WBC_ABN_MAT^LYM_NRBC	C^LMNE^SuspiciousHgbNrbcInterf	WBC abn. matrix LYM/NRBC
S^WBC^WBC_ABN_MAT^SEP_LYM_NRBC	C^LMNE^NrbcLymInterf	WBC abn. matrix LYM/NRBC
S^WBC^WBC_ABN_MAT^LYM_PHOCYTOSIS	C^LMNE^SuspiciousLym	WBC abn. matrix Lymphocytosis?
S^WBC^WBC_ABN_MAT^NEU_NOISE	C^LMNE^SmallNeu	WBC abn. matrix NEU/Noise
S^WBC^WBC_ABN_MAT^NEU_NOISE	C^LMNE^MediumBackgndNoiseExcess	WBC abn. matrix NEU/Noise

Transmitted data Alarm type, measurement type, alarm message	Transmitted data Alarm type, considered channel, alarm technical name	Description
S^WBC^WBC_ABN_MAT^NEU_EOS_NOISE	C^LMNE^HighBackgndNoiseExcess	WBC abn. matrix NEU+EOS/Noise
S^WBC^WBC_ABN_MAT^NEU_EOS_NOISE	C^SYNTHESIS_WBC^HighBackgndNoiseExcessLmne1	WBC abn. matrix NEU+EOS/Noise
S^WBC^OOR_WBC^VISIBILITY	C^SYNTHESIS_WBC^WbcLimOfVisibility	WBC out of range High visibility
S^WBC^OOR_WBC^LINEARITY	C^SYNTHESIS_WBC^WbcLimOfLinearity	WBC out of range High linearity
S^WBC^OOR_WBC^LOQ	C^SYNTHESIS_WBC^WbcLimOfQuantif	WBC out of range Low linearity
S^WBC^WBC_ABN_MAT^INTE_RF_TNC	C^SYNTHESIS_WBC^LmneLowWbc	WBC abn. matrix TNC interference
S^WBC^WBC_ABN_MAT^INTE_RF_TNC	C^SYNTHESIS_WBC^TncLowWbc	WBC abn. matrix TNC interference
S^WBC^WBC_ABN_MAT^INTE_RF_NRBC	C^SYNTHESIS_WBC^SuspiciousNrbcCount	WBC abn. matrix NRBC interference
S^WBC^WBC_ABN_MAT^INTE_RF_NRBC	C^SYNTHESIS_WBC^SuspiciousNrbcCount1	WBC abn. matrix NRBC interference
S^WBC^WBC_ABN_MAT^NRBC_PLTAGR	C^SYNTHESIS_WBC^AbnormNrbcClassification	WBC abn. matrix NRBC/PLT aggregates?

#### DIFF

Transmitted data Alarm type, measurement type, alarm message	Transmitted data Alarm type, considered channel, alarm technical name	Description
S^DIFF^WBC_ABN_MAT^INTE_RF_LYM	C^BASO^BackgndNoise	WBC abn. matrix LYM interference
S^DIFF^WBC_ABN_MAT^SEP_LYM_NEU	C^LMNE^LymNeuSep	WBC abn. matrix LYM/NEU
S^DIFF^WBC_ABN_MAT^SEP_NEU_EOS	C^LMNE^NeuEosSep	WBC abn. matrix NEU/EOS

Transmitted data Alarm type, measurement type, alarm message	Transmitted data Alarm type, considered channel, alarm technical name	Description
S^DIFF^WBC_ABN_MAT^SEP_MON_NEU	C^LMNE^MonNeuSep	WBC abn. matrix MON/NEU
S^DIFF^WBC_ABN_MAT^SEP_LYM_MON	C^LMNE^LymMonSep	WBC abn. matrix LYM/MON
S^DIFF^WBC_ABN_MAT^SEP_LYM_NEU	C^LMNE^LymNeuDensSepExcess	WBC abn. matrix LYM/NEU
S^DIFF^WBC_ABN_MAT^SEP_LYM_MON	C^LMNE^LymMonDensSepExcess	WBC abn. matrix LYM/MON
S^DIFF^WBC_ABN_MAT^SEP_LYM_MON	C^LMNE^LymMonDensSepExcess1	WBC abn. matrix LYM/MON
S^DIFF^WBC_ABN_MAT^SEP_LYM_MON	C^LMNE^LymMonDensSepExcess2	WBC abn. matrix LYM/MON
S^DIFF^WBC_ABN_MAT^SEP_NEU_EOS	C^LMNE^NeuEosDensSepExcess	WBC abn. matrix NEU/EOS
S^DIFF^WBC_ABN_MAT^SEP_MON_NEU	C^LMNE^MonNeuDensSepExcess	WBC abn. matrix MON/NEU
S^DIFF^WBC_ABN_MAT^SEP_NEU_IMG	C^LMNE^NeuImgDensSepExcess	WBC abn. matrix NEU/IMG
S^DIFF^WBC_ABN_MAT^SEP_NEU_IMG	C^LMNE^NeuImgDensSepExcess1	WBC abn. matrix NEU/IMG
S^DIFF^WBC_ABN_MAT^SEP_EOS_IMG	C^LMNE^EosImgDensSepExcess	WBC abn. matrix EOS/IMG
S^DIFF^WBC_ABN_MAT^NEU_IMG	C^LMNE^ShrinedRightNeuBox	WBC abn. matrix NEU/IMG
S^DIFF^WBC_ABN_MAT^NEU_IMG	C^LMNE^ShrinedRightNeuBox1	WBC abn. matrix NEU/IMG
S^DIFF^WBC_ABN_MAT^SEP_MON_IMM	C^LMNE^RightMon	WBC abn. matrix MON/IMM
S^DIFF^WBC_ABN_MAT^INTERF_BASO	C^SYNTHESIS_WBC^BasoLowWbc	WBC abn. matrix BASO interference
S^DIFF^WBC_ABN_MAT^BASO_IMG	C^SYNTHESIS_WBC^LicMonBasoInterf	WBC abn. matrix BASO/IMG

Transmitted data Alarm type, measurement type, alarm message	Transmitted data Alarm type, considered channel, alarm technical name	Description
S^DIFF^WBC_ABN_MAT^SEP_LYM_MON	C^SYNTHESIS_WBC^LyaMonSep	WBC abn. matrix LYM/MON
S^DIFF^WBC_ABN_MAT^SUSP_NB_LIC0	C^SYNTHESIS_WBC^NbLicExc ess0	WBC abn. matrix LIC?
S^DIFF^WBC_ABN_MAT^SUSP_NB_LIC1	C^SYNTHESIS_WBC^NbLicExc ess1	WBC abn. matrix LIC?
S^DIFF^WBC_ABN_MAT^SUSP_NB_LIC2	C^SYNTHESIS_WBC^NbLicExc ess2	WBC abn. matrix LIC?
S^DIFF^WBC_ABN_MAT^SUSP_P_LIC	C^SYNTHESIS_WBC^PLicExcess	WBC abn. matrix LIC?
S^DIFF^WBC_ABN_MAT^SUSP_NB_ALY0	C^SYNTHESIS_WBC^NbAlyExc ess0	WBC abn. matrix ALY?
S^DIFF^WBC_ABN_MAT^SUSP_NB_ALY1	C^SYNTHESIS_WBC^NbAlyExc ess1	WBC abn. matrix ALY?
S^DIFF^WBC_ABN_MAT^SUSP_NB_ALY2	C^SYNTHESIS_WBC^NbAlyExc ess2	WBC abn. matrix ALY?
S^DIFF^WBC_ABN_MAT^SUSP_NB_ALY3	C^SYNTHESIS_WBC^NbAlyExc ess3	WBC abn. matrix ALY?
S^DIFF^WBC_ABN_MAT^SUSP_NB_ALY4	C^SYNTHESIS_WBC^NbAlyExc ess4	WBC abn. matrix ALY?
S^DIFF^WBC_ABN_MAT^SUSP_P_ALY	C^SYNTHESIS_WBC^PAlyExcess	WBC abn. matrix ALY?
S^DIFF^WBC_ABN_MAT^SUSP_NB_IMM	C^SYNTHESIS_WBC^NblmmExcess	WBC abn. matrix IMM?
S^DIFF^WBC_ABN_MAT^SUSP_P_IMM	C^SYNTHESIS_WBC^PlmmExcess	WBC abn. matrix IMM?
S^DIFF^WBC_ABN_MAT^SUSP_NB_IML	C^SYNTHESIS_WBC^NblmlExc ess	WBC abn. matrix IML?
S^DIFF^WBC_ABN_MAT^SUSP_P_IML	C^SYNTHESIS_WBC^PlmlExcess	WBC abn. matrix IML?

## RBC

Transmitted data Alarm type, measurement type, alarm message	Transmitted data Alarm type, considered channel, alarm technical name	Description
S^RBC^OOR_RBC^VISIBILITY	C^SYNTHESIS_RBC_PLT^RbcLimOfVisibility	RBC out of range High visibility
S^RBC^OOR_RBC^LINEARITY	C^SYNTHESIS_RBC_PLT^RbcLimOfLinearity	RBC out of range High linearity
S^RBC^OOR_RBC^LOQ	C^SYNTHESIS_RBC_PLT^RbcLimOfQuantif	RBC out of range Low linearity
S^RBC^PLT_CONCENT	C^SYNTHESIS_RBC_PLT^PitConcentrateMode	PLT concent. mode
S^RBC^OOR_PLT^VISIBILITY	C^SYNTHESIS_RBC_PLT^PitLimOfVisibility	PLT out of range High visibility

## RBC\_OTH

Transmitted data Alarm type, measurement type, alarm message	Transmitted data Alarm type, considered channel, alarm technical name	Description
S^RBC_OTH^RBC_ABN_HIST^DBL_POP	C^RBC^RbcDoublePop	RBC abn. histogram Double population?
S^RBC_OTH^OOR_HCT^VISIBILITY	C^SYNTHESIS_RBC_PLT^HctLimOfVisibility	HCT out of range High visibility
S^RBC_OTH^OOR_HCT^LINEARITY	C^SYNTHESIS_RBC_PLT^HctLimOfLinearity	HCT out of range High linearity
S^RBC_OTH^OOR_HCT^LOQ	C^SYNTHESIS_RBC_PLT^HctLimOfQuantif	HCT out of range Low linearity
S^RBC_OTH^RBC_ABN_HIST^ABN_DISTRI	C^SYNTHESIS_RBC_PLT^AbnormRbcDistrib	RBC abn. histogram Abnormal distribution
S^RBC_OTH^OOR_HGB^VISIBILITY	C^SYNTHESIS_RBC_PLT^HgbLimOfVisibility	HGB out of range High visibility
S^RBC_OTH^OOR_HGB^LINEARITY	C^SYNTHESIS_RBC_PLT^HgbLimOfLinearity	HGB out of range High linearity
S^RBC_OTH^OOR_HGB^LOQ	C^SYNTHESIS_RBC_PLT^HgbLimOfQuantif	HGB out of range Low linearity

Transmitted data Alarm type, measurement type, alarm message	Transmitted data Alarm type, considered channel, alarm technical name	Description
S^RBC_OTH^PLT_ABN_HIST^SCH_MACRO_PLT	C^PLT^Sch	PLT abn. histogram Schistocyte / Macro PLT?
S^RBC_OTH^PLTO_ABN_MAT^RBC_DBL_POP	C^RBC_PLTO^RbcDoublePop	PLTO abn. matrix RBC Double population?

## PLT

Transmitted data Alarm type, measurement type, alarm message	Transmitted data Alarm type, considered channel, alarm technical name	Description
S^PLT^PLT_INTERF^PLTAGR	C^SYNTHESIS_WBC^PitAggr	PLT Interference PLT aggregates?
S^PLT^PLT_INTERF^PLTAGR	C^SYNTHESIS_WBC^PitAbsAggr	PLT Interference PLT aggregates?
S^PLT^PLT_INTERF^PLTAGR	C^SYNTHESIS_WBC^PitAggrLowPit	PLT Interference PLT aggregates?
S^PLT^PLT_INTERF^PLTAGR	C^SYNTHESIS_WBC^PitAggrAbnormalPit	PLT Interference PLT aggregates?
S^PLT^PLT_ABN_HIST^SEP_RBC_PLT	C^PLT^RbcPitFixedSep	PLT abn. histogram RBC/PLT
S^PLT^PLT_ABN_HIST^SEP_RBC_PLT	C^PLT^RbcPitSepOutOfRange	PLT abn. histogram RBC/PLT
S^PLT^PLT_ABN_HIST^ABN_MPV	C^PLT^AbnormMpv	PLT abn. histogram Abnormal MPV
S^PLT^PLT_ABN_HIST^ABN_PDW	C^PLT^AbnormPdw	PLT abn. histogram Abnormal PDW
S^PLT^PLT_ABN_HIST^LARGE_PLT	C^PLT^AbnormPit	PLT abn. histogram Large PLT?
S^PLT^OOR_PLT^LINEARITY	C^SYNTHESIS_RBC_PLT^PitLimOfLinearity	PLT out of range High linearity
S^PLT^OOR_PLT^LOQ	C^SYNTHESIS_RBC_PLT^PitLimOfQuantif	PLT out of range Low linearity

Transmitted data Alarm type, measurement type, alarm message	Transmitted data Alarm type, considered channel, alarm technical name	Description
S^PLT^PLT_ABN_HIST^LOW_COUNT	C^SYNTHESIS_RBC_PLT^LowPlt	PLT abn. histogram Low count
S^PLT^PLTO_ABN_MAT^SCH	C^RBC_PLTO^SchExcess	PLTO abn. matrix Schistocyte?
S^PLT^PLTO_ABN_MAT^PLTAGR_BAL_PLTO	C^SYNTHESIS_RBC_PLT^PltPltoHistoBalance	PLTO abn. matrix PLT aggr. or PLTO/PLT bal.?
S^PLT^PLT_INTERF^PLTAGR	C^SYNTHESIS_RBC_PLT^PltAggrSuspicion	PLT Interference PLT aggregates?

### PLTO

Transmitted data Alarm type, measurement type, alarm message	Transmitted data Alarm type, considered channel, alarm technical name	Description
S^PLTO^OOR_RBC^VISIBILITY	C^SYNTHESIS_RBC_PLT^RbcLimOfVisibility	RBC out of range High visibility
S^PLTO^OOR_RBC^LINEARITY	C^SYNTHESIS_RBC_PLT^RbcLimOfLinearity	RBC out of range High linearity
S^PLTO^OOR_RBC^LOQ	C^SYNTHESIS_RBC_PLT^RbcLimOfQuantif	RBC out of range Low linearity
S^PLTO^PLT_CONCENT	C^SYNTHESIS_RBC_PLT^PltConcentrateMode	PLT concent. mode
S^PLTO^OOR_PLT^VISIBILITY	C^SYNTHESIS_RBC_PLT^PltLimOfVisibility	PLT out of range High visibility
S^PLTO^PLTO_ABN_MAT^SEP_PLTOX_RBC	C^RBC_PLTO^RbcPltoFixedSep	PLTO abn. matrix PLT-Ox/RBC
S^PLTO^PLTO_ABN_MAT^LPF_RBC	C^RBC_PLTO^InvalidLpfCalc	PLTO abn. matrix LPF/RBC
S^PLTO^PLTO_ABN_MAT^SEP_LPF_RBC	C^RBC_PLTO^SchLpfInterf	PLTO abn. matrix LPF/RBC
S^PLTO^OOR_PLTO^VISIBILITY	C^SYNTHESIS_RBC_PLT^PltoLimOfVisibility	PLTO out of range High visibility

Transmitted data Alarm type, measurement type, alarm message	Transmitted data Alarm type, considered channel, alarm technical name	Description
S^PLTO^OOR_PLTO^LINEARITY	C^SYNTHESIS_RBC_PLT^PltoLimOfLinearity	PLTO out of range High linearity
S^PLTO^OOR_PLTO^LOQ	C^SYNTHESIS_RBC_PLT^PltoLimOfQuantif	PLTO out of range Low linearity
S^PLTO^PLTO_ABN_MAT^PLTAGR_BAL_PLTO	C^SYNTHESIS_RBC_PLT^PltPltoHistoBalance	PLTO abn. matrix PLT aggr. or PLTO/PLT bal.?
S^PLTO^PLT_INTERF^PLTAGR	C^SYNTHESIS_RBC_PLT^PltAggrSuspicion	PLT Interference PLT aggregates?

### RET

Transmitted data Alarm type, measurement type, alarm message	Transmitted data Alarm type, considered channel, alarm technical name	Description
S^RET^RET_ABN_MAT^SEP_RBC_PLT	C^RET^RbcPltFixedSep	RET abn. matrix RBC/PLT
S^RET^RET_ABN_MAT^SEP_RBC_PLT	C^RET^RbcPltSepOutOfRange	RET abn. matrix RBC/PLT
S^RET^RET_ABN_MAT^RET_PLT	C^RET^ReticulatedPltSepExcess	RET abn. matrix RET/PLT
S^RET^RET_ABN_MAT^SEP_RET_RBC	C^RET^RbcRetFixedSep	RET abn. matrix RET/RBC
S^RET^OOR_RET^VISIBILITY	C^SYNTHESIS_RBC_PLT^RetLimOfVisibility	RET out of range High visibility
S^RET^OOR_RET^LINEARITY	C^SYNTHESIS_RBC_PLT^RetLimOfLinearity	RET out of range High linearity
S^RET^OOR_RET^LOQ	C^SYNTHESIS_RBC_PLT^RetLimOfQuantif	RET out of range Low linearity
S^RET^OOR_RET^LINEARITY	C^SYNTHESIS_RBC_PLT^PRetLimOfLinearity	RET out of range High linearity
S^RET^RET_ABN_MAT^LOW_COUNT	C^SYNTHESIS_RBC_PLT^LowRet	RET abn. matrix Low count

## BFWBC

Transmitted data Alarm type, measurement type, alarm message	Transmitted data Alarm type, considered channel, alarm technical name	Description
S^BFWBC^BFWBC_ABN_MAT^ NOISE_PN	C^BF_LMNE^HighBackgndNois eExcess	WBC abn. matrix Noise/PN
S^BFWBC^BFWBC_ABN_MAT^ NHC	C^BF_LMNE^NucleatedCellsBf WbcInterf	WBC abn. matrix Non hemato cells?
S^BFWBC^OOR_BFWBC^VISIB ILITY	C^SYNTHESIS_BF^BfWbcLimO fVisibility	WBC out of range High visibility
S^BFWBC^OOR_BFWBC^LINE ARITY	C^SYNTHESIS_BF^BfWbcLimO fLinearity	WBC out of range High linearity
S^BFWBC^OOR_BFWBC^LOQ	C^SYNTHESIS_BF^BfWbcLimO fQuantif	WBC out of range Low linearity

## BFDIFF

Transmitted data Alarm type, measurement type, alarm message	Transmitted data Alarm type, considered channel, alarm technical name	Description
S^BFDIFF^BFWBC_ABN_MAT^ SEP_MN_PN	C^BF_LMNE^SuspiciousMnPnS ep	WBC abn. matrix MN/PN
S^BFDIFF^BFWBC_ABN_MAT^ INTERF_LEFT	C^BF_LMNE^BackgndNoise	WBC abn. matrix Left interference

## BFRBC

Transmitted data Alarm type, measurement type, alarm message	Transmitted data Alarm type, considered channel, alarm technical name	Description
S^BFRBC^BFRBC_ABN_HIST^ INTERF_RBC	C^BF_RBC^BfRbcInterf	RBC abn. histogram RBC interference?
S^BFRBC^OOR_BFRBC^VISIB ILITY	C^SYNTHESIS_BF^BfRbcLimOf Visibility	RBC out of range High visibility

Transmitted data Alarm type, measurement type, alarm message	Transmitted data Alarm type, considered channel, alarm technical name	Description
S^BFRBC^OOR_BFRBC^LINEA RITY	C^SYNTHESIS_BF^BfRbcLimOf Linearity	RBC out of range High linearity
S^BFRBC^OOR_BFRBC^LOQ	C^SYNTHESIS_BF^BfRbcLimOf Quantif	RBC out of range Low linearity

## 3.5. Curves and Matrix Transmission

### 3.5.1. General Decoding

The image data must be uncompressed using first the *base64* and secondly the *deflate* algorithms.

The data must be converted from binary to text format.

### 3.5.2. Matrices

Matrices are transmitted in the manufacturer record as follows:

- The Message Type field (14.3) is set to: MATRIX.
- The Measurement Type field (14.4) is set to: DIFF, PLTO, RET, BFDIFF.
- The Name field (14.5) is set to: LMNERESABS, RBCPLTRESABS, RETRESFLUO, MPRESABS.
- The Thresholds field (14.6) contains the encoding type and the FLOATLE-stream/deflate:base64 coded value of the matrix polygons.
- The Points field (14.7) contains the encoding type and the FLOATLE-stream/deflate:base64 coded value of the matrix points.

## Matrix Thresholds

Matrix polygons thresholds data must be in accordance with the following frame:

Number of bytes	Data	Format	Meaning
4	X display min	FLOATLE FLOATLE: IEEE 754 floating point value transmitted in Little Endian byte order (Intel)	X min value for Matrix start
4	X display max	FLOATLE	X max value for Matrix end
4	Y display min	FLOATLE	Y min value for Matrix start
4	Y display max	FLOATLE	Y max value for Matrix end
4	NumberOfList = 3	FLOATLE	One list for X coordinates of polygons, one list for Y coordinates of polygons, one list for box identifiers
4	ListLength	FLOATLE	Number of thresholds in the list
4 x ListLength	X	FLOATLE	X value of the threshold for each threshold
4 x ListLength	Y	FLOATLE	Y value of the threshold for each threshold
4 x ListLength	BoxID	FLOATLE	ID of each box

BoxID for the matrix: LMNERESABS.

BoxID	Box Name
0	LYM
1	MON
2	NEU
3	EOS
4	IMG
5	ALY
7	RN
8	RM

BoxID	Box Name
9	IMM
10	IML
11	LN
12	Background Noise Low
13	Background Noise High
15	Low Optical Correlation
16	Background Noise Bubbles
17	Erythroblasts
18	Platelet aggregates
19	MON + RM + IMM
20	LYM + ALY + IML
21	<b>5DIFF:</b> LN + NEU + RN + IMG <b>6DIFF:</b> LN + NEU
22	<b>6DIFF:</b> RN + IMG
23	Erythroblasts + Platelet aggregates

BoxID for the matrix: RBCPLTRESABS.

BoxID	Box Name
0	Small Cells
1	RBC = MIC + RBC_N + MAC
2	PLT = PLT_N + PLT_L
3	RBC_N
4	MIC
5	MAC
6	PLT_N
7	PLT_L

BoxID for the matrix: RETRESFLUO.

BoxID	Box Name
0	Background Noise
1	PLT
2	PLT_RET
3	RBC
4	RET_L
5	RET_M
6	RET_H
7	IMR
8	NRBC
9	WBC

BoxID for the matrix: MPRESABS.

BoxID	Box Name
0	Background Noise Low
1	Background Noise High
2	BFMN = LL + LYM + MON
3	LL
4	LYM
6	MON
7	BFPN = EOS + NEU
8	EOS
9	NEU
10	Immature cells = Immature polynucleated cells + Immature mononucleated cells
11	Immature polynucleated cells
12	Immature mononucleated cells

## Matrix Points

Matrix points data must be in accordance with the following frame:

Number of bytes	Data	Format	Meaning
4	X display min	FLOATLE	X min value for Matrix start
4	X display max	FLOATLE	X max value for Matrix end
4	Y display min	FLOATLE	Y min value for Matrix start
4	Y display max	FLOATLE	Y max value for Matrix end
4	X scale NB	FLOATLE	Number of X ticks in the list to display
4 x X scale NB	X scale	FLOATLE	X tick values
4	Y scale NB	FLOATLE	Number of Y ticks in the list to display
4 x Y scale NB	Y scale	FLOATLE	Y tick values
4	NumberOfList = 4	FLOATLE	One list for X points, one list for Y points, one list for the number of points for this coordinate, one list for the population to which the point belongs (X, Y coordinates)
4	ListLength	FLOATLE	Number of elements in the list
4 x ListLength	X	FLOATLE	X (Coordinate)
4 x ListLength	Y	FLOATLE	Y (Coordinate)
4 x ListLength	Qty	FLOATLE	Quantity (number of points for the (X, Y) coordinate)
4 x ListLength	Pop	FLOATLE	Population to which the point (X, Y coordinates) belongs

Population ID for the matrix: LMNERESABS.

Population ID	Population Name
0	LYM
1	MON
2	NEU
3	EOS

Population ID	Population Name
4	IMG
5	ALY
6	LL
7	RN
8	RM
9	IMM
10	IML
11	LN
12	Background Noise Low
13	Background Noise High
14	BAS
15	Low Optical Correlation
16	Background Noise Bubbles
17	Erythroblasts
18	Platelet aggregates

Population ID for the matrix: RBCPLTRESABS.

Population ID	Population Name
0	PLT_N
1	PLT_L
2	MIC
3	MAC
4	RBC_N
5	Small Cells

Population ID for the matrix: RETRESFLUO.

Population ID	Population Name
0	Background Noise
1	PLT

Population ID	Population Name
2	PLT_RET
3	RBC
4	RET_L
5	RET_M
6	RET_H
7	IMR
8	NRBC
9	WBC

Population ID for the matrix: MPRESABS.

Population ID	Population Name
0	LYM
1	MON
2	NEU
3	EOS
4	Immature polynucleated cells
6	LL
7	Immature mononucleated cells
8	Background Noise Low
9	Background Noise High

### 3.5.3. Histograms

Histograms are transmitted in the manufacturer record as follows:

- The Message Type field (14.3) is set to: HISTOGRAM.
- The Measurement Type field (14.4) is set to: DIFF, PLT, WBC, RBC, BFRBC, BFDIFF.
- The Name field (14.5) is set to: BASOALONGRES, EOSALONGABS, EOSALONGRES, LYMALONGABS, LYMALONGRES, MONALONGABS, MONALONGRES, NEUALONGABS, NEUALONGRES, PLTALONGRES,

TNCALONGRES, RBCALONGRES, RBCPLTALONGRES, MNALONGABS, MNALONGRES, PNALONGABS, PNALONGRES.

- The Thresholds field (14.6) contains the encoding type and the FLOATLE-stream/deflate:base64 coded value of the histogram thresholds.
- The Points field (14.7) contains the encoding type and the FLOATLE-stream/deflate:base64 coded value of the histogram points.

## Histogram Thresholds

Histogram thresholds data must be in accordance with the following frame:

Number of bytes	Data	Format	Meaning
4	X display min	FLOATLE	X min value for Histogram start
4	X display max	FLOATLE	X max value for Histogram end
4	Y display min	FLOATLE	Y min value for Histogram start
4	Y display max	FLOATLE	Y max value for Histogram end
4	NumberOfList = 2	FLOATLE	One list for X threshold positions and one list of threshold identifiers
4	ListLength	FLOATLE	Number of thresholds in the list
4 x ListLength	X	FLOATLE	X value of the threshold for each threshold
4 x ListLength	ThrsID	FLOATLE	ID of each threshold

ThrsID for the histogram: BASOALONGRES.

ThrsID	Threshold Name	Value
0	RBG	Setting
1	RBasoMin	Setting
2	RBasoMax	Setting

ThrsID for the histogram: RBCALONGRES.

ThrsID	Threshold Name	Value
0	VMic	Setting
1	VMac	Setting

ThrsID for the histogram: PLTALONGRES.

ThrsID	Threshold Name	Value
0	thrPitRbc	Calculated

ThrsID for the histogram: RBCPLTALONGRES.

ThrsID	Threshold Name	Value
0	thrSynthPitRbc	Calculated

ThrsID for the histogram: TNCALONGRES.

ThrsID	Threshold Name	Value
0	RTNC1	Setting
1	RTNC2	Setting
2	RTNC3	Setting

There are no thresholds (ListLength = 0) for the following histograms: EOSALONGABS, EOSALONGRES, LYMALONGABS, LYMALONGRES, MONALONGABS, MONALONGRES, NEUALONGABS, NEUALONGRES, MNALONGABS, MNALONGRES, PNALONGABS, PNALONGRES.

## Histogram Points

Histogram points data must be in accordance with the following frame:

Number of bytes	Data	Format	Meaning
4	X display min	FLOATLE	X min value for Histogram start
4	X display max	FLOATLE	X max value for Histogram end
4	Y display min	FLOATLE	Y min value for Histogram start
4	Y display max	FLOATLE	Y max value for Histogram end
4	X scale NB	FLOATLE	Number of X ticks in the list to display
4 x X scale NB	X scale	FLOATLE	X tick values
4	Y scale NB	FLOATLE	Number of Y ticks in the list to display

Number of bytes	Data	Format	Meaning
4 x Y scale NB	Y scale	FLOATLE	Y tick values
4	NumberOfList = 2	FLOATLE	One list for X positions and one list for the Y positions (quantity on each X position)
4	ListLength	FLOATLE	Number of elements in the list
4 x ListLength	X	FLOATLE	X (Coordinate)
4 x ListLength	Y	FLOATLE	Y (Coordinate)

### 3.6. Example of Data Frame

#### 3.6.1. Example of a Query: With Response

```

<- <ENQ>
-> <ACK>
<- <STX>1H|\^&||MHR1^210M2SH01011^1.7.0|||||P|LIS2-A2|
20230929091956<CR><ETX>00<CR><LF>
-> <ACK>
<- <STX>2Q|1|^2023092700000011^1^042249^1||ALL|||||O<CR><ETX>E6<CR><LF>
-> <ACK>
<- <STX>3Q|2|^2023092700000012^1^042249^2||ALL|||||O<CR><ETX>EA<CR><LF>
-> <ACK>
<- <STX>4Q|3|^2023092700000013^1^042249^3||ALL|||||O<CR><ETX>EE<CR><LF>
-> <ACK>
<- <STX>5Q|4|^2023092700000014^1^042249^4||ALL|||||O<CR><ETX>F2<CR><LF>
-> <ACK>
<- <STX>6Q|5|^2023092700000015^1^042249^5||ALL|||||O<CR><ETX>F6<CR><LF>

```

```

-> <ACK>
<- <STX>7Q|6|^2023092700000016^1^042249^6||ALL|||||O<CR><ETX>FA<CR><LF>
-> <ACK>
<- <STX>0Q|7|^2023092700000017^1^042249^7||ALL|||||O<CR><ETX>F6<CR><LF>
-> <ACK>
<- <STX>1Q|8|^2023092700000018^1^042249^8||ALL|||||O<CR><ETX>FA<CR><LF>
-> <ACK>
<- <STX>2Q|9|^2023092700000019^1^042249^9||ALL|||||O<CR><ETX>FE<CR><LF>
-> <ACK>
<- <STX>3Q|10|^2023092700000020^1^042249^10||ALL|||||O<CR><ETX>47<CR><LF>
-> <ACK>
<- <STX>4L|1|N<CR><ETX>07<CR><LF>
-> <ACK>
<- <EOT>
-> <ENQ>
<- <ACK>
-> <STX>1H|\^&||YP8K|||||P|LIS2-A2|20230929092120<CR><ETX>B4<CR><LF>
<- <ACK>
-> <STX>2P|1||00000011||PATIENT 11^TEST||19851114^37^Y|M<CR><ETX>C4<CR><LF>
<- <ACK>
-> <STX>3O|1|2023092700000011^1^042249^1||^DIF|R|20230927174534|20230927174534|||
N||||BLOOD|||||Q<CR><ETX>9B<CR><LF>
<- <ACK>
-> <STX>4P|2||00000012||PATIENT 12^TEST||19970225^26^Y|F<CR><ETX>C5<CR><LF>
<- <ACK>
-> <STX>5O|1|2023092700000012^1^042249^2||^DIF|R|20230927174536|20230927174536|||
N||||BLOOD|||||Q<CR><ETX>A3<CR><LF>
<- <ACK>
-> <STX>6P|3||00000013||PATIENT 13^TEST||19760228^47^Y|M<CR><ETX>D4<CR><LF>
<- <ACK>
-> <STX>7O|1|2023092700000013^1^042249^3||^DIF|R|20230927174538|20230927174538|||
N||||BLOOD|||||Q<CR><ETX>AB<CR><LF>

```

```

-> <ACK>
-> <STX>0P|4||00000014||PATIENT 14^TEST||19710717^52^Y|F<CR><ETX>C4<CR><LF>
<- <ACK>
-> <STX>1O|1|2023092700000014^1^042249^4||^^^DIF|R|20230927174540|20230927174540|||
N|||BLOOD|||||||Q<CR><ETX>99<CR><LF>
<- <ACK>
-> <STX>2P|5||00000015||PATIENT 15^TEST||19650903^58^Y|M<CR><ETX>D6<CR><LF>
<- <ACK>
-> <STX>3O|1|2023092700000015^1^042249^5||^^^DIF|R|20230927174542|20230927174542|||
N|||BLOOD|||||||Q<CR><ETX>A1<CR><LF>
<- <ACK>
-> <STX>4P|6||00000016||PATIENT 16^TEST||20130315^10^Y|F<CR><ETX>B6<CR><LF>
<- <ACK>
-> <STX>5O|1|2023092700000016^1^042249^6||^^^DIF|R|20230927174544|20230927174544|||
N|||BLOOD|||||||Q<CR><ETX>A9<CR><LF>
<- <ACK>
-> <STX>6P|7||00000017||PATIENT 17^TEST||19870921^36^Y|M<CR><ETX>E0<CR><LF>
<- <ACK>
-> <STX>7O|1|2023092700000017^1^042249^7||^^^DIF|R|20230927174546|20230927174546|||
N|||BLOOD|||||||Q<CR><ETX>B1<CR><LF>
<- <ACK>
-> <STX>0P|8||00000018||PATIENT 18^TEST||19360920^87^Y|F<CR><ETX>D5<CR><LF>
<- <ACK>
-> <STX>1O|1|2023092700000018^1^042249^8||^^^DIF|R|20230927174548|20230927174548|||
N|||BLOOD|||||||Q<CR><ETX>B1<CR><LF>
<- <ACK>
-> <STX>2P|9||00000019||PATIENT 19^TEST||19240128^99^Y|M<CR><ETX>E1<CR><LF>
<- <ACK>
-> <STX>3O|1|2023092700000019^1^042249^9||^^^DIF|R|20230927174550|20230927174550|||
N|||BLOOD|||||||Q<CR><ETX>A7<CR><LF>
<- <ACK>
-> <STX>4P|10||00000020||PATIENT 20^TEST||19890613^34^Y|F<CR><ETX>F3<CR><LF>
<- <ACK>

```

```

-> <STX>5O|1|2023092700000020^1^042249^10||^^^DIF|R|20230927174552|20230927174552|||
N|||BLOOD|||||||Q<CR><ETX>CD<CR><LF>
<- <ACK>
-> <STX>6L|1|N<CR><ETX>09<CR><LF>
<- <ACK>
-> <EOT>

```

### 3.6.2. Example of a Query: Known Tube Without Worklist

```

<- <ENQ>
-> <ACK>
<- <STX>1H|\^&||MHR1^210M2SH01011^1.7.0|||||P|LIS2-A2|
20230929091956<CR><ETX>00<CR><LF>
-> <ACK>
<- <STX>2Q|1|^2023092700000005^1^042249^1||ALL|||||O<CR><ETX>E9<CR><LF>
-> <ACK>
<- <STX>3L|1|N<CR><ETX>06<CR><LF>
-> <ACK>
<- <EOT>
-> <ENQ>
<- <ACK>
-> <STX>1H|\^&||YP8K|||||P|LIS2-A2|20230929092446<CR><ETX>BF<CR><LF>
<- <ACK>
-> <STX>2P|1|<CR><ETX>BB<CR><LF>
<- <ACK>
-> <STX>3O|1|2023092700000005^1^042249^1||^^^|||||N|||||Y<CR><ETX>6F<CR><LF>
<- <ACK>
-> <STX>4L|1|N<CR><ETX>07<CR><LF>
<- <ACK>
-> <EOT>

```

### 3.6.3. Example of a Query: Unknown Tube

```

-> <ENQ>
-> <ACK>
<- <STX>1H|\^&||MHR1^210M2SH01011^1.7.0||||P|LIS2-A2|
20230929091956<CR><ETX>00<CR><LF>
-> <ACK>
<- <STX>2Q|1|^2023092700000205^1^042249^1||ALL|||||O<CR><ETX>EB<CR><LF>
-> <ACK>
<- <STX>3L|1|N<CR><ETX>06<CR><LF>
-> <ACK>
<- <EOT>
-> <ENQ>
<- <ACK>
-> <STX>1H|\^&||YP8K||||P|LIS2-A2|20230929092601<CR><ETX>B8<CR><LF>
<- <ACK>
-> <STX>2P|1|<CR><ETX>BB<CR><LF>
<- <ACK>
-> <STX>3O|1|2023092700000205^1^042249^1||^|N|||||Z<CR><ETX>72<CR><LF>
<- <ACK>
-> <STX>4L|1|N<CR><ETX>07<CR><LF>
<- <ACK>
-> <EOT>

```

### 3.6.4. Example of Result With Rerun and Reflex

```

<- <ENQ>
-> <ACK>

```

```

<- <STX>1H|\^&||MHR1^210M2SH01011^1.7.0||||P|LIS2-A2|
20230302102751<CR><ETX>E3<CR><LF>
-> <ACK>
<- <STX>2P|1||00000000||PAT 2303020020^PAT 2303020020||19250709000000^97^Y|
F<CR><ETX>1F<CR><LF>
-> <ACK>
<- <STX>3O|1|202309270000020^1^041176^1||^|DIR|R||20230302082259||||BLOOD||||
AP^STANDARD(m)|||F<CR><ETX>6E<CR><LF>
-> <ACK>
<- <STX>4C|1||S^DIFF^WBC_ABN_MAT^SEP_NEU_EOS||<CR><ETX>82<CR><LF>
-> <ACK>
<- <STX>5C|2||C^LMNE^NeuEosSep||<CR><ETX>DE<CR><LF>
-> <ACK>
<- <STX>6M|1|HISTOGRAM|WBC|TNCALONGRES|FLOATLE-stream/
deflate:base64^Y2AAgXpnMNWwCEQ7ABEQNwBxgSMQQ+QYGuxBcgA=|FLOATLE-stream/
deflate:base64^xdV/
qF5zHAFwg2EYuzY2OzaLobhwcYw7m7G8N2XVyzMVzMXHbxjGEYHhvrzq9u0lzMui2tpbWWpJu
WniQtSUvSknSTCTdtLS0tf3Pp/pW<ETB>4B<CR><LF>
-> <ACK>
<- <STX>7P4Qyrdefc45z3nOeZ5zPt/3N8vSeLp5sFTW12pt5A7EfqWxtp+nY/RSpZ
+sKcvqKdETivtIKnQRTc99LKRLfRRZRvb2UE/
OxIgN9n0LBvCUIZRx0hGMYZ6GpjJEpMZirTyJnBLGbtShvzWEA7C
+mgkzJLWcZyKqzgWVbRxfO8yMt08wqvsoYe3mAt6+hiPW+zgY28wyY2s4V3eY/
36eMDtvIhVT7iYz5hG5/yGZ+znS/4k<ETB>7A<CR><LF>
-> <ACK>
<- <STX>0q/Ywdd8w7f08x3f8wM7+ZGf
+JkBfmEXv7Kb39jDXjl9cAAHchBDOJhDOJShHMBhHMEwjuQohPH0RzDCEZyLMdxPKMYzQmcy
BjGchLjqOdkxjOBBk7hVE5jlqdzBmcyibM4m3MocS7ncT6TmclFXMHuLUJiLmEal3IzjeQ0MZ1mZj
CTy7mCWVzJVzNbObQwjW0ci3XcT1t3MBcbmQe87mJm1nALdzKbbRzO3dwJwu5i0XcT<ETB>
A9<CR><LF>
-> <ACK>
<- <STX>1Qf3cC+L6eQ+7ucByizhQR5iKQ/zCI
+yjMd4nCdYzpM81fxHLvyjkTijj5woxfzPYs5tjn4bG8+/JX5/uv8zrOS5mkraXhG/
dVE8p7nxncZ4zxOIT1lfpX5LvTci+ml89MSUeJct8bw74prpHI3qS6zhrZrKWUeDmZiqut403aP
+nrt88FzXqPb9gvqanVV/lcl8a7mRw/MjB5K/ZV6MvV5miv75iGaV7tiDm
+NXNgU2dEde<ETB>E8<CR><LF>

```

```
-> <ACK>
<-
<STX>2dQZ2VSO3GqLfGuITeWzZmZ7xQGRq3345W4387Ypc3qcc2iO7SwV51JTxWdOf874+9vsL6
0E11odKrBd5bFcK60bx8yvv9Upxjfmr41nhOsVziuftr9iLxf2/07v/dvwX1/j/x
+8=<CR><ETX>10<CR><LF>
```

```
-> <ACK>
<- <STX>3M|2|HISTOGRAM|RBC|RBCALONGRES|FLOATLE-stream/
deflate:base64^Y2AAgXZnEPn2+3wXIOUAWRIODAXSYHEGhgZ7AA=|FLOA<ETB>6C<CR><LF>
```

```
-> <ACK>
<- <STX>4TLE-stream/deflate:base64^dZR5UBRnGsZRFBHHCzArHkBiSFDxABURIf6+9/
k4hEREE4kHXiErGliDB0ZcIAHBJUgQcFFGEGE8UilOBNFALFDWZSOwWsB6RROlikfwCKZYUcO
+spWt/Wenaqa7Z95v+vf08+s2M3vz2o43n20deYo3wux/vvPcVz3/
iaHBG2Tg7lwL1aC90Vwq4n3q4Xtrt8bC7DEmxcKOso6RMwRf<ETB>A5<CR><LF>
```

```
-> <ACK>
<- <STX>5KMdLL5UNpu/ESa39wgeV5GBaXwTKZcXrqX50wy0K6IZ0ukV2wFz1fL8c21vKZJjvC
+xevuyn5H2njtr7JT9xuvN6cNHn0oKqgfRaweSHwuCs35Ay0vHU6L6x2lz0vzzZwp0M6FAtwmETM
QhXqQV
+xM8swSxDzkWuNH45vfozGdc4nZyNFIIY3wDqGhS1cQc9KgneHU70gk9TkfRcxMXe2bqFO3hV5
aG4j5qWNcljNsp3aVzB<ETB>F0<CR><LF>
```

```
-> <ACK>
<- <STX>6wp9DgkIVnS6GF0vNk0r20vzDTbmopzGauvXS7KpfZ8uj6jXzmM1FT
+0FmLKQuqPMWUR1TseZ9QTVeH3NvCVUHXYamc9Qpf4b5q6giuRzzF5JZwrOM381nSq/
yBlqqLjxO85RS8fa6jnLZTps0cB5msjkcJUzXad9Hjc51y0yBv3A2Zopa/
VPnO8uZSS0csYHIJrziHO2UXLpE876jBLr2znvr2Ro7eDMLyJw7BXn/o0qdG<ETB>C7<CR><LF>
```

```
-> <ACK>
<-
<STX>7Z4ad0DnnY9wdcfZ5x6oWNCb7i7WYC7QImXJdpVX7gFWIF7wYlgHR6H9Mf40AHgjnBMP
wgPowdJTKw1uC8cTrbFvbQhcMp6C9wdTAVD0VJoh7eLh4F7RF75CNyuGomRNfBgTmFsdMT1G
29jaPM74H6R1fYumtqdYNs5Gtw1MizG4IpuLabZjAP3ji8dxqPOaQJOLhPBDuALD1fUeLmhj/
dksA9IDJqK6mB3mC+dBnYDhtXTUa<ETB>E1<CR><LF>
```

```
-> <ACK>
<-
<STX>0n3RFf0DLAniE2YhYpkL7xM08DOYFOOxJkCQkchwP5gfak3TpX74JcqX7BL0NfPRnGjPx7f
CAB7hYjW93GsbQ4etgeCHUOY2TwtcpiPe7oPwL4h1G4BTA7BaHb6COWelrktwj6PxbjltQTsIRYHL
IUxaBmuBy8HO4ng0JXIwv0xmvShYD8xP/aPyEhYhSvJYWBXMSdrNVJzPkvtQQTYW/
gXf4bk0jWoKdeDHYZPzVok1q9Dde<ETB>48<CR><LF>
```

```
-> <ACK>
```

```
<- <STX>1N6sM+g5mgYWjeisu1zsNvw6oxBrNmfUG4RC/Yc0222okIXh+/T4rh3A4ZnGbj7eCzKif/
9WfB/Xz6vFsyMivf69Fgb
+1F0XGtKMRC9F3kJQqebBZp9vvFvnhqYdz1UOSe7St9rUbK0qcx8t4PbnJ0wTRZW+khJzyfKh
+7usqBpjEyzd5Rhm8elvNBLGWvUS9EYMs7Ijvqlh3v1a4B/
5NxHheEqHqspj20Q2hnXokwtb2IE<ETB>E3<CR><LF>
```

```
-> <ACK>
<- <STX>247B8vKmc7yiaeSARPCpW1DiuwXc1Rmzv6H9BtsRg
+67Glib0Uz4qPouv8uKt9aSs6TfqQlJy1xRzca0scXtkP4Gh5KQZUy4eue57D+xQ1cfN2JwMgBSr/
FQQWFuqpHU6DuT/tQFZmtUiM9P1fB/bepVe+mqF+2paune7JU7W6jjmUq2IX5qk3z9Nl/
nnqckSuikoyqtOvs9SBHRIq4sxUZRuzXc3VG1Ru/2iVax+u+s9aog<ETB>45<CR><LF>
```

```
-> <ACK>
<-
<STX>3znA5XvVahzKz1Uq994ZRC3SpU0DFPDF9mof67RqVM7LNSKAZ3UkR870XtHO65at8Fqeyv
cTC34fv4t9NZfQ9qMBpQF16Pu/b8jveQCork/
xYO8Mqw4dBID447hY8+DmLo8D9cG5udMmFvmQrbm0mYdzgOzqUbkHkxHOFBjQc0KghNK3x
wes4M7C2YhMg1o/E8ZgR6PbDGfkr1Pn3QO/iZ7QusoVS1jVQ8sq/
Uq5WRvXXjp<ETB>5E<CR><LF>
```

```
-> <ACK>
<-
<STX>4CvVS45n91JxUmJtFIFk7tfKDX0+4C2fupNHqs86EKIC7navkMXY4dSgnEggTy3pFeXetKim/
+Sd6/8LB9PbpaFExvl0p41csafK6SsOinvPymUt7flywXHs
+V3xzPktJBU6T82WcYzKqTr8ETZUplgy3skylC18TJir0Gq+3Hce5xsc9kiA5dvlm89XS9TnkbK/
Hth0vrcSulya4nMLVkg9yel90DfeWxQCGjLad3O3lo5A<ETB>2E<CR><LF>
```

```
-> <ACK>
<- <STX>5RpU+ssW75ylPMSH/H5rOVTPytZds1cvrrcKZhR2A1qFswJr6uEd
+OqBC7N50UA01HxKWyg2J09gHR8tV+MS/RKC4s2SW
+HJlqun5OEufnGwTzdN8THisiRfr0hSkvcsSLogmC/RT6OhvRYbQUfgNea1YRdzR2XVpdITj/
9GcdqZr+rpMrcIqS+vjt+e/29/f+87naHEXpu559fBE9/bN8Youfd2/
vZlnDr4vk7RXI2O09c<ETB>B4<CR><LF>
```

```
-> <ACK>
<- <STX>6+WaXzParquvRE2Mva5OTPE+6YrTlxiQtY8MX2k/
56drsbKMW73hISzCe1f4N<CR><ETX>C0<CR><LF>
```

```
-> <ACK>
<- <STX>7M|3|HISTOGRAM|PLT|PLTALONGRES|FLOATLE-stream/
deflate:base64^Y2BgYDA2VnYCUgrJjOBaAcGhgZ7Rh5sBgA|FLOATLE-stream/
deflate:base64^jdR7UNRVFAfwUB0UdQARdBFUBGjGmzV1Nzf75xrtWnMGFg
+cjO<ETB>50<CR><LF>
```

```
-> <ACK>
```

```
<- <STX>0RAhtHB1HTMBWtRURzMSVFkzAUExHTBN1IHdzGN6RmumRARqgh
+EBhEU2o7+83g1F/tTn7uVy7/mc+9jVaDQavX6AQKmp2NARbbKs+edD
+Hus0tE6BkhaR45ksbvQDpBLSyNli32hPMRkkTEu10XZMFYqW41OjLtkm86Dhpi6U2ZmH8J8WuL
7AtVfJSSTRRDWkkE7mazGGNKZZxPyKa7ZdGZKMKRhnYwCSRuxNpOibVIYn0N<ETB>18<CR>
<LF>
```

-> <ACK>

```
<- <STX>12Vx5yWEmpV4A8NtqRcAS57ORrPYI8pWSuu4CcDnKFVSCvk
+LibiJ3HTly6pHfRUbnYxgaNlk6MmpjX6snJxi0XFbalVEnm
+t6sjPjIw3a3oyauSmsL0fbdGw1BjPq57i4ULa7wlhnDmfshR05EazX69liH87YFxudo3IHgoFxDZAEj
FdhGGGMZ
+wXxgQYUTDegjEZxIR2hZlgTlCRaYmWRrhyMWTBmw5gLLwHGfBgLYSgkQ<ETB>09<CR><LF>
```

-> <ACK>

```
<-
<STX>25CUthJGPNcjrse5TWGasTYG3CutTYaZxkHkN9rYWedBtiBXOvz1yLcBNWxEzgzUsQI5
N6OWTOTegj1/hT1sQ01ZMLJhblfxDYwcGDth5MLYBSMPRj7q3AOjAIYVxrcw9sHYD
+MAjEYRTAOWbDB
+B5GMfZ8GMYRGEdhIMA4BsMO40cYx2GcgHESxikYZ2CcxmVvwwJxjkY52FcgHERxi8wLsFww
Cih8SuMKzAqYFTCqJl<ETB>1B<CR><LF>
```

-> <ACK>

```
<-
<STX>3xFYTYxjUY1TBuwLgJ4w8YNTBqYdTBuA3jLox7MOph3IfRAKMRhgvGQxjNMB7BeAzjCe6i
BUYrDI1w5GiE1dhBmOs6CLwxode7Ca3DttIT3AXem7DYPURcXCdh0HoKvD1RF9VZ2F2dRWZm
F4F3KlxOrQgye4mmsK4Cb1LsSOgmlvh6i2ibt8D7xPeyB4weMhRC6AnjGRg+MHxg+MLwg
+EHoxeM3jB6w/CH4Q+jD4wAGAFczw6<ETB>02<CR><LF>
```

-> <ACK>

```
<- <STX>4E0RdGXxj9YOhg6GAewegPoz+MYBjBMEJE22/D//
0M8vrOYEr1kFZHnplGeb0oTzy3QDbqrsuvhQykMXMmUvngpXRnfRblxxaT+4PfyP+L8r/
lIADpg7lQVrBQ3+I5urf3uW0cYs58XYKx1zN4C4TsvlEcaFbcm08a8QJdp/5M
+8vvsjffKoZuTnXs5afm/
GAp69w8YIbTVy28AnfW9OqRtq4Zs74s5kbr7eyeaBGKL9x6Zc<ETB>31<CR><LF>
```

-> <ACK>

```
<- <STX>51onxwC5vrH3LsFRcHvn2fq1bc4YwnNSyOV3PgJCdXrqzskUv/
DE5Wbt5nv1az3N1RRl/
1nCSlw7buXtuCTfUHOWZVUVcYd7LRb328OKxu9g7L5srUrbxurxNTGc28MBuX/
DZT9ZyjNPMpm4rePmmRajxlw71TIQ7j5dvjueDPtN41vBJ7N4wgffMHM8vzx3L1yoN7Ag38Pxb15A
7nCOOhrBbUxBffr4fK+dIWuXhZqvDOKu<ETB>B8<CR><LF>
```

-> <ACK>

```
<- <STX>6wC8eP7sTe
+RoueO8xpY2vobufVxFqpex9pylgjFatvUQXavaS6aUneTW9DXOewtFXdplza9YqMh/
jRpZhasp9c1U9X5QI6E+yrZ+SEbd+8g3XW0DJ8VQRssUOUj7BkUcFmos2/
oSutFkaht1cRjhZimxNhr3rCPlzuNH9aLR83zle3cPqr1qtO8X
+XsR6pYfSffkbOst9a20bzNaytXafp4G9iePmXNRHcc7kx3hr2S8OzX<ETB>85<CR><LF>
```

-> <ACK>

```
<- <STX>76vLNdLjm9VY3E2ly5NKng6ZjSxo/6UvbenS5jz+rbx7/FasjX0c
+A9aGqn3sXfaoX67OVdq2Pt6yOk
+JtjVKKOPKmrB5plSz2m8fitO2vuR0oGrhTJ6GMvbfUOa09fH9Qp5GCecgKW37aP8/
JZTvYfvAGalXO6RYcSVOl9zSM6WRhSmSe3KsNGf3MAI3KbmXREruM6Zjbsn9MOdZ6W8=<CR>
<ETX>61<CR><LF>
```

-> <ACK>

```
<- <STX>0M|4|HISTOGRAM|RBC|RBCPLTALO<ETB>E0<CR><LF>
```

-> <ACK>

```
<- <STX>1NGRES|FLOATLE-stream/deflate:base64^Y2AAgXpnMOWg7AliGRga7BkYBJxAQgA=|
FLOATLE-stream/deflate:base64^rdV/iN9zHafwD26MzTlztuHwtQ0bw5mZYdzndI
+MHfti5nDsMAxnDsMwfGPNJeliaUnrG4traV1C19K6InRjuiRGWpektaSF1rWWPD53r9u
+vuYP5V2P3p/3+/v78369nt8kycbzTcNTO<ETB>A5<CR><LF>
```

-> <ACK>

```
<-
<STX>2jWfHBhpkhRH9pNiw8g6zfYo0ccgSWOS1JCjnpQCbXRQpJMu1lOimx566aOfAbYzyE52M0
QyP0mqGMt4aqhlMnXkmMZ0ZILPbOYyj5Q8C2imwGJaaKWNZSynnQ5WsorVFHmBNaylk5d5hVf
p4jXW8QbreZO32ECJt9nu3SzfzTA8f8CEf08sWPmErFwJuzj6jn8/5gi8Z4Cu+5hu28z0/
sINBfuQnfmYnu/iFX9nNb/zBH<ETB>B0<CR><LF>
```

-> <ACK>

```
<-
<STX>3obYyz7+JFEDh3Aoh1HFGA7nCMZyJEcxjvEcTTXHUMOxTOA4ajmeiUxiMidwldRx8mcwq
nkOI0pTGUap3MGZzKdGZf2czkHM7IPOo5n1lCwGwuZA4XMZeLuYRLmcdIXE4DKY3Mp4k8V3A
IV7GAq7mGhTRzLdexiALXcwM3spibWMLNtHALt3lbrdzOHSyljTu5i7tZxj3cy30s534e4EHaeYgVPE
wHj/Aoj7GSx3mCJ1nFU<ETB>CD<CR><LF>
```

-> <ACK>

```
<- <STX>4zzNM6zmWZ5r2p8Lw/2f9XJr1Mic+A2FuA+z4vz2Rk1uir5oid7Mengger4rcqEQeZGL/
BiMTCmVzaXImzQMZ1LDAUn691wqz6TRORefVYy8yfKINvKhPb7n5uinrCeQo2bc+
+Jaxnf9jrmHbXzH701J35h8UqzOJ32TzHvYm2eeYk2Wo0Xr1H7R4+IEJDOuspz93nPXdnW/
pDbzV/xHuu15f4kXWxFmsiHnFgVXRGrWjY<ETB>6A<CR><LF>
```

-> <ACK>

```
<- <STX>5GpuYdTcoqi/xqj5GdEvWT9lvVdTJuvTrI/3RB7siFzZEtm0MXJrQ2RbZ+Rhe  
+Rnc2RrlrFZ7mY5nN3Xocjx3tATuiP7u+Is2sr+lyqv0zi3pMJozZSfc+Xj5f7tNeW1NlpxYqaK6X/  
rMHkIMprsbxWh/fKRraurOPKcbC9/zL2f5f/YYz8pr8A<CR><ETX>57<CR><LF>
```

```
-> <ACK>
```

```
<- <STX>6M|5|HISTOGRAM|DIFF|EOSALONGABS|FLOATLE-stream/  
deflate:bas<ETB>64<CR><LF>
```

```
-> <ACK>
```

```
<- <STX>7e64^Y2AAgXpnMMXQ4AgkHEAsAA==|FLOATLE-stream/  
deflate:base64^7dV7aI5RGADw130YhrkPw9xvcx/Gvnp7de5zn/ttGOY  
+LEILkiRJ0plkSZIkSUuSJEmsJC1JS9KSJEl+H31t6fvyrz889et0Tud5znnq7T1BEI79Wb  
+GoCgzqlyQeWQ94/c8FF6jhFLKCOQkkEwqlbLJJT/  
zd81ijnGKEi5ylZuU8oAnvKCMc<ETB>66<CR><LF>
```

```
-> <ACK>
```

```
<-  
<STX>0ir4Gr7PqCCoSRzxJJBIS5JIJoVu9CKVgaSRTojRjGcy2cwih4Xksol15JFPATspplgDHOQX  
RzmCEc5xnFOcJJTnOYMZynhHOe5wEUucZkrXOUa17nBTW5xmzuUcpd73OcBD3nEY57wlGc85  
wUvecVrynjDW95Rzns  
+8JEKpVgZL3zlG9/5QeAbqEZ1aICTWtSmDnHUpR71iacBDWIEAo1pQIMSaUzZwTcSVrSmDUM0  
pR3tS<ETB>C7<CR><LF>
```

```
-> <ACK>
```

```
<-  
<STX>1aYDHelECp3pQle60Z0e9KQXvelDX1LpR38GMJBBDDGYIaQxIGMNJZwQjySBEJqPIYjRjGM  
s4xjOBiUxiMIOYySymc4MZjKL2cxhLjnMYz4LWMgiFrOEXJayOWsYCWrrWM0a1rKO9eSxgY1slp/  
NbGErBWxjOzvYyS52s4dC9rlvq/K/8D/+x78Y4Xcplv+VRV5q8JJEEXV/Fii5VU9K1adWHmR/  
VF7+Euf0fKj7YIVP9YZf<ETB>7C<CR><LF>
```

```
-> <ACK>
```

```
<- <STX>29433FdZuE7mTw==<CR><ETX>BC<CR><LF>
```

```
-> <ACK>
```

```
<- <STX>3M|6|HISTOGRAM|DIFF|EOSALONGRES|FLOATLE-stream/  
deflate:base64^Y2AAgXpnMMXwwAFigDADAA==|FLOATLE-stream/  
deflate:base64^7dT7a05xHAfw4z73ud8Z5n6b+zD22IZhmPudud+Z  
+zAsSZK0JC1JT5KWPcVpSVqSliRjKqQISZIkSZK8ZnyyNMo/sHe9+n7P53u  
+dX44vYOgP<ETB>79<CR><LF>
```

```
-> <ACK>
```

```
<- <STX>4IdTfy9BWSj4E/  
u8ynlecsVzqHxGmJJQ5fvjgyCWOBIkUkW2eRxnHwKCFNIEcWUUMpjinPGOz7xjSAICGoTQyNiaU  
lbOhJHPL3ptwLDSCSJEgmkk0Emc1jAerJYzXo2k80ucsglijMc5RjHOcFJTpHPac5wlgLOcZ4LhLn  
IJS5TyBWuCo0irnODmxRzi9vcoYS73OM  
+pTzglY94zBOe8oznvOAlryjjiNW94yzve84GPfOlzX<ETB>8D<CR><LF>
```

```
-> <ACK>
```

```
<- <STX>5/jKN77zg58E/  
oEa1KQWtalDXeoRQ30a0JBGNKYJTYmlGc1pQUta0Zo2tKUd7elARzrRmS7E0ZVudCeeHvSkf73  
pQ1/60Z8BDGQQCQxmCEMZxnBGMJJERjGaMSQxlnEkE218KaSSxgQmMol0JjOFqWQwjenMIJO  
ZzGl2c5jLPOazglUsYjFLWMoylpPFClayitWsYS3rWM8GNrKJzWxhK9vIZjs72MkudrOHveSwj/  
0clJeDH<ETB>79<CR><LF>
```

```
-> <ACK>
```

```
<- <STX>6Er90wvVqUikB6P7MClYr9qVVYwJ9tH3I2flnVpSKRy1hkN/d284cjf57+/8I+iz6ITn//  
kF<CR><ETX>29<CR><LF>
```

```
-> <ACK>
```

```
<- <STX>7M|7|HISTOGRAM|DIFF|LYMALONGABS|FLOATLE-stream/  
deflate:base64^Y2AAgXpnMMXgBKIdQCwA|FLOATLE-stream/  
deflate:base64^7dT7a45hGAfwB8OchzmfmhmfhmEOs3fmbBiGYRjmfBrmfHqTJEmSJC0tS<ETB>C  
D<CR><LF>
```

```
-> <ACK>
```

```
<- <STX>0ZIkSUvSkiRjKtaStCQtSZIkST7vPPIf+Gnf+vR0Pz/ddV/XNwhiOzXv8wkmhN  
+aRIlgGp6jGX/Okdg/SimniiAzCBJlIoUIORRQRJTjnOY8pVzJmWU85jnVFJFNZ/  
5nhm7VBDEEU9TEkikPZ1Jlpm  
+DCSFVNIYR4SjTCWbHHLJl58CClnHJooozG8HiHKEoxzjOCc4ySIoc4aznOM8FyjhlqVc4jJXuMo1  
mODm9ziNnco4<ETB>19<CR><LF>
```

```
-> <ACK>
```

```
<- <STX>1y73uE85D3jllx7zhKc84zkveEkFlbziNW+o4i3veE81H/jlJz7zha984zs/  
+MkvAjNqH7rUl476NKAh8TSiMU1oSjOa04IEWtKK1iTShra0oz0d6EgnOtOfFnQjie70oCfJ9KI3fehL  
P/ozglEMYjBDSGEowxhOKiMYySjSGM0YxjKODMaTQYTMrd  
+7kcVEJjGZKUXlGtOZQTYzmcVscpjDXOaRy3wWsJA8FRGYJeSzlGUsp<ETB>C6<CR><LF>
```

```
-> <ACK>
```

```
<- <STX>24AVrGQVhaxmDwtZx3o2sJFNbGYLWyliG9vZQTE72cVU9rCXfeznAAc5IPWvF/5X/  
vZMrF9Swr6ID/e0JNyD2LwmhO  
+dHr5LevjOPcMZi81qRbhbX8IOiHVHVdgd3pWENrCQuiVvXedhtNXeoTW1ql78B<CR><ETX>3B<C  
R><LF>
```

```
-> <ACK>
```

```
<- <STX>3M|8|HISTOGRAM|DIFF|LYMALONGRES|FLOATLE-stream/  
deflate:base64^Y2AAgXpnMMXADqldQCwA|FLOATLE<ETB>BE<CR><LF>
```

```
-> <ACK>
<- <STX>4-stream/deflate:base64^7dTpi05hGAfgg7Fvr31n7Dtj383L2Nexj33sO2Mf
+5skSZqkaZI0SZIkSZKkSZIkSZIkTZIkSZIkYXV4iv9Ayl1XT/
f5dD787I8UxbMv6+cTVQzVz0IGUSrsqcxfezL
+RjEIIbINi6IE6WSQJJtc8kxhiAKKKOYcl7hGCXd5yFNKecMHvAnJ6I0KIGNBHVpSFPSaUMHupB
Bb/ozmCQjGMMes<ETB>AC<CR><LF>
-> <ACK>
<-
<STX>5plODvPIZSkRWUseW8hnNyn2c4CDHOlwRzhKAcc4TiFFnOAKpyjmNGc4yznOc4GLXOlyV7j
KNa5zg5UuUclvb3OEu97jPAx7yiMc84SnPeM4LSnnJK17zhr84z0f
+MgnPvOFr3zjO5EMIKes5UijPBWYfuWIEpWpQlWqUZ0a1CRBLWpTh7rUoz4NaEgjGtOEpjSjOS1I
pyWtaE0b2tK09nSgl53oTBe60o3uZNCdNvSiN33oS<ETB>AD<CR><LF>
-> <ACK>
<- <STX>6z/
6M4CBDGlwQxhKJkmGMZwsRjCSUyxmDGMZx3gmMJFJTCabKUXlGtOZwUxmkcNs5jCXecxnA
QvJZRLGWcJSirGcFaxkFatZw1rWsZ4N5LGRTWxmC1vZxnby2cFOdrGpPezN
+tOL3Pi7on7JhF6pSD0RVq4+Zxwo4Xhpu6FnMf5TQsZrRjyGuf4c8j+i3BXp8Mt54d
+iPsk7pzS0FGpP7osN3Rclviz+9LDHv9r3JX/5//8+<ETB>1C<CR><LF>
-> <ACK>
<- <STX>7/MD<CR><ETX>07<CR><LF>
-> <ACK>
<- <STX>0M|9|HISTOGRAM|DIFF|MONALONGABS|FLOATLE-stream/
deflate:base64^Y2AAgXpnMMWwwwFIOIBYAA==|FLOATLE-stream/
deflate:base64^7dT5a89hAAfwj3vuuW
+GuYcv5j72nQ1zz30Nc8899803SZK0JC1JS9KStCQtSUuSJEmSJC1JkiRjkuT15aMtfwHIXa
+enqee44endxDecyDj5xCUpgcViQZBL<ETB>C2<CR><LF>
-> <ACK>
<- <STX>1FyPpf2aR+NrFFFGOYE9iSQRIUo2ueQT4wgFFFJEMSPxv+4r4y4PeUo5b/jAl/
h7RgdBdRkOryJNaUlbbkimOylIESGUll4iSSRaTyGYmc8khl
+XksY58trKTvcQ4yCEOc4SjHOM4BZzgJKco5DRnOEsR5zjPBYq5yCUuU8IVrnKNUq5zg5uUcYvb
3OEu97jPAx7yiMc84SnPeM4LynnJK17zhr84z0f+MgnPvOFr3zjO4E/U<ETB>A1<CR><LF>
-> <ACK>
<- <STX>2IWqVKM6NahJLRkoTR3qUo/
6NKAhiTSiMU1oSjOa04KwTKl1bWhLO9rTgSQ60onOJNOFrnSjOz3oSS9S6E0f+hKhH/
0ZQCoDgcRghjCUYQxnBCMZRRpR0hINBpmMYSzjyG18E5jJCYzhalkM43pzGAmS5jNH0Yyj/
kslleFLGlxuSxhKctYzgpWsoo8VrOGtaxjPRvYSD6b2MwWtrKN7exgJ7vYzR72so/9GRW98FckWtE/
l<ETB>FF<CR><LF>
-> <ACK>
```

```
<- <STX>3Tvodw/97qJl2D2Ve6eokoKwf0pDJeG8INwXCbss3mnx88rCu2KVxij4lrggFO/IP/3P//zb
+QE=<CR><ETX>D4<CR><LF>
-> <ACK>
<- <STX>4M|10|HISTOGRAM|DIFF|MONALONGRES|FLOATLE-stream/
deflate:base64^Y2AAgXpnMMUwwwFIOIBYAA==|FLOATLE-stream/
deflate:base64^7dTpi01hHMDxY98ZjH0b+9jHvnPHWAaDi7FvY98Z+85N0iRJ<ETB>20<CR><LF>
-> <ACK>
<- <STX>5kiZJkyRjmiRNkiZJkRjmiRNkiRjkiTJ53Av88ILf8D86tPznKdzufF6RsE4RzK+LUE
+enB34kEQSx+Hhv9+zoSnIFAMWUEnkkihTQiRMkhlxh5nEz//
f4CLINIEcU8oIRSynjLR76G3zMmCKpSk7okkUxzWpNCJ1LpSRoDGMllowkyyiZDOXheSwgVslJft
7GY/MQ5zhKPkCyzjnOAKpzjNGfl5yznOU8AFLnKjy1zh<ETB>C7<CR><LF>
-> <ACK>
<- <STX>6Ktco5Do3uEkRt7jNHYq5yz3u84CHPOlxJTzhKc8o5TkveEkZr3jNG97yjd84COF
+MwXvvKN7/
wg8A9UojJvQeO1qIODmtSiNnWoS3q04AkGtKlxiTThKY0ozktaEkrWtOGtrQjhfZ0oCOd6EwXupJK
N7rTg570ojd9SKMv/
ejPAAYyiMEMYSjDGM4IRjKK0URIZwwZjGUc45IAJhOZxGSymMJUphFIOjOYSTazmM0c5jKP<ET
B>07<CR><LF>
-> <ACK>
<- <STX>7+SxglYtYzBjYwMoylrOclaxiNwTYzrWs4GNbGluWxhK9vYz52sovd7GEv+9jPAQ5m/
O1CxfzfHJ0MGxl2MSXewrR4GxONTLQxJ76PlmtmtFxDY/Fe5sU7mngmUu6ecA37m+hwYh
+2ueAfYnHhNwZx5fe/Ol8xFfNnfgl=<CR><ETX>A6<CR><LF>
-> <ACK>
<- <STX>0M|11|HISTOGRAM|DIFF|NEUALONGABS|FLOATLE-stream/
deflate:base64^Y2AAgXpnMM<ETB>FD<CR><LF>
-> <ACK>
<- <STX>1XwyAllOIBYAA==|FLOATLE-stream/
deflate:base64^7dVraFZ1HMDxU2mprVx5Sc1q5qp5bdXUectnm
+aylVPnJZ22atXyOnWa5sqHkBgHMiJERgKlHmilMSJGhIwxZjIjGCJDJEZISIQMERER6XPc/
+Dzple9rD980oc5z5tz+31PFMvrf8m9Tfr7cXR/
paloHY6nFwz8tSxHaKdPqKiKMomh3xSIFNFLWkaaOQITzy<ETB>F3<CR><LF>
-> <ACK>
<- <STX>2khTba6aKbXvq4Sj
+3iZtIlaQRTYjGcN4csglj6nkU0A80ixkFLKKKeCNVRSRTU1bKaWOvZQT5qvOMDXNPANBzIEI9/
yHYc5wiG08T1NHOCeP3CSUZtZly208hM/08Yv/
Mpp2umgzN0cZznKebHi5wkV4ucZnf6CseeK5X+IOr/MlfxKOf69zJre4zR3uxu+Dd
+ABHuQhBjGyH3mEIQxIGl+SxWM8znCyeYInGcFIRjG<ETB>34<CR><LF>
```

```
-> <ACK>
<-
<STX>3apxjDWMbxNON5hmd5jhwrm8DwTyeUFXuQl8pjEZKYwIWIM52XyeYVXeY0CZjCTWRQym
znMZR7zeZ0FpCiiMBIWsoG3WEwpb7Ketyjbd5hKeUsYzkrqGAlq1jNgT5lLeuoZD0bel8q3ucDPqSa
j/YT6jhUzayic1sYSvbqGU7O9hJHbvYzWfsYS+fs496vuDLkvtd+M
+s1EDb4q5lh54ILUuHfrWFvVWFPmWF/uSH3pSGttS<ETB>2D<CR><LF>
-> <ACK>
<- <STX>4FRiRNOBZmP5n11ox57wgzfTbMeGeY7WR2L2XM9sUgmfuesN
+d8V9XaEhL6EzSmubQoaZwTodCyw6E3iX9i1sYdzHuZHkQX1fc0sKwLQjXHxc5s9H94f60hca3hP
vWkCEdVIX7m/ndSKTCNyX6B/EzSqGnlki/k5l+n/92/U3<CR><ETX>4F<CR><LF>
-> <ACK>
<- <STX>5M|12|HISTOGRAM|DIFF|NEUALONGRES|FLOATLE-stream/
deflate:base64^Y2AAgXpnM<ETB>CA<CR><LF>
-> <ACK>
<- <STX>6MVQBKIdQCwA|FLOATLE-stream/deflate:base64^7dSji05RGAfga98GY98GY98Z
+24+xr4Oxr4Nxr4Nxr5NkiRpkIRJkyRjkqRj0iRjkiRjkr4kSzlksZL4noyW/
APirafTPfd265zznI8UJWpXrchWhfGbxWlovzwnJ/+
TmWmKOOYUjEg6MomVTsJfJNrnks48CjLIGS5QRDg3uMcyj4rzklR
+JhkRRWSqSRDK1qU<ETB>21<CR><LF>
-> <ACK>
<-
<STX>78KqbSkLR1Jowd9GECMoYxkLJlKMZ3ZZJPDUIaSSx6b2U4+u9nDXvaxnwMcpIBDHOYIRz
nGcU5QyEiOcZoznOUc57nARS5xmSKucJvRFHODg9zkFre5w13ucZ8HPOQRj3nCU
+I84zkveMkrXvOGt7zjPR/4yCc
+84VID5SiNGUoSznKU4GKVKlyVUiiKtWoTjI1qEktaOHutSjPg1oSCNSaEwTmPJkM5rTgpa0jVtaE
s72t<ETB>CB<CR><LF>
-> <ACK>
<- <STX>00BjnSiM11loyvd6E4PetKL3vShL/
3ozwAGM0YgxmCBkMZRjDGcFIRjGaMYxlHOOZQCYTmcRkspjCVKYxnRnMZBazmcNc5pHNf
BawkBwWsZglGUZy1nBSlaxmjXksbj+91eTx4b2MgmNrOFrWxjOzvYmfEjF/
62SuTV75iVklfJlaNyQx4I8YI5ENeuMcXw516E/o/KfRm+9ALWeGs5oW9L9nfVeEMFobzHBO
+7x16My<ETB>14<CR><LF>
-> <ACK>
<- <STX>1X8K/
HP1+HOXgt5cCjkTU7lqLSQdfGQkUUHmwtCtibyNh7WVRjWmJBYc8n6/+T391Hsh1/mf9rL//
Ww1lc=<CR><ETX>93<CR><LF>
-> <ACK>
```

```
<- <STX>2M|13|HISTOGRAM|DIFF|BASOALONGRES|FLOATLE-stream/
deflate:base64^Y2AAgXpnMMXg4QQkHIAIcC4MjAoOCPkGuxBcgA=|FLOATLE-stream/
deflate:base64^7dR7aFZIHAFwY800nbXK+63l<ETB>5A<CR><LF>
-> <ACK>
<-
<STX>3LUvNWd6yqcfNzEptlZm3cnnJeZvTps37S0SoRAyRGCL2lHlWpmlhEhEjQkQkRCRGIaYJG
ClhESli4ue45+DL+4f/2V898OF9zu/8c97n93u
+UZSsbeV3fqKpZdHdFUdRjTzE1uf46RGIkaaiSZFURHFIBBTQSU1ZnHOPQ1kOchRTtDIKc7SRD
MtXOU6kW8qoD2FFNGZ7vSmmAEMZigjGQspcRMLmv9b9OoYcZmU8li6li<ETB>FB<CR><LF>
-> <ACK>
<- <STX>4JTXUUsdmMnzCp3zGdnbyOV9Qzy528yUN7GEv+8iywnN8zUG+5RCHOcoxfuA4J/iRn/
iZrn7hV05yitOc4TfOco7z/E4Tf3CBizRziT/5ixYuc4W/uco//
Ms1rnODm9xK5sEMtOEBHqSATjxE09rzMB3oSCGdelRHKelxHucJOtOfmSjOz3oSS9604e
+PEkxT9PG/gxgllN4msE8w7MMYSjDel7hIDCC53mBkYxiNGMYy4uM<ETB>43<CR><LF>
-> <ACK>
<- <STX>54yVKGc8EJhlziTLKmczLTOEvpvlqr/E605jODN6ggjd5i7eZyTvM4l1mM4e5zGM+7/E
+C6jkAxayIMUs4UOWUsUylrOClayimtXUs1a1fEQt61jPxp9SxgY1sYjNb2Fp
+Nxfu10pzJc2UNFdSab7EITMKwn3eHwa7V+ih88jIqkw9X3n
+hkN8x2G17znOEBLsUt9R3ppv1eFcZ4UeDwuz1jfiMYjKnyTwn857OfnInroR71hTu<ETB>F9<CR><
LF>
-> <ACK>
<- <STX>638lwfW+EPKgNeVMasinJrZaQd9mQh2kmNoR9TcjOOE
+arShvn2ZuKsqpFefUmsM5ZvNkgjvhH+Wlc2pJn9Je5Uvf3Utuv/Nr+fx8d/+v/
34IPUh63xjfbg==<CR><ETX>0E<CR><LF>
-> <ACK>
<- <EOT>
-> <ENQ>
<- <ACK>
-> <STX>1H|^&||YP8K|||||P|LIS2-A2|20230929093211<CR><ETX>B6<CR><LF>
<- <ACK>
-> <STX>2P|1||00000020||PATIENT 20^TEST||19890613^34^Y|F<CR><ETX>C1<CR><LF>
<- <ACK>
-> <STX>3O|1|2023092700000020^1^042249^10||^SLIDE_EC|R|20230927174552|
20230927174552||N|||BLOOD|||||||O<CR><ETX>4E<CR><LF>
<- <ACK>
```

```
-> <STX>4L|1|N<CR><ETX>07<CR><LF>
```

```
<- <ACK>
```

```
-> <EOT>
```

### 3.6.5. Example of Result With Nothing Left to Do

```
<- <ENQ>
```

```
-> <ACK>
```

```
<- <STX>1H|\^&||MHR1^210M2SH01011^1.7.0||||P|LIS2-A2|  
20230302102751<CR><ETX>E3<CR><LF>
```

```
-> <ACK>
```

```
<- <STX>2P|1||00000000||PAT 2303020020^PAT 2303020020||19250709000000^97^Y|  
F<CR><ETX>1F<CR><LF>
```

```
-> <ACK>
```

```
<- <STX>3O|1|2023092700000005^1^041176^1||^DIR|R||20230302082259||||BLOOD||||  
AP^STANDARD(m)|||F<CR><ETX>71<CR><LF>
```

```
-> <ACK>
```

```
<- <STX>4C|1||S^DIFF^WBC_ABN_MAT^SEP_NEU_EOS||<CR><ETX>82<CR><LF>
```

```
-> <ACK>
```

```
<- <STX>5C|2||C^LMNE^NeuEosSep||<CR><ETX>DE<CR><LF>
```

```
-> <ACK>
```

```
<- <STX>6M|1|HISTOGRAM|WBC|TNCALONGRES|FLOATLE-stream/  
deflate:base64^Y2AAgXpnMnwWCEQ7ABEQnwBxgSMQQ+QYGuxBcgA=|FLOATLE-stream/  
deflate:base64^xdV/  
qF5zHAfwg2EYu2zY2OzaLobhwcYw7rm7G8N2XVyzMVzMXHbxjGEYHhvrzq9u0lzMui2tpbWWpJu  
WniQtSuVsknSTtCTdtLS0ltf3Pp/pW<ETB>4B<CR><LF>
```

```
-> <ACK>
```

```
<- <STX>7P4Qyrdefc45z3nOeZ5zPt/3N8vSeLp5sFTW12pt5A7EfQWxtp+nY/RSpZ  
+sKcvqqKdETivtIKnQRTc99LKRLfRRZRvb2UE/  
OxlgN9n0LBvCUIZRx0hGMYZ6GpjJEpMZirTyJnBLGbTshvWEA7C  
+mgkzJLwCzYkqzgWVbRxfO8yMt08wqvsoYe3mAt6+hIPW+zgY28wyY2s4V3eY/  
36eMDtvHVT7iYz5hG5/yGZ+znS/4k<ETB>7A<CR><LF>
```

```
-> <ACK>
```

```
<- <STX>0q/Ywdd8w7f08x3f8wM7+ZGf  
+JkBfmEXv7Kb39jDXjl9cAAHchBDOJhDOJShHMbhHMEWjuQohIPH0RzDCEZyLMDxPKMYzQmcy  
BjGchLjqOdkxjOBBk7hVE5jlqdzBmcyibM4m3MocS7ncT6TmclFXMHULuJiLmEal3IzjeQ0MZ1mZj  
CTy7mCWVzJVVzNbObQwjW0ci3XcT1t3MBcbmQe87mJm1nALdzKbbRzO3dwJwu5i0XcT<ETB>  
A9<CR><LF>
```

```
-> <ACK>
```

```
<- <STX>1Qf3cC+L6eQ+7ucByizhQR5iKQ/zCI  
+yjMd4nCdYzpM81fxHLvyjkTij5woxfzPYs5tjn4bG8+/JX5/uv8zrOS5mkraXhG/  
dVE8p7nXncZ4zxOiT1lfpX5LvTci+ml89MSUeJct8bw74prpHl3qS6zhZrKWIUeDmZiqut403aP  
+nrt88FzXqPb9gvqanVV/lcl8a7mRw/MjB5K/ZV6MvV5miv75lGaV7tiDm  
+NXNgU2dEde<ETB>E8<CR><LF>
```

```
-> <ACK>
```

```
<- <STX>2dQZ2VSO3GqLfGulTEwZmZ7xQGRq3345W4387Ypc3qcc2iO7SwV51JTxWdOf874+9vsL6  
0E11odKrBd5bFcK60bx8yyv9Upxfmr41nhOsVziuftr9lXf2/07v/dvwx1/j/x  
+8=<CR><ETX>10<CR><LF>
```

```
-> <ACK>
```

```
<- <STX>3M|2|HISTOGRAM|RBC|RBCALONGRES|FLOATLE-stream/  
deflate:base64^Y2AAgXZnEPn2+3wXIOUAWRIODAxSYHEGhgZ7AA==|FLOA<ETB>6C<CR><LF>
```

```
-> <ACK>
```

```
<- <STX>4TLE-stream/deflate:base64^dZR5UBRnGsZRFBHHcZArHkBiSFDxABURIf6+9/  
k4hHEREE4kHXiErGliDB0ZciAHBjUgQcFFGEGE8UllOBNFALFDWZSOwWsb6RROlikfwCKZYUcO  
+spWt/Wenaqa7Z95v+vf08+s2M3vz2o43n20deYo3wux/vvvPcVz3/  
iaHBG2Tg7lwL1aC90Vwq4n3q4Xtrt8bC7DEmxcOso6RMwRf<ETB>A5<CR><LF>
```

```
-> <ACK>
```

```
<- <STX>5KMDLL5UNpu/ESa39wgeV5GBaXwTKZcXrqX50wy0K6lZ0ukV2wFz1fL8c21vKZJjvC  
+xevuyn5H2njtr7JT9xuvN6cNHn0oKqgfRaweSHwuCs35Ay0vHU6L6x2lz0vzzZwp0M6FAtwmETM  
QhXqQV  
+xM8swSxDzkWuNH45vfozGdc4nZyNFIIY3wDqGhS1cQc9KgneHU70gk9TkrCxmXe2bqFO3hV5  
aG4j5qWNcljNsp3aVzB<ETB>F0<CR><LF>
```

```
-> <ACK>
```

```
<- <STX>6wp9DgkIVnS6GF0OvNk0r20vzDTbmopzGauvXS7KpfZ8uj6jXzmM1FT  
+0FmLKQruqPMWUR1TseZ9QTVeH3NvCVUHXyamc9Qpf4b5q6giuRzzF5JZwrOM381nSq/  
yBlqLjxO85RS8fa6jnLZTps0cB5msjkcJUzXad9Hjc51y0yBv3A2Zopa/  
VPnO8uZSS0csYHlJrziHO2UxLpE876jBLr2znvr2Ro7eDMLjyW7BxN/oOqdG<ETB>C7<CR><LF>
```

```
-> <ACK>
```

```
<-  
<STX>7Z4ad0DnnY9wdcfZ5x6oWNcb7i7WYC7QImXJdpVX7gFWIF7wYlGHR6H9Mf40AHgjnBMP  
wgPowdjTKw1uC8cTrbFvbQhcMp6C9wdTAVD0VJoh7eLh4F7RF75CNyuGomRNfbgTmFsdMT1G  
29jaPM74H6R1fYumtqdYNs5Gtw1MizG4IpuLAbZjAP3ji8dxqPOaQJ0LhPBDuALD1fUeLmhj/  
dksA9IDJqK6mB3mC+dBnYDhtXTUa<ETB>E1<CR><LF>
```

-> <ACK>

```
<-  
<STX>0n3RFf0DLAniE2YhYpkL7xM08DOYFOOXJkCQkchwP5gfak3TpX74JcqX7BL0NfPrnGjPx7f  
CAB7hYjW93GsbQ4etgeCHUOY2TwctpiPe7oPwL4h1G4BTA7BaHb6COWelrktwJ6PxbjtQTsIRYHL  
IUxaBmuBy8HO4ng0XJIWv0xmvShYD8xP/aPyEhYhSvJYWBXMSdrNVJzPkvtQQTYW/  
gXf4bk0jWoKdeDHYZPzVok1q9Dde<ETB>48<CR><LF>
```

-> <ACK>

```
<- <STX>1N6sM+g5mgYWjeisu1zsNvw6oxBrNmfUG4RC/Yc0222oklXh+/T4rh3A4ZnGbj7eCzKif/  
9WfB/Xz6vFsyamivf69Fgb  
+1F0XGtKMRC9F3kJQqebBZp9vFvnhqYdz1UOSE7St9rUbK0qqx8t4PbnJ0wTRZW+khJzyfKh  
+7usqBpjEyzd5Rhm8elvNblGWvUS9EyMs7ljlvlh3v1a4B/  
5NxHheEqHqspj20Q2hnXokwtb2IE<ETB>E3<CR><LF>
```

-> <ACK>

```
<- <STX>247B8vKmc7yiaeSARPCpW1DiuwXc1Rmzv6H9BtsRg  
+67Glib0Uz4qPouv8uKt9aSS6TfqQlJy1xRzca0scXtkP4Gh5KQZUy4eue57D+xQ1cfN2JwMgBSr/  
FQQWFuqPHU6DuT/tQFZmtUiM9P1fb/bepVe+mqF+2paune7JU7W6jjmUq2IX5qk3z9Nl/  
nnqckSuikoyqtOvs9SBHRIq4sxUZRUzXc3VG1Ru/2iVax+u+s9aog<ETB>45<CR><LF>
```

-> <ACK>

```
<-  
<STX>3znA5XvVahzKz1Uq994ZrC3SpU0DFPDF9mof67RqVM7LNSKAZ3UkR870XtHO65at8Fqeyv  
cTC34fv4t9NZfQ9qMBpQF16Pu/b8jveQCork/  
xYO8Mqw4dBI447hY8+DmLo8D9coG5udMmFvmQrbm0mYdzgOzqUbkhkxHOFbQjC0KghNK3x  
wes4M7C2YhMg1o/E8ZGR6PbDGfkr1Pn3QO/iZ7QusoVS1jVQ8sq/  
Uq5WRvXXjp<ETB>5E<CR><LF>
```

-> <ACK>

```
<-  
<STX>4CvVS45n91JxUmJtFIFk7tfKDX0+4C2fupNHqs86EKIC7navkMXY4dSgnEgPTy3pFeXetKim/  
+Sd6/8LB9PbpaFExvl0p41csafK6SsOinvPymUt7flywXHs  
+V3xzPktJBu6T82WcYZkqTr8ETZUplgy3skyLC18TJir0Gq+3Hce5xsc9kiA5dvlm89XS9TnkbK/  
Hth0vrcSulya4nMLVkg9yel90DfeWxQCgJLad3O3lo5A<ETB>2E<CR><LF>
```

-> <ACK>

```
<- <STX>5RpU+ssW75ylPMSH/H5rOVTPytZds1cvrrcKZhr2A1qFswJr6uEd  
+OqBC7N50UA01HxKWyg2J09gHR8tV+MS/RKC4s2SW  
+HJlqun5OEufnGwTzdN8THisRfr0hSKvkcSLogmC/RT6OhvRYbQUfgNea1YRdzR2XVvpdITj/  
9GcdqZr+rpMrcIqS+vjt+e/29/f+87naHEXpu559fBE9/bN8YOufd2/  
vZlnDr4vk7RXI2O09c<ETB>B4<CR><LF>
```

-> <ACK>

```
<- <STX>6+WaXzParquvRE2Mva5OTPE+6YrTlxiQtY8MX2k/  
56drsbKMW73hISzCe1f4N<CR><ETX>C0<CR><LF>
```

-> <ACK>

```
<- <STX>7M|3|HISTOGRAM|PLT|PLTALONGRES|FLOATLE-stream/  
deflate:base64^Y2BgYDA2VnYCUgxrJjObaAcGhgZ7zRh5sBgA|FLOATLE-stream/  
deflate:base64^jdR7UNRVFAfwVUB0UdQARdBFUBGjGmzV1Nzf75xrtWnMGfG  
+cjO<ETB>50<CR><LF>
```

-> <ACK>

```
<- <STX>0RAhtHB1HTMBWtRURzMSVfKzAUExHTBN1IHdzGN6RmumRARqgh  
+EBhEU2o7+83g1F/tTNn7uVy7/mc+9jVaDQavX6AQKMP2NARbbKs+edD  
+Hus0tE6BkhaR45ksbvQDpBLSyNli32hPMRkkTEu10XZMFYqW41OjLtkm86Dhpi6U2ZmH8J8WuL  
7AtVfjSSTRRDWkkE7mazGGNKZZxPyKA7ZdGZKMKRhnYwCsruxNpOibVIYn0N<ETB>18<CR>  
<LF>
```

-> <ACK>

```
<- <STX>12Vx5yWEmvP4A8ntqRcAS57ORrPYI8pWSuu4CcDnKFVSCvk  
+LibiJ3HTly6pHfRUbnYxgaNik6MmpjX6snJxi0XFbaIVEnm  
+t6sjPjW3a3oyauSmsL0fbdGw1BjPq57i4ULa7wlhnDmfshR05EazX69IiH87Yfxudo3IHgoFxZDAEj  
FdhGGGMZ  
+wXxgQYUTDegJZxIR2hZlgTlCRAyMWRhyMWTBmw5gLlWHGfBgLYSyGkcQ<ETB>09<CR><LF>  
>
```

-> <ACK>

```
<-  
<STX>25CUthJGPNcjrse5TWGasTYG3CutTYaZxkHkN9rYWedbBtiBXOvz1yLcBNWxEzgzUsQI5  
N6OWTOtegi1/hT1sQ01ZMLJhblfxDYwcGDth5MLYBSMPRj7q3AOjAIYVxrcw9sHYD  
+MAjEiYRTAOwbDB  
+B5GMfZ8GMYRGEdhIMA4BsMO40cYx2GcgHESxikYZZCcxmVwvgJxjky52FcgHERxi8wLsFww  
CiH8SuMKzAqYFTCqJ<ETB>1B<CR><LF>
```

-> <ACK>

```
<-
<STX>3xFYTYxjUY1TBuwLgJ4w8YNTBqYdTBuA3jLox7MOph3lfrAKMRhgvGQxjNMB7BeAzjCe6i
BUYrDI1w5GiE1dhBmOs6CLWxode7Ca3DTTi3AXem7DYPURcXCdh0HoKvD1RF9VZ2F2dRWZm
F4F3KlXOrQgye4mmsK4Cb1LsSOgmlvh6i2ibt8D7xPeyB4weMHRc6AnjGRg+MHxg+MLWg
+EHoxeM3jB6w/CH4Q+jD4wAGAFcZw6<ETB>02<CR><LF>
```

-> <ACK>

```
<- <STX>4E0RdGXxj9YOhg6GAewegPoz+MYBjBMEJE22/D//
0M8vrOYEr1kFZHnplGeb0oTzy3QDbqrsuvhQykMXMmUvngpXRnfRblxxaT+4Pfyap+L8r/
lIADpg7lQVrBQ3+I5urf3uW0cYs58XYKx1zN4C4TsvlEcAFbcm08a8QJDp/5M
+8vvsjfkZuTnXs5afm/
GAp69w8YIbTVy28AnfW9OqRtq4Zs74s5kbr7eyeaBGKL9x6Zc<ETB>31<CR><LF>
```

-> <ACK>

```
<- <STX>51onxwC5vrH3LsFRcHvn2fq1bc4YwnNSyOV3PgJcDxrqzskUJ/
DE5WBt5nv1az3N1RRI/
1nCSlw7buXtuCTfUHOWZVUVcYd7LRb328OKxu9g7L5srUrbxurzNTGc28MBuX/
DZT9yjnPMPm4rePmmRajxlw71TIQ7j5dvjueDPTn41vBJ7N4wgffMHM8vzx3L1yoN7Ag38Pxb15A
7nCOOhrBbUxBffr4fk+diWuXHzqvdOKu<ETB>B8<CR><LF>
```

-> <ACK>

```
<- <STX>6wC8eP7sTe
+RoueO8xpY2vobufVxqFpex9pylgjyFatyUQXavaS6aUneTW9DXOewtFXdplza9YqMh/
jRpZhasp9c1U9X5Ql6E+yrZ+SEbd+8g3XW0DJ8VQRssUOUj7BkUcFmos2/
oSUTfKaht1cRjhzimxNhR3rCPlzuNH9aLR83zle3cPQR1qtO8X
+XsR6pYfSffkbOst9a20bzNaytXAf4G9iePmXNRHcc7kx3hR2S8OzX<ETB>85<CR><LF>
```

-> <ACK>

```
<- <STX>76vLNdLjm9VY3E2ly5NKNg6ZjSxo/6UvbenS5jz+rbx7/FasjX0c
+A9aGqn3sXfaoX67OVdq2Pt6yOk
+JtjVKKOPKmb5plSz2m8fitO2vuR0oGrhTJ6GMvbfUOa09fH9Qp5GCecgKW37aP8/
JZTVyfvAGalxO6RYcsvOI9zSM6WRhSmSe3KsNGf3MAI3KbmXREruM6ZJbsn9MOdZ6W8=<CR>
<ETX>61<CR><LF>
```

-> <ACK>

```
<- <STX>0M|4|HISTOGRAM|RBC|RBCPLTALO<ETB>E0<CR><LF>
```

-> <ACK>

```
<- <STX>1NGRES|FLOATLE-stream/deflate:base64^Y2AAgXpnMOWg7AliGRga7BkYBJxAQgA=|
FLOATLE-stream/deflate:base64^rdv/iN9zHafwD26MzTlztuHwtQ0bw5mZYdzndI
+MHfti5nDsMAxnDsMwfGPNJeliaUnrG4traV1C19K6lnRJuIRGWpektaSF1rWWPD53r9u
+vuYP5V2P3p/3+/v78369nt8kycbzTcNTO<ETB>A5<CR><LF>
```

-> <ACK>

```
<-
<STX>2jWfHBhpkhRH9pNiw8g6zfYo0ccgSWOS1JCjnpQCbcXRQpJMu1lOimx566aOfAbYzyE52M0
QyP0mqGMt4aqhlMnXkmMZ0ZILPbOYyj5Q8C2imwGJaaKWNZSynnQ5WsorVfHmBNaylk5d5hVf
p4jXW8QbreZO32ECJt9nlu3SziffZTA8f8CEf08sWpMErFwUzUz6jn8/5gi8Z4Cu+5hu28z0/
sINBfuQnfmYnu/iFX9nNb/zBH<ETB>B0<CR><LF>
```

-> <ACK>

```
<-
<STX>3obYyz7+JFEDh3Aoh1HFGA7nCMZyJEcxjvEcTTXHUMOXTOA4ajmeiUxiMidwldRx8mcwq
nkOI0pTGUap3MGZzKdGZzF2czkHM7lPOo5n1lcwGwuZA4XMZeLuYRLmcdlXE4DKY3Mp4k8V3A
IV7GAq7mGhTRzLdexiALXcwM3spibWMLNtHALt3lbrdzOHSyljTu5i7tZxj3cy30s534e4EHaeYgVPE
wHj/Aoj7GSx3mCJ1nFU<ETB>CD<CR><LF>
```

-> <ACK>

```
<- <STX>4zzNM6zmWZ5r2p8Lw/2f9XJr1Mic+A2FuA+z4vz2Rk1uir5oid7Mengger4rcqEQeZGL/
BiMTCmVzaXImzQMZ1LDAUn691wqz6TRORefVYy8yfKINvKhpB7n5uinrCeQo2bc+
+Jaxnf9jrmHbXzH701J35h8UqzOJ32TzHvYm2eeYk2Wo0Xr1H7R4+IEJD0uspz93nPXd7nW/
pDbzV/xHuu15f4kXWxFmsiHNFgVXRGrWyJ<ETB>6A<CR><LF>
```

-> <ACK>

```
<- <STX>5GpuYdTcoqi/xqj5GdEvWT9lvVdTJuvTrl/3RB7siFzZetm0MXJrQ2RbZ+Rhe
+Rnc2RrFZ7mY5nN3Xocj3tATuiP7u+Is2sr+lyqv0zi3pMJoJZSfc+Xj5f7tNeW1NlpxYqaK6X/
rMHklMprsbxWh/fKRraurOPKcbC9/zL2f5f/YyZ8pr8A<CR><ETX>57<CR><LF>
```

-> <ACK>

```
<- <STX>6M|5|HISTOGRAM|DIFF|EOSALONGABS|FLOATLE-stream/
deflate:bas<ETB>64<CR><LF>
```

-> <ACK>

```
<- <STX>7e64^Y2AAgXpnMMXQ4AgkHEAsAA==|FLOATLE-stream/
deflate:base64^7dV7al5RGADw130YhrkPw9xvcx/Gvnp7de5zn/ttGOY
+LEILkiRJ0plkSZIKSUuSJEmsJC1JS9KsJEl+H31t6fvyrz889et0Tud5znnq7T1BEI79Wb
+GoCgzqlyQeWQ94/c8FF6jhFLKCOQkkEwqlbLJIT/
zd81ijnGKEi5ylZuU8oAnvKCMc<ETB>66<CR><LF>
```

-> <ACK>

```
<-
<STX>0ir4Gr7PqCCoSrZxJJBIS5JlJoVu9CKVgaSRTojRjGcy2cwih4Xksol15JFPATspplgDHOQXx
RzmCEc5xnFOcJITnOYmZynhHOe5wEUucZkrXOUa17nBTW5xmzUcpd73OcBD3nEY57wlGc85
wUvecVrynjDW95Rzns
+8JEKPVgZL3zlG9/5QeAbqEZ1aiCTWtSmDnHUpr71iacBDWIEAo1pQIMSaUZzWtCSVrSmDUm0
pR3tS<ETB>C7<CR><LF>
```

-> <ACK>

```
<-  
<STX>1aYDHeIECp3pQle60Z0e9KQXvelDX1LpR38GMJBBDDGYIaQxIGMNJZwQjySBEJqPIYjRjGM  
s4xjOBIUxiMIOYySymc4MZjKL2cxhLjnMYz4LWMGifROEXJayjOWsYCWrwWM0a1rKO9eSxgY1slp/  
NbGErBWxjOzvYyS52s4dC9rlvq/K/8D/+x78Y4Xcplv+VRV5q8JjEEXV/Fii5VU9K1adWmHr/  
VF7+Euf0fKj7YIVP9YZf<ETB>7C<CR><LF>
```

-> <ACK>

```
<- <STX>29433FdZuE7mTw==<CR><ETX>BC<CR><LF>
```

-> <ACK>

```
<- <STX>3M|6|HISTOGRAM|DIFF|EOSALONGRES|FLOATLE-stream/  
deflate:base64^Y2AAgXpnMMXwwAFIgdADAA==|FLOATLE-stream/  
deflate:base64^7dT7a05xHAfw4z73ud8Z5n6b+zD22IzhmPudud+Z  
+zAsSZK0JC1JT5KWpCVpSvQsIIRJkqQISZIkSZK8ZnvymMo/sHe9+n7P53u  
+dX44vYOGp<ETB>79<CR><LF>
```

-> <ACK>

```
<- <STX>4IdTfy9BWSj4E/  
u8ynlecsVzqHxGmJJQ5fvjgyCWOBIIkUkW2eRxnHwKCFNIEcWUUMpjnlPGOz7xjSAICGoTQyNiaU  
lbOhJHPL3pTlWLDSCSJEgmkk0Emc1jAerJyzXo2k80ucsgljyMc5RjHOcFJtHPac5wlgLocZ4LhLn  
lJS5TyBWuco0irmODmxRzi9vcoYS73OM  
+pTzglY94zBOe8oznvOAlryjjiNW94yzve84GPfOlzX<ETB>8D<CR><LF>
```

-> <ACK>

```
<- <STX>5/jKN77zg58E/  
oEa1KQWtalDXeoRQ30a0JBGnKYJTYmlGc1pQUta0Z02tKUd7elARzrRmS7E0ZVudCeeHvSkF73  
pQ1/60Z8BDGQQCQxmCEMznBGMJJERjGaMSQxlnEkE2I8KaSSxgQmMol0JjOFqWQwjenMIJO  
ZzGI2c5jLPOazglUsYjFLWMoylpPFClayitWsYS3rWM8GNrKJzWxhK9vZjs72MkudrOHveSwj/  
0clJeDH<ETB>79<CR><LF>
```

-> <ACK>

```
<- <STX>6Er90wvVqUikB6P7MCl9r9qVVYWj9tH3I2flnVpSKRy1hkN/d284cif57+/8!+iz6ITn/  
kF<CR><ETX>29<CR><LF>
```

-> <ACK>

```
<- <STX>7M|7|HISTOGRAM|DIFF|LYMALONGABS|FLOATLE-stream/  
deflate:base64^Y2AAgXpnMMXgBKldQCwA|FLOATLE-stream/  
deflate:base64^7dT7a45hGAfwB8OchzmfmhmEos3fmbBiGYRjmfBrmfHqTJEmSJC0tS<ETB>C  
D<CR><LF>
```

-> <ACK>

```
<- <STX>0ZIkSUvSkiRjKtaStCQtSZIkST7vPPIf+Gnf+vR0Pz/ddV/XNwhiOZxV8wkmhN  
+aRllGp6jGX/Okdg/SimniiAzCBJlIoUIORRQRJTjnOY8pVzLJmWU85jnVFJFNZ/  
5nhm7VBDDEEU9TEkikPZ1Jlpm  
+DCSFVNIYR4SjTCWbHHLJl58CClnHJooZg8HiHKEoxzjOCc4ySIoc4aznOM8FyjHqVc4jJXuMo1  
rnODm9ziNco4<ETB>19<CR><LF>
```

-> <ACK>

```
<- <STX>1y73uE85D3jllx7zhKc84zkveEkFlbziNW+o4i3veE81H/jlJz7zha984zs/  
+MkvAjNQh7rUI476NKAh8TSiMU1oSjOa04IEWtKK1iTShra0oz0d6EgnOtOfmQjie70ocfJ9Kl3fehL  
P/ozglEMYjBDSGEowxhOKiMYySjSGM0YxjKODMaTQYTMrD  
+7kcVEJjGZKUxlGtOZQTYzmcVscpjDXOaRy3wWsjA8FrGYJeSziGUSp<ETB>C6<CR><LF>
```

-> <ACK>

```
<- <STX>24AVrGQVhaxmDWtZx3o2sJfNBgYLYwliG9vZQTE72cVU9rCXfeznAAc5IPWvF/5X/  
vZMrF9Swr6lD/e0JNyD2LwmhO  
+dHr5LevjOPcMZi81qRbhbX8lOiHVHvdg3pWEnRcOuiuVXedhtNXeoTW1ql78B<CR><ETX>3B<C  
R><LF>
```

-> <ACK>

```
<- <STX>3M|8|HISTOGRAM|DIFF|LYMALONGRES|FLOATLE-stream/  
deflate:base64^Y2AAgXpnMMXADqldQCwA|FLOATLE<ETB>BE<CR><LF>
```

-> <ACK>

```
<- <STX>4-stream/deflate:base64^7dTpi05hGAfgg7Fvr31n7Dtj383L2Nexj33sO2Mf  
+5skSZqkaZI0SZIkSZKkSZIkSZIkSZIkSZIkXV4iv9Ayl1XT/  
f5dD787I8UxbMv6+cTVQzvz0IGUSrsqcxfezL  
+RjEiIBiNi6IE6WSQJtC8kxhxiAKKKOYcl7hGCXd5yFNKecMHvhANj6I0KIGNBHvPSFpSaUMHupB  
Bb/ozmCQjGMMES<ETB>AC<CR><LF>
```

-> <ACK>

```
<-  
<STX>5plODvPIZSkRWUseW8hnNyn2c4CDHOlwRzhKAcc4TiFFnOAKpyjmNGc4yznOc4GLXOIyV7j  
KNa5zg5uUclvb3OEu97jPAx7yiMc84SnPeM4LSnnJK17zhre84z0f  
+MgnPvOfR3zjO5EMIKES5UijPBWYfuWIEpWpQIWqUZ0a1CRBLWpTh7rUoz4NaEgjGoEjSjOS11  
pyWtaE0b2tKO9nSgl53oTBe60o3uZNCdnvSiN33oS<ETB>AD<CR><LF>
```

-> <ACK>

```
<- <STX>6z/  
6M4CBDGllwQxhKJkmGMZwsRjCSUYxmDGMZx3gmMJFJTCabKUxlGtOzWUxmkcNs5jCXecxnA  
QvJZRGLWcJSlrGcFaxkFatZw1rWsZ4N5LGRTWxmC1vZxnby2cFOdrGbPezN  
+t0Lf3Pi7on7JhF6pSD0RVq4+Zxwo4Xhpu6FnMf5TQsZrRjyGuf4c8j+i3Bxp8Mt54d  
+iPsk7pzS0FGpP7osN3Rclviz+9LDHv9r3JX/5//8+<ETB>1C<CR><LF>
```

-> <ACK>

```
<- <STX>7/MD<CR><ETX>07<CR><LF>
```

```
-> <ACK>
<- <STX>0M|9|HISTOGRAM|DIFF|MONALONGABS|FLOATLE-stream/
deflate:base64^Y2AAgXpnMMWwwwFIOIBYAA==|FLOATLE-stream/
deflate:base64^7dT5a89hAAfwj3vuuW
+GuYcv5j72nQ1zz30Nc8899803SZK0JC1JS9KStCQtSUuSJEmSJC1JkiRJKuT15aMtfwHIXa
+enqee44endxDCEyDj5xCUpgcViQZBL<ETB>C2<CR><LF>
-> <ACK>
<- <STX>1FyPpf2aR+NrFFFGOYE9iSQRIUo2ueQT4wgFFFJEMSPv+4r4y4PeUo5b/jAl/
h7RgdBdRkoRyJNaUlbkkimOyIESGUiI4iSSRaTyGYmc8khI
+XksY58trKTvcQ4yCEOC4SjHOM4BZzgJKco5DRnOEsR5zjPBYq5yCUuU8IVrnKNUq5zg5uUcYvb
3OEU97jPax7yiMc84SnPeM4LynnJK17zhre84zOf+MgnPvOFr3zjO4E/U<ETB>A1<CR><LF>
-> <ACK>
<- <STX>2IWqVKM6NahJLRkoTR3qUo/
6NKAhITSiMU1oSjOa04KWtK1bWhLO9rTgSQ60onOJNOFrnSjOz3oSS9S6E0f+hKhH/
0ZQCoDgCrgHjCUYQxnBCMZRRpR0hINBpmMYSzjyG18E5jJCYzhalkM43pzGams5jNHOYyj/
ksllEFLGlxuSxhKctYzgpWsoo8VrOGtaxjPRvYSD6b2MwWtrKN7exgJ7vYzR72so/9GRW98FckWtE/
I<ETB>FF<CR><LF>
-> <ACK>
<- <STX>3Tvodw/97qJl2D2Ve6eokoKwf0pDJeG8INwXCbss3mnx88rCu2KVxj4lrggFO/IP/3P//zb
+QE=<CR><ETX>D4<CR><LF>
-> <ACK>
<- <STX>4M|10|HISTOGRAM|DIFF|MONALONGRES|FLOATLE-stream/
deflate:base64^Y2AAgXpnMMUwwwFIOIBYAA==|FLOATLE-stream/
deflate:base64^7dTpi01hHMDxY98ZjH0b+9jHvnPHWAaDi7FvY98Z+85N0iRj<ETB>20<CR><LF>
-> <ACK>
<- <STX>5kiZJkyRjmiRNkiZKjRjmiRNkiRjkiTJ53Av88ILf8D86tPznKdzTufF6RsE4RzK+LUE
+enB34kEQSx+Hhv9+zoSnIFAMWUEnkkihTQiRMkhlxh5nEz//
f4CLINIEcU8oIRSynjLR76G3zMMCKpSk7okkUxzWpNCJ1LpSRoDGMllowlkyyiZDOXheSwgVslJft
7GY/MQ5zhKPkcyznOAKpzjNGfl5yznOU8AFLnKJy1zh<ETB>C7<CR><LF>
-> <ACK>
<- <STX>6Ktco5Do3uEkRt7jNHYq5yz3u84CHPOIxJTzhKc8o5TkveEkZr3jNG97yjd84COF
+MwXvvKN7/
wg8A9UojJVqEo1qLOdmtSiNnWoSz3q04AkGtKlxiTThKY0ozktaEkrWtOGtrQjhfZ0oCOd6EwXupJK
N7rTg570ojd9SKMv/
ejPAAYyiMEMYSjDGM4IRjkK0URIZwwZjGUc45IAJhOzXGSymMJUphFIOjOYSTazmM0c5jKp<ET
B>07<CR><LF>
-> <ACK>
```

```
<- <STX>7+SxglYtYzBjYwMoylrOclaxiNwTyzrWs4GNbGzUwXhK9vYzg52sovd7GEv+9jPAQ5m/
O1CxfzfHJ0MGx12MSXewrR4GxONTLQxJ76PlmtmtFxDY/Fe5sU7mngmUu6ecA37m+hwYh
+2ueAFYnHhNwZx5fe/Ol8xFfNnfgl=<CR><ETX>A6<CR><LF>
-> <ACK>
<- <STX>0M|11|HISTOGRAM|DIFF|NEUALONGABS|FLOATLE-stream/
deflate:base64^Y2AAgXpnMM<ETB>FD<CR><LF>
-> <ACK>
<- <STX>1XwyAlI0IBYAA==|FLOATLE-stream/
deflate:base64^7dVraFZ1HMDxU2mprVx5Sc1q5qp5bdXUectnm
+aylVPnJZ22atXyOnWa5sqHkBgHMiJERGLhMilMSJGhIwxZjIJCJDJEZISIQMERER6XPc/
+Dzple9rD980oc5z5tz+31PFMvrf8m9TfR7cXR/
paloHY6nFwz8TsXHakKdPqKikMOMh3xSIFNFLWkaaOQITZy<ETB>F3<CR><LF>
-> <ACK>
<- <STX>2khTba6aKbXvq4Sj
+3iZtIlaQRTYjGcN4csglj6nkU0A80ixkFLKKKeCNVRSRTU1bKaWOvZQT5qvOMDXNPANBzIEI9/
yHYc5wlGO8T1NHOCeP3CSUzTzly208hM/08Yv/
Mpp2umgzN0cZZznKebHi5wkV4ucZnf6CseeK5X+IOr/MlFXKO6f9zJre4zR3uxu+Dd
+ABHuQhBjGyH3mEIQxIGl+SxWM8znCyeYInGcFIRjG<ETB>34<CR><LF>
-> <ACK>
<-
<STX>3apxjDWMbxNON5hmd5jhwm8DwTyeUFXuQl8pjEZKYwIWIM52XyeYVXeY0CZjCTWRQym
znMZR7zeZ0FpCimBIWsog3WEwpb7KEtyjbd5hKeUsYzkrqGAlq1jNGt5iLeuoZD0bel8q3ucDPqSa
j/iYT6jhUzayic1sYSvbqGU7O9hJHbvYzWfsYS+fs496vuDLkvtd+M
+s1EDb4q5lh54ILUuHfrWfVvWFPmWF/uSH3pSGttS<ETB>2D<CR><LF>
-> <ACK>
<- <STX>4FRiRNOBZmP5n11ox57wgzftbMeGeY7WR2L2XM9sUgmfuesN
+d8V9XaEhL6EzSmubQoaZwTodCyw6E3iX9i1sYdzHuZHkQX1fc0sKwLQjXHXc5s9H94f60hca3hP
vWkCEdVIX7m/ndSKTCNyX6B/EzSqcGnlki/k5+n/92/U3<CR><ETX>4F<CR><LF>
-> <ACK>
<- <STX>5M|12|HISTOGRAM|DIFF|NEUALONGRES|FLOATLE-stream/
deflate:base64^Y2AAgXpnM<ETB>CA<CR><LF>
-> <ACK>
<- <STX>6MVQBkIdQCwA|FLOATLE-stream/deflate:base64^7dSJi05RGAfga98GY98GY98Z
+24+xr4Oxr4Nxr4Nxr5NkiRpkirJkyRjKqRJOiRjkiRjkr4kSZikSZL40yW/
APirafTPfd265zrn8UJWpXxrchWhfGbxWlovzwnJ/+/
TmWmKQOYUjEg6MomVTSiJFJNrnks48CjLIGS5QRDG3uMcyj4rzklR
+JhKRRWSqSRDK1qU<ETB>21<CR><LF>
```

```
-> <ACK>
<-
<STX>78KqbSkLR1Jowd9GECMoYxkLJlKMZ3ZZJPDUIaSSx6b2U4+u9nDXvaxnwMcpIBDHOYIRz
nGcU5QyEIocZoznOUc57nARS5xmSKucJVrFHODG9zkFre5w13ucZ8HPOQRj3nCU
+I84zkveMkrXvOGt7zjPR/4yCc
+84VID5SiNGUoSznKU4GKVKlyVUiiKtWoTjI1qEktaIOHutSjPg1oSCNSaEwTmPJkM5rTgpa0jVtaE
s72t<ETB>CB<CR><LF>
```

```
-> <ACK>
<- <STX>00BjnSiM11Ioyvd6E4PetKL3vShL/
3ozwAGMh0YgxmCBkMZrjDGcFIRjGaMYxIHOOZQCYTmcRkspjCVKYxnRnMZBazmcNc5pHNf
BawkBwWsZglGUZy1nBSlaxmjXksbj+91eTx4b2MgmNrOFrWxjOzvYmfEjF/
62SuTV75IVklfJlaNyQx4I8iYI5ENeuMcXw516E/o/KfRm+9ALWeGs5oW9L9nfVeEMFobzHBO
+7x16My<ETB>14<CR><LF>
```

```
-> <ACK>
<- <STX>1X8K/
HP1+HOXgt5cCjkTU7lqLSQdfGQkUUHmWtCtibyNh7WVRjWmJBYc8n6/+T391Hsh1/mf9rL//
Wv1lc=<CR><ETX>93<CR><LF>
```

```
-> <ACK>
<- <STX>2M|13|HISTOGRAM|DIFF|BASOALONGRES|FLOATLE-stream/
deflate:base64^Y2AAgXpnmMXg4QQkHIAIiCc4MjAoOCPkGuxBcgA=FLOATLE-stream/
deflate:base64^7dR7afZIHafwY800nbXK+63l<ETB>5A<CR><LF>
```

```
-> <ACK>
<-
<STX>3LUvNWd6yqcfNzEptlZm3cnnJeZvTps37S0SoRAyRGCL2lhIWpmlhEhEjQkQkRCRGiAyJG
CihESli4ue45+DL+4f/2V898OF9zu/8c97n93u
+UZSsbeV3fqKpZdHdFUdRjtQzE1uf46RGlkaaiSZFURHFIBBTQSU1ZNHOPQ1kOchRTtDIKc7SRD
MtXOU6kW8qoD2FFNGZ7vSmmAEMZigljGQspcRMLmv9b9OoYCazmU8li6li<ETB>FB<CR><LF>
```

```
-> <ACK>
<- <STX>4JTXUUsdmMnzCp3zGdnbyOV9Qzy528yUN7GEv+8iynwN8zUG+5RCHOcoxfuA4J/iRn/
iZRn7hV05yitOco4TfOco7z/E4Tf3CBizRziT/5ixYuc4W/uco//
Ms1rnODm9xK5sEMtOEBHqSAtxEO9rZMB3oSCGdelRHKelxHucJOtOfrnSjOz3oSS9604e
+PEkxT9GP/gxglIN4msE8w7MMYSjDel7hIDCC53mBkYxiNGMYy4uM<ETB>43<CR><LF>
```

```
-> <ACK>
<- <STX>54yVKGc8EJhlzTLKmczLTOEvpIqr/E605jODN6ggjd5i7eZyTvM4I1mM4e5zGM+7/E
+C6jkAxayiMUs4UOWUsUylrOClayimtXUsla1fEQt61jPxsXgY1sYjNb2Fp
+Nxfu1OpzJc2UNFdSab7EITMKwn3eHwA7V+iH88jlqkw9X3n
+hkN8x2G17znOEbLsUt9R3ppv1eFcZ4UeDwuz1jfmYjKnyTwn857OfnInroR71hTu<ETB>F9<CR><LF>
```

```
-> <ACK>
<- <STX>638lwfw+EPKgNeVMasinJrZaQd9mQh2kmNoR9TcjOOE
+arSVhn2ZuKsqpFefUmsM5ZvNkgjvhH+Wlc2pJn9Je5Uvf3Utuv/Nr+fx8d/+v/
34IPUh63xfBg==<CR><ETX>0E<CR><LF>
```


```
-> <ACK>
<- <EOT>
-> <ENQ>
<- <ACK>
-> <STX>1H|^&|||YP8K|||||P|LIS2-A2|20230929093529<CR><ETX>C2<CR><LF>
```

```
<- <ACK>
-> <STX>2P|1|<CR><ETX>BB<CR><LF>
<- <ACK>
-> <STX>3O|1|2023092700000005^1^041176^1||^|^|N|||||Y<CR><ETX>6D<CR><LF>
```

```
<- <ACK>
-> <STX>4L|1|N<CR><ETX>07<CR><LF>
<- <ACK>
-> <EOT>
```

### 3.6.6. Example of Data Frame for Matrices

#### DIFF



4M|15|MATRIX|DIFF|LMNERESABS|FLOATLE-stream/  
deflate:base64^nZShbsMwElat4k4aDwnYA0wKrmY7pLCwYA9QULA3WEDAwMBQNFb  
NgXuAgsHCPMPQYOFG4e6Wu+bsnLe2lj7FsX2//  
ztbNgbB2mOX2uNuYLvrYMB34P92ve8AZOyB/  
uNN2rDmBzmlsCzH88xUjMD5kBF4FhN8VP6/4T  
+PenmtBBnwc3eJsKnTlpx43Xvse1/tBl/OXMVm0Xupwbba+9zwR/  
mWfa4Q53RB1FM5sHjYp+CxS/uUZaX7+imNfvFarJ+ewje7BHnx  
+A2sY69yYtXJncLzJx5yDjWxaZqdC/U4fgYU4nyk1p68oPeZ2LOD/  
xZYKT5Te0sPrMt5NqlOqCn1tdpsgJ0f4nbCY1EONZL5MBxbiPITvcc  
+ND8r8t4lvGfkmVXPhrQ25PfFhbQ2TeVCFgmswsMJ4Nun8eTGVDdBnhwvJ/  
0FtLuQ6QXWn82VD5Jv7+17aNK27nNczeT+bHw==|FLOATLE-stream/  
deflate:base64^7Z0PaF13lt+VxJ5oEifRJE6imXVmPFnNjLtNu  
+rWbd3Wba7O827V1my9WxW8i2DV1ksFq1K1FdS0LvsQxggjDDCGCPEwvhjDDCGC  
OMMQ8hjDHGGGOMCSGIYIIJphgBhNM6Dujz9fnvOsn20lm//YJftx7f/f3O3++5/zO7/  
c7996njg7/+NKR9PR/6pFx9Sf/W5Hx4B1dOxulA8aZYTz0UbZ3iibbfX  
+UKN0NURRKHsoOxpIG/fHuOf9IK3pVE2QW8P1zsbZWQjDDdKL/  
Uf0L6g7ECO7RxH6FdwdJ67GmUrtHuSLPvpN5pk205f6TDM353IsDvRGaHvENfebzDx3  
ZHa7KBfAd8PUt1u+vVzvpX+m6jbVZJ3M/gMQ2cv9bu5J1m2gt8uZO  
+HTw99imS7fmhupG0veu2iz1bw+QD9hHGv+9tp0wOtTbTRvT3lv5I  
+m7nuhldBu2FKP7R3QLMXDLaRSfyGob2Vspkim1fh05/  
sPEiBUDqMUL852WRLajuCnH5/DFyk6w7qqxy74bfbwg+E7Qh6aoxk3/  
S6deC7nXvCuQv9R5M9tlv4oPTxuk4w6Upy7018+pFvo4Vv9XLcDC3pP1ay/  
47Upxt5sn8OcD5G22Habbfw903oPsh1L+3dhvvQdy+6KSYUyKvxJt3J2wHLMZoP/  
e66bPFIL5003+MdhuhJztsQf8uZFaR3SX7JnTdAm6KwYotu2k7AB/  
FsE3J7sKux2lcFxbjdK9FTO1KeuwFK9FX3NwOD8XMvdhoX+LfAy2NyUHqpeNO  
+nRTNwoPxd9u2ji9/Zz3QG8HvNcha6bTzvh2GMRK4doo/7D8B2xGGvjyX6uT6eFX/  
ckzBWvN1jMTQUYyQc/sJh7RpFp2ML/19FnH/lo7snOW7k/  
Qt8q96vos4e2HyB3P7zkM7LtfvqPJ5kUr/  
dAv9MilirOy0793O9G332UPdDZYDFjFNGLOLtlLw6E867ub810dAYUZzZRN0W  
+krfLou4kOdWYdJp4a+bLMA51yVDCcON1I9gs33ItA8a+zIK9t30HaL/dmjsBa9NtO  
+BhsbnMDzkj7uhvZfSbxGfhLPi  
+j7kFk77OISDZ4dFHJQP7QW3LtpV4aE5ZGvqq7jaazEneLsbjXK8UU42ykKjXGqUeqPM  
NcpN6s5TN9Mo9xrlidKM8apSPuefXS41yjusvqPNysVE  
+bZRZaF9oIKO0XaHNce45zxo8TjXKR9Bd4v5Z+jjfK9A  
+Rf9bIDvUH6PdSfqc495VdHF5bnN+HvrXaXeN+iVozcLH252gpbH6WeR6AK/  
71C2ntmehd5Lr82B5E6wW0ftL6s9yXUOWs+C2RJ  
+L2OUWtM6i1+VGmUevea7ryH4WW1/CHjegfQNCpOw  
+NQ4LoKJ63O3Uc7ASzhehN817DCLHler4QK2uMFxAdstcD1H34+Q7SJ9TiP3Avwvgt8  
5aH2KLKeQ4ww6nuD6DnXz6KC+R6Evf17geh75l+h3Ouk5b+ELM2BRp81NrufR4xZ1p

+B5Ab28yHeXOD8Bf/nbLdp/  
lnC6ijxz6PcZdr0CT40r2egqfXFyWS3FXTU2PKYOY8hk2P6ObyWoHWBtvKvC8nWly38/  
AK0TiP7Avdvl8ssss/A4za2qGF7x  
+wT9LiPTMv0u4I0NfqfQ687YPYxdQ9tdQzXLMbHCFs5mu6dBLsFeAvD6+h6CR0vlpNig  
Ou8CL+5JONI+M/SdhneDyz8QvFDsq3PrKlvdegfxkM/P6X1J1FPtXLRhorc  
+A0h50fWMSBS9ybQlfpRjlgq/Op8zrDtc85VQvfnE42SgHaSneh211/  
pnEjnPoN2kRb/aD/yT3p2zVr6bpMwXNaeQaTzY7Rv8jtjreDyJbDaz3l/cM/S8g6yEL/  
5dvT3F+kvT2Oog8pwGJ  
+e7G55TIFkLn55Gxil4CM8LntPOcxxjrKA7FPUz9JnGr4u21iS6xy05xOtw8g1i5wztNHcc  
RY+E7SZoP2ExfpDvn8o0ZynznlfRddpZDrK/RlozaPPBLQuwH8ZeuP0OWkR  
+07Dx2nsSTgcp+04NE6j9yJtDyY5XS/tPUSjhpwnoLOQ5F3iel/S0WnJpw6gyTXfu+MhQ  
+K/  
kTC7zQ6HUY2jdm94H0E2ceQQ2NICprH6ScbHoDWUeTwdkMcZddD3DuT2vl4HYRG1c  
KXziDzLHpehOf+hNUsNosWflfDpifQxe9fhd9BZNSYOYDOB  
+BxmL7HkVvjtsb5rIVPyq8Og90cWExj772J5jHq+kn/5sCl8MWY0e+KP87Rd2JZNdJ  
+nvdWeoPlpvG/kGL9WeVdppzq8jgeo5ajBe/dwu8FCvY/PYYAQue5Ftr8V6eyrpgcz5/  
ah9oNcrFsxQP4U+sxYx4Qzyjyc7H4OP1ucTCeujHEmybCPNIMWay+tlyahl36TCV/F  
+wMJ10vQnqY/UeSZtxg3m7ZYr07bjG2ZGvFxd3oUoOHn2v9pTE0ZbH/  
GMd2Y7SvYYtxi/3YLLlvQGOCviexmcb0CLiNwv1YY7aP6z2P9rbyw/  
GEnjVXFg0nanEL630CngxZxfAkai7Q7Dr1ZdNyxMD1gMRYVR+QXI  
+Cp2HUEelXuU4yQbHuk9Y9Pcr1EXgtQuvERdwTvqeQbdFifqyizx5o78OOlyz8RjHn0k  
6+dcpiDpCc+6DTbzFXaO+3CL1LFvHqHO2Oco8sNBryhrXR/  
GDt4O633ZfO94HMcmZS7yvFsGkwx0H4zYK95m/5rdrdgl/TVe5uAlryO81LByzigNNQ/B  
+G3ynqteZZ4P4xdPPjLtoloyUwEb6nLPzgmMUcJ9nHLPbzC6nPjIV/  
Kx6fB9/51Ebj0ftd7d4HYbGPHhWaX8afsMwa43j6HnQYI00Y  
+EP2qtMwEvri7Mwa1OtXfq5nsfW01zvtZgPNBeeTO1mLNaUiscTFnOK/  
FPr1YNgst9iLphAhhmL+v/rAs03y4mX1roT4DFpsYZxvsqdyRenLOZ6jX3F8pFk1yGLTY/  
kmUp2PcD5VXgNWfjKP0m01H7X9fBx/RCaqs1tvue1kZaD+6l7ZzF3It  
+eJB2V5BBc4pyCqJ7HHucp+4s56fq/wT20TprGPryK60jJ9NxiDajFusL7X29fszChzS  
+Xb5dFuNo6+Fpi/  
Wy1kdat2hu07ztctYt5sY96LPHYh08YhF3teeUtBwf0v7tJG0VK2TfMcpNi/  
3LsMU4kW9oT3bQYnz0WuwnTyYbyPeOWqw3DlvkpmvIKPtpXJ2yWOfPWayzRGPWYv  
2n2DBtsU+sWaxrxzRrSvI7fu07tEcO2cRJ7Vnnk0IIPic/Jhiz1e1WLu1xx2gbYu92WLfa/  
G2HykfY6e7DIGYt8zhz09mCPPC60R1BM2Z/4KE+jdY5yW1r/7rNYU  
+R5UXtjSxtmwbA39teQU4fp6MW+4DTyFPFBsq/iYePe+VQFLvPQUP  
+oTin9dVBizhbo53qHQM9i3C617i312Jca/+e18maz+R/  
2nflhiforxzLCDhqbtlN4Rear3ZbzNHZ6cs4r/Gkvacfl/5F8XkvB49ZzE/  
7oDuPo7Kp2geumSx3hxMGOX91gmLOXiY/sNcKw67XpqbNR9rDjnf  
+QQyL1uskXtSjvMVY3WuRk9qPbbVuOmexTzoDTfeDi

+jhONy2WINJlXmw0DyiWHuc8zz3n7f17ZwEf63TjyZ7LijE2Gnk1Z5NuU3tU7Xfn6L/  
eYu5q2bh+4pdyuNdRM9xC//RPKr1uGhq3arnMq5DHF20/9RcMID+qlcMVA7J  
+Y9ZxPF5i32YnoPir5THu2wR66Yt8n9a44yjj9tA+zbtvw+lv6zRfSYsVc1JoT/  
EJS8vfbSWjMpijheK9+/xWLuOJ7ua0+gWODH5eoVY/  
ZTP8x1Ddm0p1WuTfmigxZrnX0W7wNMW6wDhhNWk+g/aDFvas2xIHjqcQuC5/  
TPvAQOipGKTc3YJF30TpP8+2MxXNN7Vdr2POExbMb5Wzz/  
mqPRb51Ev2I73GutZ5UHN0eXnF1yZrHqR+1rh22yEUdgtZli/3XEH2dtvb4S  
+B30Wlu0Vw0bpEX0zpHOWfv62NL68Uj0L9okfsesNij1CzmTuW7LisQxTXNB/  
pOeVRA573tQZR3irvA6vg6f3qCUPtdYwxtxm0eO6I52vD8NXaQPvVHD/  
HoaF4oxgzgk2dp3IPW8La62NxHelPuPj1sozXeTedvifho9ywcq3nLDIHxY5Lbka3pul8f  
beeyl5zFan50GqzPc32+Rc9oBNnuobWaxtKIRXxWTmqXxTyqnPSoxbhWTnXIYv9wAPK  
u08dpDVise5Xfu2QRh7QeGAXXwaSv5jHFQfn1Xov1ueyofYHT3IYc2mMcxw6y3wGL5yb  
aWyn  
+ulwaV9rjam951CKWH7bYw40gv2L8ELg5jxsWe17t3YWd8uraP6iv1qiz0Na8rD1N3WJPC  
w4cLkH3HOc16rV/0tpZexGt/ZRv07pL62+3w26LZw+z8FS+1OXVs9ntFu  
+ZjjEGM0PiuFjFu8DKMerPfxR9Fcs1pVY0N7kcP01/pla7a6hQ9qjSK/0LM/5dUXLFZ0hxlV  
+dBZi2eGe9BVeyTNZzMW+d0xizyH9vNL8FKO9yZFOxOKZ4MctU7QXIXrceULF2h7zSL/  
orWV9gaKdZMW6wTNKVqXnYGG8viHkPcsMlywLmKB8pZKLZpXnSZTqKnc1Vi3noH  
O2VhzlpsR/XsxZRONGa7T8tskSntU  
+bh7Z5y2erSgGam2IMbnlcQ98DiDPkkVOXOvuG9xXDU22RW5a6zHtReSHyxa5Ha29Ziz  
e49Qa9Rr6LlrkK+R3y5R+i2e6yIMvWlzpQYtc/VWLVNGf8JHfay0qnl9a8/5K88cZMNf  
+RfbdnTAdB+udFs  
+KFXtPjtsptzplMafesvBVzU96bjFs8dxK2A5Z7B81Jnda7Mm1xjuPva5YrB8HLOaVgxZrul  
sWubyaHQ/uAXfiRs9YxNn9yKd84VmL+fkA9tSeqW4xLvUcQfjr  
+UAdPsaseU8zhuyDFmfv68jr9y9b5G5GLd7IVs60ZvGsVvu1Axbv4un5wlcSS89d5S/  
KTwxvbK9yAnzOJxkvWswNLv8S9h2zyJ  
+cRt9byKA8906L53DKNSvGan5X3JyJvZ55HgO/  
wxa5ea1n91g8h9H68hz3z6CL8pTKiW21WCfreYrWsf684SoyDIqs1bXeIY0WLD6pHQP/  
YYvcofiDyh9oDaWxp7yV3gW7zrFOW6dzyWLNvWixR8rPNg4lJlctnnMoFz1u8XzyAm3HL  
J7fKqYqYwuwQWYntqr6S9w47EU3sG7QmGOe+3el65YLH  
+0ZyoOXUZmhrLyu0P0EfrKM1rBfSPwNvtrOcCkb5Qj1b0TpxE3GNf  
+V6dlo8M9XzW3Hslp+cSPSGLPa9WteOjt4LFs/Z5Sdary/  
6bHrpsUaZ8JiPzK8e51IWv5x7JF7nLQImft9JQvrWI7zW0noXUO3PQs5Bx8jnHUuko  
+sNdiHVu12OP8Ejt+RntPkF/zud4f0zuLn9nq+27KUyt  
+Od0Vi3zkssV7BO4bHldfrsS7aF6WwEFxfBm+s/TVXK356Ab0vK/ySJ  
+is7edhuZIZNezRD2nk+0PW7wLeA3ZlyzxdMe6ILFvHDOlr91zWJ/Vad  
+Dv3r0Lpu8ZxNzyT1DsLIRFv7MffBLzk/b5FLGrFYp1205j2vn9+31ff3LkJK3kflp5YpFvkBr/  
MsWMcf5av6com4ZnG5Y+OcVi/  
2C4r9yxcswMWMMeHpcs8jXzFu9gK9fhGD6yWEt6n3sWe0k9q1mi/


12LHKI8Ues3H3evNnZr9UZ5r1E2NMrbjflmx583Sjf13+P4YqN0Ncq7jfKLymr/  
HzfK32mUtxrlp/R5mT5vQv/7jfiD+r2W6PxGo/xT  
+v6Afi81Sk8lcmU53zppsfSvNC6Vvm36dROuZapUrtJi  
+fxh0p0dJxKNPWsUe9LTKR2epdAvPNxskRf/CST8IF6Z0B8jyY5xVP5/  
kOpXjJ5P48v6R0Vf692O3vMcrjz8fUu+/7mHIAH2/zCXVfcf1CZfX+Vxbv+H8JbR  
+v92nn53ds9f1159eJ3V/Cpj+qrl4Tn198j7zJVueulb7V8/t9q+sRj4Wes/A9sub/  
ae6vczwx3qjdNhquZ/OF/pW6T3sW22nugXqOI1W6m9Db4OtHtdRV+tb/  
X7Cj3e5XmyU641yuS/anOS+071PXb0v6HYmWS8jj9PwPO1GW/3ew/  
XvsfgGayO4bLb43sqPWY2+c9vGsYfzldpvpa/Tf4gMPbTdYvHtYg  
+8xGcjdR3USY4Oi29uNqdjp8W3Xh9YfK  
+jNj1gKpkGaNdr8c3dB4meY9APja0W36J1WXw/  
KN2I1yDyraPNTtp1UzKtztchTphZF0n8dtDeCp7AZgYa330EbyTpk4aPdiZ5s25nw2Wjx3Z/  
TG6bPVu6pbzf9VO8+tQHdtyQ79lh8K9RDXT/3tyY/  
0LGfthnXLclGG0u2kWydqX4bvITXVmxQlmcNLZTLx7Lb6zHbL4trJImA5bfP+1w  
+lbyR3c7yphUVh8/7rN4vtU2bXTYhy5bPKXLouX02Xxraj8tBP53Nd20GcjdETL9d  
+Z7LzRwj83W3wnudHie3Lx7YV3R5JbPrsuYS49u+C1JeG8wcJXntGn0+i7LvXZYDF  
+hMsHFmO9k7YdFt/P9tJ3l8WY0JjZaeHvO+i32eLb0e1cu008NhYJi57EU7quEJNX  
+sJ/uhl/xbBdFnFQY0njdiv3P0DfgYRVj4Ufy667Lb453mnnxnaPHb9lCMmhckU4XFjHvYV/  
EKI1TyZz9Y0fSvZFBvqjx2WpXpZ/i567Ef7vF98yK4RqLhYUvaS4bRlCbi  
+9VZbPsnwPYSepI2w7RvN6Wseuy5hile7Kd8x62GGeK84plm9BmUVcks9tqg  
+agxfSA9y7jy8pGggfTHbsEevxfz3QkTjU3FNY2PbdY8fmSrbeme89liEvs6kO1+X3wXr  
DjsdcPwUvxweit9MX7yukXjSfDc5XihOa3TosYIL+QD2xBbs0tmy2+Sx+k38O  
+aNdNvWLLJffzY9dCzC6/TM48dFrFgJ9dDFt/  
Lb7CljwPIOMLxewvrrLL6V15y8O7Xvtogjkm5751ov9nPeD22tg2RjrcGGLL6xvo8NdnGu4u  
uzMXxkgP7Ke4xYzOEaOwMWz+a9zZiF3bda/P7IPovfIVCew9tMW/xWhRftXURjG/  
313MPp6TI01/  
jbQf1myzWgspzKdfp9eMWucQPLH6TQDy3W7yfsAkMdP9A4q14V4CN7Dac9Bqz  
+A65oK/iz2aL35vYxD2nr7lzxJrn6l6L7+xIT/ntTmt+92UQflvBvSfhoHisHJCeMSg3IXPN/RZ  
+LhtL4MWv7uh8TxAP9ljm0Wef5i6nUkOzTda52gdNpL4eLsJZB6H3qjFHOF6nbX4zQa9  
A7OHorjdb+HfPRbfi+9KOG5iPHdb/EaJ5s4ca3ZZjD3FkyGLGLAdXfJ6XbR2WLwbrPit  
+FG1yB/KJ+Tzwxa/k7EdG0jeXoucjtYqoxa/d6H5bjTJuQeZJpN9hOnRxLsXvOR/  
uymHLfYDGIPOY9aafw/HMdn+f8Qi5uv+gDWvG/  
stxtK4RTwcpE54a12vmDtjsY6ZvFbMr3oKpn8KF+YQBfZWHsy2eyQxe/  
V7LCIjWpJlop92xMmm8Fkm8VYUuweS/YcxOZ7LX4jQT6/M9m82+Kdj20WY7jfwl  
+0Lc6XfvWfRRzHZenpfw3Nt1zI3uDYqrUtn3oXu2mvPwwlbnPq2mu/w  
+ydp77mSa7RxGp77uEu/m+j1QiW+2fY826LF+4BflYPTuYwsV201T3lo8fP  
+nt9bttV8mbfxPKDn1r5G3kv093aeM9Nzy1ca5+vJ1dyEnt5T8Dyo+3nNlm/  
ofV1ml8vzNJ5/XELmW+j4yJrfU/S8nJ4jvkh+yPHwfOIZ9LlnWMzY/  
GM7SwyuyfWXYD7nJetvjthpvwroO9cpu3wPA4+DgO97Gf6/8QXK5bPOc9idw3oXcFzJV

b9DI3IxicRub7Ft/Cu5x3wG8RGe6By/uVeDZ/CR1+Wonvsc9ThOfnyO00PQ/  
3k8oq3o79zVDXOyiXsM8I2t  
+zeOfRNnj7+QZ85BGY69nCRfTy4wPaXrN4b8HnzxfJud6weE/qHXR6gA7L0Dtm8f30Tezp  
+Mzbam74ChidQgbXdw6bKM/uvjHN+cUJUP0KXL8HE  
+X9u8Sszslj8gqd83G3cCd97yOR17s/uV19jr/Xkqg/Dw+3rOWk/vsnxFDTuc/  
4SeemP8ZNZsDqLDJL1IPT83i301nsEC/jaJvzBsXCfdl/OHLi+FZiClmPdAY8vwo0B2Ljfus/  
qXVAfd2P4mxcfnzPfwitwn1Xie5Xr+NgVZFhHDHL8FhsHnxcZ/e5H9L  
+tsVvNHxmEZ82kVO+A256nuO+9Bb0/TiDLHeg  
+6XFOWY3qZuHdjXh92Y13pmsY585MH2JvPIDi/ecX63E70RcRW4v+6H/  
ERi4vp9bPFu8iA2dvz+DcGwcox8Rp+6A3w34bKjEtxcr4K+c+m/xLEPP5s5YxMRXK/  
He1dvE8jsW30i4TlqtXjeNFTZW4rdNZJ/  
14OsYvkKb9yvzezbXLX47ppuxsoic62knm1+28IHD2MPP91u8O/8CPvwQH/  
zY4ptRf5ZzOtnzhsV7EYsJ  
+2sW74YoRtyBnuPj40zvfvk9HwP67utVZB2Fnuzhh9MgnEX49V5nbb4/  
YePwfsWgJ7E3hfgrTHXCQ/F0wmLuUTztNvFsXsTG7s/TVm8pzUNRrPw/  
U3sXEXXGWW8Ab1XK0FT71j+gDnhdcbOy5X4jml9Jd6D+j6+spJ8bJF  
+88isZ4uu0wV8wmPp1xbvG7uPbKrEs1aPnQfo95A2N/  
EDx9bHxucW8UtzuPvIMjwelt8CODtO99BX78lcpa+Pt45KxAa/52NS75PKv37OM8AjFr  
+Bsgxer1TiW5wO8HVcprDZF9j3PeRYAJsOYoiexk3FN1zOj9DxdWT  
+EeunerDPW5X4tsn1faMSv7OksX4B+72NnTUvnsWefryE7RQLZpFjFsvWwAGnT7Cb2/  
M8flLoq8btNep/Sj/eR3X31IPJpTfaQ8sNKvLfgPB8h/yS4OybXse89ZL  
+P3PLdw2B7iz4vMn+fxU9uQeNz+Mxjs3vwcV7epwYWL+NHL1fiebiPI+/j430OHJ331/  
ip6+k+6naXPfXeha9vfGz6fL4Z2g8s1sVHLZ77u4/ftPjGxH2shv/  
UaLceHPVU9CcW79nN0f8uej/Etl+Cm9ZhWkvW0dPLj8Bkl/  
lqNv4CG7ncM8jsfttZie8GTid5Z5DTsfEx+SW0nEc3fn0CzO6C00+InS9V4t2PzKp8q+19fl  
+icXbQ4rdqXqvEb1QtW6zV7yOP0xzH/o7v9yrN71mcg8d9sF2P/7j81yzmhh5wv4HN/Oh  
+6rFrFtw11t/F1p9afM+t9ZXLuoF5w9eaDy3mzRns+shi3XYGWy5arLMuoYven1Hs  
+Qqbvo2fac2hGHbT4I2OZfp/jU1d1nvYRb75iPqP0fcTi3W/xpP73M+JAQ/ov74S75XOI  
+M1aPg4uoMtP0KXewk7b/Mz4sYr2OKChQ  
+7z3jMcdRHKw8MnLTX1dzaxbj241fUv4sPvY89P8Uuejfnl2B9CWzeT3a/avH9i9bezu8Y  
+H8JVu9V4vcD3bd/TKzWuvvoONI22eOdsXSV+E8yxXwFP5Qy+tvhdg33Q2oc  
+d8BwymJ/VQcPvUt80+ldcq/z8e5xzOPoOvDUXugjsHUsXqzEb7HTR+ZXif  
+vMJam4HUSbDyGOC3Noa7ry8zH6yqxdniLmDNhsRZcQLdxcP0cnH1f6+1vg4v331JZX  
ZMpnrvzFm8t6Q1+yOL7x61PniBGHbPYk/5EjHJ8esGE2H3RiVigY/  
dA2DsuH5Gfr0z9kol9hZV5HIAHJ/  
Q57TFd8+up9tC3wK8B0Z3Lb4dcUwWwd9t4lh9ZbGuVf7oLnyWoXOQ9ifR0/v8/Up8J+A  
++yk+9kWymd57Ogdej7CFYvmGSvwm4E2LmaP538f9GeQ6ii+5LD4WPGYrByS/  
kT4z8HkHH3JsvqadcgNd4HMM3bS39HGsfNtDMFbc/9Jit/Rj/O/  
VSvXepfa3N5DrpsW3Zfvxlw3EttdZZ2hdeweel1ci9+o6+Vh+jXXDOYv82nrGj/

vUJeymtWo38UK2vArdH8DL5VbexjGZs/htJfLx  
+I5cPEXrPzQXfBxnTorkQNwfFyPdyvxWx3O6wFyicc7jLVXiTefWOQ/  
9I3TmMVVVGp9XwfvL2h3F7w09tyGyg9oP3Ha4I3P71XiHa5b  
+Os65laNHb0jfnJivVOD95uV2L/P0P6yxfvK3v6I5+K+1c411w2ih32Wext3cd8nt9v0f8+/  
ZyX4/mjSvwGmvJ0rv9Fi3yV6+Fx4oeMeeUfvd9n2OC6xTx1GPw/t/idjC  
+4f9Vi3XUPX1FM83HTzXroe4zfVYqRR/gavI9YrM3uW3xT4/S2ME  
+9X4n10TrscxW9nX8X8foVfM3fyZSsXfjQrEW+4H3GjvKpB7n3Gn7yHmuBDuTV/  
ld7plu0rVPvOPUwXucs3s31Mac18OfY8HXmghr6T1vkdI9ifjwNjnon1/sdBsvN8TXnD/  
AtH1Pue7+0yMt9RP20xfd60/il1g0e166j33rwc/0VW/  
O4c5y9OHYbsafbbA5d3c5ukxvYvncG8jq40Xvf3rbZYv3a99jrnWsFYtucq31ztdc/  
xybuG2OINc88r+Z7PY5tBz/eXzvHPZ2X9HYWwBf77+MbyqvAladvO9DI/  
XFD4mR7n81fPb1SsSF9/  
DP97Cnx7Gb9FWOcjM6dxNr3mbevE2bMXB7GX27iJPnwVHj0/  
XzWPBiJWKL1g6z6PY2PncNbO+Cm8cXxVC1PW3x3OfTsuw9JzjAXjege9L  
+Nj30fMLfNn99ozFb2524tdVi1yB9r/H8M/1+NULHN/A565i9xn0fhvf1d70ETb1/se4fpO4/  
WOOGzg6b+Wmj+DjP2PsuCzvMCYmaNdDDHuFOOP93Oc1jtw/  
Xq1Evu8OerzLGL9N302VyJO9Ull52PU8AO6fg9NV+CpvMIX+ryHLJ/  
ig20sx0PEYtRjHL1VinTiPXhct/NLHuGKN20t5q1vlcAQ/dRxp4Hfe5iNs7f77BuOgoxL  
+eQV+2kees/i9ko1g28VYPYVpvcF4/D40HVO/voGdtFc7j95LFr+V8yJY6lmWY/ga4/8l/  
qpc3l2L70iuQUtzgeYmH49fwmjZL1n8TvuGpOfwe8ydbcsnlt+BIYX0Unz2rLfb4FesHg  
+or3yDYv16hKyf2qx1rpGu8/Bx8fBhMW+w3FwW21CjvWV+L1t5au1x  
+tgrH6Bxb4En0787KhFLu6uNef83Q576XuNPI3Y8JDF3v2hRS7qDPgsWKwhNMcpX6lcg  
t/TnlK2+dRijOp5qfYnDxKuE2DxTiX2Rb+A/yOLb31ercTvuvlYcX2+tvht+E  
+4nqHPD4hV3vaszZredaha  
+Kv71Pcq8Vs26yqRF1SOWrk953XQ4nsY9zPP4SIHvnpY6zKMYsMTFR/5s2jxG/  
eO65v47jlyud+9h7weT25YfgfItDxWav10Hn9z/r3Erw786B3o+nkXPncV+/6sEr9v/hZ+4z5/  
HLs6n0n84gl+63a+Q7sli2f9wtFjInbLmpfvgnUXf/AyBeX6KO14Rz8JtDxAdc/JA6/  
QUzrrERMeJHYdBceb4HXQ/zN44THBo/fPI6cl88t72Ffx/FVcPsN/HUOH3Ef8xj/  
00rkDZWjf4n48BvE7evo/Tpxyfv3uZUjv4GMv8QWmltfJk4dBuV50TX84Cf4+G9W4vvq/  
djGaXTT/ivkVr6toxK/  
P30Cn3WsfqQGLoOPIxZy27LEeRUvu49sHHcerCB1vhe5xifByOX9TV8yv3Zx+p1i  
+dvr1fi/264LI6RY/h3K/FNzhFoLfnMPZoHfhPMNoKZnsHou6Cvoefrr1P46yFkeWSRz/  
F4sI9+P6jEvvK8xVz4U3zRMf0Eu02A  
+zv4oP0rc71skevVMxqndxPMf4vYKN8XVsr3HYWn89Oz9Rfwlx+nB8BYc4j7uM/  
VbnfHdAZ/m4HuXuoOWXyTfpix4HK7v/oapwZOWnM7Di9UYj/  
zFmPgDn6mdzlm8c8T6Pt9xo7mz3GwexddHWPIBK5Z/  
Jb8L9K4eQn9v4dtx5HJ45Hb73Xk1Lh1enr/4X34+P/qUxnl+HuN8qeNsrVR/k+j/  
ltGqTTK/Opt/6yymcdqazK9I8b5T9Xvt+/cv/cRTu3475G2cn17zfKHvc5rv9Bo+xolN9uID  
+hzvv/R2j+I9r8c2Tx/e4/qaz69m/T/j9Bc6BR/h3Xv8c9l+13Kqtx3u/3u60b5X80yj/j6Hr



RET



```
3M|14|MATRIX|RET|RETRESFLUO|FLOATLE-stream/deflate:base64^hY/  
NDYAgDIU7gmN49MjBxAKDwSiO4iiMok15pFQTX/KF/  
rwCJRKVTONkJoTDV2shKe2ZwVx6Ugsni0p0sccZkPvYSYYr6395XLWDv6NHDt4j80h7  
CII3cYQaulNWcF96DUHfMC/CzWeoThTjzEM/zB5aiO07FEyw0=|FLOATLE-stream/  
deflate:base64^7X0PiGXVmefttjppY6svTq86iYkvM51JNjGTN4kmTibZPt  
+9ZuLEqDUZtZsnLVmJvj7LjwBhzWZYV9SAiFiDyaRgoRuUgjhUhTNCJF08hFpHk00hQi  
UjSNPGRoCpHmISKFNLLv+/3q++71yq729lmZpbXcHjv3XvOd74/v+/POffc6izDv/  
9TZLVP/OulrLrtx5Nrf9syy7LOuLXGLY3bME1+z/  
P7PK8v2qRvi9fb48+FcVsetypNPvNjNhfuz3FMYp8W6fTHbSbPsgH7ILzeZeuQ9ix/  
z3L8kG0p0O3xfotjOub8t0m7y3st/s5IT+PmwvyJ49r8PkdaJfvOk398ryj/POfMwvvgO  
+ZLeRE/3V9IHc0Au6S/j/vl+z/hd/lqW+lQc26YMmq9N/  
bTyCf2S9BPprFGWDudf5n3ZZEibgt76DUgPco/GbZ1zSe99zil  
+B4E34Qh9h9TbIvsnOV62XCKdFXMcgs4658B3ySlby549ju9y/  
oxzJXMZhAfzsqRMA15fzP  
+bc8emMNtnnwXOU1GXnaBz8DvHMYucZy3l1qdsyxf5bX5oGPpwwzyzJNeRXmIM2Eqy  
91+KxzTo35XOL901aMMbdKcoxy6XprjsscmzC/R/sL  
+MvnuU997c9fpHHGyxngSbTHMsUhZu415lc88eW+FOUOrO2zWPKdF/Jc/  
1+UTmUbDXMNH52Ry7S/wtmQa83jX3K7U5cWzDB28gLvqkU5rjS/YZUQfwRcWcfrC/
```

```
eBjSTuvUiTDTCXaVPXVdsu+jn68GXXXJ30bf8edN1Cn6yiZr5AP6rijT6vj7+  
+P20bj9cpwnToLGuKXx9yvHbXH8/ei4HR  
+3mfHvI8afb43b7vH3lfHnaNyeH7cXef21cTvB8afGnw  
+OP4+NPy8ff54Zf77APmjnOPerpP9mPqF5RzG5vs5r2fj3znH7STGhfTif0BqStzj  
+yLidHrdd477v8h5kvG38+4nx53v55Do+K947Rpk/R57XkLPY75/ziU7eGLEz4/  
YbZvEq53+Z/Y+yzx9ynlMcC34/pPvxjz8Xxu0dynamcelsmHw  
+R9m9RV39Ge5zIXMc49k3y3ufny0H3x0jrQD6ZD/dkXofuL6MNwdOfj79/  
ppjwfowynyHvi7yGT+j+Q/J6hP3e4tzPU0/ic0cx4fE06bxAOuF8PMn2i3G/  
xzgGfD5NnbxEGd/gfEu01RnOtUjdXzsef3MxoQ9ZD+UTvZ8IOABed7HPc8QNdp0+9fAC  
+TpFeY5Th+DlUVH91ynrTvrBcfIH/I/hulq4Ocg+L3KuHdTrTcUEH/  
AX2PDRcXuGeDIJ3QIPa9T506S7t5iME8/v0p7rIO116gn  
+cZy6PE0ZRuH3Wepnkbyep7+RzGhJzmL8e  
+reV9+D9p9zvs8534yn2D7OOkqFnytmMj/HvUwyN1Xldv3C8cCbH57MekL2WdJ/  
y3K9D77vUU+IO8DxMgax4DuAuUdkNc9xUT/sO0H4/ZgPtH7iJ  
+r5O0UbXSCNGGDW4vJ3OD5KY4/SlnZ6d5HzHqDspzIDxcMf59kPcxj+IZPv  
+EcfAN6vut0G+V14XDdyn/IY7/38TA07nHgSsYT05Q/  
69QH89QD2+w35C0j3LMnbQ5bAG+gO3XKNuA892eT/wd8aekXUH3I/Z5n/I+S15X  
+fsV6v1F4kf+8R3GhNeoyz9mbK6IFWDqGvKMOLnAsZ8lb8/Q3sfYXso9/  
skPXqf84O056vkw59tBWcFvUbhNrtTJcfJwjinO9Qj1gjbv4/  
j8WEzrvch7lwAH53Ekab9BG6lu4epa66JPeq+RpgX1kN/Cxr/  
A4cJz9D1Kf0MmzI0tp8lmRxxcoD3AAv/  
g19Qzswp9+n3HoMHIS3jpiORXL1ijLU9Rhn9eXaXfE2nXa8jXSXqUejlAv  
+D3Pcc9znPzINc414PwVZQXfj/  
P6UdoQsRlylmeDBrCIWuSzbFawLXMH0Pq7gzvVW5XTQLbCyPvUD9PcK7TtPuLufvT4  
dx9FHY6QFqK7x9wzDKx8si4/SXrFdUNsOtPiYm3yc8RfkLHZe61g  
+qGJynnKnWzSNtjngeKCV3cP0s5nyFOQA94/CJ9CveQMw+R1w9pU8TFy+iHH  
+aO4eeozXOCYwsEQNHeP2qkCdO08Ynab8Z+sabIA  
+4+vn42n9k3lSnHfnPkPy9QxmG5FNyO8NP  
+BRshxyRFV7DxNj6Fuf8kLy8QJpvkdevEzOgjTrgHK9jz9lrfc2x63QDkvUMfjrk5+Xc893i7lj  
dIXjdtGvjuaOwavlS3DXYi3ySO65/hXq5i3ygTHlz/KtB4jbjkjo6Q7t9qfCab418oT985XF  
+nmRf1Q1P07all9DDDtWdJ4kWPe5vgF8vYyx2L+AVtJWuDtNertJfY5RBpv5O4nRn  
+Gfldpr9OkC2zuJS/g68eFx2Ho7duMARh3eeE5TeuMjDXyCvVzmLzso4Ud1G3V7x/  
N3kXzRPkXfEYPpFor4+oA9n8XO5xY0Q8HKIcb+dey5/m5zL1+zbl  
+zzrq6Occ4Zx7BDteoS0X8+9Lnm94X3l6mfNznHilTB32/  
ly12MZ5DpSergMtZGFXX7EunNMH7eSHzAtgfi91uc7wPOh3sYh3h5JfsflnZuYJ2i2HiWc  
82w9n0n9zwGOqNg/1ni6gR1ejdz9jnK16zZQGcXc9sR0oHvfpEynsi93niRNDD2BWLmo/  
9eQZ+TH2CO47nXnQelyf9cTHg  
+nXsc31N4fQX6j1HWs9Q39HJL4TXG3cTULuaLGcZkYG5AO4PfXj7xmQ+pC8Vj1C4/  
ISbfj0wD2LZH1Jfe6hX5Dj415OU9x7ytUz+3yMNxNn/Xkx4PUudCRfQB
```

+LJ6YDDPawDXqb+IftN1OFc7jkK8x+mTS8jFoX9E9TX67TvCdL6lPeaCfO/  
wvGK9VeSzx3Up9ade0gfgJvWme9k/s+wpDz76S  
+FQMRI5HPDhKrixWeH0H7XuL9Kfb5HPX4Lsc/RP1cRV28R14foY5g+2toZ  
+kf369ivl4h3x/InmuflyvYeA9yP05a0MezvPcu+TpHHSmGwd4Vf3+D  
+Rz9kVeuofxv5h7HVDPfwbh4gLh7nDlska8PKf889Ynv2KsAjH/2/TNf869TnqedM  
+ST9SPWGO+RhnRR3tQ4PErTNEHAR/g9RbG83eoj/F1hON20levpY6uJF/  
3sAYAnfulST18lvNh/j8jJqGnP6KdYOev0d8zYgt7Y79L4vVeWsxIOxjsDbn  
+P1Z4vYMMYdrl/J7J/2/  
R7wpPyfCL3cRV0DT4dlCzzCbpAFOr0v9zXMO7nvpcAm8MdV2mORfEGWPynCfXapl8  
RDyPcMbfEEbTbKfW9O9Rfaw8TqiHY+yM8hefw15T9GDDyZ+1rvUY5V7H4t93p3J  
+3+InUGXSE3az8Oekad/Qe012760M2U9W/Hn1+m7F+nDyFmA0+PUZYniXPg6gTtAZ0/  
TT5ezz2flaYgL+Ee4i5q2t8rJjSA/5fJM8Z180ms0D4I7114/  
kXGamDqHwrfS4MPQh6sW4DXQe5r+gXa5XL2Uf0JLOS8NiTmr2M/4AixBzb+C  
+ldmMA6+wv0i5OU/2vEOHT6Hfr7KfrRC8QsaGnPR7U+fOmvSHcHY  
+labQtbGu2yTD38BWPBG9QlXkPGW2i354j1mcJz9RnKfjX1skgsdlsJP  
+ALNcYT5BXzvUicSPZriKuMcQrte9QZcHoj7ao1/QfkryB/T5KPPQ8SK8ozyBmx1K  
+0Cv3mc9r2aMkEexc0rGX+kP+xFPUHMzRauO+2Pol/y9/2MU29Svm9QBx+xAU/  
aA9Fa4i+JH8V/tG8VEx0ipiNXHSROpICeVLvcyHi2k7xgDsXvTjFz75OvZ8I9h5h/  
1PUF3ztCHV+N33wvHn9ZRpn/37F4yB/4E6gXz/ddx+VEzy52flww8KXwvPc47/  
QrtBnt/Q/srR4B9xDRHhm+8WXu9iftU1j7G9SixDZ78af36VGAZt+P0vC68HYfPKt/  
lcS7g4gFeB3892u431Nsdpa8+F4gvxCrsT2b+3pVa2jY8wpi+EHGD+05QS8/  
K3xf/hvE/OPs+136+jukqTXZ/Yxj6L  
+fulWOFH7eJo4RJ09RB4gdkBUYOEledtHGh0n7byg7fGKVnolrqmPE+yni  
+yDn7vHee8TD09THfn7/PcYxPYNpF57/PkOebqXd7mJ82UHdXEk/gk9eThrA12O55/  
eniO03iZfvF76Ggp0f4nfQzTj2dtL8FrEL3ruUFzqFnyu2lw  
+22Edr2Q9zz3XAgZEensMgV6A2+AP6DewF2yMOAn+IG9pT+E  
+F568+sZoy8gpwPX75AV4ga7nigk/s/SfZ2iPhxkLzlG3j/  
MTtZDqLQjhm4XrHvOucI5ni1jtL326w4QzrulMch0NeXTMzn4JnR5beG1AupJ0Eae  
+G2ORRyC//9T4fs8/3P8/e+Lydrn6sJ5/iZ1BJ3Bx  
+QbV1Ae6GBv4fvj8MWbC9+j2Uk6lxMHD9KPtJb/  
Fn0CFEGWw9Q5bHIGGLmXfTDHX5I2MAS  
+tYf4TO77D1dxznO0pdawNxEpElcPo96FXHI57Qd5L6Gtr2fPAArj9KmmOPOwvdBkefu  
oy7W6Fo/6t8om/EmppyyWz0fiP18RXi9HLS6lInkHU/fellbX+cOFUuBO60LwC  
+sF79lekjPyA+fpV+u5M2BGaQP4Er5JzfkSZxHL6BmAVb/y/q+tf054eoB  
+gH8QT7kNcVvu5H/vgy7XcVfeNo7rleGFohj6pjYLN7aGvw/DB5upY22kia93l  
+0DxI3SOOlP49GuRBRFA9Ccz8gjxiHGRBRP4ufQ254Oe09S7iCTa/  
iXJat9CHnl9DP3rGBFPd2gDjXiUPqAm0T4U5EYseJg3Mc5ZYBP5hE8QM6OQG  
+jDmwH7VAxwL2x6hzuCnmG8/  
ax1E3CiuP8/5gC3YF3H3MepNzwGA4RPUH2qMhxDj0NM6ZQAWsB77NnWSCo8BoA

+9IbZDZvfo21cQZ9FTN1FvAB8HaM8XaP8/  
oo8dyD03Q2bw8wXqCbEDsQf6OExsfZD7s5hfco7P0IdA6++oV  
+jtHzj3KdoGdglWEIuwThtS7z8g5rvUA+zzY15HPvWJ+/  
648PgFfD1EW95DP7iSvV33pIF18G302/9WTGI5dJMoD+ykZwRnyV+buNF  
+I85NIXbcXXgcsygvZe8wDalzydz399GLIPu7+b166k76BUxAXEJ/o  
+8fQWx9C3iCDEK9TLul3ZfRt+RrcHrrZwX2EW+xbqloo2g47+m3T5Ln52lzXFPz7fuLHzP  
+jTt+QXqCfbFvNANxiNuvE9/2VN4HkUu+wpX+hxlqyID4jvix735JKZoTo/2epT6fzp3+/  
yIGP6IfOD3L+kXu+hjPySu2vTLXcSM9nuwPt5N+XYSS8gJiN/lu18INhG/  
4dvAP3wPeeF71Pfr5Bd2+R3Kid/7Cq8fgb0vES+ly3O0zTL5fom83Ue8/  
BbbE9Tvo5ThPuIV+L6j8PUccLWj8L2ez5BHyPIP5BN6up+Y  
+jx1Awwi9uFaThmRI26gXqBUL2Liz2hjnNtHn1P8Rh48xE/EauSNO2l/5BnFhO8QA6B/  
BfWTkffvU3ctj13hvW+SR9RtmA++ijoX/no5ZQbf+QXcRg2QzwoqlD3OQ77ID/  
IOPgx5oV9ztC3dTbrSX7P2ZDHriEeHiYGYO+naL9bSRPrTgRe/ilm9rbPESMQT  
+lrXqO9iViAPZKpHVXuAZbo2YDP6/  
QRvATxODbOSf4Rr65hf6w7j4KXENn7yLejpMviraDrYgJpEXbiRujhOLiAW35RNMl48cp  
W4Psv208HijZ4nw5x8Rr8AteEGsgG99nphR3YxYdwt5hM4x5tu02xP0c/  
hmmfueDrCNWAlfuf2gCywP3jCb+Rj6O3Pqb/HeR81xm3k3B/  
qB3NeSR2Dv5dyfBzKfd//rzgv5sS9FcyKQP7gH2B9dPE7Bep13Vi7HLOAxvo  
+Z72by8njrRWP537Wk3PfN5n/zOkt4f6BL33iaWM2UesULxZi2dnzpnWrsKf/  
4IWTd2NdLWvuo/Fp5vryv8ec11pAVc6BkHYtGOWs  
+wZYXld0gTugsRMr9nAb6wL4zhZ9HEyb17EDPza6iXiCn9sR+l3bA3KgXED+gq6vIr/  
bUcR8xR2cEDIO/  
wJtqSu3VzxBXOqu4kzyt0g7nSOsy2llnjs6R3zeJb51v1L7z31BG0HuGssI2mLNFfd5Y  
+HO5jHxhDNbtf1q4LHrOPUMcfJkY/h5xgHod+zFYFyAH4N2H3v5JrdLnWWGdge3w/  
OzGeZvx58/yyVnfPfnkDK7O3GfmZ7+H/  
L1xvjInzVzqvQqdYV8OZ3ab7zm0da6YZ8er5GeyN855J3/vYGT  
+foPGLnLmivkZdL0noPcInS8Y81xxx/y9isr5751frxK/v4Dzr9DtjL5ewQ6Lx/  
fUdB57MTzzjpD3jc/i73A+2U4U93ivJtn8pPrqTR/RyXqWgFdDx67zXHJ6u9O9HluWnOJ/  
4q0Wjz3LJ2m5O8zyAY3EA8Zdag3jqz3Q39dV6/S7rADsbrHQXJrn1bsBy4Hc+jB/  
yumxTJdensFQl5b0VgUc6Gx7lZqDflvmZ9g7tA3otvPJGJ3j13n4khgess8c  
+VoK43CufJF98XsgPfn9HyPzl+bv1nTM32XQOy2LxvlJ1vm7/  
CkgleNT2J5FHSAMcLicujTS86/3juYY1/IPR/4LINfpeTvusiv5UvCocZsxITcsTZvjvPd+QQ/  
Xfn3K2TbFvW2Yh5nMI5bDWOGXIF0IWzOUs/CEcboXZ  
+MuO4FbLUDrR55bvE6cAuZ5K/DMH7l+232nWe/HvU4oJzL4bvehamS9ymTY0zve/  
Sox4H5+yKJPC8FO65xjN7NKclXt56xwZz6z0M6XQ3bYXoXegcCtRb2c7D20jNv5ArkCT  
xbxecsP7+eT9Yu+I5cMmfuHx3aYiQsUP96n0OxTu8gtc3ff1JOMGU/  
xa8OxyyY555F83j1d1b0zIHHPEbr/bcO5RWG5bvL1FOyt6/  
ireJ3N9DUPip3irVW1WOU3vtaphzJPEcoJvXDvIk8IR8r5iXz9w2FP/  
VV7GuZv6um7+JPvjdV/o6N3vOJuSAL/


VrmuV6xqR9+y2bKOfG9sTguM3+vTrlrX0cbRHq90E/Yr5JjRO  
+BLQS6KXxXTuiax2fhoGX19+Xa5u9LiV/FFtEozf27G75LX9JV2zy/  
aYxi1ZJ5TCoDLfHYDrRk70RbKaag76p5jIj9ulbHV9+8JkM875nHd  
+G9Ob9s3DX3n9l8rivnt/k75pV5yqYqrynviDviDeaT3xKNxg3MI9dc  
+ZY1DuL4jvGBWG7NK8bcW2d10qrvz8oHoUj2TfGHcVKYUTxJMYD2WQh8JDM39GV7r  
vmfjtn9fcC9U6e4pPqC+k1meNo0LCr/HbJPD4p9rT42TGPby3elx8I+/  
KzKlwT7S7HxXpeehb2ZOvSPD5mwQ4d83diO+Y  
+3w3f9amYIsxj6EL4Ljwh6PAv2ruFX6uBZsrDirn9szfuYxxPwXe  
+ubxUTFoYpV3jT02RfN8ID3gPnxVdTTsIPwPwvhV0lwl74rvmWxtXvO0zGOH8CW9qr6U  
TygnxfCcxVfEu2Li7ICby/eyMNe6uc9s9Es+p3DcN4/  
JVfKYrT6KFXPm9exs0LPmki7IG6XVsd3t4Fwnpn7FD6zvL6OiHm4ZXW/a5nXFJn5O/  
0Rz8KLYorym3Qm+vK1dfMY2zKPM13zd5hjhJPWovJz2U/  
zdM2rxnkT7J11zzGL5rnmirgRbx0gs6rwJ900uP4OfN30TEen/  
Gpo9feZZdMu5VbMlo0li3Aivcnuqj+Vc4bmNXwZfsuPy  
+TFV4PtIszwwjx5G5njps0+8+YxbVvk2zv195CXzWKI  
+0kfPPB5l7gH1sUJe8FmlOr777Cv9KaYpBqyZr+uEkZF5LJV/  
yY9nzesarTuGwd9WyJfWTLKpYI3HPNcl33Pm+VV1pHSvHK/+ikNal82Zx6lB518zX5/J/  
oqnWjocqnlRB14rBwublfl3VMq8Re+Z/L0BYXTX3k1XzOljrTck5F+ywEuaW/  
hQTFNokf1zXmjDqSzZQXF0zr4vXzN9dHwQeKvO/uyD/  
k036gVY7jM1oP8WfYjxmd8M46UW2UaxRbTfrjqGZsMbv2k+smb  
+twLkm71AT3NXoV8/2C6Z7wloPTo0j  
+PCknxOdiXcEn3F8FXzGKcafTMG5HWMylbJHFMlgVfVqJJN+Uxy4Try9wrnGgY6irVLVI  
+TLpivERXLZ8Jabs58r2KY/G+5gN7lvA5T3JQ8wt+KeY5R/JWdM/  
O9BulzaJ6bZBdcVz4XDfnnvNX/toVitPAP3UHXq+a19gl/  
Ze918xpgxTy3prAn0bfwtzCS21ExVTipkteoiumqqzKORcxxOCD/iiU94mTY6Lcc  
+rbMY1PJOUfmda7ynXgQRoWdlWDFkoyTXbrChHltV5rXCp1wbZY8q64cprDva/  
WaSraWDG3zOmMYajT5UAp21b6T1hO95Pve/YZsc  
+axQHoahDIXrL4um8k95ikW98zXUYpfkgd9sXcq3MmvhuZ1VMxVG/pl3FoiTfC/O/  
e80zeP5cPfkfXcEOoHsa+Y1rmw  
+R16H5n4qOl3zNQTuXK6esmMcP2axtXp8qnquOqILNFZu159k19xXJvGr  
+t3QUa7pBNtWTK7Jxcn3LXzvm68uR+dpCthya1wQ35J7XyhTWiMISH/  
Flrlb2qnVK9NlXmQ6youlaYM69VNnyDPrkY5INsapEWxiufrAcczJm8TU27Y  
+IntZkQ95TjBHul8zx2zKvy+fNc7Fqu8zqz5Qgt2JXRr2CZi953TY0ryNK8/w/  
fYp39Dq9YBwPTTfExc2F81jYivwUfG+YlzPHIvyJcxVmu8hJ/OaVzZu577HpDWr8o/  
Wv7KhsDI0x0ab9+J6V3WM+i5Rr6rj182flcn/  
1sz3FrRGHprncbXhmdx1lp2PyMvl6liEDHtyz+HLnGfZPNeqpbReM/NnN  
+JHcUz3FMe7gZ9F89pf9b3spDVMZV5jL5nXqOBROUn+L/  
3OmV9dp4rzKWaOzHMD7Kj4o9zVDVhQTZMljzOyp/  
C5S14XzLGXzGvyWsxUfZfqzPbQbbNOJHcbzJzv1dNt8Z+0jVyW4r5+kk5TvpSPBB

+FAsq81wlnsVjlbweko2kw8y8Tlni+JT8+fPAvDZaMY8IHfLT41pC8Xc50NN6YM782XTL/  
BIRmer7SH3zfUvFS1XnbXl8b3kvhtzXld8qwbvnuNwxfzZ6qJ5jawaqDLfN1wxr4E24rp5r  
QHes9z3M1LydbBkKZPXOUvsU2os  
+RCWFK9XyXcZcNMY30MrzfO08tK6ef5QHfVc0GcvYDEzX/P0rL62XjZ/  
vjpn9edurSCrfFB5VDX8xrok7IH5r18aL6vN2DTXoT2offlXstGH2ibrzsUzzBG8SIFPMqXIY  
c3Yn8Kz4Npm675mgGxQLqcM9+fURyc4Vzgd0/u+w  
+95JhdM8fzaviu2l4xe4mbJfPcLz6kG83ZM89tWodk5vunivvr56itYRqolgHjIKdl8lcmz6RP  
F8q90jetvleh+abCTkPNPfmPpf6CYOyr/Dbon5FS364bp4j5bOr5ljfk3u8Fj5uCDIhzup  
+pZioprq+a76GGYX7+3Jfr+3JXX+VeZzQWmA26G/Jfl2j+DgKOOOpbeL6YXA  
+yd0nZFs3PEiyY5+sl0hZP0t2q+Z7VOvuofkzEeiLWY1B/rIu9zh6FT/  
Fb0naL5muOkXk9LUwpluzOfV91ZL42qYK918zj1Q25x/Zl0pc/  
rpnvfSe2VfM1eMe89tUe1bx5zlvj9Y2ck9drdcUF3F/nfeUO5Y+9uccOXLs  
+D3sZ5nVrSZ3M5J4rkvmzJJ0JUTySLw+CzOC5a17bqP  
+6uS9luecj5RedsVihLtfNc458MgvXNa5LnFS8H/OOsKm5latBP0ueS/  
ukJRxs5ljqDv9UC3fM88/  
A6meR1qnbLPf1x8b5vOS1VZW89IGfjbhCfhVzwKNyeJnCVhbHDFN9fd8IX/Ip6Uj3N/  
Jicv9fNq9vNnwi9xikPtKjz9gQDbAn9xjfDp+VuQxLVj8Dpr3KuWD/jzhMTOrs0qz5Pta  
+3DGueLpAWopTl6tjBbSBpYH5+qkd7JKZ55Qh9SI/6Af+bsh9rTUKGlxxemheP23k9bx  
+fg4866yOxsimyrFaEykmrHFO8a89SN3rB/sMzWOU6s8+9QJ96tzjZ2dPE  
+q1s049vq8f05QuVqxuWse25STNmOR81ppZF6PLF93QDMKF6sU65I83XRsjl+  
+4FG23wfpZPfQ1a+2Zd7DO4xFu7O/ZxDxRiwO/c98g5xc33udaLwE2Prari3wHv41zWPk/  
J76b1n/oxkRJuX7Ks8ObL6GUzInDZprZivj/c0/  
GLVvOaOsXXOfN9pleBpYPW6f9k8dq9xDuiult6ZNU8niyYr3cWzfcUdD1LvmbpBJ5UF5  
ZBF4k4alW83zGPY1XyWDXPPaZUX0spj0qPpfn  
+YdfqeXnWvNZbpbqzCQyv3cwgZY0LNPbITqeV+55Fn/  
jcnfuaSXFNUVYyKLBpC96RuX8l8rpS+bnHzYluTO6+dRkXfW1X3Jb/  
Uzt33z3LdKunquvm6+xy1M4Jw43q/  
AeUCcG8S5fryrgncw8H4Nzg3inRG8f4z3M38IPcKf2b  
+TT9772HivJ5+8GwMaeJ8d7x3cIE/OHOL9Nbx38ybyv3vq+hvImDMAj  
+XyNvDn9P8zvE653uE84XE5zjAdLrc/zTgc8XKNfn5AtnJfh+Ad45maMceG  
+tJB2MPcl7fcqBuX6Vt85b9knrRdlDX/eR7IF+/yHpo89h9r  
+HvOpvNN1EOQ9TpgPko8N7D5PPRN4GpD1Pej/jPPq7l4/k/  
vdRbiedN9n3AdKcJ83bSotX/L5A2Q6SR4w/zu+/Dt8fYZ+H2B/XZmkX0H  
+MNryJ8+pvauldmysL/9sreJ/sCr638Tm+F3MV35HA+yN4H  
+m3+W4E3rHYx3cq8l4f3m96qPjk/ysCZ1w/VUu+16UaJ+PvjDUd+rUUL/k943Xd7yXfU  
+kwJqK1zPemMI7T3FWq89EL9zd44u8q+X4K7mX8ngLv4kXjUuNTfOq35ojzZ2xtq/  
fLkvOawjWNFZ8p6KBlfi8L99q8rjqPiyWxwJd8zNHwQrFPEWaG/SCN9Tg8dEXXWC/  
crkc0ZepbPhdVqSza7rKkx7zD5Pk4WmsbH/tJxFviJyZUm6jU1+WnwVW7Bd9YU/  
u9hW23alEO2TXFubaYr/wE2IkYsx2+lh57je9NH9qO16jj7e5vX89W/Dd95BN126S//



8uj6jXzmM1FT  
+0FmLKQuqPMWUR1TseZ9QTVeH3NvCVUHXyamc9Qpf4b5q6giuRzzF5JZwrOM381n  
Sq/yBlqqLjxO85RS8fa6jnLZTps0cB5msjkcJUzXad9Hjc51y0yBv3A2Zopa/  
VPnO8uZSS0csYHlJrziHO2UXLpE876jBLr2znvr2Ro7eDMLyjW7BXn/  
o0qdGZ4ad0DnnY9wdcfZ5x6oWNCb7i7WYC7QlMxJdpVX7gFWIF7wYIghR6H9Mf40AH  
gjnBMPwgPowdjTKw1uC8cTrbFvbQhcMp6C9wdTAVD0VJoh7eLh4F7RF75CNyuGomR  
NfbgTmFsdMT1G29jaPM74H6R1fYumtqdYNS5Gtw1MizG4lpuLabZjAP3ji8dxqPOaQJOL  
hPBDuALD1fUeLmhj/dksA9IDJqK6mB3mC  
+dBnYDhtXTUan3RFf0DLAniE2YhYpkL7xM08DOYFOOxJkCQkchwP5gfak3TpX74JcqX  
7BL0NfPRnGjPx7fCAB7hYjW93GsbQ4etgeCHUOY2TwctpiPe7oPwL4h1G4BTA7BaHb6  
COwelrktw6PxbjltQTsIRYHLIUxaBmuBy8HO4ng0JXIWv0xmvShYD8xp/  
aPyEhYhSvJYWBXMSdrNVJzPkvtQQTYW/gXf4bk0jWoKdeDHYZPzVok1q9DdeN6sM  
+g5mgYWjeisu1zsNvw6oxBrNmfUG4RC/Yc0222okIXh+/T4rh3A4ZnGbj7eCzKif/9WfB/  
Xz6vFsyaMivf69Fgb  
+1FOXGtKMRC9F3kJQqebBZp9vFvnhqYdz1UOSe7St9rUbK0qxx8t4PbnJ0wTRZW  
+khJzyfKh+7usqBpjEyzd5Rhm8elvNbLGWvUS9EyMs7Ijvlqih3v1a4B/  
5NxHheEqHqspj20Q2hnXokwtb2IE47B8vKmc7yiaeSARPCpW1DiuwXc1Rmzv6H9BtsRg  
+67Glib0Uz4qPouv8uKt9aSs6TfqQlJy1xRzca0scXtkP4Gh5KQZUY4eue57D  
+xQ1cfN2JwMgBSr/FQQWFuqpHU6DuT/tQFZmtUiM9P1fB/bepVe+mqF  
+2paune7JU7W6jjimUq2IX5qk3z9NI/nnqckSuikoyqtOvs9SBHRlq4sxUZRuzXc3VG1Ru/  
2iVax+u  
+s9aogznA5XvVahzKz1Uq994ZRc3SpU0DFPDF9mof67RqVM7LNSKAZ3Ukr870XtHO6  
5at8FqeyvcTC34fv4t9NZfQ9qMBpQF16Pu/b8jveQCork/  
xYO8Mqw4dBID447hY8+DmLo8D9coG5udMmFvmQrbm0mYdzgOzqUbkHkxHOFbQjC  
0KghNK3xwes4M7C2YhMg1o/E8ZgR6PbDGfkr1Pn3QO/iZ7QusoVS1jVQ8sq/  
Uq5WRvXXjpCvVS45n91JxUmJtFfK7fKDX0+4C2fupNHqs86EKIC7navkMXY4dSgnEgp  
Ty3pFeXetKim/+Sd6/8LB9PbpaFExvl0p41csafK6SsOinvPymUt7flywXHs  
+V3zxpKtJBU6T82WcYZkqTr8ETZUplgy3skyLC18TJir0Gq  
+3Hce5xsc9kiA5dvlm89XS9TnkbK/  
Hth0vrcSulya4nMLVkgH9yel90DfeWxQCGjLad3O3lo5ARpU+ssW75yIPMSh/  
H5rOVTpytZds1cvrcKZhr2A1qFsw0Jr6uEd+OqBC7N50UA01HxKWyg2J09gHR8tV  
+MS/RKC4s2SW+HJlqun5OEufnGwTzdN8THisIRfr0hSKvkcSLogmC/  
RT6OhvRYbQUfgNea1YRdzR2XVvpdlTj/9GcdqZr+rpMrclqS+vjt+e/29f  
+87naHExpu559fBE9/bN8Y0ufd2/vZlnDr4vk7RXI2O09c+WaXzParquyVrE2MVA5OTPTe  
+6YrTlxiQtY8MX2k/56drsbKMW73hISzCe1f4N


### PLT



0M|3|HISTOGRAM|PLT|PLTALONGRES|FLOATLE-stream/  
deflate:base64^Y2BgYDA2VnYCUgxrJjOBaAcGhgZ7zRh5sBgA|FLOATLE-stream/  
deflate:base64^jdR7UNRVFAfvUB0UdQARdBFUBGjGmzV1Nzf75xrtWnMGFG  
+cjORAhtHB1HTMBWtRURzMSVfKzAUExHTBN1IHdzGN6RmumRARqgh  
+EBhEU2o7+83g1F/tTNn7uVy7/mc+9jVaDQavX6AQKmp2NARbbKs+edD  
+Hus0tE6BkhaR45ksbvQDpBLSyNli32hPMRkkTEu10XZMFYqW41OjLtkm86Dhpi6U2Zm  
H8J8WuL7AtVFjSSTRRDWkkE7mazGGNKZZxPyKA7ZdGZKMKRhnYwCSRuxNpOibVIY  
n0N2Vx5yWEmvP4A8NtqRcAS57ORrPY18pWSuu4CcDnKFVSCvk  
+LibiJ3HTly6pHfRUbnYxgaNik6MmpjX6snJxi0XFbalVEnm  
+t6sjPJlw3a3oyauSmsL0fbdGw1BjPq57i4ULa7wlhnDmfshR05EazX69IiH87YFxudo3IHg  
oFxZDAEjFdhGGGMZ  
+wXxgQYUTDegjEZxlR2hZlgTlcRAyMWRrhyMWTBmw5gLwHGfBgLYSyGkcQ5CUthJG  
PNCjgrse5TWGgasTYG3CutTYaZxkHkN9rYWedbBtiBXOvz1yLcBNWxEzgzUsQl5N6OWT  
OTegj1/hT1sQ01ZMLJhblfxDYwcGDth5MLYBSMPRj7q3AOjAIYVxrcw9sHYD  
+MAjEIYRTAOwbDB  
+B5GMfz8GMYRGEhdhIMA4BsMO40cYx2GcgHESxikY2ZCcxmVwvGjXjkY52FcgHERxi  
8wLsFwwCiH8SuMKzAqYFTCqlJxYYTjxjUY1TBuwLgJ4w8YNTBqYdTBUA3jLox7MOph  
3IfRAKMRhgvGQxjNMB7BeAzjCe6iBUYrDI1w5GiE1dhBmOs6CLwxode7Ca3DTTIT3AXe  
m7DYPURcXCdh0HoKvD1RF9VZ2F2dRWZmF4F3KlxOrQgye4mmsK4Cb1LsSOgmlvh6i  
2ibt8D7xPeyB4weMHRc6AnjGRg+MHxg+MLwg+EHoxeM3jB6w/CH4Q  
+jD4wAGAFczw6E0RdGxXj9YOhg6GAewegPoz+MYBjBMEJE22/D//  
0M8vrOYEr1kFZHnplGeb0oTzy3QDbqrsuvhQykMXMmUvngpXRnfrBlxxaT+4Pfyap+L8r/  
lIADpg7lQvRBQ3+I5urf3uW0cYs58XYKx1zN4C4TsvlEcaFbcm08a8QJdP/5M  
+8vvsJlfKozuTnXs5afm/  
GAp69w8YIbTVy28AnfW9OqRtq4Zs74s5kbr7eyeaBGKL9x6Zc1onxwC5vrH3LsFRcHvn2  
fq1bc4YwnNSyOV3PgJcDxRqzkskUV/DE5WBt5nv1az3N1RRI/  
1nCSlw7buXtuCTfUHOWZVUVcYd7LRb328OKxu9g7L5srUrbxurzNTGc28MBuX/  
DZT9ZyjNMPm4rePmmRajxlw71TIQ7j5dvjueDPtN41vBJ7N4wgffMHM8vzx3L1yoN7Ag  
38Pxb15A7nCOOhrBbUxBffr4fK+dIWuXHzqvdOKuwC8eP7sTe


```
+RoueO8xpY2vobufVxFqpex9pylgjFatvUQXavaS6aUneTW9DXOewtFXdplza9YqMh/
jRpZhasp9c1U9X5QI6E+yrZ+SEbd+8g3XW0DJ8VQRssUOuj7BkUcFmos2/
oSuTfKaht1cRjhzimxNhrR3rCPlzuNH9aLR83zle3cPQr1qtO8X
+XsR6pYfSfkbOst9a20bzNaytXAf4G9iePmXNRHcc7kx3hR2S8OzX6vLNdljm9VY3E2ly
5NKng6ZjSxo/6UvbenS5jz+rbx7/FasjX0c+A9aGqn3sXfaoX67OVdq2Pt6yOk
+JtjVKKOPKmr5pISz2m8fitO2vuR0oGrhTJ6GMvbfUoa09fH9Qp5GCecgKW37aP8/
JZTvYfvAGalxO6RYcsvgO19zSM6WRhSmSe3KsNGf3MAI3KbmXREruM6ZJbsn9MOdZ6
W8=
```

**BASO**



```
2M|13|HISTOGRAM|DIFF|BASOALONGRES|FLOATLE-stream/
deflate:base64^Y2AAgXpnMMXg4QQkHIAliCc4MjAoOCPkGuxBcgA=|FLOATLE-stream/
deflate:base64^7dR7aFZIHafwY800nbXK
+63lLUvNWd6yqcfNzEptlZm3cnnJeZvTps37S0SoRAyRGCL2lhIWpmlhEhEjQkQkRCRGi
AyJGCIhESli4ue45+DL+4f/2V898OF9zu/8c97n93u
+UZSsbeV3fqKpZdHdFUdRjtQzE1uf46RGlkaaiSZFURHFIBBTQSU1ZNhOPQ1kOchRTt
DIKc7SRDMtXOU6kKW8qoD2FFNGZ7vSmmAEMZigljGQspcRMLmv9b9OoYCaZmU8li6li
JTXUUsdmMnzCp3zGdnbyOV9Qzy528yUN7GEv+8iynwN8zUG+5RCHOCoxfuA4J/iRn/
iZRN7hV05yitOc4TfOco7z/E4Tf3CBizRziT/5ixYuc4W/uco//
Ms1rnODm9xK5sEMtOEBHqSAtxEO9rzMB3oSCGdelRHKelxHucJOtOfrnSjOz3oSS960
4e+PEkxT9GP/
gXgllN4msE8w7MMYSjDel7hIDCC53mBkYxiNGMYy4uM4yVKGc8EJhlziTLKmczLTOEV
pvlqr/E605jODN6ggjd5i7eZyTvM4l1mM4e5zGM+7/E
+C6jkAxayiMUs4UOWUsUylrOClayimtXUsla1fEQt61jPx9SxgY1sYjNb2Fp
+Nxfu10pzJc2UNFdSab7EITMKwn3eHwA7V+iH88jlqkw9X3n
+hkN8x2G17znOEbLsUt9R3ppv1eFcZ4UeDwuz1jfMYjKnyTwn857OfnInroR71hTu38lwfW
+EPKgnEvMasinJrZaQd9mQh2kmNoR9TcjOOE+arSVhn2ZuKsqpFefUmsM5ZvNkgjhvH
+Wlc2pJn9Je5Uvf3Utuv/Nr+fX8d/+v/34IPUH63xfBg==
```

**WBC**



```
6M|1|HISTOGRAM|WBC|TNCALONGRES|FLOATLE-stream/
deflate:base64^Y2AAgXpnMNWwCEQ7ABEQNwBxgSMQQ+QYGuxBcgA=|FLOATLE-stream/
deflate:base64^xdV/
qF5zHAfwg2EYU2zY2OzaLobhwcYw7rm7G8N2XVyzMVzMXHbxjGEYHhvrzq9u0lzMui2t
pbWWpJuWniQtSuVsknSTtCTdtLS0lft3Pp/pWP4Qyrdefc45z3nOeZ5zPt/
3N8vSeLp5sFTW12pt5A7EfQWxtp+nY/RSpZ
+sKcvqqKdETivtlKnQRTc99LKRLfRRZRvb2UE/
OxlgN9n0LBvCUIZRx0hGMYZ6GpjlJEpMZirTyJnBLGbTShvzWEA7C
+mgkzJLWcZyKqzgWVbRxfO8yMt08wqvsoYe3mAt6+hIPW+zgY28wyY2s4V3eY/
36eMDtvlhVT7iYz5hG5/yGZ+znS/4kq/Ywdd8w7f08x3f8wM7+ZGf
+JkbfmEXv7Kb39jDXjl9cAAHchBDOJhDOJShHMbhHMEWjuQohIPH0RzDCEZyLMdxPK
MYzQmcyBjGchLjqOdkxjOBBk7hVE5jlqdzBmcyibM4m3MocS7ncT6TmclFXMhULuJiLm
Eal3IzjeQ0MZ1mZjCTy7mCWVzJVVzNbObQwjW0ci3XcT1t3MBcbmQe87mJm1nALdzK
bbRzO3dwJwu5i0XcTQf3cC+L6eQ+7ucByizhQR5iKQ/zCl
+yjMd4nCdyzpm81fxHLvyjktIjj5woxfzPYs5tjn4bG8+/JX5/uv8zrOS5mkraXhG/
dVE8p7nXncZ4zxOiT1lfpX5LvTci
+m189MSUeJct8bw74prpHI3qS6zhrZrKWlUeDmZiqut403aP+nrt88FzXqPb9gvqanVV/
lcl8a7mRw/MjB5K/ZV6MvV5miv75lGaV7tiDm
+NXNgU2dEdedQZ2VSO3GqLfGulTEwZmZ7xQGRq3345W4387Ypc3qcc2iO7SwV51JT
xWdOf874+9vsL60E11odKrBd5bFcK60bx8yyv9Upxjfmr41nhOsVziuftr9iLxf2/07v/
dvwX1/j/x+8=
```

### 3.6.8. Example of Data Frame for Statistics

```

<- <ENQ>
-> <ACK>
<- <STX>1H|\^&||MHR1^210M2SH01010^1.7.1|||||P|LIS2-A2|
20230306080334<CR><ETX>E8<CR><LF>
-> <ACK>
<- <STX>2M|1|STATS|SAMPLE_ORDER^20221130202301^20230306080334^UPTIME|
SamplingNumber^^11431\AnalyticalSamplingStatistics^^12961\SlideSamplingStatistics^^1286<
CR><ETX>5C<CR><LF>
-> <ACK>
<- <STX>3M|2|STATS|RACK^20221130202301^20230306080334^UPTIME|
RackLoadedFromLoader^^1407\RackEjectedOnTray^^1411\RackLoadedFromConveyor^^0\Rac
kEjectedOnConveyor^^0<CR><ETX>FD<CR><LF>
-> <ACK>
<- <STX>4M|3|STATS|PATIENT^20221130202301^20230306080334^UPTIME|
RunReportPatientCBC^^5932\RunReportPatientDIF^^5715\RunReportPatientRET^^11\RunRepor
tPatientCBR^^299\RunReportPatientDIR^^249\RunReportPatientCBF^^17\RunReportPatientDIF_
LV^^0\RunReport<ETB>82<CR><LF>
-> <ACK>
<-
<STX>5PatientRBC_PLTO^^200\RunReportPatientSLIDE_E^^3\RunReportPatientSLIDE_EC^^12
82\RunReportPatientSLIDE_C^^1<CR><ETX>68<CR><LF>
-> <ACK>
<- <STX>6M|4|STATS|CONTROL^20221130202301^20230306080334^UPTIME|
RunReportControlCBC^^0\RunReportControlDIF^^160\RunReportControlRET^^157\RunReportC
ontrolCBR^^0\RunReportControlDIR^^0\RunReportControlCBF^^94\RunReportControlDIF_LV^^
0\RunReportControl<ETB>45<CR><LF>
-> <ACK>
<- <STX>7RBC_PLTO^^157<CR><ETX>15<CR><LF>
-> <ACK>
<- <STX>0M|5|STATS|RERUN_REFLEX^20221130202301^20230306080334^UPTIME|
RerunPatient^^26\ReflexPatient^^191\RerunControl^^0\ReflexControl^^157<CR><ETX>79<CR>
<LF>
-> <ACK>

```

```

<- <STX>1M|6|STATS|REPEAT_CALI^20221130202301^20230306080334^UPTIME|
RunRepeatabilityDIF^^5\RunRepeatabilityDIR^^25\RunCalibrationDIF^^0<CR><ETX>C1<CR><L
F>
-> <ACK>
<- <STX>2M|7|STATS|INVALID_RUN^20221130202301^20230306080334^UPTIME|
BadIdentification^^0\BadExtraction^^3\BadSampling^^66\BadCondition^^4\BadCalculation^^1\
StopRequested^^0\BadSmearStain^^22<CR><ETX>AC<CR><LF>
-> <ACK>
<- <STX>3M|8|STATS|FAILURE^20221130202301^20230306080334^UPTIME|
FailureNumber^INST^24\FailureNumber^CIM^0\FailureNumber^SPS^30<CR><ETX>37<CR><LF
>
-> <ACK>
<- <STX>4M|9|STATS|REAGENTS^20221130202301^20230306080334^UPTIME|
Volume^BASOLYSE^56753700\Volume^BUFFER^9808215\Volume^CLEANER^12607750\Volum
e^DILUENT^432352008\Volume^FLUOCYTE^9310200\Volume^LYSEBIO^7029900\Volume^MET
HANOL^10749060\Volume^NUCEDIFF^<ETB>57<CR><LF>
-> <ACK>
<-
<STX>514087800\Volume^STAIN1^6730781\Volume^STAIN2^938164<CR><ETX>11<CR><LF>
-> <ACK>
<- <STX>6M|10|STATS|TECHNICAL_CYCLES^20221130202301^20230306080334^UPTIME|
AnalyzerStartup^^53\AnalyzerShutdown^^38\AnalyzerRecovery^^21\AnalyzerBgNoiseCheck^^1
69\AnalyzerCleaner^^148\AnalyzerMinoclaire^^6\CIMStartup^^0\CIMShutdown^^0\CIMRecovery
^^0\S<ETB>36<CR><LF>
-> <ACK>
<-
<STX>7PSSStartupDaily^^42\SPSSStartupWeekly^^2\SPSSShutdownDaily^^34\SPSSShutdownWeek
ly^^1\SPSRecovery^^11<CR><ETX>25<CR><LF>
-> <ACK>
<- <STX>0L|1|N<CR><ETX>03<CR><LF>
-> <ACK>
<- <EOT>

```

### 3.6.9. Example of Data Frame for Remote Command

```
-> <ENQ>
<- <ACK>
-> <STX>1H|\^&||YP8K|||||P|LIS2-A2|20230306081037<CR><ETX>AE<CR><LF>
<- <ACK>
-> <STX>2M|1|EXECUTE|QC^EXTQC_CORRECT<CR><ETX>2F<CR><LF>
<- <ACK>
-> <STX>3L|1|N<CR><ETX>06<CR><LF>
<- <ACK>
-> <EOT>
<- <ENQ>
-> <ACK>
<- <STX>1H|\^&||MHR1^210M2SH01010^1.7.1|||||P|LIS2-A2|
20230306081028<CR><ETX>E9<CR><LF>
-> <ACK>
<- <STX>2M|1|EXECUTE|QC^EXTQC_CORRECT||ACCEPTED<CR><ETX>60<CR><LF>
-> <ACK>
<- <STX>3L|1|N<CR><ETX>06<CR><LF>
-> <ACK>
<- <EOT>
<- <ENQ>
-> <ACK>
<- <STX>1H|\^&||MHR1^210M2SH01010^1.7.1|||||P|LIS2-A2|
20230306081028<CR><ETX>E9<CR><LF>
-> <ACK>
<- <STX>2M|1|EXECUTE|QC^EXTQC_CORRECT||DONE<CR><ETX>4D<CR><LF>
-> <ACK>
<- <STX>3L|1|N<CR><ETX>06<CR><LF>
-> <ACK>
<- <EOT>
```