Yumizen C LDL

- Yumizen C230
- Yumizen C240

REF 1300148018

REAGENT 1 38 mL

REAGENT 2 14 mL

IVD  2797

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Yumizen C LDL

Referência

1. Centers for Disease Control/National Institutes of Health Manual, "Biosafety in Microbiological and Biomedical Laboratories", 1988. I have also seen this as: Richardson JH and Barkley WE. eds. Biosafety in Microbiological and Biomedical Laboratories, U.S. Dept. of Health and Human Services, Public Health Service, HHS Publication No. (CDC) 84-8395, Washington, DC (1984).
2. National Committee for Clinical Laboratory Standards, Preparation and Testing of Reagent Water in the Clinical Laboratory - Third Edition; Approved Guideline NCCLS Document C3-A3 (1997).
3. Gotto AM. Lipoprotein metabolism and the etiology of hyperlipidemia, Hospital Practice (1988) **23** (Suppl. 1): 4-13.
4. Crouse JR, Parks JS, Schey HM, Kahl FR. Studies of low density lipoprotein molecular weight in human beings with coronary artery disease. J. Lipid Res. (1985) **26** (5): 566-574.
5. Badimon JJ, Badimon L, Fuester V. Regression of Atherosclerotic Lesions by High Density Lipoprotein Plasma Fraction in the Cholesterol-Fed Rabbit. Journal of Clinical Investigation (1990) **85**: 1234-1241.
6. Castelli WP, Doyle JT, Gordon T, Hames CG, Hjortland MC, Hulley SB, Kagan A, Zukel WJ. HDL Cholesterol and other lipids in coronary heart disease. Circulation (1977) **55**: 767-772.
7. Barr DP, Russ EM, Eder HA. Protein-lipid relationships in human plasma. Am. J. Med. (1951) **11**: 480.
8. Gordon T, Castelli WP, Hjortland MC, Kannel WB, Dawber TR. High density lipoprotein as a protective factor against coronary heart disease. Am. J. Med. (1977) **62**: 707-714.
9. Guder WG, Zawta B. The Quality of Diagnostics Samples. 1st ed. Darmstadt: GIT Verlag (2001): 22-23.
10. Bachorik PS, Ross JW. National Cholesterol Education Program Recommendations for Measurement of Low-Density Lipoprotein Cholesterol: Executive Summary, Clin. Chem. (1995) **41** (10): 1414-1420.
11. Council Directive (2000/54/EC). Official Journal of the European Communities. No. L262 from October 17, 2000: 21-45.
12. Evaluation of detection capability for clinical laboratory measurement procedures. Approved Guideline, 2nd ed., CLSI (NCCLS) document EP17-A2 (2012) **32** (8).
13. Evaluation of Precision of Quantitative Measurement Procedures. Approved Guideline, CLSI (NCCLS) document EP05-A3 (2014) **24** (25).
14. Evaluation of Linearity of Quantitative Measurement Procedures. 2nd Edition, CLSI (NCCLS) guideline EP06-Ed2 (2020) **40** (16).
15. Measurement Procedure Comparison and Bias Estimation Using Patient Samples. Approved Guideline, 3rd ed., CLSI (NCCLS) document EP09c (2018) **38** (12).
16. Passing H, Bablok W. A new biometrical procedure for testing the equality of measurements from two different analytical methods. J. Clin. Chem. Clin. Biochem. (1983) **21**: 709-720.
17. Young DS. Effects of Drugs on Clinical Laboratory Tests. 5th Edition, Washington, DC, AACC Press (2000).
18. Young DS. Effects of Preanalytical Variables on Clinical Laboratory Tests. 2nd Edition, Washington, DC, AACC Press (1997) **3**: 120-132.