

# MINOTROL CRP

## Hematology Reference Control

9755 Research Drive  
 Irvine, California 92618  
 Telephone: (949) 453-0500  
 (888) 903-5001  
 Fax: (949) 753-0533

Lot No.: **MC412**

Exp. Date: **2019-10-05** (Y/M/D)

	Parameter/Units	Low Abnormal Lot # <b>MC412-1</b>			Normal Lot # <b>MC412-2</b>			High Abnormal Lot # <b>MC412-3</b>		
		Assay	Range	Limits	Assay	Range	Limits	Assay	Range	Limits
MICROS CRP 200	WBC K/ $\mu$ L	<b>1.9</b>	0.4	1.5 - 2.3	<b>7.4</b>	1.0	6.4 - 8.4	<b>19.7</b>	2.2	17.5 - 21.9
	RBC M/ $\mu$ L	<b>2.31</b>	0.15	2.16 - 2.46	<b>4.59</b>	0.18	4.41 - 4.77	<b>5.55</b>	0.25	5.30 - 5.80
	HGB g/dL	<b>6.0</b>	0.5	5.5 - 6.5	<b>13.3</b>	0.6	12.7 - 13.9	<b>17.6</b>	0.7	16.9 - 18.3
	HCT %	<b>16.2</b>	2.0	14.2 - 18.2	<b>36.7</b>	2.5	34.2 - 39.2	<b>49.4</b>	3.0	46.4 - 52.4
	MCV fL	<b>70</b>	4	66 - 74	<b>80</b>	4	76 - 84	<b>89</b>	4	85 - 93
	MCH pg	<b>26.0</b>	2.4	23.6 - 28.4	<b>29.0</b>	2.4	26.6 - 31.4	<b>31.7</b>	2.8	28.9 - 34.5
	MCHC g/dL	<b>37.1</b>	3.0	34.1 - 40.1	<b>36.2</b>	3.0	33.2 - 39.2	<b>35.6</b>	3.0	32.6 - 38.6
	RDW %	<b>14.4</b>	4.0	10.4 - 18.4	<b>13.8</b>	4.0	9.8 - 17.8	<b>12.8</b>	4.0	8.8 - 16.8
	PLT K/ $\mu$ L	<b>76</b>	20	56 - 96	<b>258</b>	40	218 - 298	<b>510</b>	60	450 - 570
	MPV fL	<b>8.1</b>	2.0	6.1 - 10.1	<b>7.9</b>	2.0	5.9 - 9.9	<b>7.7</b>	2.0	5.7 - 9.7
	Lym %	<b>66.0</b>	7.0	59.0 - 73.0	<b>37.5</b>	6.0	31.5 - 43.5	<b>17.5</b>	6.0	11.5 - 23.5
	Lym K/ $\mu$ L	<b>1.3</b>	0.2	1.1 - 1.5	<b>2.8</b>	0.5	2.3 - 3.3	<b>3.4</b>	1.3	2.1 - 4.7
	Mon %	<b>9.5</b>	6.0	3.5 - 15.5	<b>7.5</b>	5.0	2.5 - 12.5	<b>6.0</b>	4.0	2.0 - 10.0
	Mon K/ $\mu$ L	<b>0.2</b>	0.2	0.0 - 0.4	<b>0.6</b>	0.4	0.2 - 1.0	<b>1.2</b>	0.8	0.4 - 2.0
	Gran %	<b>24.5</b>	7.0	17.5 - 31.5	<b>55.0</b>	6.0	49.0 - 61.0	<b>76.5</b>	7.0	69.5 - 83.5
	Gran K/ $\mu$ L	<b>0.5</b>	0.2	0.3 - 0.7	<b>4.1</b>	0.5	3.6 - 4.6	<b>15.1</b>	1.5	13.6 - 16.6
	CRP mg/dL	<b>0.48</b>	0.30	0.18 - 0.78	<b>2.10</b>	0.70	1.40 - 2.80	<b>4.12</b>	1.00	3.12 - 5.12

For technical assistance, call the Customer Support Center at 1-888-903-5001.