

# MINOTROL CRP

## Hematology Reference Control

9755 Research Drive  
Irvine, California 92618  
Telephone: (949) 453-0500  
(888) 903-5001  
Fax: (949) 753-0533

Lot No.: **MC413**

Exp. Date: **2020-01-05** (Y/M/D)

	Parameter/Units	Low Abnormal Lot # <b>MC413-1</b>			Normal Lot # <b>MC413-2</b>			High Abnormal Lot # <b>MC413-3</b>		
		Assay	Range	Limits	Assay	Range	Limits	Assay	Range	Limits
MICROS CRP 200	WBC K/ $\mu$ L	<b>2.0</b>	0.4	1.6 - 2.4	<b>7.1</b>	1.0	6.1 - 8.1	<b>19.1</b>	2.2	16.9 - 21.3
	RBC M/ $\mu$ L	<b>2.31</b>	0.15	2.16 - 2.46	<b>4.50</b>	0.18	4.32 - 4.68	<b>5.62</b>	0.25	5.37 - 5.87
	HGB g/dL	<b>6.1</b>	0.5	5.6 - 6.6	<b>13.3</b>	0.6	12.7 - 13.9	<b>17.5</b>	0.7	16.8 - 18.2
	HCT %	<b>16.4</b>	2.0	14.4 - 18.4	<b>36.5</b>	2.5	34.0 - 39.0	<b>49.5</b>	3.0	46.5 - 52.5
	MCV fL	<b>71</b>	4	67 - 75	<b>81</b>	4	77 - 85	<b>88</b>	4	84 - 92
	MCH pg	<b>26.4</b>	2.4	24.0 - 28.8	<b>29.6</b>	2.4	27.2 - 32.0	<b>31.1</b>	2.8	28.3 - 33.9
	MCHC g/dL	<b>37.2</b>	3.0	34.2 - 40.2	<b>36.5</b>	3.0	33.5 - 39.5	<b>35.4</b>	3.0	32.4 - 38.4
	RDW %	<b>14.0</b>	4.0	10.0 - 18.0	<b>13.5</b>	4.0	9.5 - 17.5	<b>13.0</b>	4.0	9.0 - 17.0
	PLT K/ $\mu$ L	<b>77</b>	20	57 - 97	<b>250</b>	40	210 - 290	<b>490</b>	60	430 - 550
	MPV fL	<b>7.9</b>	2.0	5.9 - 9.9	<b>7.5</b>	2.0	5.5 - 9.5	<b>7.3</b>	2.0	5.3 - 9.3
	Lym %	<b>68.5</b>	7.0	61.5 - 75.5	<b>35.0</b>	6.0	29.0 - 41.0	<b>18.5</b>	6.0	12.5 - 24.5
	Lym K/ $\mu$ L	<b>1.4</b>	0.2	1.2 - 1.6	<b>2.5</b>	0.5	2.0 - 3.0	<b>3.5</b>	1.3	2.2 - 4.8
	Mon %	<b>8.0</b>	6.0	2.0 - 14.0	<b>9.0</b>	5.0	4.0 - 14.0	<b>5.5</b>	4.0	1.5 - 9.5
	Mon K/ $\mu$ L	<b>0.2</b>	0.2	0.0 - 0.4	<b>0.6</b>	0.4	0.2 - 1.0	<b>1.1</b>	0.8	0.3 - 1.9
	Gran %	<b>23.5</b>	7.0	16.5 - 30.5	<b>56.0</b>	6.0	50.0 - 62.0	<b>76.0</b>	7.0	69.0 - 83.0
	Gran K/ $\mu$ L	<b>0.5</b>	0.2	0.3 - 0.7	<b>4.0</b>	0.5	3.5 - 4.5	<b>14.5</b>	1.5	13.0 - 16.0
	CRP mg/dL	<b>0.39</b>	0.30	0.09 - 0.69	<b>2.07</b>	0.70	1.37 - 2.77	<b>3.95</b>	1.00	2.95 - 4.95

For technical assistance, call the Customer Support Center at 1-888-903-5001.