

MINOTROL CRP

Hematology Reference Control

9755 Research Drive
 Irvine, California 92618
 Telephone: (949) 453-0500
 (888) 903-5001
 Fax: (949) 753-0533

Lot No.: **MC414**
 Exp. Date: **2020-04-05** (Y/M/D)

	Parameter/Units	Low Abnormal Lot # MC414-1			Normal Lot # MC414-2			High Abnormal Lot # MC414-3		
		Assay	Range	Limits	Assay	Range	Limits	Assay	Range	Limits
MICROS CRP 200	WBC K/ μ L	1.8	0.4	1.4 - 2.2	7.3	1.0	6.3 - 8.3	19.4	2.2	17.2 - 21.6
	RBC M/ μ L	2.40	0.15	2.25 - 2.55	4.55	0.18	4.37 - 4.73	5.60	0.25	5.35 - 5.85
	HGB g/dL	6.1	0.5	5.6 - 6.6	13.3	0.6	12.7 - 13.9	17.5	0.7	16.8 - 18.2
	HCT %	16.8	2.0	14.8 - 18.8	36.9	2.5	34.4 - 39.4	49.8	3.0	46.8 - 52.8
	MCV fL	70	4	66 - 74	81	4	77 - 85	89	4	85 - 93
	MCH pg	25.4	2.4	23.0 - 27.8	29.2	2.4	26.8 - 31.6	31.3	2.8	28.5 - 34.1
	MCHC g/dL	36.3	3.0	33.3 - 39.3	36.1	3.0	33.1 - 39.1	35.1	3.0	32.1 - 38.1
	RDW %	14.5	4.0	10.5 - 18.5	14.0	4.0	10.0 - 18.0	13.5	4.0	9.5 - 17.5
	PLT K/ μ L	77	20	57 - 97	263	40	223 - 303	550	60	490 - 610
	MPV fL	8.3	2.0	6.3 - 10.3	7.9	2.0	5.9 - 9.9	7.9	2.0	5.9 - 9.9
	Lym %	66.0	8.0	58.0 - 74.0	34.0	6.0	28.0 - 40.0	16.5	6.0	10.5 - 22.5
	Lym K/ μ L	1.2	0.2	1.0 - 1.4	2.5	0.5	2.0 - 3.0	3.2	1.3	1.9 - 4.5
	Mon %	9.5	7.0	2.5 - 16.5	7.0	5.0	2.0 - 12.0	6.0	4.0	2.0 - 10.0
	Mon K/ μ L	0.2	0.2	0.0 - 0.4	0.5	0.4	0.1 - 0.9	1.2	0.8	0.4 - 2.0
	Gran %	24.5	7.0	17.5 - 31.5	59.0	6.0	53.0 - 65.0	77.5	7.0	70.5 - 84.5
	Gran K/ μ L	0.4	0.2	0.2 - 0.6	4.3	0.5	3.8 - 4.8	15.0	1.5	13.5 - 16.5
	CRP mg/dL	0.50	0.30	0.20 - 0.80	2.62	0.70	1.92 - 3.32	4.09	1.00	3.09 - 5.09

For technical assistance, call the Customer Support Center at 1-888-903-5001.