

MINOTROL CRP

Hematology Reference Control

9755 Research Drive
Irvine, California 92618
Telephone: (949) 453-0500
(888) 903-5001
Fax: (949) 753-0533

Lot No.: **MC415**

Exp. Date: **2020-07-05** (Y/M/D)

	Parameter/Units	Low Abnormal Lot # MC415-1			Normal Lot # MC415-2			High Abnormal Lot # MC415-3		
		Assay	Range	Limits	Assay	Range	Limits	Assay	Range	Limits
MICROS CRP 200	WBC K/ μ L	2.0	0.4	1.6 - 2.4	7.5	1.0	6.5 - 8.5	20.4	2.2	18.2 - 22.6
	RBC M/ μ L	2.43	0.15	2.28 - 2.58	4.61	0.18	4.43 - 4.79	5.62	0.25	5.37 - 5.87
	HGB g/dL	6.2	0.5	5.7 - 6.7	13.3	0.6	12.7 - 13.9	18.0	0.7	17.3 - 18.7
	HCT %	16.8	2.0	14.8 - 18.8	36.9	2.5	34.4 - 39.4	50.6	3.0	47.6 - 53.6
	MCV fL	69	4	65 - 73	80	4	76 - 84	90	4	86 - 94
	MCH pg	25.5	2.4	23.1 - 27.9	28.9	2.4	26.5 - 31.3	32.0	2.8	29.2 - 34.8
	MCHC g/dL	37.0	3.0	34.0 - 40.0	36.1	3.0	33.1 - 39.1	35.6	3.0	32.6 - 38.6
	RDW %	13.7	4.0	9.7 - 17.7	14.6	4.0	10.6 - 18.6	13.5	4.0	9.5 - 17.5
	PLT K/ μ L	65	20	45 - 85	245	40	205 - 285	500	60	440 - 560
	MPV fL	8.8	2.0	6.8 - 10.8	8.1	2.0	6.1 - 10.1	8.0	2.0	6.0 - 10.0
	Lym %	64.5	8.0	56.5 - 72.5	34.0	6.0	28.0 - 40.0	17.0	6.0	11.0 - 23.0
	Lym K/ μ L	1.3	0.2	1.1 - 1.5	2.6	0.5	2.1 - 3.1	3.5	1.3	2.2 - 4.8
	Mon %	10.5	7.0	3.5 - 17.5	7.5	5.0	2.5 - 12.5	6.0	4.0	2.0 - 10.0
	Mon K/ μ L	0.2	0.2	0.0 - 0.4	0.6	0.4	0.2 - 1.0	1.2	0.8	0.4 - 2.0
	Gran %	25.0	7.0	18.0 - 32.0	58.5	6.0	52.5 - 64.5	77.0	7.0	70.0 - 84.0
	Gran K/ μ L	0.5	0.2	0.3 - 0.7	4.4	0.5	3.9 - 4.9	15.7	1.5	14.2 - 17.2
	CRP mg/dL	0.59	0.30	0.29 - 0.89	2.11	0.70	1.41 - 2.81	4.68	1.00	3.68 - 5.68

For technical assistance, call the Customer Support Center at 1-888-903-5001.