

MINOTROL CRP

Hematology Reference Control

9755 Research Drive
 Irvine, California 92618
 Telephone: (949) 453-0500
 (888) 903-5001
 Fax: (949) 753-0533

Lot No.: **MC416**
 Exp. Date: **2020-10-05** (Y/M/D)

	Parameter/Units	Low Abnormal Lot # MC416-1			Normal Lot # MC416-2			High Abnormal Lot # MC416-3		
		Assay	Range	Limits	Assay	Range	Limits	Assay	Range	Limits
MICROS CRP 200	WBC K/ μ L	1.9	0.4	1.5 - 2.3	7.6	1.0	6.6 - 8.6	20.0	2.2	17.8 - 22.2
	RBC M/ μ L	2.32	0.15	2.17 - 2.47	4.50	0.18	4.32 - 4.68	5.58	0.25	5.33 - 5.83
	HGB g/dL	6.0	0.5	5.5 - 6.5	13.2	0.6	12.6 - 13.8	17.7	0.7	17.0 - 18.4
	HCT %	16.2	2.0	14.2 - 18.2	36.9	2.5	34.4 - 39.4	49.7	3.0	46.7 - 52.7
	MCV fL	70	4	66 - 74	82	4	78 - 86	89	4	85 - 93
	MCH pg	25.9	2.4	23.5 - 28.3	29.3	2.4	26.9 - 31.7	31.7	2.8	28.9 - 34.5
	MCHC g/dL	36.9	3.0	33.9 - 39.9	35.8	3.0	32.8 - 38.8	35.6	3.0	32.6 - 38.6
	RDW %	14.0	4.0	10.0 - 18.0	13.5	4.0	9.5 - 17.5	13.0	4.0	9.0 - 17.0
	PLT K/ μ L	72	20	52 - 92	265	40	225 - 305	505	60	445 - 565
	MPV fL	8.0	2.0	6.0 - 10.0	7.5	2.0	5.5 - 9.5	7.3	2.0	5.3 - 9.3
	Lym %	63.5	8.0	55.5 - 71.5	34.0	6.0	28.0 - 40.0	17.0	6.0	11.0 - 23.0
	Lym K/ μ L	1.2	0.2	1.0 - 1.4	2.6	0.5	2.1 - 3.1	3.4	1.3	2.1 - 4.7
	Mon %	10.5	7.0	3.5 - 17.5	7.5	5.0	2.5 - 12.5	5.0	4.0	1.0 - 9.0
	Mon K/ μ L	0.2	0.2	0.0 - 0.4	0.6	0.4	0.2 - 1.0	1.0	0.8	0.2 - 1.8
	Gran %	26.0	7.0	19.0 - 33.0	58.5	6.0	52.5 - 64.5	78.0	7.0	71.0 - 85.0
	Gran K/ μ L	0.5	0.2	0.3 - 0.7	4.4	0.5	3.9 - 4.9	15.6	1.5	14.1 - 17.1
	CRP mg/dL	0.42	0.30	0.12 - 0.72	2.05	0.70	1.35 - 2.75	3.55	1.00	2.55 - 4.55

For technical assistance, call the Customer Support Center at 1-888-903-5001.