

MINOTROL CRP

Hematology Reference Control

9755 Research Drive
Irvine, California 92618
Telephone: (949) 453-0500
(888) 903-5001
Fax: (949) 753-0533

Lot No.: **MC417**

Exp. Date: **2021-01-05** (Y/M/D)

	Parameter/Units	Low Abnormal Lot # MC417-1			Normal Lot # MC417-2			High Abnormal Lot # MC417-3		
		Assay	Range	Limits	Assay	Range	Limits	Assay	Range	Limits
MICROS CRP 200	WBC K/ μ L	2.0	0.4	1.6 - 2.4	7.6	1.0	6.6 - 8.6	19.8	2.2	17.6 - 22.0
	RBC M/ μ L	2.42	0.15	2.27 - 2.57	4.56	0.18	4.38 - 4.74	5.50	0.25	5.25 - 5.75
	HGB g/dL	6.1	0.5	5.6 - 6.6	13.5	0.6	12.9 - 14.1	17.7	0.7	17.0 - 18.4
	HCT %	16.9	2.0	14.9 - 18.9	37.8	2.5	35.3 - 40.3	50.1	3.0	47.1 - 53.1
	MCV fL	70	4	66 - 74	83	4	79 - 87	91	4	87 - 95
	MCH pg	25.2	2.4	22.8 - 27.6	29.6	2.4	27.2 - 32.0	32.2	2.8	29.4 - 35.0
	MCHC g/dL	36.0	3.0	33.0 - 39.0	35.7	3.0	32.7 - 38.7	35.4	3.0	32.4 - 38.4
	RDW %	14.0	4.0	10.0 - 18.0	13.5	4.0	9.5 - 17.5	12.5	4.0	8.5 - 16.5
	PLT K/ μ L	74	20	54 - 94	252	40	212 - 292	507	60	447 - 567
	MPV fL	8.9	2.0	6.9 - 10.9	7.9	2.0	5.9 - 9.9	7.3	2.0	5.3 - 9.3
	Lym %	66.0	8.0	58.0 - 74.0	35.5	6.0	29.5 - 41.5	17.0	6.0	11.0 - 23.0
	Lym K/ μ L	1.3	0.2	1.1 - 1.5	2.7	0.5	2.2 - 3.2	3.4	1.2	2.2 - 4.6
	Mon %	10.5	7.0	3.5 - 17.5	7.8	5.0	2.8 - 12.8	5.0	4.5	0.5 - 9.5
	Mon K/ μ L	0.2	0.2	0.0 - 0.4	0.6	0.4	0.2 - 1.0	1.0	0.9	0.1 - 1.9
	Gran %	23.5	7.0	16.5 - 30.5	56.7	6.0	50.7 - 62.7	78.0	7.0	71.0 - 85.0
	Gran K/ μ L	0.5	0.2	0.3 - 0.7	4.3	0.5	3.8 - 4.8	15.4	1.5	13.9 - 16.9
	CRP mg/dL	0.46	0.30	0.16 - 0.76	1.92	0.70	1.22 - 2.62	3.94	1.00	2.94 - 4.94

For technical assistance, call the Customer Support Center at 1-888-903-5001.