

MINOTROL CRP

Hematology Reference Control

9755 Research Drive
 Irvine, California 92618
 Telephone: (949) 453-0500
 (888) 903-5001
 Fax: (949) 753-0533

Lot No.: **MC418**

Exp. Date: **2021-04-05** (Y/M/D)

	Parameter/Units	Low Abnormal Lot # MC418-1			Normal Lot # MC418-2			High Abnormal Lot # MC418-3		
		Assay	Range	Limits	Assay	Range	Limits	Assay	Range	Limits
MICROS CRP 200	WBC K/ μ L	2.1	0.4	1.7 - 2.5	7.6	1.0	6.6 - 8.6	19.9	2.2	17.7 - 22.1
	RBC M/ μ L	2.37	0.15	2.22 - 2.52	4.66	0.18	4.48 - 4.84	5.55	0.25	5.30 - 5.80
	HGB g/dL	6.1	0.5	5.6 - 6.6	13.7	0.6	13.1 - 14.3	17.6	0.7	16.9 - 18.3
	HCT %	16.8	2.0	14.8 - 18.8	38.2	2.5	35.7 - 40.7	50.0	3.0	47.0 - 53.0
	MCV fL	71	4	67 - 75	82	4	78 - 86	90	4	86 - 94
	MCH pg	25.7	2.4	23.3 - 28.1	29.4	2.4	27.0 - 31.8	31.7	2.8	28.9 - 34.5
	MCHC g/dL	36.3	3.0	33.3 - 39.3	35.9	3.0	32.9 - 38.9	35.2	3.0	32.2 - 38.2
	RDW %	14.0	4.0	10.0 - 18.0	13.5	4.0	9.5 - 17.5	13.0	4.0	9.0 - 17.0
	PLT K/ μ L	82	20	62 - 102	274	40	234 - 314	530	60	470 - 590
	MPV fL	8.2	2.0	6.2 - 10.2	8.1	2.0	6.1 - 10.1	7.8	2.0	5.8 - 9.8
	Lym %	66.0	8.0	58.0 - 74.0	37.0	6.0	31.0 - 43.0	17.0	6.0	11.0 - 23.0
	Lym K/ μ L	1.4	0.3	1.1 - 1.7	2.8	0.5	2.3 - 3.3	3.4	1.3	2.1 - 4.7
	Mon %	9.5	7.0	2.5 - 16.5	8.0	5.0	3.0 - 13.0	5.5	5.0	0.5 - 10.5
	Mon K/ μ L	0.2	0.2	0.0 - 0.4	0.6	0.4	0.2 - 1.0	1.1	1.0	0.1 - 2.1
	Gran %	24.5	7.0	17.5 - 31.5	55.0	6.0	49.0 - 61.0	77.5	7.0	70.5 - 84.5
	Gran K/ μ L	0.5	0.2	0.3 - 0.7	4.2	0.6	3.6 - 4.8	15.4	1.5	13.9 - 16.9
	CRP mg/dL	0.46	0.30	0.16 - 0.76	2.09	0.70	1.39 - 2.79	3.89	1.00	2.89 - 4.89

For technical assistance, call the Customer Support Center at 1-888-903-5001.