

# MINOTROL CRP

## Hematology Reference Control

9755 Research Drive  
Irvine, California 92618  
Telephone: (949) 453-0500  
(888) 903-5001  
Fax: (949) 753-0533

Lot No.: **MC419**  
Exp. Date: **2021-07-05** (Y / M / D)

	Parameter/Units	Low Abnormal Lot # <b>MC419-1</b>			Normal Lot # <b>MC419-2</b>			High Abnormal Lot # <b>MC419-3</b>		
		Assay	Range	Limits	Assay	Range	Limits	Assay	Range	Limits
MICROS CRP 200	WBC K/ $\mu$ L	<b>2.0</b>	0.4	1.6 - 2.4	<b>7.5</b>	1.0	6.5 - 8.5	<b>20.2</b>	2.2	18.0 - 22.4
	RBC M/ $\mu$ L	<b>2.32</b>	0.15	2.17 - 2.47	<b>4.51</b>	0.18	4.33 - 4.69	<b>5.51</b>	0.25	5.26 - 5.76
	HGB g/dL	<b>5.9</b>	0.5	5.4 - 6.4	<b>13.1</b>	0.6	12.5 - 13.7	<b>17.8</b>	0.7	17.1 - 18.5
	HCT %	<b>16.2</b>	2.0	14.2 - 18.2	<b>36.5</b>	2.5	34.0 - 39.0	<b>50.1</b>	3.0	47.1 - 53.1
	MCV fL	<b>70</b>	4	66 - 74	<b>81</b>	4	77 - 85	<b>91</b>	4	87 - 95
	MCH pg	<b>25.4</b>	2.4	23.0 - 27.8	<b>29.0</b>	2.4	26.6 - 31.4	<b>32.3</b>	2.8	29.5 - 35.1
	MCHC g/dL	<b>36.3</b>	3.0	33.3 - 39.3	<b>35.9</b>	3.0	32.9 - 38.9	<b>35.5</b>	3.0	32.5 - 38.5
	RDW %	<b>13.6</b>	4.0	9.6 - 17.6	<b>13.8</b>	4.0	9.8 - 17.8	<b>13.0</b>	4.0	9.0 - 17.0
	PLT K/ $\mu$ L	<b>67</b>	20	47 - 87	<b>254</b>	40	214 - 294	<b>510</b>	60	450 - 570
	MPV fL	<b>8.1</b>	2.0	6.1 - 10.1	<b>7.4</b>	2.0	5.4 - 9.4	<b>7.2</b>	2.0	5.2 - 9.2
	Lym %	<b>69.5</b>	8.0	61.5 - 77.5	<b>35.5</b>	6.0	29.5 - 41.5	<b>18.5</b>	5.0	13.5 - 23.5
	Lym K/ $\mu$ L	<b>1.4</b>	0.2	1.2 - 1.6	<b>2.7</b>	0.5	2.2 - 3.2	<b>3.7</b>	1.1	2.6 - 4.8
	Mon %	<b>7.0</b>	7.0	0.0 - 14.0	<b>7.5</b>	5.0	2.5 - 12.5	<b>5.0</b>	5.0	0.0 - 10.0
	Mon K/ $\mu$ L	<b>0.2</b>	0.2	0.0 - 0.4	<b>0.6</b>	0.4	0.2 - 1.0	<b>1.0</b>	1.0	0.0 - 2.0
	Gran %	<b>23.5</b>	7.0	16.5 - 30.5	<b>57.0</b>	6.0	51.0 - 63.0	<b>76.5</b>	7.0	69.5 - 83.5
	Gran K/ $\mu$ L	<b>0.5</b>	0.2	0.3 - 0.7	<b>4.3</b>	0.5	3.8 - 4.8	<b>15.5</b>	1.5	14.0 - 17.0
	CRP mg/dL	<b>0.55</b>	0.30	0.25 - 0.85	<b>2.20</b>	0.70	1.50 - 2.90	<b>4.10</b>	1.00	3.10 - 5.10

For technical assistance, call the Customer Support Center at 1-888-903-5001.