

# MINOTROL CRP

## Hematology Reference Control

9755 Research Drive  
Irvine, California 92618  
Telephone: (949) 453-0500  
(888) 903-5001  
Fax: (949) 753-0533

Lot No.: **MC420**

Exp. Date: **2021-10-05** (Y/M/D)

	Parameter/Units	Low Abnormal Lot # <b>MC420-1</b>			Normal Lot # <b>MC420-2</b>			High Abnormal Lot # <b>MC420-3</b>		
		Assay	Range	Limits	Assay	Range	Limits	Assay	Range	Limits
MICROS CRP 200	WBC K/ $\mu$ L	<b>2.0</b>	0.4	1.6 - 2.4	<b>7.6</b>	1.0	6.6 - 8.6	<b>19.8</b>	2.2	17.6 - 22.0
	RBC M/ $\mu$ L	<b>2.33</b>	0.15	2.18 - 2.48	<b>4.55</b>	0.18	4.37 - 4.73	<b>5.48</b>	0.25	5.23 - 5.73
	HGB g/dL	<b>6.0</b>	0.5	5.5 - 6.5	<b>13.4</b>	0.6	12.8 - 14.0	<b>17.6</b>	0.7	16.9 - 18.3
	HCT %	<b>16.5</b>	2.0	14.5 - 18.5	<b>37.8</b>	2.5	35.3 - 40.3	<b>49.9</b>	3.0	46.9 - 52.9
	MCV fL	<b>71</b>	4	67 - 75	<b>83</b>	4	79 - 87	<b>91</b>	4	87 - 95
	MCH pg	<b>25.8</b>	2.4	23.4 - 28.2	<b>29.5</b>	2.4	27.1 - 31.9	<b>32.1</b>	2.8	29.3 - 34.9
	MCHC g/dL	<b>36.3</b>	3.0	33.3 - 39.3	<b>35.5</b>	3.0	32.5 - 38.5	<b>35.3</b>	3.0	32.3 - 38.3
	RDW %	<b>14.0</b>	4.0	10.0 - 18.0	<b>13.5</b>	4.0	9.5 - 17.5	<b>13.0</b>	4.0	9.0 - 17.0
	PLT K/ $\mu$ L	<b>71</b>	20	51 - 91	<b>256</b>	40	216 - 296	<b>481</b>	60	421 - 541
	MPV fL	<b>9.3</b>	2.0	7.3 - 11.3	<b>8.5</b>	2.0	6.5 - 10.5	<b>8.1</b>	2.0	6.1 - 10.1
	Lym %	<b>69.5</b>	8.0	61.5 - 77.5	<b>38.5</b>	6.0	32.5 - 44.5	<b>17.0</b>	5.0	12.0 - 22.0
	Lym K/ $\mu$ L	<b>1.4</b>	0.2	1.2 - 1.6	<b>2.9</b>	0.5	2.4 - 3.4	<b>3.4</b>	1.1	2.3 - 4.5
	Mon %	<b>7.5</b>	7.0	0.5 - 14.5	<b>6.5</b>	5.0	1.5 - 11.5	<b>6.0</b>	5.0	1.0 - 11.0
	Mon K/ $\mu$ L	<b>0.2</b>	0.2	0.0 - 0.4	<b>0.5</b>	0.4	0.1 - 0.9	<b>1.2</b>	1.0	0.2 - 2.2
	Gran %	<b>23.0</b>	7.0	16.0 - 30.0	<b>55.0</b>	6.0	49.0 - 61.0	<b>77.0</b>	7.0	70.0 - 84.0
	Gran K/ $\mu$ L	<b>0.5</b>	0.2	0.3 - 0.7	<b>4.2</b>	0.5	3.7 - 4.7	<b>15.2</b>	1.5	13.7 - 16.7
	CRP mg/dL	<b>0.62</b>	0.30	0.32 - 0.92	<b>2.38</b>	0.70	1.68 - 3.08	<b>3.95</b>	1.00	2.95 - 4.95

For technical assistance, call the Customer Support Center at 1-888-903-5001.