

# MINOTROL CRP

## Hematology Reference Control

9755 Research Drive  
Irvine, California 92618  
Telephone: (949) 453-0500  
(888) 903-5001  
Fax: (949) 753-0533

Lot No.: **MC421**

Exp. Date: **2022-01-05** (Y/M/D)

	Parameter/Units	Low Abnormal Lot # <b>MC421-1</b>			Normal Lot # <b>MC421-2</b>			High Abnormal Lot # <b>MC421-3</b>		
		Assay	Range	Limits	Assay	Range	Limits	Assay	Range	Limits
MICROS CRP 200	WBC K/ $\mu$ L	<b>2.2</b>	0.4	1.8 - 2.6	<b>7.3</b>	1.0	6.3 - 8.3	<b>19.2</b>	2.2	17.0 - 21.4
	RBC M/ $\mu$ L	<b>2.29</b>	0.15	2.14 - 2.44	<b>4.46</b>	0.18	4.28 - 4.64	<b>5.58</b>	0.25	5.33 - 5.83
	HGB g/dL	<b>6.0</b>	0.5	5.5 - 6.5	<b>13.3</b>	0.6	12.7 - 13.9	<b>17.8</b>	0.7	17.1 - 18.5
	HCT %	<b>16.5</b>	2.0	14.5 - 18.5	<b>37.0</b>	2.5	34.5 - 39.5	<b>50.8</b>	3.0	47.8 - 53.8
	MCV fL	<b>72</b>	4	68 - 76	<b>83</b>	4	79 - 87	<b>91</b>	4	87 - 95
	MCH pg	<b>26.2</b>	2.4	23.8 - 28.6	<b>29.8</b>	2.4	27.4 - 32.2	<b>31.9</b>	2.8	29.1 - 34.7
	MCHC g/dL	<b>36.4</b>	3.0	33.4 - 39.4	<b>35.9</b>	3.0	32.9 - 38.9	<b>35.1</b>	3.0	32.1 - 38.1
	RDW %	<b>13.5</b>	4.0	9.5 - 17.5	<b>14.0</b>	4.0	10.0 - 18.0	<b>13.5</b>	4.0	9.5 - 17.5
	PLT K/ $\mu$ L	<b>74</b>	20	54 - 94	<b>256</b>	40	216 - 296	<b>510</b>	60	450 - 570
	MPV fL	<b>8.8</b>	2.0	6.8 - 10.8	<b>8.0</b>	2.0	6.0 - 10.0	<b>7.8</b>	2.0	5.8 - 9.8
	Lym %	<b>64.5</b>	8.0	56.5 - 72.5	<b>33.5</b>	6.0	27.5 - 39.5	<b>16.0</b>	5.0	11.0 - 21.0
	Lym K/ $\mu$ L	<b>1.4</b>	0.2	1.2 - 1.6	<b>2.4</b>	0.5	1.9 - 2.9	<b>3.1</b>	1.1	2.0 - 4.2
	Mon %	<b>10.0</b>	7.0	3.0 - 17.0	<b>8.5</b>	5.0	3.5 - 13.5	<b>5.0</b>	5.0	0.0 - 10.0
	Mon K/ $\mu$ L	<b>0.2</b>	0.2	0.0 - 0.4	<b>0.6</b>	0.4	0.2 - 1.0	<b>1.0</b>	1.0	0.0 - 2.0
	Gran %	<b>25.5</b>	7.0	18.5 - 32.5	<b>58.0</b>	6.0	52.0 - 64.0	<b>79.0</b>	7.0	72.0 - 86.0
	Gran K/ $\mu$ L	<b>0.6</b>	0.2	0.4 - 0.8	<b>4.2</b>	0.5	3.7 - 4.7	<b>15.2</b>	1.5	13.7 - 16.7
	CRP mg/dL	<b>0.56</b>	0.30	0.26 - 0.86	<b>2.15</b>	0.70	1.45 - 2.85	<b>3.40</b>	1.00	2.40 - 4.40

For technical assistance, call the Customer Support Center at 1-888-903-5001.