

MINOTROL CRP

Hematology Reference Control

9755 Research Drive
 Irvine, California 92618
 Telephone: (949) 453-0500
 (888) 903-5001
 Fax: (949) 753-0533

Lot No.: **MC422**

Exp. Date: **2022-04-05** (Y/M/D)

	Parameter/Units	Low Abnormal Lot # MC422-1			Normal Lot # MC422-2			High Abnormal Lot # MC422-3		
		Assay	Range	Limits	Assay	Range	Limits	Assay	Range	Limits
MICROS CRP 200	WBC K/ μ L	2.2	0.4	1.8 - 2.6	8.6	1.0	7.6 - 9.6	21.0	2.2	18.8 - 23.2
	RBC M/ μ L	2.33	0.15	2.18 - 2.48	4.53	0.18	4.35 - 4.71	5.56	0.25	5.31 - 5.81
	HGB g/dL	6.2	0.5	5.7 - 6.7	13.4	0.6	12.8 - 14.0	17.4	0.7	16.7 - 18.1
	HCT %	16.8	2.0	14.8 - 18.8	37.1	2.5	34.6 - 39.6	49.5	3.0	46.5 - 52.5
	MCV fL	72	4	68 - 76	82	4	78 - 86	89	4	85 - 93
	MCH pg	26.6	2.4	24.2 - 29.0	29.6	2.4	27.2 - 32.0	31.3	2.8	28.5 - 34.1
	MCHC g/dL	37.0	3.0	34.0 - 40.0	36.1	3.0	33.1 - 39.1	35.2	3.0	32.2 - 38.2
	RDW %	13.5	4.0	9.5 - 17.5	13.0	4.0	9.0 - 17.0	12.0	4.0	8.0 - 16.0
	PLT K/ μ L	79	20	59 - 99	268	40	228 - 308	499	60	439 - 559
	MPV fL	8.5	2.0	6.5 - 10.5	8.0	2.0	6.0 - 10.0	8.1	2.0	6.1 - 10.1
	Lym %	66.5	8.0	58.5 - 74.5	33.5	6.0	27.5 - 39.5	17.0	5.0	12.0 - 22.0
	Lym K/ μ L	1.5	0.2	1.3 - 1.7	2.9	0.5	2.4 - 3.4	3.6	1.1	2.5 - 4.7
	Mon %	10.0	7.0	3.0 - 17.0	8.0	5.0	3.0 - 13.0	5.0	4.5	0.5 - 9.5
	Mon K/ μ L	0.2	0.2	0.0 - 0.4	0.7	0.4	0.3 - 1.1	1.1	1.0	0.1 - 2.1
	Gran %	23.5	7.0	16.5 - 30.5	58.5	6.0	52.5 - 64.5	78.0	7.0	71.0 - 85.0
	Gran K/ μ L	0.5	0.2	0.3 - 0.7	5.0	0.5	4.5 - 5.5	16.4	1.5	14.9 - 17.9
	CRP mg/dL	0.51	0.30	0.21 - 0.81	2.23	0.70	1.53 - 2.93	4.23	1.00	3.23 - 5.23

For technical assistance, call the Customer Support Center at 1-888-903-5001.