

MINOTROL CRP

Hematology Reference Control

9755 Research Drive
Irvine, California 92618
Telephone: (949) 453-0500
(888) 903-5001
Fax: (949) 753-0533

Lot No.: **MC424**

Exp. Date: **2022-10-05** (Y/M/D)

	Parameter/Units	Low Abnormal Lot # MC424-1			Normal Lot # MC424-2			High Abnormal Lot # MC424-3		
		Assay	Range	Limits	Assay	Range	Limits	Assay	Range	Limits
MICROS CRP 200	WBC K/ μ L	2.1	0.4	1.7 - 2.5	7.3	1.0	6.3 - 8.3	19.0	2.2	16.8 - 21.2
	RBC M/ μ L	2.26	0.15	2.11 - 2.41	4.47	0.18	4.29 - 4.65	5.59	0.25	5.34 - 5.84
	HGB g/dL	6.0	0.5	5.5 - 6.5	13.6	0.6	13.0 - 14.2	17.9	0.7	17.2 - 18.6
	HCT %	16.3	2.0	14.3 - 18.3	37.5	2.5	35.0 - 40.0	50.9	3.0	47.9 - 53.9
	MCV fL	72	4	68 - 76	84	4	80 - 88	91	4	87 - 95
	MCH pg	26.5	2.4	24.1 - 28.9	30.4	2.4	28.0 - 32.8	32.0	2.8	29.2 - 34.8
	MCHC g/dL	36.9	3.0	33.9 - 39.9	36.2	3.0	33.2 - 39.2	35.2	3.0	32.2 - 38.2
	RDW %	13.5	4.0	9.5 - 17.5	13.0	4.0	9.0 - 17.0	12.5	4.0	8.5 - 16.5
	PLT K/ μ L	75	20	55 - 95	262	40	222 - 302	489	60	429 - 549
	MPV fL	8.5	2.0	6.5 - 10.5	8.0	2.0	6.0 - 10.0	7.7	2.0	5.7 - 9.7
	Lym %	64.0	8.0	56.0 - 72.0	36.0	6.0	30.0 - 42.0	17.0	5.0	12.0 - 22.0
	Lym K/ μ L	1.3	0.2	1.1 - 1.5	2.6	0.5	2.1 - 3.1	3.2	1.1	2.1 - 4.3
	Mon %	10.5	7.0	3.5 - 17.5	8.5	5.0	3.5 - 13.5	6.0	5.0	1.0 - 11.0
	Mon K/ μ L	0.2	0.2	0.0 - 0.4	0.6	0.4	0.2 - 1.0	1.1	1.0	0.1 - 2.1
	Gran %	25.5	7.0	18.5 - 32.5	55.5	6.0	49.5 - 61.5	77.0	7.0	70.0 - 84.0
	Gran K/ μ L	0.5	0.2	0.3 - 0.7	4.1	0.5	3.6 - 4.6	14.6	1.5	13.1 - 16.1
	CRP mg/dL	0.45	0.30	0.15 - 0.75	2.03	0.70	1.33 - 2.73	3.50	1.00	2.50 - 4.50

For technical assistance, call the Customer Support Center at 1-888-903-5001.