

MINOTROL CRP

Hematology Reference Control

9755 Research Drive
 Irvine, California 92618
 Telephone: (949) 453-0500
 (888) 903-5001
 Fax: (949) 753-0533

Lot No.: **MC425**
 Exp. Date: **2023-01-05** (Y/M/D)

	Parameter/Units	Low Abnormal Lot # MC425-1			Normal Lot # MC425-2			High Abnormal Lot # MC425-3		
		Assay	Range	Limits	Assay	Range	Limits	Assay	Range	Limits
MICROS CRP 200	WBC K/ μ L	2.1	0.4	1.7 - 2.5	7.6	1.0	6.6 - 8.6	20.4	2.2	18.2 - 22.6
	RBC M/ μ L	2.24	0.15	2.09 - 2.39	4.57	0.18	4.39 - 4.75	5.52	0.25	5.27 - 5.77
	HGB g/dL	6.1	0.5	5.6 - 6.6	13.5	0.6	12.9 - 14.1	18.0	0.7	17.3 - 18.7
	HCT %	16.4	2.0	14.4 - 18.4	37.0	2.5	34.5 - 39.5	50.8	3.0	47.8 - 53.8
	MCV fL	73	4	69 - 77	81	4	77 - 85	92	4	88 - 96
	MCH pg	27.2	2.4	24.8 - 29.6	29.5	2.4	27.1 - 31.9	32.6	2.8	29.8 - 35.4
	MCHC g/dL	37.3	3.0	34.3 - 40.3	36.5	3.0	33.5 - 39.5	35.4	3.0	32.4 - 38.4
	RDW %	13.5	4.0	9.5 - 17.5	13.5	4.0	9.5 - 17.5	13.0	4.0	9.0 - 17.0
	PLT K/ μ L	82	20	62 - 102	271	40	231 - 311	522	60	462 - 582
	MPV fL	8.4	2.0	6.4 - 10.4	7.8	2.0	5.8 - 9.8	8.0	2.0	6.0 - 10.0
	Lym %	65.5	8.0	57.5 - 73.5	35.5	6.0	29.5 - 41.5	17.5	5.0	12.5 - 22.5
	Lym K/ μ L	1.4	0.2	1.2 - 1.6	2.7	0.5	2.2 - 3.2	3.6	1.1	2.5 - 4.7
	Mon %	9.5	7.0	2.5 - 16.5	7.0	5.0	2.0 - 12.0	7.0	5.0	2.0 - 12.0
	Mon K/ μ L	0.2	0.2	0.0 - 0.4	0.5	0.4	0.1 - 0.9	1.4	1.0	0.4 - 2.4
	Gran %	25.0	7.0	18.0 - 32.0	57.5	6.0	51.5 - 63.5	75.5	7.0	68.5 - 82.5
	Gran K/ μ L	0.5	0.2	0.3 - 0.7	4.4	0.5	3.9 - 4.9	15.4	1.5	13.9 - 16.9
	CRP mg/dL	0.50	0.30	0.20 - 0.80	2.48	0.70	1.78 - 3.18	4.01	1.00	3.01 - 5.01

For technical assistance, call the Customer Support Center at 1-888-903-5001.