

MINOTROL CRP

Hematology Reference Control

9755 Research Drive
Irvine, California 92618
Telephone: (949) 453-0500
(888) 903-5001
Fax: (949) 753-0533

Lot No.: **MC426**

Exp. Date: **2023-04-05** (Y/M/D)

	Parameter/Units	Low Abnormal Lot # MC426-1			Normal Lot # MC426-2			High Abnormal Lot # MC426-3		
		Assay	Range	Limits	Assay	Range	Limits	Assay	Range	Limits
MICROS CRP 200	WBC K/ μ L	2.0	0.4	1.6 - 2.4	7.4	1.0	6.4 - 8.4	19.7	2.2	17.5 - 21.9
	RBC M/ μ L	2.34	0.15	2.19 - 2.49	4.48	0.18	4.30 - 4.66	5.48	0.25	5.23 - 5.73
	HGB g/dL	6.2	0.5	5.7 - 6.7	13.2	0.6	12.6 - 13.8	17.4	0.7	16.7 - 18.1
	HCT %	16.8	2.0	14.8 - 18.8	36.3	2.5	33.8 - 38.8	49.3	3.0	46.3 - 52.3
	MCV fL	72	4	68 - 76	81	4	77 - 85	90	4	86 - 94
	MCH pg	26.5	2.4	24.1 - 28.9	29.5	2.4	27.1 - 31.9	31.8	2.8	29.0 - 34.6
	MCHC g/dL	36.8	3.0	33.8 - 39.8	36.4	3.0	33.4 - 39.4	35.3	3.0	32.3 - 38.3
	RDW %	14.0	4.0	10.0 - 18.0	13.5	4.0	9.5 - 17.5	13.0	4.0	9.0 - 17.0
	PLT K/ μ L	77	20	57 - 97	270	40	230 - 310	520	60	460 - 580
	MPV fL	8.8	2.0	6.8 - 10.8	8.4	2.0	6.4 - 10.4	8.3	2.0	6.3 - 10.3
	Lym %	66.0	8.0	58.0 - 74.0	34.0	6.0	28.0 - 40.0	17.0	5.0	12.0 - 22.0
	Lym K/ μ L	1.3	0.2	1.1 - 1.5	2.5	0.5	2.0 - 3.0	3.3	1.1	2.2 - 4.4
	Mon %	9.0	7.0	2.0 - 16.0	7.5	5.0	2.5 - 12.5	6.0	5.0	1.0 - 11.0
	Mon K/ μ L	0.2	0.2	0.0 - 0.4	0.6	0.4	0.2 - 1.0	1.2	1.0	0.2 - 2.2
	Gran %	25.0	7.0	18.0 - 32.0	58.5	6.0	52.5 - 64.5	77.0	7.0	70.0 - 84.0
	Gran K/ μ L	0.5	0.2	0.3 - 0.7	4.3	0.5	3.8 - 4.8	15.2	1.5	13.7 - 16.7
	CRP mg/dL	0.40	0.30	0.10 - 0.70	1.93	0.70	1.23 - 2.63	4.38	1.00	3.38 - 5.38

For technical assistance, call the Customer Support Center at 1-888-903-5001.