

# ABX Minotrol CRP



**LOT**

**MC 425**

**CONTROL**



**(Exp.) 2023-01-05**  
(YYYY-MM-DD)

Rev 1

PARAMETRES PARAMETERS		UNITES UNITS		ABX Lysebio													
				CONTROL			1	CONTROL			2	CONTROL			3	TOLERANCES TOLERANCE	
				MICROS		PENTRA	TOLERANCES TOLERANCE	MICROS		PENTRA	MICROS		PENTRA				
				EMI CRP	LC-767G	MS CRP		EMI CRP	LC-767G	MS CRP	EMI CRP	LC-767G	MS CRP				
GB WBC	10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l	2.1	2.1	2.2	± 0.4	7.9		7.9	8.2	± 1.0	20.8	21.6	21.8	± 2.2			
GR RBC	10 <sup>6</sup> /mm <sup>3</sup> ;10 <sup>12</sup> /l	2.27	2.26	2.24	± 0.15	4.59	4.59	4.60	± 0.18	5.39	5.43	5.49	± 0.25				
HB HGB	g/dl	6.2	6.2	5.9	± 0.4	13.7	13.6	13.6	± 0.5	18.0	18.1	18.3	± 0.7				
	g/l	62	62	59	± 4	137	136	136	± 5	180	181	183	± 7				
	mmol/l	3.85	3.85	3.66	± 0.25	8.51	8.45	8.45	± 0.31	11.18	11.24	11.36	± 0.44				
HT HCT	%	16.3	16.5	16.4	± 2.0	37.6	37.2	36.8	± 2.5	49.6	49.4	48.9	± 3.0				
	l/l	0.163	0.165	0.164	± 0.020	0.376	0.372	0.368	± 0.025	0.496	0.494	0.489	± 0.030				
VGM MCV	µm <sup>3</sup> ; fl	72	73	73.0	± 4	82	81	80.0	± 4	92	91	89.0	± 4				
TGMH MCH	pg	27.3	27.4	26.3	± 2.4	29.8	29.6	29.6	± 2.4	33.4	33.3	33.3	± 2.8				
	fmol	1.70	1.70	1.64	± 0.15	1.85	1.84	1.84	± 0.15	2.07	2.07	2.07	± 0.18				
CCMH MCHC	g/dl	37.9	37.6	36.1	± 3.0	36.4	36.6	37.0	± 3.0	36.3	36.6	37.5	± 3.0				
	g/l	379	376	361	± 30	364	366	370	± 30	363	366	375	± 30				
	mmol/l	23.56	23.34	22.41	± 1.86	22.60	22.72	22.95	± 1.86	22.54	22.75	23.26	± 1.86				
IDR RDW	%	13.5	13.0	14.5	± 3.0	13.0	12.0	13.0	± 3.0	12.5	12.0	12.5	± 3.0				
PLAQ. PLTS	10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l	76	79	68	± 20	271	271	263	± 40	511	514	506	± 60				
VPM MPV	µm <sup>3</sup> ; fl	8.1	8.9	9.4	± 2.0	8.4	8.2	8.9	± 2.0	8.7	8.6	8.9	± 2.0				
LYMPHO	#	1.3	1.3	N/A	± 0.2	2.7	2.7	N/A	± 0.5	3.9	3.8	N/A	± 1.1				
	%	64.0	63.0	N/A	± 8.0	34.5	34.0	N/A	± 6.0	18.5	17.5	N/A	± 5.0				
MONO	#	0.2	0.2	N/A	± 0.2	0.6	0.6	N/A	± 0.4	1.4	1.4	N/A	± 1.0				
	%	8.5	8.5	N/A	± 7.0	7.0	7.0	N/A	± 5.0	6.5	6.5	N/A	± 5.0				
GRANULO	#	0.6	0.6	N/A	± 0.2	4.6	4.7	N/A	± 0.5	15.7	16.4	N/A	± 1.5				
	%	27.5	28.5	N/A	± 7.0	58.5	59.0	N/A	± 6.0	75.0	76.0	N/A	± 7.0				
CRP	mg/dl	0.54	0.57	0.54	± 0.40	2.77	2.95	2.48	± 0.80	4.52	4.67	4.63	± 1.50				
	mg/l	5.4	5.7	5.4	± 4.0	27.7	29.5	24.8	± 8.0	45.2	46.7	46.3	± 15.0				

Ref: TEMP-1070 Rev.42 BACKVERSO 1300092170