

# MINOTROL CRP

## Hematology Reference Control

9755 Research Drive  
 Irvine, California 92618  
 Telephone: (949) 453-0500  
 (888) 903-5001  
 Fax: (949) 753-0533

Lot No.: **MC427**  
 Exp. Date: **2023-07-05** (Y / M / D)

	Parameter/Units	Low Abnormal Lot # <b>MC427-1</b>			Normal Lot # <b>MC427-2</b>			High Abnormal Lot # <b>MC427-3</b>		
		Assay	Range	Limits	Assay	Range	Limits	Assay	Range	Limits
MICROS CRP 200	WBC K/ $\mu$ L	<b>2.0</b>	0.4	1.6 - 2.4	<b>7.9</b>	1.0	6.9 - 8.9	<b>20.0</b>	2.2	17.8 - 22.2
	RBC M/ $\mu$ L	<b>2.34</b>	0.15	2.19 - 2.49	<b>4.47</b>	0.18	4.29 - 4.65	<b>5.67</b>	0.25	5.42 - 5.92
	HGB g/dL	<b>6.1</b>	0.5	5.6 - 6.6	<b>13.4</b>	0.6	12.8 - 14.0	<b>18.3</b>	0.7	17.6 - 19.0
	HCT %	<b>16.6</b>	2.0	14.6 - 18.6	<b>36.7</b>	2.5	34.2 - 39.2	<b>51.6</b>	3.0	48.6 - 54.6
	MCV fL	<b>71</b>	4	67 - 75	<b>82</b>	4	78 - 86	<b>91</b>	4	87 - 95
	MCH pg	<b>26.1</b>	2.4	23.7 - 28.5	<b>30.0</b>	2.4	27.6 - 32.4	<b>32.3</b>	2.8	29.5 - 35.1
	MCHC g/dL	<b>36.7</b>	3.0	33.7 - 39.7	<b>36.6</b>	3.0	33.6 - 39.6	<b>35.5</b>	3.0	32.5 - 38.5
	RDW %	<b>14.0</b>	4.0	10.0 - 18.0	<b>13.5</b>	4.0	9.5 - 17.5	<b>13.0</b>	4.0	9.0 - 17.0
	PLT K/ $\mu$ L	<b>76</b>	20	56 - 96	<b>256</b>	40	216 - 296	<b>509</b>	60	449 - 569
	MPV fL	<b>8.7</b>	2.0	6.7 - 10.7	<b>8.2</b>	2.0	6.2 - 10.2	<b>7.9</b>	2.0	5.9 - 9.9
	Lym %	<b>68.0</b>	8.0	60.0 - 76.0	<b>33.5</b>	6.0	27.5 - 39.5	<b>17.0</b>	5.0	12.0 - 22.0
	Lym K/ $\mu$ L	<b>1.4</b>	0.2	1.2 - 1.6	<b>2.6</b>	0.5	2.1 - 3.1	<b>3.4</b>	1.1	2.3 - 4.5
	Mon %	<b>9.5</b>	7.0	2.5 - 16.5	<b>8.0</b>	5.0	3.0 - 13.0	<b>5.5</b>	5.0	0.5 - 10.5
	Mon K/ $\mu$ L	<b>0.2</b>	0.2	0.0 - 0.4	<b>0.6</b>	0.4	0.2 - 1.0	<b>1.1</b>	1.0	0.1 - 2.1
	Gran %	<b>22.5</b>	7.0	15.5 - 29.5	<b>58.5</b>	6.0	52.5 - 64.5	<b>77.5</b>	7.0	70.5 - 84.5
	Gran K/ $\mu$ L	<b>0.5</b>	0.2	0.3 - 0.7	<b>4.6</b>	0.5	4.1 - 5.1	<b>15.5</b>	1.5	14.0 - 17.0
	CRP mg/dL	<b>0.30</b>	0.30	0.00 - 0.60	<b>2.08</b>	0.70	1.38 - 2.78	<b>4.08</b>	1.00	3.08 - 5.08

For technical assistance, call the Customer Support Center at 1-888-903-5001.