

ABX Minotrol CRP



LOT

MC 428

CONTROL



(Exp.) 2023-10-05

(YYYY-MM-DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	ABX Lysebio																	
		CONTROL					1	CONTROL					2	CONTROL					3
		MICROS			PENTRA		TOLERANCES TOLERANCE	MICROS			PENTRA		TOLERANCES TOLERANCE	MICROS			PENTRA		TOLERANCES TOLERANCE
		EMI CRP	LC-767G		MS CRP			EMI CRP	LC-767G		MS CRP			EMI CRP	LC-767G		MS CRP		
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	2.0	2.0		2.1		± 0.4	7.5	7.4		7.8		± 1.0	19.9	20.2		20.8		± 2.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	2.32	2.30		2.28		± 0.15	4.66	4.63		4.66		± 0.18	5.50	5.54		5.62		± 0.25
HB HGB	g/dl	6.0	6.0		5.7		± 0.4	13.6	13.5		13.6		± 0.5	17.6	17.7		17.9		± 0.7
	g/l	60	60		57		± 4	136	135		136		± 5	176	177		179		± 7
	mmol/l	3.73	3.73		3.54		± 0.25	8.45	8.38		8.45		± 0.31	10.93	10.99		11.12		± 0.44
HT HCT	%	16.2	16.3		16.0		± 2.0	37.3	37.0		36.8		± 2.5	49.0	48.8		48.8		± 3.0
	l/l	0.162	0.163		0.160		± 0.020	0.373	0.370		0.368		± 0.025	0.490	0.488		0.488		± 0.030
VGM MCV	µm ³ ; fl	70	71		70.0		± 4	80	80		79.0		± 4	89	88		86.0		± 4
TGMH MCH	pg	25.9	26.1		25.0		± 2.4	29.2	29.2		29.2		± 2.4	32.0	31.9		31.9		± 2.8
	fmol	1.61	1.62		1.55		± 0.15	1.81	1.81		1.81		± 0.15	1.99	1.98		1.98		± 0.18
CCMH MCHC	g/dl	36.9	36.7		35.7		± 3.0	36.5	36.4		36.9		± 3.0	36.0	36.3		37.0		± 3.0
	g/l	369	367		357		± 30	365	364		369		± 30	360	363		370		± 30
	mmol/l	22.94	22.82		22.18		± 1.86	22.65	22.63		22.94		± 1.86	22.33	22.55		23.00		± 1.86
IDR RDW	%	14.0	13.0		15.5		± 3.0	12.5	13.0		13.0		± 3.0	12.0	12.0		12.0		± 3.0
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	78	75		68		± 20	262	256		259		± 40	532	538		554		± 60
VPM MPV	µm ³ ; fl	7.6	8.5		8.9		± 2.0	8.1	8.2		8.7		± 2.0	8.0	8.1		8.4		± 2.0
LYMPHO	#	1.3	1.3		N/A		± 0.2	2.7	2.6		N/A		± 0.5	3.5	3.4		N/A		± 1.1
	%	64.0	63.5		N/A		± 8.0	36.0	35.0		N/A		± 6.0	17.5	17.0		N/A		± 5.0
MONO	#	0.2	0.2		N/A		± 0.2	0.6	0.6		N/A		± 0.2	1.3	1.3		N/A		± 1.0
	%	8.0	8.0		N/A		± 7.0	7.5	8.0		N/A		± 5.0	6.5	6.5		N/A		± 5.0
GRANULO	#	0.6	0.6		N/A		± 0.2	4.2	4.2		N/A		± 0.5	15.1	15.5		N/A		± 1.5
	%	28.0	28.5		N/A		± 7.0	56.5	57.0		N/A		± 6.0	76.0	76.5		N/A		± 7.0
CRP	mg/dl	0.48	0.51		0.44		± 0.40	2.49	2.57		2.29		± 0.80	4.62	4.95		4.46		± 1.50
	mg/l	4.8	5.1		4.4		± 4.0	24.9	25.7		22.9		± 8.0	46.2	49.5		44.6		± 15.0

Ref: TEMP-1070 Rev.42 BACKVERSO 1300092170