

MINOTROL CRP

Hematology Reference Control

9755 Research Drive
 Irvine, California 92618
 Telephone: (949) 453-0500
 (888) 903-5001
 Fax: (949) 753-0533

Lot No.: **MC430**

Exp. Date: **2024-04-05** (Y/M/D)

	Parameter/Units	Low Abnormal Lot # MC430-1			Normal Lot # MC430-2			High Abnormal Lot # MC430-3		
		Assay	Range	Limits	Assay	Range	Limits	Assay	Range	Limits
MICROS CRP 200	WBC K/ μ L	2.1	0.4	1.7 - 2.5	7.3	1.0	6.3 - 8.3	19.8	2.2	17.6 - 22.0
	RBC M/ μ L	2.21	0.15	2.06 - 2.36	4.56	0.18	4.38 - 4.74	5.68	0.25	5.43 - 5.93
	HGB g/dL	5.9	0.5	5.4 - 6.4	13.5	0.6	12.9 - 14.1	18.0	0.7	17.3 - 18.7
	HCT %	15.9	2.0	13.9 - 17.9	37.4	2.5	34.9 - 39.9	51.1	3.0	48.1 - 54.1
	MCV fL	72	4	68 - 76	82	4	78 - 86	90	4	86 - 94
	MCH pg	26.7	2.4	24.3 - 29.1	29.6	2.4	27.2 - 32.0	31.7	2.8	28.9 - 34.5
	MCHC g/dL	37.1	3.0	34.1 - 40.1	36.1	3.0	33.1 - 39.1	35.2	3.0	32.2 - 38.2
	RDW %	13.5	4.0	9.5 - 17.5	13.5	4.0	9.5 - 17.5	13.5	4.0	9.5 - 17.5
	PLT K/ μ L	72	20	52 - 92	264	40	224 - 304	515	60	455 - 575
	MPV fL	8.2	2.0	6.2 - 10.2	7.4	2.0	5.4 - 9.4	7.8	2.0	5.8 - 9.8
	Lym %	68.0	8.0	60.0 - 76.0	37.0	6.0	31.0 - 43.0	15.5	5.0	10.5 - 20.5
	Lym K/ μ L	1.4	0.2	1.2 - 1.6	2.7	0.5	2.2 - 3.2	3.1	1.1	2.0 - 4.2
	Mon %	11.0	7.0	4.0 - 18.0	8.0	5.0	3.0 - 13.0	9.0	5.0	4.0 - 14.0
	Mon K/ μ L	0.2	0.2	0.0 - 0.4	0.6	0.4	0.2 - 1.0	1.8	1.0	0.8 - 2.8
	Gran %	21.0	7.0	14.0 - 28.0	55.0	6.0	49.0 - 61.0	75.5	7.0	68.5 - 82.5
	Gran K/ μ L	0.4	0.2	0.2 - 0.6	4.0	0.5	3.5 - 4.5	14.9	1.5	13.4 - 16.4
	CRP mg/dL	0.39	0.30	0.09 - 0.69	2.19	0.70	1.49 - 2.89	4.02	1.00	3.02 - 5.02

For technical assistance, call the Customer Support Center at 1-888-903-5001.