

# ABX Minocal



**LOT**

**CX 486**

Rev **2**

**CAL**



**(Exp.) 2023-10-05**  
(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	ABX Lysebio								TOLERANCES TOLERANCE
		YUMIZEN		YUMIZEN						
		H1500 H2500 V1.0 to V1.4	H1500 H2500 Since V1.7							
GB WBC	10 <sup>9</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l	<b>8.70</b>	<b>8.72</b>							± 0.20
GR RBC	10 <sup>6</sup> /mm <sup>3</sup> ;10 <sup>12</sup> /l	<b>4.51</b>	<b>4.50</b>							± 0.06
HB HGB	g/dl	<b>13.2</b>	<b>13.2</b>							± 0.2
	g/l	<b>132</b>	<b>132</b>							± 2
	mmol/l	<b>8.20</b>	<b>8.20</b>							± 0.12
PLA PLT	10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l	<b>241</b>	<b>241</b>							± 10
VGM MCV	µm <sup>3</sup> ; fl	<b>82.5</b>	<b>83.0</b>							± 4.0

1300116192  
Ref: TEMP-0387 Rev.50