

ABX Minocal



LOT

CX 477

CAL



(Exp.) 2023-01-05
(YYYY-MM-DD)

Rev **1**

PARAMETRES PARAMETERS	UNITES UNITS	ABX Lysebio								TOLERANCES TOLERANCE
		YUMIZEN								
		H1500 H2500								
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.75								± 0.20
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.57								± 0.06
HB HGB	g/dl	14.0								± 0.2
	g/l	140								± 2
	mmol/l	8.69								± 0.12
HT HCT	%	N/A								N/A
	l/l	N/A								N/A
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	264								± 10
VMP MPV	µm ³ ; fl	N/A								N/A
VGM MCV	µm ³ ; fl	86.0								± 4.0

1300060884
Ref: TEMP-0387 Rev.48