

ABX Minocal



LOT

CX 485

CAL



(Exp.) 2023-09-05
(YYYY-MM-DD)

Rev **1**

| PARAMETRES PARAMETERS | UNITES UNITS | ABX Lysebio | | | | | | | | TOLERANCES TOLERANCE |
|--------------------------|---|----------------|--|--|--|--|--|--|--|-------------------------|
| | | YUMIZEN | | | | | | | | |
| | | H1500 H2500 | | | | | | | | |
| GB WBC | 10 ³ /mm ³ ;10 ⁹ /l | 8.75 | | | | | | | | ± 0.20 |
| GR RBC | 10 ⁶ /mm ³ ;10 ¹² /l | 4.55 | | | | | | | | ± 0.06 |
| HB HGB | g/dl | 13.4 | | | | | | | | ± 0.2 |
| | g/l | 134 | | | | | | | | ± 2 |
| | mmol/l | 8.32 | | | | | | | | ± 0.12 |
| HT HCT | % | N/A | | | | | | | | N/A |
| | l/l | N/A | | | | | | | | N/A |
| PLA PLT | 10 ³ /mm ³ ;10 ⁹ /l | 241 | | | | | | | | ± 10 |
| VMP MPV | µm ³ ; fl | N/A | | | | | | | | N/A |
| VGM MCV | µm ³ ; fl | 82.0 | | | | | | | | ± 4.0 |

1300060884
Ref: TEMP-0387 Rev.49