

ABX Minocal



LOT

CX 490

CAL



(Exp.) 2024-02-05
(YYYY - MM - DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	ABX Lysebio								TOLERANCES TOLERANCE
		YUMIZEN		YUMIZEN						
		H1500 H2500 HELO 1.x	H1500 H2500 HELO 2.0							
GB WBC	10 ⁹ /mm ³ ;10 ⁹ /l	8.75	8.75							± 0.20
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.48	4.51							± 0.06
HB HGB	g/dl	13.2	13.2							± 0.2
	g/l	132	132							± 2
	mmol/l	8.20	8.20							± 0.12
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	246	243							± 10
VGM MCV	µm ³ ; fl	83.0	83.7							± 4.0

1300116192
Ref: TEMP-0387 Rev.51