

ABX Minotrol Retic



LOT RX 437
Rev 1

CONTROL

(Exp.) 2022-11-05
(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	ABX											
		CONTROL 1			TOLERANCES TOLERANCE	CONTROL 2			TOLERANCES TOLERANCE	CONTROL 3			TOLERANCES TOLERANCE
		YUMIZEN				YUMIZEN				YUMIZEN			
Méthode Manuelle / Manual RETIC	%	1.3			± 0.7	5.4			± 2.2	10.7			± 3.2
RETIC	%	1.81			± 0.91	5.66			± 1.42	11.38			± 2.28
RETIC	10 ⁶ /mm ³	0.063			± 0.032	0.199			± 0.050	0.403			± 0.081
	10 ⁹ /l	63			± 32	199			± 50	403			± 81
MFI MRI *	%	30.2			± 8.0	38.7			± 10.0	41.6			± 12.0
	-	0.302			± 0.080	0.387			± 0.100	0.416			± 0.120
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	3.49			± 0.30	3.51			± 0.30	3.54			± 0.30
PIC RANGE	-	13			± 8	13			± 8	13			± 8

1300102875
Ref: TEMP-0382 Rev.42

* For research use only in Canada and USA