

ABX Minotrol Retic



LOT RX 438
Rev 1

CONTROL

(Exp.) 2023-01-05
(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	ABX											
		CONTROL 1			TOLERANCES TOLERANCE	CONTROL 2			TOLERANCES TOLERANCE	CONTROL 3			TOLERANCES TOLERANCE
		YUMIZEN				YUMIZEN				YUMIZEN			
Méthode Manuelle / Manual RETIC	%	1.5			± 0.8	5.8			± 2.4	11.2			± 3.4
RETIC	%	1.70			± 0.85	5.58			± 1.40	12.17			± 2.43
RETIC	10 ⁶ /mm ³	0.060			± 0.030	0.198			± 0.050	0.432			± 0.086
	10 ⁹ /l	60			± 30	198			± 50	432			± 86
MFI MRI *	%	29.3			± 8.0	42.2			± 10.0	46.5			± 12.0
	-	0.293			± 0.080	0.422			± 0.100	0.465			± 0.120
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	3.54			± 0.30	3.54			± 0.30	3.55			± 0.30
PIC RANGE	-	13			± 8	13			± 8	13			± 8

* For research use only in Canada and USA

1300102875
Ref: TEMP-0382 Rev.42