

ABX Minotrol Retic



LOT RX 440
Rev 1

CONTROL

(Exp.) 2023-05-05
(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	ABX											
		CONTROL 1			TOLERANCES TOLERANCE	CONTROL 2			TOLERANCES TOLERANCE	CONTROL 3			TOLERANCES TOLERANCE
		YUMIZEN				YUMIZEN				YUMIZEN			
Méthode Manuelle / Manual RETIC	%	1.3			± 0.7	5.5			± 2.2	11.9			± 3.6
RETIC	%	1.60			± 0.80	5.27			± 1.32	10.71			± 2.14
RETIC	10 ⁶ /mm ³	0.056			± 0.028	0.182			± 0.046	0.367			± 0.073
	10 ⁹ /l	56			± 28	182			± 46	367			± 73
MFI MRI *	%	29.8			± 8.0	41.6			± 10.0	44.6			± 12.0
	-	0.298			± 0.080	0.416			± 0.100	0.446			± 0.120
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	3.49			± 0.30	3.46			± 0.30	3.43			± 0.30
PIC RANGE	-	13			± 8	13			± 8	13			± 8

* For research use only in Canada and USA

1300102875
Ref: TEMP-0382 Rev.42