

ABX Minocal



LOT

CX 491

CAL



(Exp.) 2024-03-05
(YYYY - MM - DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	ABX Lysebio								TOLERANCES TOLERANCE
		YUMIZEN		YUMIZEN						
		H1500	H1500	H2500	H2500					
		HELO 1.x	HELO 2.0							
GB WBC	10 ⁹ /mm ³ ;10 ⁹ /l	8.65	8.69							± 0.20
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.54	4.61							± 0.06
HB HGB	g/dl	13.3	13.4							± 0.2
	g/l	133	134							± 2
	mmol/l	8.26	8.32							± 0.12
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	245	241							± 10
VGM MCV	μm ³ ; fl	82.5	82.5							± 4.0

1300116192
Ref: TEMP-0387 Rev.51