

ABX Minocal



LOT

CX 492
Rev 1

CAL

(Exp.) **2024-04-05**
(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	ABX Lysebio									TOLERANCES TOLERANCE
		MICROS	MICROS	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	
		EMI CRP	LC-767G	MS CRP	XLR	DX NEXUS DF NEXUS	MS60		60 60 C+ ES60	XL80	
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.2	9.1	9.3	9.2	9.4	9.1		9.2	9.2	± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.47	4.46	4.48	4.50	4.51	4.48		4.50	4.50	± 0.06
HB HGB	g/dl	13.2	13.2	13.4	13.4	13.3	13.3		13.4	13.4	± 0.2
	g/l	132	132	134	134	133	133		134	134	± 2
	mmol/l	8.20	8.20	8.32	8.32	8.26	8.26		8.32	8.32	± 0.12
HT HCT	%	37.5	37.2	36.5	36.7	36.5	36.3		36.7	36.7	± 1.0
	l/l	0.375	0.372	0.365	0.367	0.365	0.363		0.367	0.367	± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	233	239	222	226	235	227		229	226	± 10
VMP MPV	µm ³ ; fl	8.9	8.9	N/A	N/A	N/A	N/A		N/A	N/A	± 0.5

Ref: TEMP-0387 Rev.51 BACK / VERSO 1300090308