

ABX Minocal



LOT

CX 492

CAL



(Exp.) 2024-04-05
(YYYY - MM - DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	ABX Lysebio								TOLERANCES TOLERANCE
		YUMIZEN		YUMIZEN						
		H1500	H1500	H2500	H2500					
		HELO 1.x	HELO 2.0							
GB WBC	10 ⁹ /mm ³ ;10 ⁹ /l	8.75	8.88							± 0.20
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.52	4.61							± 0.06
HB HGB	g/dl	13.4	13.5							± 0.2
	g/l	134	135							± 2
	mmol/l	8.32	8.38							± 0.12
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	243	240							± 10
VGM MCV	µm ³ ; fl	83.5	84.8							± 4.0

1300116192
Ref: TEMP-0387 Rev.51